

JUN 05 2002

## 2.4 510(k) Summary

### 510(k) Summary

#### U-Systems Ultrasound System

#### U-Systems INC.

Prepared November 21, 2001

Product Name: USI-2000 Ultrasound System with Needle Guide

Manufacturer: U-Systems Inc.

Generic Name Diagnostic Ultrasound System accessory

Classification Name: Ultrasound Imaging System and Transducers (Class II); Classification codes:  
IYO 892.1560 System, Imaging Pulsed Echo, Ultrasonic  
IYN 892.1550 System, Imaging, Pulsed Doppler, Ultrasonic  
ITX 892.1570 Transducer, Ultrasonic, Diagnostic

Contact Person: Sheila W. Pickering Ph.D.  
2081 Longden Circle  
Los Altos, California 94024  
Telephone/Fax 650 969 6114  
e-mail: swpraqa@aol.com

#### A. Legally Marketed Predicate Device

The modification to the USI-2000 is substantially equivalent to the UltraGuide 1000 System.

#### B. Device Description

##### Indications For Use:

The modified USI-2000 includes an accessory for use with the existing diagnostic ultrasound system. The accessory displays graphics depicting the position and future path of a rigid interventional instrument, such as a biopsy needle or an aspiration needle, on a computer monitor screen that also shows the ultrasound image of the target organs.

#### C. Intended Use

The accessory is intended to be used in clinical interventions and for anatomical structures where ultrasound is currently used for visualizing such procedures.

#### D. Substantial Equivalence

The USI-2000 is substantially equivalent to the UltraGuide 1000 System, which is currently in commercial distribution.

CONFIDENTIAL  
U-Systems Inc. 510(k) Notification  
Device Modification

**E. Performance Data**

The USI-2000 Needle Guide performance has been validated according to the company's quality assurance procedures.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

JUN 05 2002

U-Systems, Inc.  
% Sheila W. Pickering, Ph.D.  
2081 Longden Circle  
LOS ALTOS CA 94024

Re: K013902

Trade Name: U-Systems USI-2000 Diagnostic Ultrasound System  
Addition of Needle Guide Accessory

Regulation Number: 21 CFR 892.1550

Regulation Name: Ultrasonic pulsed doppler imaging system

Regulation Number: 21 CFR 892.1560

Regulation Name: Ultrasonic pulsed echo imaging system

Regulation Number: 21 CFR 892.1570

Regulation Name: Diagnostic ultrasonic transducer

Regulatory Class: II

Product Code: 90 IYN, IYO, and ITX

Dated: May 1, 2002

Received: May 6, 2002

Dear Dr. Pickering:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

This determination of substantial equivalence applies to the following transducers intended for use with the U-Systems USI-2000 Diagnostic Ultrasound System, as described in your premarket notification:

Transducer Model Number

7.5 MHz Transducer

10 MHz Transducer

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This determination of substantial equivalence is granted on the condition that prior to shipping the first device, you submit a postclearance special report. This report should contain complete information, including acoustic output measurements based on production line devices, requested in Appendix G, (enclosed) of the Center's September 30, 1997 "Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Transducers." If the special report is incomplete or contains unacceptable values (e.g., acoustic output greater than approved levels), then the 510(k) clearance may not apply to the production units which as a result may be considered adulterated or misbranded.

The special report should reference the manufacturer's 510(k) number. It should be clearly and prominently marked "ADD-TO-FILE" and should be submitted in duplicate to:

Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center (HFZ-401)  
9200 Corporate Boulevard  
Rockville, Maryland 20850

This letter will allow you to begin marketing your device as described in your premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus permits your device to proceed to market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801, please contact the Office of Compliance at (301) 594-4591. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597 or at its Internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Page 3 – Dr. Pickering

If you have any questions regarding the content of this letter, please contact Rodrigo C. Perez at (301) 594-1212.

Sincerely yours,

A handwritten signature in black ink that reads "Nancy C. Brogdon". The signature is written in a cursive style with a large, prominent 'N' and 'B'.

Nancy C. Brogdon  
Director, Division of Reproductive,  
Abdominal and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure(s)

10013902

FDA Submission Cover Sheet

510(k) Number (if known): ~~K003479~~ K013902

Device Name: Modified U-Systems USI-2000 Diagnostic Ultrasound System

Indications For Use:

The modified USI-2000 includes an accessory for use with the existing diagnostic ultrasound system. The accessory displays graphics depicting the position and future path of a rigid interventional instrument, such as a biopsy needle or an aspiration needle, on a computer monitor screen that also shows the ultrasound image of the target organs.

The accessory is intended to be used in clinical interventions and for anatomical structures where ultrasound is currently used for visualizing such procedures.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence Of CDRH, Office Of Device Evaluation (ODE)

Prescription Use

OR

Over-The-Counter Use

(Per 21CFR 801)

Nancy C Brogdon  
(Division Sign-Off)  
Division of Reproductive, Abdominal,  
and Radiological Devices  
510(k) Number K013902

RA  
001

100

**Diagnostic Ultrasound Indications for Use**

510(k) Number(s):

Device Name: USI-2000  
 Diagnostic Ultrasound Pulsed Echo System  
 Diagnostic Ultrasound Pulsed Doppler Imaging System

Intended Use:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)*	Other (specify)
Ophthalmic										
Fetal										
Abdominal		N	N	N		N	N		Note 1	
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (breast, thyroid, testes)		N	N	N		N	N		Note 1	
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Laposcopic										
Peripheral Vascular		N	N	N		N	N		Note 1	
Musculo-skeletal Conventional		N	N	N		N	N		Note 1	
Musculo-skeletal Superficial										

Note 1: Combined includes: B/M; B/PWD; B/Color Doppler; B/Amplitude Doppler; B/Color Doppler/PWD and B/Amplitude Doppler/PWD

The USI Needle Guide accessory is intended for small parts use for breast biopsy.

N = new indication  
 P = previously cleared by FDA

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)  
 Concurrence of CDRH, Office of Device Evaluation (ODE)

  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 510(k) Number           K013902

**Diagnostic Ultrasound Indications for Use**

510(k) Number:

Device Name: 7.5 MHz Transducer  
Diagnostic Ultrasound Transducer

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:


Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)*	Other (specify)
Ophthalmic										
Fetal										
Abdominal	N	N	N			N	N	Note 1		
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (breast, thyroid, testes)	N	N	N			N	N	Note 1		
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Traneseophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Laposcopic										
Peripheral Vascular	N	N	N			N	N	Note 1		
Musculo-skeletal Conventional	N	N	N			N	N	Note 1		
Musculo-skeletal Superficial										

Note 1: Combined includes: B/M; B/PWD; B/Color Doppler; B/Amplitude Doppler; B/Color Doppler/PWD and B/Amplitude Doppler/PWD

The USI Needle Guide accessory is intended for small parts use for breast biopsy.

N = new indication  
P = previously cleared by FDA

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)  
Concurrence of CDRH, Office of Device Evaluation (ODE)

  
 (Division Sign-Off)  
 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 510(k) Number           K013902

**Diagnostic Ultrasound Indications for Use**

510(k) Number:

Device Name: 10 MHz Transducer  
Diagnostic Ultrasound Transducer

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application	Mode of Operation									
	A	B	M	PWD	CWD	Color Doppler	Amplitude Doppler	Color Velocity Imaging	Combined (specify)*	Other (specify)
Ophthalmic										
Fetal										
Abdominal		N	N	N		N	N		Note 1	
Intraoperative (specify)										
Intraoperative Neurological										
Pediatric										
Small Organ (breast, thyroid, testes)		N	N	N		N	N		Note 1	
Neonatal Cephalic										
Adult Cephalic										
Cardiac										
Transesophageal										
Transrectal										
Transvaginal										
Transurethral										
Intravascular										
Laposcopic										
Peripheral Vascular		N	N	N		N	N		Note 1	
Musculo-skeletal Conventional		N	N	N		N	N		Note 1	
Musculo-skeletal Superficial										

Note 1: Combined includes: B/M; B/PWD; B/Color Doppler; B/Amplitude Doppler; B/Color Doppler/PWD and B/Amplitude Doppler/PWD

The USI Needle Guide accessory is intended for small parts use for breast biopsy.

N = new indication  
P = previously cleared by FDA

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)  
Concurrence of CDRH, Office of Device Evaluation (ODE)

  
 Division Sign-Off  
 Division of Reproductive, Abdominal,  
 and Radiological Devices  
 510(k) Number K013902