

K903840

NOV 20 1996



Diagnostics

## 510(k) Summary

---

**Introduction** According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence.

---

**1. Submitter name, address, contact** Boehringer Mannheim Corporation  
2400 Bisso Lane  
P.O. Box 4117  
Concord, CA 94524-4117  
(510) 674 - 0690, extension 8413

Contact Person: Yvette Lloyd

Date Prepared: September 20, 1996

---

**2. Device name** Proprietary name: CEDIA® Phenytoin II Assay  
Common name: Homogeneous enzyme immunoassay for the determination of phenytoin.

Classification name: Enzyme immunoassay, Diphenylhydantoin

---

**3. Predicate device** The Boehringer Mannheim CEDIA® Phenytoin II Assay is substantially equivalent to other products in commercial distribution intended for similar use. Most notably it is substantially equivalent to the currently marketed CEDIA® Phenytoin Assay (K905689).

---



**4.  
Device  
Description**

The CEDIA<sup>®</sup> Phenytoin II Assay is based on the bacterial enzyme  $\beta$ -galactosidase, which has been genetically engineered into two inactive fragments. These fragments spontaneously reassociate to form fully active enzyme that, in the assay format, cleaves a substrate, generating a color change that can be measured spectrophotometrically.

In the assay, phenytoin in the sample competes with analyte conjugated to one inactive fragment of  $\beta$ -galactosidase for antibody binding site. If analyte is present in the sample, it binds to antibody, leaving the inactive enzyme fragments free to form active enzyme. If analyte is not present in the sample, antibody binds to analyte conjugated on the inactive fragment, inhibiting the reassociation of inactive  $\beta$ -galactosidase fragments, and no active enzyme is formed.

---

**5.  
Intended use**

Immunoassay for the in vitro quantitative determination of phenytoin in human serum and plasma.

---

**6.  
Comparison  
to predicate  
device**

The Boehringer Mannheim CEDIA<sup>®</sup> Phenytoin II Assay is substantially equivalent to other products in commercial distribution intended for similar use. Most notably it is substantially equivalent to the currently marketed CEDIA<sup>®</sup> Phenytoin Assay (K905689).

The following table compares the CEDIA<sup>®</sup> Phenytoin II Assay with the predicate device, CEDIA<sup>®</sup> Phenytoin Assay. Specific data on the performance of the test have been incorporated into the draft labeling in attachment 5. Labeling for the predicate device is provided in attachment 6.

**Similarities:**

- Intended Use: Immunoassay for the in vitro quantitative determination of phenytoin
  - Sample type: Serum and plasma
  - Assay range: 0 - 40  $\mu\text{g/mL}$
  - Same Antibody
  - Same Conjugate
- 

*Continued on next page*



510(k) Summary, Continued

6. Comparison to predicate device cont.

Differences:

Feature	CEDIA® Phenytoin II	CEDIA® Phenytoin
Reaction test principle	Spectrophotometric 570 nm	Spectrophotometric 415 nm
Instrument required	Hitachi 911	Hitachi 704
Enzyme Substrate	CPRG (Chlorophenol red-β-D-Galactopyranoside)	m-CNPG (m-Cyano-p-nitrophenol-β-D-galactopyranoside)

Performance Characteristics:

Feature	CEDIA® Phenytoin II			CEDIA® Phenytoin		
Precision	Modified NCCLS (µg/mL):			Modified NCCLS (µg/mL):		
Level	<u>Sample</u>	<u>Control 2</u>	<u>Control 3</u>	<u>Low</u>	<u>Mid</u>	<u>High</u>
N	120	120	120	120	120	119
Within-Run %CV	6.30	14.80	26.75	10.7	17.7	35.4
Total %CV	3.2	2.0	1.3	3.6	2.0	1.5
	6.30	14.80	26.75	10.7	17.7	35.4
	5.1	3.1	2.3	4.7	3.8	3.3

Continued on next page



510(k) Summary, Continued

6. Comparison to predicate device, (cont.)

Performance Characteristics:



Feature	CEDIA® Phenytoin II	TDx Phenytoin
Lower Detection Limit	0.6 µg/mL	0.6 µg/mL
Linearity	0.6 - 40.0 µg/mL	0.6 - 40.0 µg/mL
Method Comparison	Vs CEDIA® Phenytoin <u>Least Squares</u> $y = 1.00x - 0.17$ $r = 0.998$ $SEE = 0.55$ $N = 114$  <u>Deming's</u> $y = 1.00x - 0.19$ $r = 0.998$ $SEE = 0.39$ $N = 114$	Vs Abbott TDx® Phenytoin <u>Least Squares</u> $y = 0.984x - 0.10$ $r = 0.994$ $N = 92$

Continued on next page



510(k) Summary, Continued

6. Comparison to predicate device, (cont.)

Performance Characteristics:

Feature	CEDIA® Phenytoin II	CEDIA Phenytoin
Interfering substances	No interference at:	No interference at:
Bilirubin	66 mg/dL	50 mg/dL
Hemoglobin	1000 mg/dL	1000 mg/dL
Lipemia	1000 mg/dL	1000 mg/dL
Total Protein	12.0 g/dL	13 g/dL
Rheumatoid Factor	100 IU/mL	180 IU/mL
Specificity	% Cross-reactivity	% Cross-reactivity
HPPH	1.8%	1.4%
5-MPPH	5.7%	4.8%