

NOV 20 1997

K971329

10. 510(K) SUMMARY

Summary Of Safety And Effectiveness

10.1. General Information

Classification: Class II
Magnetoencephalograph

Common/Usual Name: Magnetoencephalographic (MEG) Device

Proprietary Name: CTF "Whole-Cortex MEG System"

Establishment Registration: to be submitted

Manufacturer:

CTF Systems Inc.
15 - 1750 McLean Ave.
Port Coquitlam, BC
Canada V3C 1M9
Phone: (604) 941-8561
Fax: (604) 941-8565

Performance Standards: no applicable performance standard have been issued under section 514 of the Food, Drug and Cosmetic Act.

10.2. Intended Uses

The CTF "Whole-Cortex MEG System" system is intended for use as a magnetoencephalographic (MEG) device which non-invasively detects and displays biomagnetic signals produced by electrically active nerve tissue in the brain. When interpreted by a trained clinician, the data enhances the diagnostic capability by providing useful information about the location relative to brain anatomy of active nerve tissue responsible for critical brain functions.

10.3. Device Description

The CTF "Whole-Cortex MEG System" integrates up to 200 dc-SQUID axial gradiometers with workstation computers and data acquisition software in order to measure the magnetic signals generated by the intercellular dendritic currents. These detectors positioned in a helmet shaped array which gives the user the ability to record the electrical activity of the entire surface of the brain simultaneously without having to move the position of the probe.

10.4. Safety and Effectiveness

The CTF "Whole-Cortex MEG System" is substantially equivalent to both the Neuromag-122 (K962764) and the Biomagnetic Technologies Magnes Single (K901215A) in safety and effectiveness. The following chart has been compiled to demonstrate that the CTF "Whole-Cortex MEG System" is substantial equivalence to these devices.

Substantial Equivalence Chart:

Parameter	CTF "Whole-Cortex MEG System"	Neuromag-122 (K962764)	Biomagnetic Technologies Magnes Single (K901215A)
No. of SQUID detectors/channels for MEG data:	64 to 200	122	37
Operating Principle	superconducting flux transformer coupled with dc-SQUID controlled by digital flux-locked loop	superconducting flux transformer coupled with dc-SQUID controlled by analog flux-locked loop	superconducting flux transformer coupled with dc-SQUID controlled by analog flux-locked loop
No. of auxiliary channels for other types of data:	88	166	51
Gradiometer:	1 axial first order gradiometer per location	2 orthogonal planar first order gradiometers per location	1 axial first order gradiometer per location
Intersensor spacing:	32 mm (150 sensor configuration)	43-44 mm	20 mm
Gradiometer placement:	64 to 200 locations distributed across the helmet shaped lower tip of a dewar (Optional Caucasian or Oriental head shape).	61 locations distributed across the helmet shaped lower tip of a dewar.	37 locations positioned in a circular array over a concave spherical surface.
Cryogen used:	Liquid Helium	Liquid Helium	Liquid Helium
Coverage:	One acquisition to cover entire head.	One acquisition to cover entire head.	Six to ten acquisitions to cover entire head.
Gantry:	Floor mounted, standard gantry is fixed. Optional gantry tilts to 90°.	Floor mounted, standard gantry tilts up to 30°. Optional gantry tilts to 45°.	Suspended from ceiling, gantry can tilt up to 45°.
Patient Position:	Seated, or lying on back with optional bed.	Seated or supine. Optional chair insert for children	Seated, or lying on back or side.
Head position indicator:	Included	Available	Available
Computer:	HP workstation with UNIX environment	HP workstation with UNIX environment	SUN workstation with UNIX environment
Networking capabilities:	Ethernet connections to other workstations included	Ethernet connections to other workstations available	Ethernet connections to other workstations available

Magnetic shielded room accessories:	Interior DC lights, video camera and monitor and two-way intercom for monitoring patients	Video monitor and two-way intercom for monitoring patients	Interior DC lights, video cameras and two-way intercom for monitoring patients
-------------------------------------	---	--	--

Intended use comparison:

The CTF "Whole-Cortex MEG System" non-invasively measures the magnetoencephalographic (MEG) signals produced by the active tissue of the brain. These signals are displayed and may be interpreted by trained physicians to help localize these active areas. The locations may be correlated to anatomical information of the brain.

Based on the product literature: The Neuromag-122 system is intended for use as a magnetoencephalographic (MEG) device which non-invasively detects and displays biomagnetic signals produced by electrically active nerve tissues in the brain. When interpreted by a trained clinician, the data enhances the diagnostic capability by providing useful information about the location relative to brain anatomy of active nerve tissue responsible for critical brain functions.

Based on the product literature: The Magnes Single non-invasively detects small biomagnetic signals produced by brain and provides information about the location of electrically active nerve tissue responsible for producing these signals. The data is presented to the physicians in an MEG image, from which they may draw information about the location of critical brain function relative to brain anatomy.



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room - WO66-G609
Silver Spring, MD 20993-0002

Stephen E. Robinson, Ph.D.
Senior Scientist
CTF Systems, Inc.
15 - 1750 McLean Avenue
Port Coquitlam, British Columbia
CANADA V3C 1M9

APR - 9 2012

Re: K971329

Trade/Device Name: Whole-Cortex MEG System
Regulation Number: 21 CFR 882.1400
Regulation Name: Electroencephalograph
Regulatory Class: II
Product Code: OLY, GWQ
Dated (Date on orig SE ltr): September 5, 1997
Received (Date on orig SE ltr): September 8, 1997

Dear Mr. Robinson:

This letter corrects our substantially equivalent letter of November 20, 1997.

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Malvina B. Eydelman". The signature is written in a cursive style and is positioned to the left of the typed name.

Malvina B. Eydelman, M.D.

Director

Division of Ophthalmic, Neurological,
and Ear, Nose and Throat Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

510(k) Number (if known): K971329

Device Name: Whole-cortex MEG System (with optional EEG)

Indications For Use:

The CTF Systems Inc. "Whole-Cortex MEG System" non-invasively measures the magnetoencephalographic (MEG) signals (and, optionally, electroencephalographic EEG signals) produced by electrically active tissue of the brain. These signals are recorded by a computerized data acquisition system, displayed, and may then be interpreted by trained physicians to help localize these active areas. The locations may then be correlated with anatomical information of the brain. MEG is routinely used to identify the locations of visual, auditory, somatosensory, and motor cortex in the brain.* MEG is also used to non-invasively locate regions of epileptic activity within the brain. The localization information provided by MEG may be used, in conjunction with other diagnostic data, in neurosurgical planning.

*when routinely used in conjunction with evoked response averaging devices

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-Off)
Division of Cardiovascular, Respiratory,
and Neurological Devices

510(k) Number K971329

Prescription Use X
(Per 21 CFR 801.109)

OR

Over-The-Counter Use _____

(Optional Format 1-2-96)