

510(k) SUMMARY
as required per 807.92(c)1. Submitters Name, Address:

Siemens Medical Systems, Inc.
Electromedical Systems Group, PCS
Danvers, MA 01923
Tel: (508) 750-7500
Fax: (508) 777-3398

Official Correspondent:
David Simard, Director
Quality Assurance & Regulatory Affairs

Contact person for this submission:
Wulf R. Trepte
Quality Assurance
Date submission was prepared: 16/10/97

2. Trade Name, Common Name and Classification Name:A. Trade Name:

KION Anesthesia System

B. Common Name, Classification Name, Class and Regulation Number:

Common Name	Classification Number	Class	Regulation Number
Gas-machine, Anesthesia	73BSZ	II	21 CFR 868.5160
Gas-machine, Analgesia	73ELI	II	21 CFR 868.5160
Cardiac monitor	74DRT	II	21 CFR 870.2300
Pulse rate monitor	74BWS	II	21 CFR 870.2300
Pulse oximeter	74DQA	II	21 CFR 870.2700
Breathing Frequency Monitor	73BZQ	II	21 CFR 868.2375
Clinical Electronic Thermometer	80BWX	II	21 CFR 880.2910
Indwelling Blood Pressure Monitor	74CAA	II	21 CFR 870.1110
Heart Rate Monitor, Neonatal	74FLO	II	21 CFR 870.2300
Ventilatory Effort Monitor (Apnea Detector)	73FLS	II	21 CFR 868.2375
Monitor Blood Pressure, Neonatal, Invasive	74FLP	II	21 CFR 870.1110
Arrhythmia detector & Alarm	74DSI	III	21 CFR 870.1025
Medical Cathode-Ray Tube Display	74DXJ	II	21 CFR 870.2450
ST Segment Monitor with Alarm	74MLD	III	21 CFR 870.1025
Non-indwelling Blood Pressure Monitor	74DXN	II	21 CFR 870.1130
End-tidal Carbon-Dioxide Monitor	73CCK	II	21 CFR 868.1400

510(k) Notification
Siemens KION

3. Predicate Device Identification:

<i>Legally marketed devices to which equivalence is being claimed</i>	<i>510(k) #</i>
Siemens Servo Ventilator 900C	K841529
Servo Ventilator 900C/D Electronic gas supply unit	K871354
Servo Ventilator 900C/D CAV 1.0 Master Program Flow & Prs, CO2 Lung Monitor	K861707
Servo Ventilator 900C/D P.C. Interface	K852771
Servo Ventilator 900C/D CAV 4.0 Software Servo Graphics	K923444/A
Servo Ventilator 900C/D Man. Vent. Acc.	K854036
Servo Ventilator 900C/D/E CAV 3.0 Software	K890932
Siemens Servo Anesthesia Circle 985	K893786
Siemens Vaporizers 950/951/952	K841157
Ohmeda Tec5 Vaporizer	K942091
Siemens Servo Ventilator 300	K960010
Siemens SC9000/SC9015 Bedside Monitoring System	K946306

4. Device Description:

KION is a modular anesthesia system. It is available as a set of hardware and software options, which can be added to take the system from being a low-cost, general purpose anesthesia machine, to a complete state-of-the-art anesthesia workstation. This gives the opportunity to specify a system and then upgrade and modify the system as needs change in the future. The KION anesthesia system is designed to enable the operator to work with the basic components of an anesthesia system in the most suitable way for each patient case. As such, the system is mounted on wheels to allow it to be easily moved into an optimum position, where many of the components can also be individually rotated. The KION system is intended for use on adult and pediatric patients

Physically, the KION consists of three main sections:

1. The Trolley which contains the main pneumatic inlets, power supply and power inlets, battery, communication ports, and the gas delivery system.
2. The Column consists of a revolving Vaporizer Holder, Servo Bellows Absorber and Patient Cassette.
3. The Integrated User Interface consists of the Control Panel and the Monitor.

The KION offers the following standard ventilation capabilities:

- Volume controlled ventilation
- Manual ventilation

5. Intended Use:

The intended use for the KION is as a modular anesthesia system for use on Adult and Infant patient populations. The unit is designed to be used in the operating room or any other hospital location where administration of anesthesia is necessary.

6. Table of Device Similarities and differences to predicate device

Manufacturer	Siemens-Elema AB, Electromedical Systems Division	Siemens-Elema AB, Electromedical Systems Division	-
510(k) number	Predicate devices <ul style="list-style-type: none"> • Siemens Servo ventilator 900C - K841529 • Siemens Servo Ventilator 300 - K960010 • Servo Anesthesia Circle 985 -K893786. • Siemens Vaporisers 950, 951, 952 - K841157. • Ohmeda Tec 5 vaporiser - K942091. • Siemens Servo Screen - K960168 	Modified Device Kion Anesthesia System K number to be decided	-

<i>Parameters General</i>	<i>Predicate devices: All</i>	<i>Modified Device</i> Kion Anesthesia System	<i>Explanation of the differences compared to the Predicate devices</i>
Applications and Intended Use	Infants to Adults	Same	-
Gas Delivery System	Pneumatic Powered Servo Controlled "Bag in Bottle" system.	Same	-
Battery Back up	No	Yes	Battery backup for the 900C is only available with an external system.

<i>Parameters Gas Dosage</i>	<i>Predicate devices</i> Siemens Servo ventilator 900C - K841529 Servo Anesthesia Circle 985 - K893786.	<i>Modified Device</i> Kion Anesthesia System	<i>Explanation of the differences compared to the Predicate devices</i>
Fresh gas dosage Range	Range:500 ml/min - 40 l/min	Range:100 ml/min - 18lpm	Enhanced functionality on Kion for low flow applications
Accuracy	+ /- 10%	Same	-
Oxygen Percent	21% / 25% - 100% (Air / N ₂ O)	21% / 28% - 100% (Air / N ₂ O)	Enhanced anti-hypoxic measures
Gas Mixture Setting	O ₂ /Air or O ₂ /N ₂ O	Same	-
Anti-hypoxia device (Limits O₂/N₂O fresh gas ratio)	Yes	Same	-
Ventilator Integral	Yes	Same	-

<i>Anesthetic Vapour Delivery Modules</i>	<i>Predicate devices</i> Siemens Vaporizers 950/951/952 and Ohmeda Tec5	<i>Modified Device</i> Kion Anesthesia System Vaporizers	<i>Explanation of the differences compared to the Predicate devices</i>
Operating Principle	High pressure precision injector	Same	-
Vaporizer mounting	Single vaporizer mounting	Magazine for three vaporizers mounting	Enhanced functionality.

510(k) Notification
Siemens KION

Breathing Systems	<u>Predicate devices</u> SV 900 and Servo Anesthetic Circle 985	<u>Modified Device</u> Kion Anesthesia System	<i>Explanation of the differences compared to the Predicate devices</i>
Circle System	Yes	Same	-
Auxiliary Fresh Gas Outlet	No	Yes	Enhanced functionality to support external breathing systems (e.g. Bain, T-piece)
Adjustable Pressure Regulating (APL) valve	No	Adjustable to 90 cm H ₂ O	Improved functionality

Ventilation Modes	<u>Predicate devices</u> SV 900C	<u>Modified Device</u> Kion Anesthesia System	<i>Explanation of the differences compared to the Predicate devices</i>
Manual	Yes	Same	-
Volume Controlled Ventilation	Yes	Same	-
Pressure Controlled Ventilation	Yes	No	Planned for future enhancement in Kion
Pressure Support	Yes	No	Planned for future enhancement in Kion

Ventilatory Control Settings	<u>Predicate devices</u> SV 300	<u>Modified Device</u> Kion Anesthesia System	<i>Explanation of the differences compared to the Predicate devices</i>
Preset Minute Volume	0.5 - 45LPM	Same	-
Tidal Volume	2 - 1500 ml	20 - 1500 ml	Altered functionality, no neonatal capability on Kion.
CMV Rate	0.5 - 120	6 - 99	Enhanced functionality due to specific clinical anesthesia requirements
I:E Ratio	1:4 - 4:1 in seven steps.	1:3 - 4:1 in five steps.	Enhanced functionality due to specific clinical anesthesia requirements
Pause Time	0, 5, 10, 20, 30% of cycle	10% of cycle in Volume Control Mode (except for I:E - 4:1 = 0%)	Enhanced functionality due to specific clinical anesthesia requirements
Inspiratory Gas Flow	Maximum 180 lpm	Same	-

Monitored Ventilatory & Respiratory parameters	<u>Predicate Devices</u> SV 900 and Servo Screen 390	<u>Modified Device</u> Kion Anesthesia System	<i>Explanation of the differences compared to the Predicate devices</i>
Inspired Tidal Volume	Yes	Same	-
Expired Tidal Volume	Yes	Same	-
Rate	Yes	Same	-
Peak airway Pressure	Yes	Same	-
Mean Airway Pressure	Yes	Same	-
Pause Pressure	Yes	Same	-
PEEP/CPAP	Yes	Same	-
Airway pressure / Analogue Display	Yes	Same	-

Ventilatory Alarms	<u>Predicate Device</u> SV900	<u>Modified Device</u> Kion Anesthesia System	<i>Explanation of the differences compared to the Predicate devices</i>
	This device has the following alarms: High Airway Pressure, High Continuous Airway Pressure, Disconnect, Low Minute Volume, Low O ₂ Concentration	Same	Enhanced functionality. In addition Kion also have the following alarms: Low Air Supply Wall and Tank, Low O ₂ Supply Wall and Tank, Low N ₂ O Supply Wall and Tank, High Air Supply Wall and Tank, High O ₂ Supply Wall and Tank, High N ₂ O Supply Wall and Tank, Low Battery Capacity

7. Assessment of non-clinical performance data for equivalence:

See the testresults under section In-house Testing, Section T, in this submission.

8. Assessment of clinical performance data for equivalence:

Clinical trials ongoing.

9. Biocompatibility:

Not applicable

10. Sterilization:

Not applicable

11. Standards and Guidances:

KION complies to the following standards:

- IEC 601-1 (EN 60601-1) Medical Electrical Equipment, part 1: General requirements for safety.
- IEC 601-1-1 (EN 60601-1-1) Medical Electrical Equipment, part 1.1: Safety Requirements for medical electrical systems.
- IEC 601-1-2 (EN 60601-1-2) Medical Electric Equipment, part 1.2: Electromagnetic Compatibility- Requirements and Tests
- ISO/DIS 8835-1.2 (Forthcoming IEC 60601-2-13, 2nd Ed.) Anesthesia Workstations and Their Components.
- prEN 794-1 Lung ventilators-Part 1: Particular Requirements for Critical Care Ventilators.
- ISO 7767 Oxygen monitors for monitoring patient breathing mixtures- Safety Requirements.



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

SEP 13 1999

Mr. David Simard
Simens Medical Systems, Inc.
16 Electronics Avenue
Danvers, MA 01923

Re: K973971
KION Anesthesia System
Regulatory Class: III (three)
Product Code: 73 BSZ and 74 DSI
Dated: July 20, 1999
Received: July 21, 1999

Dear Mr. Simard:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Pre-market Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Mr. David Simard

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4648. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597, or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,



Thomas J. Callahan, Ph.D.
Director
Division of Cardiovascular,
Respiratory, and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if known): K973971

Device Name: KION Anesthesia System

Indications for Use:

Use of the KION Anesthesia System is indicated for adult or infant populations in an environment where patient care is provided by Healthcare Professionals, trained in the administration of anesthesia, when the professional determines that a device is required to assist the breathing of the patient. This device can be used to administer anesthesia while controlling the entire ventilation for patients without any ability to breathe, as well as for supporting patients with reduced ability to breathe.

MRI Compatibility Statement:

Siemens KION is not compatible for use in a MRI magnetic field.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

_____ Concurrency of CDRH, Office of Device Evaluation (ODE) _____

Prescription Use
(Per 21 CFR 801.109)

J. A. Weir

OR

Over-The-Counter Use _____

(Optional Format 1-2-96)

(Division Sign-Off)
Division of Cardiovascular, Respiratory,
and Neurological Devices
510(k) Number K973971