

**Elecsys® PSA on Elecsys 1010**

K974189

- 
- Introduction** According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence.
- 
- 1. Submitter name, address, contact** Boehringer Mannheim Corporation  
4300 Hacienda Drive  
P.O. Box 9002  
Pleasanton, CA 945660900  
(510) 730-8215
- Contact Person: Patricia M. Klimley
- Date Prepared: November 6, 1997
- 
- 2. Device name** Proprietary name: Elecsys® PSA Assay
- Common name: Electrochemiluminescence assay for the determination of Prostate-Specific Antigen (PSA).
- Classification name: System, Test , Prostate-Specific antigen
- 
- 3. Predicate device** We claim substantial equivalence to the Elecsys® PSA Assay on Elecsys 2010.
- 
- 4. Device Description** The Elecsys® test principle is based on sandwich principle. Total duration of assay: 18 minutes (37° C).
- 1st incubation (9 minutes): Sample (40 µL), a biotinylated monoclonal PSA-specific antibody (60 µL), and a monoclonal PSA-specific antibody labeled with a ruthenium complex (60 µL) react to form a sandwich complex.
  - 2nd incubation (9 minutes): After addition of streptavidin-coated microparticles (40 µL), the entire complex is bound to the solid phase via interaction of biotin and streptavidin.

*Continued on next page*

## 510(k) Summary, Continued

---

**4.  
Device  
Description**

•The reaction mixture is aspirated into the measuring cell where the microparticles are magnetically captured onto the surface of the electrode. Unbound substances are then removed with ProCell. Application of a voltage to the electrode then induces chemiluminescent emission which is measured by a photomultiplier (0.4 second read frame).  
•Results are determined via a calibration curve which is instrument-specifically generated by 2-point calibration and a master curve provided via the reagent bar code.

---

**5.  
Intended use**

Immunoassay for the in vitro quantitative determination of Prostate-Specific Antigen in human serum and plasma to aid in the management of prostate cancer patients.

---

**6.  
Comparison  
to predicate  
device**

The Boehringer Mannheim Elecsys® PSA Assay has been approved for use on the Elecsys 2010 immunoassay analyzer (K964351/S1). The application of the Elecsys® PSA Assay on the Elecsys 1010 immunoassay analyzer is substantially equivalent to the same assay (Elecsys PSA Assay) on the Elecsys 2010.

The following table compares the Elecsys® PSA on the Elecsys 1010 with the predicate device, Elecsys® PSA on the Elecsys 2010. Specific data on the performance of this test for both the Elecsys 1010 and 2010 have been incorporated into the draft labeling in attachment 5. Labeling for the predicate device provided in attachment 6 will be replaced upon the approval of this premarket approval submission with the combined Elecsys 2010 and 1010 insert (attachment 5).

**Similarities:**

- Intended Use: Immunoassay for the in vitro quantitative determination of prostate-specific antigen. The assay is further indicated for serial measurement of PSA to aid in the management of cancer patients.
- Assay range: 0-100 ng/ml
- Assay methodology: Sandwich immunoassay
- Cross-Reactivity: 0% to PAP
- Kit (cat. no. 1731262) approved for use on the Elecsys 2010 (K964351/S1)
- Sample and reagent volumes
- Reaction temperature and incubation times
- Package insert
- Performance specifications

---

*Continued on next page*

## 510(k) Summary, Continued

6.  
Comparison  
to predicate  
device cont.

**Differences:**

Feature	Elecsys 1010	Elecsys 2010
Instrument required	Elecsys 1010	Elecsys 2010
Instrument System	Batch	Random access
Reagent Storage Temp (C)	20-25C	37C

**Performance Characteristics:**

Feature	Elecsys 1010	Elecsys 2010
Lower Detection Limit	0.006 ng/mL Functional: 0.07 ng/mL	0.002 ng/mL Functional: 0.03 ng/mL
Linearity	0.01 - 100 ng/mL (with a deviation from a linear line of $\pm 10\%$ )	0.01 - 100 ng/mL (with a deviation from a linear line of $\pm 10\%$ )
Method Comparison	Vs Elecsys 2010  <u>Least Squares</u> $y=0.9978x + 0.0040$ $r=0.9915$ $N=91$  <u>Passing Bablock</u> $y=0.9825x + 0.0935$ $r=0.9915$ $N=91$	Not shown
Hook Effect	No Hook Effect up to 15,000 ng/ml PSA	No Hook Effect up to 15,000 ng/ml PSA

*Continued on next page*



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

JAN - 6 1998

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

Ms. Patricia M. Klimley  
Manager, Regulatory Affairs  
Boehringer Mannheim Corporation  
Laboratory Diagnostics  
4300 Hacienda Drive  
Pleasanton, California 94566-0900

Re: K974189  
Trade Name: Elecsys® PSA Assay  
Regulatory Class: II  
Product Code: LTJ  
Dated: November 6, 1997  
Received: November 7, 1997

Dear Ms. Klimley:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal Laws or Regulations.


Page 2

Under the Clinical Laboratory Improvement Amendments of 1988 (CLIA-88), this device may require a CLIA complexity categorization. To determine if it does, you should contact the Centers for Disease Control and Prevention (CDC) at (770)488-7655.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll free number (800) 638-2041 or at (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>"

Sincerely yours,



Steven I. Gutman, M.D., M.B.A.  
Director  
Division of Clinical  
Laboratory Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

