


1 **Product description:**

2 The Siemens SONOCUR Basic system (see figure below) is a mobile ESWT
3 unit, which consists of the following main system components:



4

5

6 (1) Shock head (or pressure wave generator) with coupling bellows;

7 (2) transportable trolley with charging unit and water conditioning system;

8 (3) control console; and

9 (4) articulating arm to position the shock head.

10

11 The SONOCUR Basic system features a finely adjustable electromagnetic shock
12 wave generator. The shock head, or pressure wave generator, is suspended by
13 an articulating arm which allows flexible movement of the tube in three different
14 planes. The flexible movement enables precise positioning of the shock tube for
15 the pain therapy application. The focus of the therapeutic shock wave can be
16 adjusted for variable penetration depths (0-50mm).

17 From the control console, the following parameters can be selected and/or dis-
18 played:

- 19
- number of impulses per treatment
 - number of impulses currently applied during session
 - energy level indication (the shock wave energy can be selected between
20 densities ranging from 0.04 mJ/mm² to 0.5 mJ/mm² in eight energy levels)
 - warning and error messages

21

22

23

24

Intended Use – Indications:

The Siemens SONOCUR Basic is a non-surgical alternative for the treatment of chronic lateral epicondylitis (commonly referred to as tennis elbow) for patients with symptoms of chronic lateral epicondylitis for 6 months or more and a history of unsuccessful conservative treatments.

Contraindications:

There are no known contraindications to ESWT with the SONOCUR Basic for treatment of chronic lateral epicondylitis.

Warnings:

The following warnings pertain to the use of the SONOCUR Basic for treatment of chronic lateral epicondylitis (tennis elbow).

- Operators of the SONOCUR Basic should be aware of the proper use of the device in delivering the correct number of shocks and in localizing the proper area to be treated.
- ESWT with the SONOCUR Basic should be prescribed by and performed under the supervision of a physician trained and experienced in the care of patients with lateral epicondylitis.
- If the patient moves after correct positioning, re-perform localization if necessary. Failure to maintain correct positioning could result in misdirection of the shockwave and injury to adjacent nerves or blood vessels.
- If patients experience severe pain/discomfort at the application site during treatment, the system operator should decrease the penetration depth of the therapeutic shock wave focus (0-50mm) by increasing the water level in the coupling bellows.
- If patients experience a vaso-vagal reaction during treatment, the patient should be reclined to a supine position until symptoms disappear.
- Patients currently undergoing systemic anticoagulation therapy (example— coumadin, heparin) should consult their physicians regarding temporary discontinuation of such medications before ESWT to prevent potential ecchymosis/bruising.
- Patients on daily aspirin therapy should temporarily discontinue aspirin intake 1 week before ESWT therapy.

67 Precautions:

68

69 The following precautions pertain to the use of the SONOCUR Basic for
70 treatment of chronic lateral epicondylitis (tennis elbow).

71

- 72 • Electromagnetic compatibility (EMC):

73 If electromagnetic interference between the extracorporeal shock wave
74 system and nearby electronic equipment is suspected (as evidenced by
75 erratic behavior with either device), it is recommended that their distance be
76 increased until proper operation resumes. If it is necessary to operate an
77 electronic device in close proximity to the ESWT system during treatment, the
78 device and the ESWT system should be tested for proper simultaneous
79 operation prior to clinical use.

- 80 • Never remove any of the cabinet covers to the system's electronics. The high
81 voltage power supply circuits utilized by extracorporeal shock wave systems
82 use voltages that are capable of causing serious injury or death from electric
83 shock.

- 84 • If the device malfunctions during treatment or the treatment is discontinued,
85 the therapeutic effects may not be as noticeable.

86

87 The safety and effectiveness of the SONOCUR Basic has not been established
88 for:

- 89 • Pregnant women;
- 90 • Patients younger than 18 years of age;
- 91 • Patients with a coagulation abnormality, thrombopathy, infection, tumor,
92 cervical compression syndrome, cervical or upper extremity arthritis, local
93 arthrosis, neurologic abnormality, or radial nerve entrapment;
- 94 • Patients who have had previous surgery for lateral epicondylitis;
- 95 • Patients who suffer from severe systemic diseases that may lead to sensory
96 changes or neuropathic pain. For example, this may include diseases such
97 as gout, diabetes mellitus, rheumatoid arthritis;
- 98 • Patients with cardiac pacemaker;
- 99 • Patients who received physical or occupational therapy less than four (4)
100 weeks prior to ESWT;
- 101 • Patients who received a local steroid injection less than six (6) weeks prior to
102 ESWT; and
- 103 • Patients with tennis elbow affecting both arms or who have had previous
104 surgery for this condition.

105

106 **Adverse Events:**

107

108 Adverse events observed during a clinical study of 114 patients that were
 109 associated with extracorporeal shock wave therapy (ESWT) include those listed
 110 below, categorized by frequency:

111

112 Adverse events reported in >20% of patients:

113 ♦ pain at, or surrounding the treatment site

114

115 Adverse events reported in <20% of patients

- ♦ nausea
- ♦ myalgia
- ♦ joint disorder
- ♦ pallor
- ♦ dizziness
- ♦ hypertonia
- ♦ hypesthesia
- ♦ paresthesia
- ♦ tremor
- ♦ vasodilation
- ♦ application site reaction
- ♦ sweating

116

117

118 The number and frequency of each reported event is summarized in Table 1 below:

119

120 **Table 1: Device Related Adverse Events at 12 Week Follow-up**

	Active			Placebo		
	Number of Patients [1]	Number of Occurrences	% of Patients [2]	Number of Patients [1]	Number of Occurrences	% of Patients [2]
Pain	28	60	50%	13	32	22.4
Nausea	10	10	17.9%	0	0	0%
Application Site Reaction	6	8	10.7%	5	5	8.6%
Sweating	5	5	8.9%	0	0	0%
Dizziness	4	4	7.1%	0	0	0%
Hypertonia	3	5	5.4%	3	3	5.2%
Hypesthesia	3	5	5.4%	1	2	1.7%
Paresthesia	3	4	5.4%	8	12	13.8%
Joint Stiffness	2	2	3.6%	0	0	0%
Myalgia	2	2	3.6%	0	0	0%
Tremor	2	2	3.6%	0	0	0%
Vasodilation	2	2	3.6%	0	0	0%
Pallor	1	1	1.8%	0	0	0%
Accidental Injury	0	0	0%	2	3	3.4%
Headache	0	0	0%	2	7	3.4%
Peripheral Edema	0	0	0%	1	1	1.7%
Twitching	0	0	0%	1	1	1.7%
Sinusitis	0	0	0%	1	2	1.7%

121 1. Number of patients experiencing at least one adverse event

122 2. Based on the total number of patients in each treatment group: active treatment group=56 patients, placebo group=58 patients

123 During the study, three patients exhibited EKG changes that were determined by
124 the investigators and cardiologists not to be treatment related.

125

126 **Potential Adverse Events:**

127

128 Potential adverse events not seen during the study include:

129

- Neuropathy
- Tendon rupture
- Local hematoma
- Misdirection of energy

130

131

132

133

134

135 ***Clinical Study:***

136

137 **Study Design:**

138 The Siemens Sonocur Basic multi-center pivotal trial was a randomized, double-
139 blind (patients and evaluators), placebo-controlled, parallel treatment study. A
140 total of 114 patients were enrolled in the study at 3 investigational sites.

141

142 Patients with chronic tennis elbow were examined and randomized to one of two
143 treatment groups (active, placebo). Each patient was scheduled to receive three
144 treatments: once a week for a three-week period. For all completed treatments,
145 a maximum of 2100 impulses per treatment session was delivered for a total
146 energy delivery of 9.27J for all three sessions. The procedure for the active and
147 placebo treatments was performed identically except that for patients receiving
148 the placebo treatment, a sound-reflecting pad was placed between the treatment
149 site and the shock wave head. No local anesthetic injection or analgesic was
150 allowed during treatment.

151

152 During the study, assessments of pain level and functional activity were
153 performed. At each visit, the pain intensity was evaluated using the Thomsen
154 provocation test (resisted wrist extension). The patient was asked to record the
155 level of pain that he/she was experiencing on a visual analog scale (VAS), which
156 was a 100mm scale with 0 for no pain and 100 for intolerable pain. In addition,
157 functional improvement was also examined using an Upper Extremity Functional
158 Scale, or UEFS, test. For this test, patients in the study were asked to score
159 their ability to perform specific daily chores (such as opening jars/doors, washing
160 dishes) on a scale from 1 to 10, with a score of 1 meaning that the patient had no
161 problem at all and a score of 10 meaning that the patient could not perform the
162 activity. Additional measures of efficacy were examined including the patient's
163 overall impression, grip strength, activity evaluation (ability to perform activities
164 that were limited by his/her tennis elbow condition), and pain medication
165 consumption. Safety assessments included an assessment of adverse events,
166 physical examination, X-rays, vital signs, 12-lead EKG, clinical labs, and
167 proportion of patients who couldn't tolerate treatment.

168

169 Patients were scheduled for follow-up evaluations, occurring at 1-, 4-, 8-, 12-
170 weeks, 6-months, and 12-months post-treatment. The primary analysis of the
171 safety and efficacy data was performed after all patients were enrolled, treated,
172 and completed their 12-week follow-up requirements.

173

174 **Inclusion/ Exclusion Criteria:**

175

176 The *inclusion criteria* included:

- 177 • history of lateral epicondylitis for at least 6-months;
- 178 • pain that is unresponsive to two of three conventional therapy programs
179 (local steroid injections, physical/occupational therapy (PT/OT), non-
180 steroidal anti-inflammatories (NSAIDs));
- 181 • pain by palpation of the lateral epicondyle;
- 182 • baseline pain that was ≥ 40 during resisted wrist extension (“Thomsen
183 provocation test”) on a 100mm visual analog scale (VAS); and
- 184 • signed informed consent.

185

186 The *exclusion criteria* included:

- 187 • < 18 years of age;
- 188 • received local injections within 6 weeks, physical/ occupational therapies
189 within 4 weeks, or non-steroidal anti-inflammatories within 1 week prior to
190 randomization;
- 191 • received systemic therapeutic anticoagulants;
- 192 • active bilateral epicondylitis;
- 193 • history *and/or* physical findings of cervical compression syndrome,
194 cervical or upper extremity arthritis, local arthrosis or neurologic
195 abnormality, rheumatoid disease, or radial nerve entrapment;
- 196 • arthrosis of the elbow, as confirmed by X-ray diagnosis (AP, lateral views)
- 197 • previous surgery for lateral epicondylitis;
- 198 • participated in a Workman’s Compensation Program or planned to apply
199 for the Program;
- 200 • thrombopathy, infection, tumor, or other severe systemic diseases;
- 201 • pregnancy;
- 202 • participated in a study with any experimental therapy within the last 30
203 days.

204

205

206 **Efficacy and Safety Endpoints:**

207

- 208 ♦ The primary efficacy endpoint was at least a 50% reduction from baseline to
209 Week 12 post-treatment in the pain visual analog scale (VAS) during resisted
210 wrist extension. The VAS is a 100-mm straight line with 0 for no pain and 100
211 for intolerable pain.

212

213

214
215
216
217
218
219
220
221
222
223
224
225
226
227
228
229
230
231
232
233
234
235
236
237
238
239
240
241
242
243
244
245
246
247
248
249
250
251
252
253
254

- ◆ The secondary efficacy endpoint was an improvement from baseline to Week 12 post-treatment in the patient's mean upper extremity function score. Function was assessed using the Upper Extremity Function Scale (UEFS) (Pransky et al.). The UEFS assesses limitations on 8 activities on a scale of 1-10 (1=no problem, 10=major problem/can't do).

Additional measures of efficacy were examined including the patient's overall impression (measured on 100mm VAS), grip strength testing (measured using a standard dynamometer), activity evaluation (ability to perform activities that were limited by his/her tennis elbow condition- measured similar to UEFS), and pain medication consumption. The patients grip strength was measured. Safety assessments included an assessment of adverse events, physical examination, X-rays, vital signs, 12-lead EKG, clinical labs, and proportion of patients who couldn't tolerate treatment.

Study Population:

Of the 114 patients enrolled in the study and included in the intent-to-treat (ITT) cohort, 56 patients were assigned to the active treatment group and 58 patients were assigned to the placebo group. Two (3.6%) active treatment group patients could not tolerate treatment and discontinued from the study prior to completing all three scheduled treatments. A third active treatment group patient was discontinued due to a low platelet count, which was found to be a pre-existing condition prior to study participation. Of the 58 placebo patients, 3 (5.2%) patients discontinued prior to the 12- week follow-up period to seek alternative therapy.

Patient demographics and treatment history are summarized in Table 2, below. The mean age for the active treatment group was 47 years (ranging from 35-71 years), and the mean age for the placebo group was 47 years (ranging from 35-60 years). There were 27 male (48.2%) and 29 female (51.8%) patients in the active treatment group and 27 male (46.6%) and 31 female (53.4%) patients in the placebo group. The mean height was 171 cm and the mean weight was 77kg. Physical exam and medical histories at baseline were also similar between the treatment groups.

255
256**Table 2: Patient Demographics and Treatment History**

Characteristic	Active Treatment Patients (N=56)	Placebo Treatment Patients (N=58)
Age (years)		
Mean	47	47.3
Range	35-71	35-60
Gender		
Male	27 (48.2%)	27 (46.6%)
Female	29 (51.8%)	31 (53.4%)
Height (cm)		
Mean	170.9	171.8
Range	152.4-188.0	149.9-190.5
Weight (kg)		
Mean	75.9	78.9
Range	50.9-120.0	53.0-120.2
Affected Arm		
Right	35 (62.5%)	41 (70.7%)
Left	21 (37.5%)	17 (29.3%)
Prior Therapies*		
All three	41 (73.2%)	43 (74.1%)
Steroid Injections & PT/OT	4 (7.1%)	6 (10.3%)
Steroid Injections & NSAIDs	6 (10.7%)	5 (8.6%)
PT/OT & NSAIDs	5 (8.9%)	4 (6.9%)
Symptom Duration (months)**		
Mean	21.3	20.8
Range	6.0-178.0	6.0-176.0

* PT/OT= physical and occupational therapy, NSAIDs= non-steroidal anti-inflammatories

** from date of initial diagnosis by a physician to enrollment into the study

257
258
259
260

Each treatment group (active, placebo) had lateral epicondylitis for an average of 21 months (12-month median) prior to randomization. In total, seventy-six (66.7%) patients had their right arm affected and 38 (33.3%) patients had their left arm affected. More than 70% of the patients in each treatment group had all three types of therapies (injections, PT/OT, NSAIDs) prior to enrollment, and 54 (93.1%) placebo patients and 51 (91.1%) active treatment patients had steroid injections.

268
269

Efficacy and Safety Results:

270
271

Efficacy Results:

272
273
274
275
276
277
278

For the ITT population (refer to Table 3, below), the placebo and active treatment groups had comparable pain scores at the baseline evaluation. The average pain score for patients who received the active treatment was 74 at baseline and 37.6 at 12 weeks. The average score for the placebo patients was 75.6 at baseline and 51.3 at 12 weeks.

279
280**Table 3: Summary of Patient Pain Assessments from Baseline to Week 12**

Treatment Group		Baseline	Week 1	Week 4	Week 8	Week 12
Active	N	56	56	56	56	56
	Mean	73.98	55.55	49.09	40.77	37.59
	SD	15.79	25.18	26.79	28.67	28.68
Placebo	N	58	58	58	58	58
	Mean	75.57	63.97	60.57	54.81	51.33
	SD	16.00	23.19	25.48	25.12	29.65

281
282
283
284

N=number of patients
SD= standard deviation

285 The primary efficacy endpoint was at least a 50% reduction from baseline to 12-
286 weeks post-treatment in the pain visual analog scale (VAS) during resisted wrist
287 extension. For the intent-to-treat cohort, the results show that the active
288 treatment group had 34/56 (60.7%) of the patients and the placebo group had
289 17/58 (29.3%) of the patients achieving at least a 50% reduction in pain during
290 provocation at Week 12 compared with baseline. There was a statistically
291 significant ($p=0.001$) between group difference.

292

293 The secondary efficacy endpoint was an improvement from baseline to Week 12
294 post-treatment in the patient's mean upper extremity function scale (UEFS)
295 score. For the ITT population (refer to Table 4, below), the placebo and active
296 treatment groups had comparable mean upper extremity function scores (UEFS)
297 at the baseline evaluation. The mean UEFS score for the active treatment group
298 was 4.68 at baseline (SD=1.78) and the mean UEFS score for the placebo group
299 was 4.63 (SD=1.8) at baseline. At Week 12, there was a statistically significant
300 ($p=0.01$) difference between groups in the mean UEFS scores, compared with
301 baseline.

302

303

304

Table 4: Summary of UEFS Assessments from Baseline to Week 12

Treatment Group		Baseline	Week 1	Week 4	Week 8	Week 12
Active	N	56	53	51	52	53
	Mean	4.68	3.23	2.80	2.54	2.25
	SD	1.78	1.89	1.70	1.52	1.57
Placebo	N	58	57	57	55	54
	Mean	4.63	3.71	3.79	3.54	3.23
	SD	1.80	1.77	1.98	2.12	2.09

305

306

307

308

N=number of patients
SD= standard deviation

309 The percent improvement in the average efficacy scores (pain, UEFS, patient's
 310 overall impression, activity level, and grip strength testing) at 12-weeks
 311 compared with baseline is summarized in Table 5, below.
 312

313 **Table 5: Percent Improvement in Average Efficacy Scores at 12 Weeks, Compared**
 314 **with Baseline**

	SONOCUR Basic			Placebo (Mock)		
	average score at beginning of study	average score at 12 weeks	% improvement	average score at beginning of study	average score at 12 weeks	% improvement
Pain *	74	37.6	49%	75.6	51.3	32%
UEFS *	4.7	2.3	51%	4.6	3.2	30%
Activity Evaluation *	7.7	3.5	55%	7.4	5	32%
Overall Impression *	70.3	32.8	53%	66	46.2	30%
Grip Strength Testing	71	87.1	23%	72.5	81.5	12%

315 * statistically significant ($p < 0.05$) between group difference. p-value is calculated using one way
 316 ANOVA

317

318

319

320

321 Safety Results:

322

323 In general, the nature, severity, frequency, duration and resolution of adverse
 324 events were similar in the active and placebo group, with the exception of certain
 325 vasovagal responses (i.e. nausea, sweating, dizziness, hypesthesia) and reports
 326 of pain during treatment for the active group.

327

328

329 Table 1 (page 7: Adverse Events) summarizes the number and frequency of
 330 adverse events that were categorized as being possibly or probably related to the
 331 study treatment.

332

333 The following table (Table 6) shows the occurrence of adverse events for the
 334 active treatment group that were judged to be possibly or probably treatment
 335 related over the course of the 12-month study period.

336

337

338
339

Table 6. Active Treatment Group: Adverse Events Through 12 Months of Follow-up

Adverse Events	Treatment Period				Follow-up Period [1]					
	Day of Treatment (Treatments 1, 2, 3) [N=56]		Between Treatments 1-2 [N=56]	Between Treatments 2-3 [N=54]	<=1 Week [N=53]	>1-4 Weeks [N=51]	>4-8 Weeks [N=52]	>8-12 Weeks [N=53]	>12 Weeks- 6 Months [N=48]	>6-12 Months [N=46]
	# of patients [2]	# of occurrences	# of patients [2]	# of patients [2]	# of patients [2]	# of patients [2]	# of patients [2]	# of patients [2]	# of patients [2]	# of patients [2]
Pain	24	46	6	5	4	5	3	1	0	0
Nausea	10	10	0	0	0	0	0	0	0	0
Application Site Reaction	4	6	2	2	3	2	2	1	1	*
Sweating	5	5	0	0	0	0	0	0	0	0
Dizziness	4	4	0	0	0	0	0	0	0	0
Hypertonia	2	3	1	0	2	1	1	1	0	0
Hypesthesia	3	5	0	0	0	0	0	0	0	0
Paresthesia	3	4	0	0	0	0	0	0	0	0
Joint Stiffness	1	1	1	2	1	1	1	1	1	*
Myalgia	1	1	1	0	0	0	0	0	0	0
Tremor	2	2	0	0	0	0	0	0	0	0
Vasodilation	2	2	0	0	0	0	0	0	0	0
Pallor	1	1	0	0	0	0	0	0	0	0

340

341 [1] relative to the last treatment, using protocol defined windows. No new adverse events judged
 342 as being possibly or probably related to treatment were reported after the 6-month follow-up
 343 period.

344 [2] number of patients experiencing at least one occurrence within each time interval

345 * patient discontinued from study during the 6-month follow-up visit and subsequently had surgery
 346 for tennis elbow. No additional study data is available for this patient after the 6-month follow-up
 347 period.

348

349

350

347
348
349

351 At the time of the 12-week follow-up visit, all device related adverse events had
352 resolved, except for one patient who had moderate elbow stiffness and mild
353 swelling that was still ongoing at the time of the 6-month follow-up visit. This
354 patient, who had an x-ray with normal findings at baseline and 12 weeks, was
355 unresponsive to treatment and terminated study participation soon after the 6-
356 month follow-up visit for surgery. There were no new device related adverse
357 effects reported during the long-term (3-12month) follow-up period.

358

359 In addition to adverse events, lab values, physical exam results, X-rays, vital
360 signs, and EKGs were assessed. No significant between group differences were
361 observed.

362

363 ***Training:***

364

365 **The Physician**

366 The physician prescribing ESWT with SONOCUR should be trained and
367 experienced in the care of patients with tendinopathies and /or lateral
368 epicondylitis. This includes conventional therapies for this indication, as well as
369 knowledge of exclusion criteria and precautions as described in this manual. The
370 physician should identify the region to be treated and supervise the treatment
371 sessions performed by a trained operator.

372

373 **The Operator**

374 The operator should be trained by an authorized* application specialist. The
375 training includes:

376

- 377 • Study of the operating manual
- 378 • SONOCUR functionality of all keys and menu options
- 379 • SONOCUR transportation and storage
- 380 • System set-up and daily maintenance procedures
- 381 • System malfunctions, error codes and proceedings
- 382 • Patient seating and positioning of the patient's arm.
- 383 • Positioning of the shockwave head and treatment protocol
- 384 • Knowledge of potential adverse events (AE) and procedures if AE's occur.
- 385 • Cleaning and disinfecting of coupling bellows

386

387 The operator should perform at least two to four (2-4) procedures under the
388 guidance of the application specialist. The name of the trainee and the trainer
389 should be kept on file with the clinic along with the training date.
390

391 * The application specialist will be provided by the manufacturer or distributor.
392

393 *System start-up*

394 Establish the line connection to the SONOCUR Basic

395 *Switching on the therapy system*

396

397 Press this button to turn on the power switches on the front panel of the
398 SONOCUR Basic.

399

400 Run the daily function and safety tests before beginning therapy.
401

402 *Daily tests*

403 Before beginning therapy

404 As the user, you must make sure that all safety-related equipment is functioning
405 properly and the system is ready for operation.

406 Check the function of the unit foot brakes

407 Inspect the articulating arm travel and counterbalance

408 Inspect all readings on control display

409 *Switching off the therapy system*

410 Press this button to turn off power via the switch on the front panel of the
411 SONOCUR Basic

412

413

414

415