

1 Patient Labeling (Date)

2 **Synvisc-One™**

3

4 Be sure to read the following important information carefully. This information does not take the
5 place of your doctor's advice. If you do not understand this information or want to know more, ask
6 your doctor.

7

8 **Glossary of Terms**

- 9 Hyaluronan (pronounced hy-al-u-ROE-nan): is a natural substance that is present in very high
10 amounts in joints. It acts like a lubricant and a shock absorber in the joint and is needed for the
11 joint to work properly.
- 12 Non-steroidal anti-inflammatory drugs: also known as "NSAIDs"; medication used to treat pain or
13 swelling. There are many examples of NSAIDs, including (but not limited to) aspirin and ibuprofen.
14 Some of these are over-the-counter drugs, and some can only be obtained by prescription.
- 15 Osteoarthritis (pronounced OS-te-o-arth-RI-tis): (OA) is a type of arthritis that involves the
16 wearing down of cartilage (the protective covering on the ends of your bones) and loss of
17 cushioning fluid in the joint

18 **Table of Contents**

19

20 **Synvisc-One™**

21 Glossary of Terms

22 Table of Contents

23 What is Synvisc-One™?

24 How is Synvisc-One™ used? (Indications)

25 Are there any reasons why I should not receive Synvisc-One™? (Contraindications)

26 What should my doctor warn me about?

27 What are the risks of using Synvisc-One™?

28 What are the benefits of Synvisc-One™?

29 How is Synvisc-One™ given?

30 What do I need to do after I get Synvisc-One™?

31 What other treatments are available for OA?

32 Non-drug treatments

33 Drug therapy
34 When should I call my doctor? (Troubleshooting)
35 What adverse events were observed in the clinical study?
36 How do I get more information about Synvisc-One? (User Assistance)
37

38 **What is Synvisc-One™?**

39 Synvisc-One is a gel-like mixture that comes in a syringe containing 6 mL (1 ½ teaspoon) and is
40 injected into your knee. It is made up of hylan A fluid, hylan B gel, and salt water. Hylan A and
41 hylan B are made from a substance called hyaluronan (pronounced hy-al-u-ROE-nan), also
42 known as sodium hyaluronate that comes from chicken combs. Hyaluronan is a natural
43 substance found in the body and is present in very high amounts in joints. The body's own
44 hyaluronan acts like a lubricant and a shock absorber in the joint and is needed for the joint to
45 work properly.

46

47 **How is Synvisc-One™ used? (Indications)**

48 The FDA-approved indication for Synvisc-One is:

49 Synvisc-One is indicated for the treatment of pain in osteoarthritis (OA) of the knee in patients
50 who have failed to respond adequately to conservative non-pharmacologic therapy and simple
51 analgesics, e.g., acetaminophen.

52 **Are there any reasons why I should not receive Synvisc-** 53 **One™? (Contraindications)**

54 Your doctor will determine if there is any reason why you are not an appropriate candidate for
55 Synvisc-One. You should be aware that Synvisc-One:

- 56 • Should not be used in patients who have had any prior allergic reactions to Synvisc,
57 Synvisc-One or any hyaluronan-based products. Signs of an allergic reaction may
58 include swelling of your face, tongue, or throat; difficulty breathing or swallowing;
59 shortness of breath; wheezing; chest pain; a tightness in your throat; sleepiness; rash;
60 itching; hives; flushing; and/or fever.

- 61 • Should not be used in patients with a knee joint infection, skin disease or infection around
62 the area where the injection will be given, or circulatory problems in the legs.

63 **What should my doctor warn me about?**

64 The following are important treatment considerations for you to discuss with your doctor and
65 understand in order to help avoid unsatisfactory results and complications:

- 66 • Synvisc-One is only for injection into the knee, performed by a doctor or other qualified
67 health care professional. Synvisc-One has not been tested to show pain relief in joints
68 other than the knee.
- 69 • Synvisc-One has not been tested to show better pain relief when combined with other
70 injected medicines.
- 71 • Tell your doctor if you are allergic to products from birds such as feathers, eggs, and
72 poultry.
- 73 • Tell your doctor if you have significant swelling or blood clots in the leg.
- 74 • Synvisc-One has not been tested in pregnant women, or women who are nursing. You
75 should tell your doctor if you think you are pregnant, or if you are nursing a child.
- 76 • Synvisc-One has not been tested in children (\leq 21 years of age).

77 **What are the risks of using Synvisc-One™?**

78 The side effects (also called reactions) sometimes seen after any injection into the knee,
79 including Synvisc-One, include: pain, swelling, heat, redness, and/or fluid build-up around the
80 knee. These reactions are generally mild and do not last long. Reactions are generally treated
81 by resting and applying ice to the injected knee. Sometimes it is necessary to give pain relievers
82 by mouth such as acetaminophen or NSAIDs, or to give injections of steroids, or to remove fluid
83 from the knee joint. Patients rarely undergo arthroscopy (a surgical inspection of the knee joint)
84 or other medical procedures related to these reactions.

85 Other side effects seen with Synvisc or Synvisc-One are: rashes, hives, itching, muscle
86 pain/cramps, flushing and/or swelling of your face, fast heartbeat, nausea (or feeling sick to your

87 stomach), dizziness, fever, chills, headache, difficulty breathing, swelling in your arms and/or legs,
88 prickly feeling of your skin, and in rare cases a low number of platelets in the blood (platelets are
89 a type of blood cell that are needed to help clot your blood when you are cut or injured). Rare
90 cases of knee joint infection have been reported. If any of the above side effects or symptoms
91 appear after you are given Synvisc-One, or if you have any other problems, you should call your
92 doctor.

93

94 **What are the benefits of Synvisc-One™?**

95 As shown in a medical study of 253 patients with osteoarthritis (OA) of the knee, where
96 approximately half received either a single injection of Synvisc-One or an injection of the same
97 volume of salt water (a "Saline Control" injection), the major benefits of Synvisc-One are pain
98 relief and improvement in other symptoms related to OA of the knee.

99

100 **How is Synvisc-One™ given?**

101 Your doctor will inject Synvisc-One into your knee.

102 **What do I need to do after I get Synvisc-One™?**

103 It is recommended you avoid strenuous activities (for example, high-impact sports such as tennis
104 or jogging) or prolonged weight-bearing activities for approximately 48 hours following the
105 injection. You should consult your doctor regarding the appropriate time to resume such activities.

106 **What other treatments are available for OA?**

107 If you have OA, there are other things you can do besides getting Synvisc-One. These include:

108 ***Non-drug treatments***

- 109 ▪ Avoiding activities that cause knee pain
- 110 ▪ Exercise or physical therapy
- 111 ▪ Weight loss
- 112 ▪ Removal of excess fluid from your knee

113 **Drug therapy**

- 114 ▪ Pain relievers such as acetaminophen and narcotics
- 115 ▪ Drugs that reduce inflammation (signs of inflammation are swelling, pain or redness), such as
- 116 aspirin and other non-steroidal anti-inflammatory drugs (NSAIDs, for example ibuprofen and
- 117 naproxen)
- 118 ▪ Steroids that are injected directly into your knee

119

120 **When should I call my doctor? (Troubleshooting)**

121 If any of the side effects or symptoms described above appear after you are given Synvisc-One,
122 or if you have any other problems, you should call your doctor.

123

124 **What did the clinical studies show?**

125 A study was conducted in 6 countries outside the United States with 21 physicians. The patients
126 in the study had mild to moderate knee OA, moderate to severe pain, and did not have sufficient
127 relief of their pain and symptoms with medications taken by mouth.

128 A total of 253 patients in the study were assigned by chance to receive either a single injection of
129 Synvisc-One (n=123 patients), or an injection of the same volume of salt water (a "Saline Control"
130 injection) (n=130 patients). Neither the patients nor the doctors evaluating them knew which
131 treatment they received. Any fluid that was present in the patient's knee was removed before the
132 injection. The patients were seen by their doctor at standard times over 6 months. Information
133 was collected about how much pain they were experiencing doing various types of activities, how
134 much they were limited in their daily activities by their OA, and on their overall condition. Their
135 doctor also provided an overall rating of their OA.

136 The main measure of the study was how much pain the subjects had doing five common types of
137 activities over the 6 months duration of the study. Daily activity limitations and overall evaluations
138 were also compared between the group of patients receiving Synvisc-One injection and the group
139 receiving salt water injection.

140 The study showed that patients receiving Synvisc-One had significantly less pain over 6 months,
141 and felt significantly better than the patients who received the salt water injections. The
142 difference in pain score reduction from baseline to 6 months between the Synvisc-One and salt
143 water control injection was 0.15 out of a 5 point scale for the measurement of OA pain in the knee.

144 **What adverse events were observed in the clinical** 145 **study?**

146 The following are the most common adverse events that occurred during the clinical trial of
147 Synvisc-One:

- 148 • Pain in the knee or at the injection site
- 149 • Stiffness, swelling or warmth in or around the knee
- 150 • Changes in the way that you walk (e.g. limping)

151

152 Severe adverse events were not observed in the Synvisc-One trial. Joint infections did not occur
153 in the injected knee in the Synvisc-One clinical trial. The most commonly occurring adverse
154 events outside of the injected knee were headache, back pain, sore throat, the flu and one
155 episode of feeling faint.

156

157 **How do I get more information about Synvisc-One?** 158 **(User Assistance)**

159 If you have any questions or would like to find out more about Synvisc-One, you may call

160 Genzyme Biosurgery at 1-888-3-SYNVISC (1-888-379-6847) or visit www.synvisc.com. Synvisc
161 and GENZYME are registered trademarks of Genzyme Corporation. Synvisc-One is a trademarks
162 of Genzyme Corporation.

163

164 **Manufactured and Distributed by:** 165 Genzyme Biosurgery

166 A division of Genzyme Corporation

167 1125 Pleasant View Terrace

1 **SYNVISC-ONE™**

2
3 Package contents provided sterile.

Genzyme Biosurgery, a division of Genzyme Corporation
1125 Pleasant View Terrace
Ridgefield, New Jersey 07657
Telephone: 1-888-3-SYNVISC (1-888-379-6847)
www.synvisc.com

4
5
6

Information for Prescribers

7 Caution: Federal law restricts this device to sale by or on the order of a physician (or properly
8 licensed practitioner).

9
10 **DESCRIPTION**

11 Synvisc-One™ (hylan G-F 20) is an elastoviscous high molecular weight fluid containing hylan A
12 and hylan B polymers produced from chicken combs. Hylans are derivatives of hyaluronan
13 (sodium hyaluronate). Hylan G-F 20 is unique in that the hyaluronan is chemically crosslinked.
14 Hyaluronan is a long-chain polymer containing repeating disaccharide units of Na-glucuronate-N-
15 acetylglucosamine.

16
17 **INDICATIONS FOR USE**

18 Synvisc-One is indicated for the treatment of pain in osteoarthritis (OA) of the knee in patients
19 who have failed to respond adequately to conservative nonpharmacologic therapy and simple
20 analgesics, e.g., acetaminophen.

21
22 **CONTRAINDICATIONS**

- 23 • Do not administer to patients with known hypersensitivity (allergy) to hyaluronan (sodium
24 hyaluronate) preparations.
- 25 • Do not inject Synvisc-One in the knees of patients having knee joint infections or skin diseases
26 or infections in the area of the injection site.

28 **WARNINGS**

- 29 • Do not concomitantly use disinfectants containing quaternary ammonium salts for skin
30 preparation because hyaluronan can precipitate in their presence.
- 31 • Do not inject Synvisc-One extra-articularly or into the synovial tissues and capsule.
- 32 • Intravascular injections of Synvisc-One may cause systemic adverse events.

33

34 **PRECAUTIONS**

35 **General**

- 36 • The safety and efficacy of Synvisc-One in locations other than the knee and for conditions other
37 than osteoarthritis have not been established.
- 38 • The safety and effectiveness of the use of Synvisc-One concomitantly with other intra-articular
39 injectables have not been established.
- 40 • Use caution when injecting Synvisc-One into patients who are allergic to avian proteins,
41 feathers or egg products.
- 42 • The safety and efficacy of Synvisc-One in severely inflamed knee joints have not been
43 established.
- 44 • Strict aseptic administration technique must be followed.
- 45 • **STERILE CONTENTS.** The syringe is intended for single use. The contents of the syringe must
46 be used immediately after its packaging is opened. Discard any unused Synvisc-One.
- 47 • Do not use Synvisc-One if package is opened or damaged. Store in original packaging
48 (protected from light) at room temperature below 86° F (30° C). **DO NOT FREEZE.**
- 49 • Remove any synovial fluid or effusion before injecting Synvisc-One.
- 50 • Synvisc-One should be used with caution when there is evidence of lymphatic or venous stasis in
51 the leg to be injected.

52 **Information for Patients**

- 53 • Provide patients with a copy of the Patient Labeling prior to use.

- 54 • Mild to moderate pain, swelling and/or effusion of the injected knee have been reported in
55 clinical trials that were related to intra-articular injection of Synvisc-One. These events were
56 typically transient and usually resolved on their own or with conservative treatment.
- 57 • As with any invasive joint procedure, it is recommended that the patient avoid strenuous
58 activities (for example, high impact sports such as soccer, tennis or jogging) or prolonged weight-
59 bearing activities for approximately 48 hours following the intra-articular injection. The patient
60 should consult his or her physician regarding the appropriate time to resume such activities.

61 **Use in Specific Populations**

- 62 • **Pregnancy:** The safety and effectiveness of Synvisc-One have not been established in
63 pregnant women.
- 64 • **Nursing mothers:** It is not known if Synvisc-One is excreted in human milk. The safety and
65 effectiveness of Synvisc-One have not been established in lactating women.
- 66 • **Pediatrics:** The safety and effectiveness of Synvisc-One have not been established in pediatric
67 patients. Pediatric patients are defined as patients ≤ 21 years of age.

68

69 **POTENTIAL ADVERSE EFFECTS OF THE DEVICE ON HEALTH**

70 **Reported Device-Related Adverse Events**

71 The most commonly reported adverse events associated with Synvisc-One are the following:

72

- 73 • Arthralgia
- 74 • Arthritis
- 75 • Arthropathy
- 76 • Injection site pain
- 77 • Joint effusion

78

79 A complete list of the frequency and rate of adverse events identified in the clinical study are
80 provided in the Safety section (Table 3).

81

82 **Potential Adverse Events**

83 The following adverse events are among those that may occur in association with intra-articular
84 injections, including Synvisc-One:

85

- 86 • Arthralgia
- 87 • Joint stiffness
- 88 • Joint effusion
- 89 • Joint swelling
- 90 • Joint warmth
- 91 • Injection site pain
- 92 • Arthritis
- 93 • Arthropathy
- 94 • Gait disturbance

95
96 A complete list of the frequency and rate of adverse events identified in the clinical study are
97 provided in the Safety section (Table 2).

98
99 **Post-marketing Experience**

100 Synvisc® (3 injection regimen) post-marketing experience has identified the following systemic
101 events to occur rarely with administration: rash, hives, itching, fever, nausea, headache, dizziness,
102 chills, muscle cramps, paresthesia, peripheral edema, malaise, respiratory difficulties, flushing
103 and facial swelling. There have been rare reports of thrombocytopenia coincident with Synvisc (3
104 injection regimen) injection.

105
106 ***Pivotal Clinical Trial***

107 Study Design: To determine the safety and effectiveness of a single injection regimen of Synvisc-
108 One in the reduction of the pain score in osteoarthritis of the knee, a prospective, randomized,
109 double-blind, 2-arm (parallel group) clinical trial in 21 centers in six European countries was
110 conducted. A total of 253 patients were randomly assigned to study treatment; 123 received 6
111 mL of Synvisc-One and 130 received 6 mL of Phosphate-Buffered Saline. Neither the patients
112 nor the clinical observers knew the patients' treatment allocations. The outcome measures
113 collected included the Western Ontario and McMaster Universities Osteoarthritis Index (WOMAC;
114 Likert 3.1 A version), patient global assessment (PTGA), clinical observer global assessment
115 (COGA), and use of rescue analgesic (see Treatment and Evaluation Schedule). The intent-to-
116 treat (ITT) population (all patients randomized) was used for the primary analysis. The primary
117 efficacy analysis was a comparison over 26 weeks between the two treatment groups of change
118 from baseline in the WOMAC A (Pain) Subscale (see Patient Population and Demographics),
119 performed by analysis of covariance (ANCOVA).

120 **Patient Population and Demographics**

121 Study patients had primary osteoarthritis of the knee per American College of Rheumatology
122 criteria and were at least 40 years old. The diagnosis was confirmed via recent radiograph
123 showing at least one osteophyte in the target knee. Study patients had continued target knee
124 pain despite use of conservative treatment and analgesics/non-steroidal anti-inflammatory drugs
125 (NSAIDs). Patients with severe disease (Grade IV) per Kellgren-Lawrence criteria, or who had
126 prior arthroplasty in the target knee, were excluded. At the beginning of the study, subjects had
127 moderate or severe target knee pain when walking on a flat surface (on a 5-point Likert scale
128 where 0 = none, 1 = mild, 2 = moderate, 3 = severe, 4 = extreme), and an average score of 1.5 to
129 3.5 on the five questions of the WOMAC A (Pain) Subscale. The WOMAC A Subscale asks
130 study subjects to rate their degree of pain when:

- 131 1. Walking on a flat surface
- 132 2. Going up and down stairs
- 133 3. Resting during the night
- 134 4. Sitting or lying
- 135 5. Standing upright

136 Table 1 summarizes the demographics and baseline characteristics. There were no clinically
137 meaningful differences between treatment groups in any baseline parameter.

Table 1: Summary of Demographic and Baseline Characteristics

Parameter/Category	Synvisc-One (N=124)*	Saline Control (N=129)*	Total (N=253)
Age, n *	124	129	253
Mean (SD)	63.6 (9.6)	62.5 (9.2)	63.0 (9.4)
Range	42, 83	43, 84	42, 84
Sex, n *	124	129	253
Female, n (%)	92 (74%)	88 (68%)	180 (71%)
Race, n *	124	129	253
Caucasian, n (%)	118 (95%)	125 (97%)	243 (96%)
Non-Caucasian, n (%)	6 (5%)	4 (3%)	10 (4%)
Body Mass Index (kg/m ²), n *	123	129	252
Mean (SD)	29.1 (4.8)	29.8 (5.7)	29.4 (5.3)
Range	20.7, 46.0	19.5, 52.4	19.5, 52.4
Prior Corticosteroids In Target Knee, n **	123	130	253
Yes – n (%)	40 (32%)	31 (24%)	71 (28%)
Prior Arthroscopy In Target Knee, n **	123	130	253
Yes – n (%)	26 (21%)	28 (22%)	54 (21%)
Tibio-Femoral Joint Modified Kellgren-Lawrence Numerical Grading System **			
Grade II	63 (51%)	51 (39%)	114 (45%)
Grade III	60 (49%)	78 (60%)	138 (55%)
Grade IV	0	1 (1%)	1 (0%)
Total WOMAC Score (0-96); Mean (SD) *	55.1 (10.5)	54.8 (9.4)	
WOMAC A Score (0-4); Mean (SD) *	2.30 (0.43)	2.25 (0.41)	
PTGA – Mean (SD) (0-4) *	2.57 (0.67)	2.50 (0.64)	
COGA – Mean (SD) (0-4) *	2.44 (0.76)	2.49 (0.75)	

* ITT Population

** Safety Population

139
140
141
142

143 **Treatment and Evaluation Schedule**

144 *Initial Treatment Phase*

145 Patients were followed for 26 weeks. Study visits were scheduled for screening, baseline, and
146 weeks 1, 4, 8, 12, 18 and 26. Injections were performed aseptically at the baseline visit after
147 arthrocentesis to withdraw any effusion or synovial fluid present. Patients were not permitted to
148 take long-acting NSAIDs (including cyclo-oxygenase II inhibitors), opioid analgesics or
149 corticosteroids (by any route) during the study, but were permitted to take up to 4 g per day of
150 acetaminophen as needed for "rescue" of injected knee pain. "Rescue" medication was not
151 permitted within 48 hours of any study visit. Injected knee assessment, patient and clinician
152 global assessments (PTGA & COGA), WOMAC and safety evaluations were performed at each
153 study visit.

154 *Repeat Treatment Phase*

155 If patients in either blinded treatment group had at least mild pain in the injected knee at the week
156 26 visit (and did not experience any significant clinical concerns after the first treatment
157 administration), they were offered an injection of (open-label) Synvisc-One. Those who chose to
158 receive the second injection were followed for 4 weeks for safety only.

159

160 **Adverse Event Summary:**

161 The frequency and type of adverse events (AEs) were similar between the group of patients that
162 received Synvisc-One and the group that received saline control.

163 **Initial Treatment Phase:** The overall proportions of patients with Treatment-Emergent AEs
164 regardless of device relatedness (Synvisc-One: n=70, 56.9%; Saline Control: n=79, 60.8%) and
165 with injected knee AEs regardless of device relatedness (Synvisc-One: n=44, 35.8%; Saline
166 Control: n=44, 33.8%) were comparable between the two treatment groups (See Table 2).

167 Table 3 lists the incidences of AEs in the injected knee that were assessed by the investigator to
168 be device-related, defined as related to either the study injection or the study treatment.

169

170
171

Table 2: Patients with Adverse Events in the Injected Knee Regardless of Relatedness

MedDRA Preferred Term	Synvisc-One N=123 n (%)	Saline Control N=130 n (%)
Any Treatment-Emergent Adverse Event	44 (35.8%)	44 (33.8%)
Arthralgia	31 (25.2%)	28 (21.5%)
Joint stiffness	10 (8.1%)	13 (10.0%)
Joint effusion	7 (5.7%)	7 (5.4%)
Joint swelling	5 (4.1%)	7 (5.4%)
Joint warmth	2 (1.6%)	5 (3.8%)
Post-traumatic pain	0	3 (2.3%)
Injection site pain	1 (0.8%)	1 (0.8%)
Synovial cyst	0	2 (1.5%)
Arthritis	1 (0.8%)	0
Arthropathy	1 (0.8%)	0
Gait disturbance	1 (0.8%)	0
Joint range of motion decreased	0	1 (0.8%)
Osteoarthritis	0	1 (0.8%)

172 Note: Patients are counted once for each unique AE regardless of device relatedness, and may have had
173 more than one unique AE
174

175 **Table 3: Patients with Device-Related Adverse Events in the Injected Knee**

MedDRA Preferred Term	Synvisc-One N=123 n (%)	Saline Control N=130 n (%)
Any Device-Related Adverse Event	7 (5.7%)	4 (3.1%)
Arthralgia	2 (1.6%)	3 (2.3%)
Arthritis	1 (0.8%)	0
Arthropathy	1 (0.8%)	0
Injection site pain	1 (0.8%)	1 (0.8%)
Joint effusion	2 (1.6%)	0

176 Note: Patients are counted once for each unique AE, and may have had more than one unique AE
177

178 Device-related AEs involving the injected knee were mild or moderate in nature and were treated
179 symptomatically. There were no serious AEs in the injected knee in either the Synvisc-One or the
180 saline control group.

181 **Repeat Treatment Phase:** The repeat treatment phase evaluated the safety profile of the initial
182 phase of patients receiving a second injection of Synvisc-One. One hundred and sixty patients

			Poor, Very poor)	
COGA	Over 26 weeks	0.71*	The odds [probability (Worse) / Probability (Better)] for Synvisc-One for over 26 weeks and at 26 weeks is approximately 71%, and 56%, respectively, to the odds for control. COGA: Clinical Observer Global Assessment has 5 scales (Very well, Well, Fair, Poor, Very poor)	Blinded clinical observers were 1.41 times more likely to assess patients treated with Synvisc-One as showing overall improvement in disease status compared to those patients treated with saline control over 26 weeks and 1.79 times more likely to assess patients treated with Synvisc-One as showing overall improvement in disease status compared to those patients treated with saline control at 26 weeks.
	At week 26	0.56*		
OMERACT-OARSI Responder	Over 26 weeks	0.66	This responder analysis did not reach statistical significance between the treatment groups.	
	At week 26	0.69		
Estimate of Treatment Difference (Analysis of Covariance)				
WOMAC C	Over 26 weeks	-0.18	The study did not show a statistically significant difference in functional improvement between the treatment groups.	
	At week 26	-0.11		

226 * Statistically significant at the 5% significance level; not adjusted for multiplicity

227 ¹Odds ratio = Odds for Synvisc-One/Odds for Control

228 Odds ratio = [Probability (Worse) / Probability (Better) for Synvisc-One] / [Probability (Worse) /
229 Probability (Better) for Control]

230 If odds ratio < 1, then in favor of Synvisc-One


231

232 The WOMAC A1 responder rate (where response was defined as a 1-or-more category
233 improvement from baseline and the patient did not withdraw from the study) was significantly
234 higher in the Synvisc-One group than in the saline control group. Seventy-one percent (71%) of
235 the patients were responders at week 18 in the Synvisc-One group (versus 54% in the saline
236 control group). At week 26, 64% of patients in the Synvisc-One group were responders, while
237 only 50% of patients in the saline control group were responders.
238

222

223

Figure 1: Plot for Categorical Secondary Endpoints – ITT Population



224

225

Table 5: Clinical Meaning of Secondary Efficacy Endpoints

	Odds Ratio		Definition	Explanation
Generalized Estimating Equation for categorical data.				
WOMAC A1	Over 26 weeks	0.64*	The odds [probability (Worse) / Probability (Better)] for Synvisc-One for over 26 weeks and at 26 weeks is approximately 64%, and 56%, respectively, to the odds for control.	Synvisc-One patients were 1.56 times more likely to self-report pain relief while walking on a flat surface compared to those patients treated with saline control over 26 weeks and 1.79 times more likely to self-report pain relief while walking on a flat surface compared to those patients treated with saline control at 26 weeks.
	At week 26	0.56*		
PTGA	Over 26 weeks	0.69*	The odds [probability (Worse) / Probability (Better)] for Synvisc-One for over 26 weeks and at 26 weeks is approximately 69%, and 51%, respectively, to the odds for control. PTGA: Patient Global Assessment has 5 scales (Very well, Well, Fair,	Synvisc-One patients were 1.45 times more likely to self-report improvement in overall health status compared to those patients treated with saline control over 26 weeks and 1.96 times more likely to self-report improvement in overall health status compared to those patients treated with saline control at 26 weeks.
	At week 26	0.51*		

284 Precaution: Do not over tighten or apply excessive leverage when attaching the needle or
285 removing the needle guard, as this may break the syringe tip.

286 • Inject the full 6 ml in one knee only.

287

288 **MANUFACTURED AND DISTRIBUTED BY:**

289 Genzyme Biosurgery a division of Genzyme Corporation

290 1125 Pleasant View Terrace

291 Ridgefield, New Jersey 07657

292 Telephone: 1-888-3-SYNVISC (1-888-379-6847)

293 Covered by U.S. patents #4,636,524, #4,713,448, #5,099,013, #5,143,724.

294 Synvisc-One™ is a Trademark of Genzyme Corporation. SYNVISC and GENZYME are
295 registered trademarks of Genzyme Corporation.

296

- 258 • Sodium chloride 51 mg
- 259 • Disodium hydrogen phosphate 0.96 mg
- 260 • Sodium dihydrogen phosphate monohydrate 0.24 mg
- 261 • Water for injection q.s. to 6.0 ml

262

263 **HOW SUPPLIED**

264 Synvisc-One is supplied in a 10-ml glass syringe containing 3 doses (48 mg) of hylan G-F 20.

265 The contents of the syringe are sterile and non-pyrogenic.

266

267 **DIRECTIONS FOR USE**

268 Precaution: Do not use Synvisc-One if the package has been opened or damaged. Store in the
269 original packaging (protected from light) at room temperature below 86°F (30°C). DO NOT
270 FREEZE.

271 Precaution: The syringe containing Synvisc-One is intended for single use. The contents of the
272 syringe must be used immediately after the syringe has been removed from its packaging.

273 Precaution: Do not concomitantly use disinfectants containing quaternary ammonium salts for
274 skin preparation because hyaluronan can precipitate in their presence.

275 Synvisc-One is administered as a single intra-articular injection. Strict aseptic administration
276 technique must be followed.

277 • Using an 18 to 20 gauge needle, remove synovial fluid or effusion before injecting
278 Synvisc-One.


279 • Do not use the same syringe for removing synovial fluid and for injecting Synvisc-One;
280 however, the same 18 to 20 gauge needle should be used.

281 • Twist the tip cap before pulling it off, as this will minimize product leakage.

282 • To ensure a tight seal and prevent leakage during administration, secure the needle
283 tightly while firmly holding the luer hub.

239

240 **Figure 2: Patient Responder Rate on WOMAC A1 (Walking Pain) -ITT Population**



241

242

243 Note: Analyzed using generalized estimating equation (GEE) for binary outcomes

244 * Statistically significant at the 5% significance level; not adjusted for multiplicity

245

246

247 **DETAILED DEVICE DESCRIPTION**

248 Synvisc-One combines the three doses of Synvisc (hylan G-F 20) which consists of hylan A

249 (average molecular weight 6,000,000 daltons) and hylan B hydrated gel in a buffered

250 physiological sodium chloride solution, pH 7.2. Synvisc-One has an elasticity (storage modulus

251 G') at 2.5 Hz of 111 ± 13 Pascals (Pa) and a viscosity (loss modulus G'') of 25 ± 2 Pa (elasticity

252 and viscosity of knee synovial fluid of 18 to 27-year-old humans measured with a comparable

253 method at 2.5 Hz: $G' = 117 \pm 13$ Pa; $G'' = 45 \pm 8$ Pa.)

254

255 Each 10 mL syringe of Synvisc-One combines the three 2 mL doses (16 mg each) of a complete

256 Synvisc treatment regimen (48 mg). Each Synvisc-One 10 mL syringe contains:

- 257
- Hylan polymers (hylan A + hylan B) 48 mg

211

212

213

214

Table 4: Primary Efficacy Results: WOMAC A (Pain) Score Overall Change from Baseline Over 26 Weeks – ITT Population

	Baseline Mean (SE) (0-4 Scale)	Mean Post-treatment (SE) (0-4 Scale)	Estimated Change (SE)	Estimated Difference from Saline Control (95% CI)	p-value (ANCOVA)
Synvisc-One (n=124)	2.30 (0.04)	1.43 (0.06)	-0.84 (0.06)	-0.15 (-0.302, -0.002)	0.047
Saline Control (n=129)	2.25 (0.04)	1.59 (0.06)	-0.69 (0.06)		

215

WOMAC A scale using 5 point Likert scale, where 0=no pain and 4 =extreme pain

216

Repeated measures Analysis of Covariance was used for the WOMAC A pain score change from the baseline.

217

218

Synvisc-One also demonstrated superiority to saline control in multiple pre-defined secondary

219

outcome measures, which included PTGA over and at 26 weeks, COGA over and at 26 weeks,

220

and pain while walking on a flat surface (WOMAC A1) over and at 26 weeks (see Figure 1 &

221

Table 5).

183 were treated during this phase of the study, of which 77 patients received a second injection of
184 Synvisc-One. Of these 77 patients, 4 (5.2%) experienced five device-related AEs in the injected
185 knee. All such events were mild to moderate and were treated symptomatically. These events
186 were arthralgia (n=2), arthritis (n=1), injection site hematoma (n=1) and injection site pain (n=1).
187 Patients who developed injected knee AEs during the initial phase of the study, and who
188 subsequently received repeat treatment, did not experience injected knee AEs upon repeat
189 exposure to Synvisc-One.

190 **Overall Injected Knee Safety Summary:** The safety profile of Synvisc-One is similar to the
191 Clinical and Post-marketing experience seen with Synvisc (3 injection regimen) where pain,
192 swelling and effusion were the most frequently occurring AEs in the injected knee. There have
193 been post-marketing reports for Synvisc indicating that in some cases the joint effusion may be
194 large and can cause pronounced pain; it is important to remove and to analyze the fluid to rule
195 out infection or crystalline arthropathies. These types of severe AEs were not observed in either
196 the initial or repeat treatment phase of the Synvisc-One trial. Joint infections did not occur in any
197 of the clinical trials of Synvisc or Synvisc-One and have been reported only rarely during clinical
198 use of Synvisc.

199

200 **Adverse Events Outside of the Injected Knee**

201 Overall 101 patients (Synvisc-One: n=47, 38.2%; Saline Control: n=54, 41.5%) experienced at
202 least one AE outside the injected knee regardless of device relatedness. The most commonly
203 occurring (5% or greater in either group) AEs outside the injected knee were headache, back pain,
204 nasopharyngitis and influenza. In the Synvisc-One group there was one AE of syncope
205 considered device-related.

206 No new systemic AEs were identified during this study as compared to Synvisc.¹

207

208 **Primary Efficacy Endpoint:**

209 The primary endpoint for the study, the difference between the treatment groups in change from
210 baseline over 26 Weeks in the WOMAC A Pain Score (Table 4) was met.