

MAY 29 2001

K010014

Alfa Scientific Designs, Inc.

11494 Sorrento Valley Road, Suite M
San Diego, CA 92121

510(K) Summary

In accordance with the Safe Medical Devices Act of 1990, a 510(K) summary is provided as outlined in 21 CFR 807.92.

Submitter

Name: Alfa Scientific Designs, Inc.

Address: 11494 Sorrento Valley Road, Suite M
San Diego, CA 92121

Telephone: (858) 350-9798

Fax: (858) 350-9709

Email: asdi@worldnet.att.net

Device Name

Trade Name: *Instant-View™ Methadone Urine Dip Strip Test*

Common Name: Methadone Test

Classification: 21 CFR 862.3620, Class II

Predicate Device

QuickScreen One Step Methadone Screening Test, Phamatech, Inc., 510(K) #K982938.

Device Description

A one-step lateral flow chromatographic immunoassay.

Intended Use

This *Methadone Urine Test* is intended for use in drug rehabilitation clinics, physician offices and reference laboratories. It provides qualitative screening results for Methadone in human urine at a cutoff concentration of 300ng/ml.

Summary of the Similarities to the Predicate Device

- Intended Use:
Professional use only.
- Interpretation of results:
The appearance of only a C line, indicates a positive result, and the Methadone level in the urine specimen tested is at or higher than 300ng/ml; or, the appearance of two lines, both C line and T line, indicates a negative result, and the Methadone level in the urine specimen tested is below 300ng/ml.
- Technological Characteristics:
Both devices are one step, qualitative, competitive binding immunoassay test.

Accuracy studies

The accuracy of the device was evaluated with 100 human urine specimens with GC/MS data at three POL sites and one reference laboratory. The correlation between the results from GC/MS and that from the *Instant-View™ Methadone Urine Dip Strip Test* is 97%.

Discussion and Conclusion

- The correlation of results from the *Instant-View™ Methadone Urine Dip Strip Test* and that from GC/MS is 97%, indicating that the *Instant-View™ Methadone Urine Dip Strip Test* is substantially equivalent to the existing legally marketed product.
 - The test results from three physician's offices laboratories and one clinical laboratory conducted by personnel with diverse educational backgrounds and working experience agreed 99%, indicating that the *Instant-View™ Methadone Urine Dip Strip Test* is suitable for use by health care professionals with diverse educational backgrounds and work experiences.
-



MAY 29 2001

Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

Naishu Wang, MD, Ph.D.
President
Alfa Scientific Designs, Inc.
11494 Sorrento Valley Road
Suite M
San Diego, CA 92121

Re: 510(k) Number: K010014
Trade/Device Name: Instant-View™ Methadone (MTD) Urine Dip Strip Test
Regulation Number: 862.3620
Regulatory Class: II
Product Code: DJR
Dated: December 26, 2000
Received: January 2, 2001

Dear Dr. Wang:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Pre-market Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895.

A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,



Steven I. Gutman, M.D., M.B.A.
Director
Division of Clinical Laboratory Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

510 (K) NUMBER (IF KNOWN): K010014

DEVICE NAME: Instant-View™ Methadone (MTD) Urine Dip Strip Test

INDICATIONS FOR USE:

It is for health care professional, in-vitro diagnostic use only.

This test is a qualitative one step lateral flow immunoassay intended for use in drug rehabilitation clinics, physician offices and reference laboratories. It provides qualitative screening results for Methadone in human urine at a cutoff concentration of 300ng/ml.

This assay provides only a preliminary result. Clinical consideration and professional judgment must be applied to any drug of abuse test result, particularly in evaluating a preliminary positive result. In order to obtain a confirmed analytical result, a more specific alternate chemical method is needed. Gas Chromatography/Mass Spectroscopy (GC/MS) is the preferred confirmation method.

Fred Loey

(Division Sign-Off)
Division of Clinical Laboratory Devices
510(k) Number K010014

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED.)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

OR

Over-The-Counter Use (Optional Format 1-2-96)