

AUG 01 2002

K012693

	<b>510(k) Summary MISTRAL</b> rev1	1 of 6
---	---------------------------------------	--------

This 510(k) summary of safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

**1. Submitter's name, address, telephone number, contact person, and date summary prepared:**

a. ACUTRONIC Medical Systems AG  
Fabrik im Schiffli  
CH 8816 Hirzel  
Tel. 0041 1 729 70 80  
Fax. 0041 1 729 70 81  
e-mail: [p.cavelti@acutronic-medical.ch](mailto:p.cavelti@acutronic-medical.ch)

b. Contact Person: Pius Cavelti  
Director Sales / Marketing

c. Date Summary Prepared: 21.06.2002

**2. Name of device, including trade name and classification name:**

a. Trade/Proprietary Name: MISTRAL

b. Classification Name: Ventilator Continuous

**3. Identification of the predicate device or legally marketed device or devices to which substantial equivalence is being claimed:**

Company: ACUTRONIC Medical Systems AG  
Device: VS-200s Universal Jet Ventilator  
510(k) : K900745C  
Date Cleared: 04/29/91

**4. A description of the device that is the subject of the 510(k), including explanation of how the device functions, basic scientific concepts, significant physical and performance characteristics (design, material, physical properties):**

**Airway pressure measurement:**

One of the key features of the MISTRAL is the safety concept to avoid having excessive pressures occurring in the patient's airways. To reach this goal, a high pressure resistant sensor with excellent resolution in the working range of 1 to 200 mbar is used.

**Display:**

A medium sized screen displays the main parameters of ventilation and allows quick information of the anaesthetist about the status of the ventilation. A graphic indication of the airway pressure makes quantitative reading possible and allows quick access for parameter changes. Special functions are activated in menus and displayed on the screen.

0074

File: 510(k) Summary MISTRAL rev1	Date: 15.02.2001	Index: 01
-----------------------------------	------------------	-----------

**Air / Oxygen blender:**

The MISTRAL has a built in blender, which allows concentrations between 21 and 100 % oxygen. The access for changes is via a knob on the front panel. This parameter is accessible without going into menus, to allow a rapid change in case of emergency. The MISTRAL uses an oxygen sensor, which is calibrated each time the unit, is switched on. If the ventilator can't reach the selected FiO<sub>2</sub>, an alarm message appears on the screen, informing about a problem in the gas supply.

**Bypass flow:**

The MISTRAL has an additional outlet for gas, which is mechanically adjusted to 15 LPM and a concentration between 21 % and 100 %. This feature allows mask ventilation of the patient for induction in anaesthesia.

**Microprocessor and Memory:**

H8 532 (manufactured by Hitachi)

Clock frequency 10MHz

Memory devices:

27C512 (several types available)

64kB EPROM

DS1386 (manufactured by Dallas)

32kB RAM timekeeper external watchdog

**Sensors:**

200 cmH<sub>2</sub>O pp-pressure sensor (Honeywell)

monitoring pressure condition in Jet tube sampling rate depending on jet frequency max 2.5 samples per second

72.5 PSI driving-pressure sensor (Honeywell)

monitoring driving pressure sampling rate maximum 200 ms

**Fault-Detection:**

**PP sensor is periodically checked for disconnection, shortcut and defective membrane (period: depending on frequency setting) *failure results in error message***

Driving pressure sensor is checked for disconnection & shortcut  
(period: depending on frequency setting) *failure results in error message*

Oxygen sensor is periodically checked for disconnection, shortcut and capacity  
(period: depending on actual regulated difference capacity only in auto calibration mode (startup routine)) *failure results in error message*

Full time gas input detection *failure results in alarm message*

fully regulated oxygen concentration with timeout function *failure results in alarm message*

0075

## Alarm functions

### General information

In case of an alarm, the Mistral will display the alarm message on the screen and an acoustic alarm sounds. The acoustic alarm is ceased as soon as the alarm condition is cleared. However, the screen message remains until it is confirmed by pressing the alarm reset button. If you press alarm reset button while the alarm conditions is not cleared, there is an alarm mute for 1 minute activated, which depresses the acoustic alarm for this period. The screen inverses to make the user aware that the alarm mute is active.

**The Mistral has an alarm concept to allow a safe use of the equipment for patient and operator. The following section shows the different alarm types and their sources:**

### PP too high:

This alarm is activated if the pressure in the Jet line remains above the set limit for the PP pressure. This alarm is a safety feature for the use of Jet Ventilation with single lumen catheters. It allows a security shut-off of the Jet Ventilator in case of an airway obstruction.

## System alarms

### General information

These alarms signal the user a state of improper function of the device. There is no immediate danger to the patient but the device should be replaced and serviced by authorised personal (exception: low level of pressure input).

### Specification and characteristics of the *ventilation* alarm

#### Application:

The ventilation device makes sure that max. oxygen level in the Jet ventilator is below a level of concern. The ventilation is periodically checked for ether disconnect (alarm level 1) or overload (alarm level 2) of the fan. This is an alarm which can not be disabled.

### Specification and characteristics of the *FIO2 not adjustable* alarm

#### Application:

The FIO2 not adjustable alarm is evoked as soon as the O2 level can not be adjusted because of any other reason than missing ether O2 or airway pressure. There could be a mechanical problem or a system malfunction of the O2 regulation. If the O2 sensor could not be calibrated this alarm is also set active.

### Specification and characteristics of the low O2 pressure alarm

#### Application:

The O2 pressure is continuously watched by the controller and goes on if there is not enough pressure on the O2 inlet. This is an alarm which can not be disabled. The alarm resets itself if pressure is sufficient.

Specification and characteristics of the low air pressure alarm

Application:

The air pressure is continuously watched by the controller and goes on if there is not enough pressure on the air inlet. This is an alarm which can not be disabled. The alarm resets itself if pressure is sufficient.

Specification and characteristics of low bat alarm

Application:

Battery level is measured and in case of a worn out battery the message "low bat" is displayed instead of the time and date. There is no alarm going off because the battery is not used unless the power failure fault. This battery has not a data backup function. In case of "low bat" a replacement of the battery is sufficient.

Specification and characteristics of Power failure alarm

Application:

Main line is less than 95 VAC the beeper on the back panel sounds with interval supplied by battery power. This is an alarm, which can not be disabled.

Specification and characteristics of MAINTENANCE REQUIRED alarm

Application:

Mistral has an integrated timer function for servicing intervals. After one year, on the display the message "MAINTENANCE REQUIRED" appears. The unit may still be used, however, a service engineer must be called for servicing of the equipment to assure a proper function.

**Error functions**

Overview

The error message signals a system malfunction. The error messages can be reset by the alarm button but if the malfunction is still detected the message immediately returns to the screen. In any of these cases the Monsoon needs a service.

Type of error	characteristic
PP sensor	the PP sensor is out of range (probably disconnected or defect)
DP sensor	the DP sensor is out of range (probably disconnected or defect)
O2 sensor is weak	O2 signal is getting to weak replace the sensor

**0077**

**5. Statement of intended use:**

The MISTRAL Universal Jet Ventilator is designed to use for short term application in bronchoscopy and laryngoscopy.

**Indication for Use**

Jet ventilation applied with the Mistral Universal Jet Ventilator is useful in airway surgery, as it is performed in thoracic surgery units and ENT surgery. Jet ventilation is the optimal ventilation technique during the application of LASER light, where the presence of an ETT bears the risk for ignition, airway fire and burn injuries. Jet ventilation is useful for the removal of foreign bodies from the airway (e.g. after accidental aspiration of foreign bodies) via rigid bronchoscopes, for supplemental oxygenation during lung surgery (when the operated lung must be recruited for additional oxygenation, but may not be ventilated conventionally for surgical reasons), for surgery of the lower trachea close to the carina (when during the reconstruction phase, no tight airway sealing can be achieved), for radiation therapy of lung metastasis (when the tidal movements of the chest during spontaneous respiration or conventional ventilation would preclude the focusing of the radiation beam onto the target).

0078

**6. Statement of how the technological characteristics of the device compare to those of the predicate or legally marketed device.**

Technical feature	VS-200S	MISTRAL	Comment
EEP Limit End Expiratory Pressure Limit	Yes, measured in between the pulses of Jet insufflations. Valve opens during pause to direct pressure on pressure sensor	Yes, measured with a high pressure resistant pressure sensor. Resolution 1 mbar	Due to availability of high pressure resistant sensor, quicker response of the pressure limit = increased patient safety
Display of EEP pressure level	No	Yes, displayed after each jet insufflation	Allows trend information on EEP pressure.
Bypass flow	No	Yes, adjusted at factory upon customer request. Allows driving of nebulizer or for continues flow for anaesthesia induction	No need for conventional anaesthesia machine for induction
Built-in air-oxygen blender	No. Use of Bird blender	Yes, electronic blender	Compact size without need for external tubings
Oxygen measurement	No. Need of separate oxygen monitor	Yes, with automatic calibration of sensor during start-up procedure	Ease of use.
Humidification	Yes, injection of distilled water into jet nozzle	No	Experiences in the past have shown, that injection of water may create some lesions in the trachea if drop size is too big.
In-service information	No	Yes, message appears on display when maintenance is needed	Facilitates the in hospital control of equipment servicing conditions
Flexible software for future improvements	No, digital electronic, no microprocessor	Yes, allows use of different language for operator without change of front label.	Menu driven operators interface
Inlet pressure controlled by servo valve	No	Yes, if pressure limit is exceeded for more then 3 s, the inlet pressure valves close for patient safety	In case of a jet valve failure, the inlet pressure is shut off automatically to avoid overpressure in the patients lung



AUG 01 2002

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Acutronic Medical Systems AG  
C/O Terry Torzala, R.A.C.  
Official Correspondent  
Medical Equipment Development, Incorporated  
P.O. Box 85820  
Tucson, Arizona 85754-5820

Re: K012693  
Trade/Device Name: Mistral Critical Care Jet Ventilator and Accessories  
Regulation Number: 21 CFR 868.5895  
Regulation Name: Continuous Ventilator  
Regulatory Class: II  
Product Code: CBK  
Dated: June 26, 2002  
Received: June 26, 2002

Dear Mr. Torzala:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Page 2 – Mr. Torzala

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4646. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely,



Timothy A. Ulatowski  
Director

Division of Anesthesiology, General Hospital,  
Infection Control and Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

INDICATIONS FOR USE STATEMENT

510(k) Number (if known): K012693


Device Name: Mistral Universal Jet Ventilator

Indications For Use:

The Mistral Universal Jet Ventilator is intended for use during bronchoscopy or laryngoscopy for applications up to 45 minutes duration on adults in a hospital or other clinical setting.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)  
Division of Dental, Infection Control,  
and General Hospital Devices

510(k) Number K012693

Prescription Use   
(Per 21 CFR 801.109)

OR

Over-The-Counter Use