

NOV 07 2001

ellman Surgitron 120 IEC

510 (k) Summary

K013255

**1. Submitter name and address:**

Frank Lin, Ph.D.  
Director of R&D engineering  
ellman international  
1135 Railroad Avenue  
Hewlett, New York 11557

**2. Device name and classification:**

- 2.1 Device Name: Surgitron 120 IEC (Also known as Surgitron 4.0 Dual RF)  
2.2 Classification: Class 2 device, 21 CFR 878.4400

**3. Description of the device:**

The ellman Surgitron 120 IEC enhanced capability Electrosurgery Generator described herein is a compact source of high power RF energy to be employed for a variety of radiosurgery procedures. This action is achieved by front panel selection of waveforms and power level. All selection is effected through push buttons and lamps which give the operator feedback of status. Power level for each mode is indicated by front panel digital displays which also show the status of self-test and monitoring. This display is interlocked with the controls to prevent operation when FAIL is displayed. The final output power control is made through foot and/or hand switches. Both Monopolar and Bipolar electrodes are provided. It is designed to comply with international safety standards.

**4. The intended use/indication for use of the device:**

**4.1 Cutting**

Snoring, Submucosal palatal shrinkage, traditional uvulopalatoplasty (RAUP), myringotomy with effective hemorrhage control, epistaxis treatment, and turbinate shrinkage, Skin Incisions, Biopsy, Cysts, Abscesses, Tumors, Cosmetic Repairs, Development of Skin Flaps, Skin Tags, Blepharoplasty,

**4.2 Blended Cutting and Coagulation**

Snoring, Submucosal palatal shrinkage, traditional uvulopalatoplasty (RAUP), myringotomy with effective hemorrhage control, epistaxis treatment, and turbinate shrinkage, Skin Tags, Papilloma Keloids, Keratosis, Verrucae, Basal Cell Carcinoma, Nevi, Fistulas, Epithelioma, Cosmetic Repairs, Cysts, Abscesses, Development of Skin Flaps,

**4.3 Hemostasis**

Control of Bleeding, Epilation, Telangiectasia

#### 4.4 Fulguration

Basal Cell Carcinoma, Papilloma, Cyst Destruction, Tumors, Verrucae, Hemostasis.

#### 4.5 Bipolar

Pinpoint, Precise Coagulation, Pinpoint Hemostasis in any field (Wet or Dry), Snoring, Submucosal palatal shrinkage, traditional uvulopalatoplasty (RAUP), myringotomy with effective hemorrhage control, epistaxis treatment, and turbinate shrinkage.

### 5. Identification to predicate devices

5.1. Sugitron IEC II with general use indication K001253

5.2 ERBOTOM ICC 200 with general use indication K933157

### 6. Summary of the technological characteristics of the new device in comparison to the predicate devices.

Substantial equivalence  
In Technological Characteristics comparison

FEATURE	<b>ollman</b> SURGITRON 120 IEC (New Application Device)	<b>ollman</b> SURGITRON IEC II WITH GENERAL USE INDICATION K001253 PREDICATE	<b>ERBOTOM</b> ICC200 WITH INDICATION K933157 PREDICATE
Indications For Use	Snoring, submucosal palatal shrinkage, traditional uvulopalatoplasty (RAUP), myringotomy with effective hemorrhage control, epistaxis treat, and turbinate shrinkage. Also, see page 5 for detail.	Same As New Device	Same As New Device
Design Specification	UL544 and IEC601-2-2	Same As New Device	Same As New Device
Output Energy	120 Watt	100 Watt	200 Watt
Output Waveform (s)	4.0 MHz Sin-wave CW, Fully Rectified, Partially Rectified, and 1.7 MHz for Fulgurating Spark-Gap	Same As New Device	350KHz Sine-Shaped and 1MHz Pulse-modulated
Standards Met	UL2601 and IEC 601-1 601-2-2, BS15724:2.2	UL2601 and IEC 601-1 601-2-2, BS15724:2.2	IEC 601-1
Delivery system and configuration	Monopolar and Bipolar	Monopolar and Bipolar	Same As New Device
Biocompatibility Test	Electrodes identical to predicate device	Same As New Device	Same As New Device
Sterilization Method(s)	Refer to page 6, Note I	Same As New Device	Not Indicated



DEC 20 2001

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Mr. Frank Lin  
Director of Engineering  
Research and Development Department  
Ellman International, Inc.  
1135 Railroad Avenue  
Hewlett, New York 11557-2316

Re: K013255  
Trade/Device Name: Surgitron 120 IEC  
Regulation Number: 878.4400  
Regulation Name: Electrosurgical cutting and coagulation device and accessories  
Regulatory Class: II  
Product Code: GEI  
Dated: September 27, 2001  
Received: September 28, 2001

Dear Mr. Lin:

This letter corrects our substantially equivalent letter of November 7, 2001 regarding the file number.

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 - Mr. Frank Lin

This letter will allow you to continue marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at their toll free number (800) 638-2041 or at (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Celia M. Witten". The signature is written in a cursive style with a large initial "C" and "W".

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative  
and Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

NOV 07 2001

510K Notification  
Surgitron - General & Plastic Surgery Use

page 5

ollman international

510(k) Number (if known): K013255

Device Name: SURGITRON 120 IFC

Indication For Use: is identical to the Surgitron as a preamendment device such as:

- \* Cutting  
Snoring, Submucosal palatal shrinkage, traditional uvulopalatoplasty (RAUP), myringotomy with effective hemorrhage control, epistaxis treatment, and turbinate shrinkage, Skin Incisions, Biopsy, Cysts, Abscesses, Tumors, Cosmetic Repairs, Development of Skin Flaps, Skin Tags, Blepharoplasty,
- \* Blended Cutting and Coagulation  
Snoring, Submucosal palatal shrinkage, traditional uvulopalatoplasty (RAUP), myringotomy with effective hemorrhage control, epistaxis treatment, and turbinate shrinkage, Skin Tags, Papilloma Keloids, Keratosis, Verrucae, Basal Cell Carcinoma, Nevi, Fistulas, Epithelioma, Cosmetic Repairs, Cysts, Abscesses, Development of Skin Flaps,
- \* Hemostasis  
Control of Bleeding, Epilation, Telangiectasia
- \* Fulguration  
Basal Cell Carcinoma, Papilloma, Cyst Destruction, Tumors, Verrucae, Hemostasis.
- \* Bipolar  
Pinpoint, Precise Coagulation, Pinpoint Hemostasis in any field (Wet or Dry), Snoring, Submucosal palatal shrinkage, traditional uvulopalatoplasty (RAUP), myringotomy with effective hemorrhage control, epistaxis treatment, and turbinate shrinkage

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Juan Walker  
(Division Sign-Off)  
Division of General, Restorative  
and Neurological Devices

Prescription Use   
(Per 21 CFR 801.109)

Over-The-Counter Use

(Optional Format 1-2-96)

510(k) Number   
K013255