

**Allergen ImmunoCAP™  
510(k) Submission  
Section 7. Summary of Safety and Effectiveness**

---

**7. SUMMARY OF SAFETY AND EFFECTIVENESS**

This summary of safety and effectiveness information is being submitted in accordance with the requirements of The Safety Medical Devices Act of 1990 (SMDA 1990) and 21 CFR Part 807.92.

**Premarket Notification 510(k) Number:** K993388

**Date of Summary Preparation:** October 5, 1999

**Distributor:** Pharmacia & Upjohn  
Diagnostics Division, US Operation  
7425-248-1  
7000 Portage Road  
Kalamazoo, MI 49001

**Manufacturer:** Pharmacia & Upjohn, Diagnostics AB  
S-751 82 Uppsala, Sweden  
and  
MIAB  
Dragarbrunnsgatan 65  
S-75320 Uppsala

**Company Contact Person:** Karen Matis  
Manager, Regulatory Affairs and Quality Management  
Diagnostics Division  
US Operation  
7000 Portage Road  
7425-248-01  
Kalamazoo, MI 49001  
(614) 794-3324 (Phone)  
(614) 794-0266 (Fax)

**Device Name:** Allergen ImmunoCAP™- c5, c6, c73, e201, e213, f50, f51, f54, f55, f56, f57, f58, f59, f60, f61, f205, f213, f218, f225, f231, f232, f233, f244, f254, f299, fx7, fx13, fx14, fx15, fx18, fx20, g201, gx3, gx4, i8, i73, i201, m207, o1, p4, t208, tx6, tx9, wx3

**Allergen ImmunoCAP™**  
**510(k) Submission**  
**Section 7. Summary of Safety and Effectiveness**

---

**Common Name:** Solid phase components of immunological test system to measure allergen specific IgE antibodies.

**Classification:**

<b>Product Name</b>	<b>Product Code</b>	<b>Class</b>	<b>CFR</b>
Allergen ImmunoCAP™ c5, c6, c73, e201, e213, f50, f51, f54, f55, f56, f57, f58, f59, f60, f61, f205, f213, f218, f225, f231, f232, f233, f244, f254, f299, fx7, fx13, fx14, fx15, fx18, fx20, g201, gx3, gx4, i8, i73, i201, m207, o1, p4, t208, tx6, tx9, wx3	82DHB	II	866.5750

**Predicate Test Systems For The Measurement of Specific IgE**

Pharmacia CAP System™ RAST FEIA	K894190, K911903
UniCAP™ Specific IgE Assay	K962274

**Intended Use Statement :**

**Allergen ImmunoCAP™** is the solid phase component of the Pharmacia & Upjohn *in vitro* immunodiagnostic systems which measure specific IgE to the respective allergen bound to the ImmunoCAP™. Allergen ImmunoCAP™ are intended to be used with Pharmacia CAP System™ RAST FEIA and UniCAP® Specific IgE *in vitro* diagnostic assays.

**Pharmacia CAP System RAST® FEIA and UniCAP® Specific IgE** are intended for *in vitro* diagnostic use as an aid in the clinical diagnosis of IgE mediated allergic disorders in conjunction with other findings, and are to be used in clinical laboratories, as well as, physician office laboratories.

**Allergen ImmunoCAP™**  
**510(k) Submission**  
**Section 7. Summary of Safety and Effectiveness**

---

**General Description**

**Allergen ImmunoCAP™**

Allergen ImmunoCAP™ consists of a cellulose sponge matrix to which allergenic components are covalently coupled. The matrix is encased in a small round plastic capsule. This capsule is at the same time a holder of the matrix for convenient automation and a reaction chamber.

The sponge matrix is manufactured from activated cellulose derivative to which allergen extract solution is added under defined optimized conditions for the allergen coupling. This solid phase is an excellent carrier of allergens and provides favorable reaction conditions.

**UniCAP®/Pharmacia CAP System™ RAST FEIA Specific IgE Test Principle**

The allergen of interest, covalently coupled to ImmunoCAP, reacts with the specific IgE in the patient serum specimen. After washing away non-specific IgE, enzyme labeled antibodies against IgE are added to form a complex. After incubation, unbound enzyme-anti-IgE is washed away and the bound complex is then incubated with a developing agent. After stopping the reaction, the fluorescence of the eluate is measured. The higher the response value, the more specific IgE is present in the specimen. To evaluate the test results, the response for the patient samples is compared directly to the response for the calibrators.

**Performance Characteristics Of Allergen ImmunoCAP™**

The safety and effectiveness of the test systems Pharmacia CAP System™ RAST FEIA and UniCAP™ Specific IgE for the determination of specific IgE antibodies have been established in previous 510(k) submissions. This 510(k) submission includes data to add 44 additional Allergen ImmunoCAP™ to the Pharmacia CAP System™ and UniCAP™ test systems for the measurement of specific IgE.

RAST inhibition verifies the immunological specificity of IgE binding for each allergen. The function of Allergen ImmunoCAP™ is further verified by testing clinical serum samples, with a history or indication of allergy to the specific allergen, and established negative samples. The analysis was performed in both Pharmacia CAP System™ and UniCAP™ test systems and results show an outstanding agreement of outcome concerning positive and negative samples in both systems.

The importance of each allergen is demonstrated with relevant literature references covering frequency, clinical use and description of related allergens. Reproducibility between production lots and stability studies complete the picture by showing the constant quality of Allergen ImmunoCAP™.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

DEC - 1 1999

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

Ms. Karen E. Matis  
Manager, Regulatory Affairs and Quality Management  
Pharmacia & Upjohn  
7000 Portage Road  
Kalamazoo, Michigan 49001-0199

Re: K993388  
Trade Name: Allergen ImmunoCAP™  
Regulatory Class: II  
Product Code: DHB  
Dated: October 8, 1999  
Received: October 8, 1999

Dear Ms. Matis:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.


Page 2

Under the Clinical Laboratory Improvement Amendments of 1988 (CLIA-88), this device may require a CLIA complexity categorization. To determine if it does, you should contact the Centers for Disease Control and Prevention (CDC) at (770)488-7655.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll free number (800) 638-2041 or at (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>"

Sincerely yours,



Steven I. Gutman, M.D., M.B.A.  
Director  
Division of Clinical Laboratory Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

**Allergen ImmunoCAP™**  
**510(k) Submission**  
**Section 1. Indications For Use Statement**

510(k) Number K993388  
 Device Name: Allergen ImmunoCAP™

Product number	Code	Name
14-4450-01	c5	Ampicilloyl
14-4451-01	c6	Amoxicilloyl
14-4347-01	c73	Insulin human
14-4824-01	e201	Canary bird feathers
14-4832-01	e213	Parrot feathers
14-4454-01	f50	Chub mackerel
14-4455-01	f51	Bamboo shoot
14-4458-01	f54	Sweet potato
14-4459-01	f55	Common millet
14-4460-01	f56	Foxtail millet
14-4461-01	f57	Japanese millet
14-4462-01	f58	Pacific squid
14-4463-01	f59	Octopus
14-4464-01	f60	Jack mackerel, Scad
14-4465-01	f61	Sardine, Pilchard
14-4837-01	f205	Herring
14-4838-01	f213	Rabbit
14-4816-01	f218	Paprika, Sweet pepper
14-4433-01	f225	Pumpkin
14-4806-01	f231	Milk, boiled
14-4804-01	f232	Ovalbumin
14-4805-01	f233	Ovomucoid
14-4839-01	f244	Cucumber
14-4434-01	f254	Plaice
14-4843-01	f299	Sweet chestnut
14-4844-01	fx7	Food Mix (f25, f45, f47, f48, f85)
14-4848-01	fx13	Food Mix (f12, f15, f31, f35)
14-4849-01	fx14	Food Mix (f25, f214, f216, f218)
14-4850-01	fx15	Food Mix (f33, f49, f92, f95)
14-4852-01	fx18	Food Mix (f12, f13, f14)
14-4853-01	fx20	Food Mix (f4, f5, f6, f9)
14-4821-01	g201	Barley
14-4193-01	gx3	Grass Pollen Mix (g1, g5, g6, g12, g13)
14-4194-01	gx4	Grass Pollen Mix (g1, g5, g7, g12, g13)
14-4474-01	i8	Moth
14-4328-01	i73	Blood worm
14-4833-01	i201	Horse bot fly
14-4831-01	m207	Aspergillus niger
14-4476-01	o1	Cotton, crude fibers
14-4475-01	p4	Anisakis*
14-4823-01	t208	Linden
14-4202-01	tx6	Tree Pollen Mix (t1, t3, t5, t7, t10)
14-4274-01	tx9	Tree Pollen Mix (t2, t3, t4, t7, t12)
14-4196-01	wx3	Weed Pollen Mix (w6, w9, w10, w12, w20)

**Allergen ImmunoCAP™** is the solid phase component of the Pharmacia & Upjohn in vitro immunodiagnostic systems which measure specific IgE to the respective allergen bound to the ImmunoCAP. Allergen ImmunoCAP are intended to be used with Pharmacia CAP System RAST FEIA and UniCAP Specific IgE in vitro diagnostic assays.

**Pharmacia CAP System RAST® FEIA and UniCAP® Specific IgE** are intended for in vitro diagnostic use as an aid in the clinical diagnosis of IgE mediated allergic disorders in conjunction with other findings, and are to be used in clinical laboratories, as well as, physician office laboratories.

Concurrence of CDRH, Office of Device Evaluation (ODE)

*Peter E. Maffei*

(Division Sign-Off)

Division of Clinical Laboratory Devices

510(k) Number

K993388

Prescription Use

OR

Over-The-Counter Use