

K100152

**510(k) Summary**

OCT 22 2010

for

**Sirona Dental Systems**

**Sirona Dental CAD/CAM System**

**1 Sponsor**

Sirona Dental Systems GmbH

Fabrikstrasse 31

D-64625 Bensheim

Germany

Contact Person: Fritz Kolle

Telephone: 49 6251 16 32 94

Date Prepared: October 20, 2010

**2 Device Name**

Proprietary Name: Sirona Dental CAD/CAM-System

Common/Usual Name: Abutment, implant, dental, endosseous

Classification Names: Endosseous dental implant abutment

**3 Predicate Devices**

Replace® NP, K091756, Brånemark®, K091756, Tissue level NN, K081005, OsseoSpeed™, K081666, Frialit® / Xive®, K032158, Osseotite K072642, Tapered Screw-Vent®, K060880.

Zirconium oxide sleeve for Camlog titanium bases (K032448).

Lava Software (K062493).

#### 4 Intended Use

---

The Sirona Dental CAD/CAM System is intended for use in partially or fully edentulous mandibles and maxillae in support of single or multiple-unit cement retained restorations. The system consists of three major parts: TiBase, InCoris mesostructure, and CAD/CAM software. Specifically, the InCoris mesostructure and TiBase components make up a two-piece abutment which is used in conjunction with endosseous dental implants to restore the function and aesthetics in the oral cavity. The InCoris mesostructure may also be used in conjunction with the Camlog Titanium base CAD/CAM (types K2244.xxxx) (K083496) in the Camlog Implant System. The CAD/CAM software is intended to design and fabricate the InCoris mesostructure. The InCoris mesostructure and TiBase two-piece abutment is compatible with the following implants systems:

- *Nobel Biocare Replace (K020646)*
- *Nobel Biocare Branemark (K022562)*
- *Friadent Xive (K013867)*
- *Biomet 3i Osseotite (K980549)*
- *Astra Tech Osseospeed (K091239)*
- *Zimmer Tapered Screw-Vent (K061410)*
- *Straumann SynOcta (K061176)*

#### 5 Device Description

The Sirona Dental CAD/CAM-System takes optical impressions and records the topographical characteristics of teeth, dental impressions, or stone models. Dental restorative prosthetic devices are manufactured using computer aided design and fabrication. The system also features the processing of mesostructures, a dental restorative prosthetic device used in conjunction with endosseous dental implant abutments.

The system that features the processing of mesostructures comprises

- Titanium bases TiBase and Camlog
- inCoris ZI meso blocks
- Sirona Dental CAD/CAM Design and fabricating devices

Titanium bases are used as an implant prosthetic titanium base for adhesion to mesostructures to restore function and aesthetics in the oral cavity.

inCoris ZI meso blocks are used in manufacturing individually designed mesostructures, which are glued to a fitting titanium base after milling and sintering.

---

Sirona Dental CAD/CAM design and fabricating devices feature the processing of mesostructures, a dental restorative prosthetic device used in conjunction with endosseous dental implant abutments, i.e. it is an accessory to it. This component consists of the devices CEREC3, CEREC AC, inEos, inEos Blue, CEREC MCXL and inLab MCXL.

## **5.1 TiBase**

### **5.1.1 Device Function**

The Sirona TiBase is a premanufactured prosthetic component directly connected to dedicated endosseous dental implants and is intended for use as an aid in prosthetic rehabilitation.

The Sirona offering consists of the titanium base TiBase, the abutment screw and the scanbody. The parts are marketed non-steril and for single use only.

The Sirona TiBase is bonded to an individually designed mesostructure, a ceramic prosthetic/restoration, that supports the final restoration.. The mesostructure is milled from an inCoris ZI meso block with Sirona CAD/CAM milling machines CEREC or inLab, and sintered afterwards.

The two piece abutment is mounted onto the implant and fixed with a screw.

The scope of delivery contains a scanbody (ABS plastic) which is mounted on a TiBase in order to acquire the topographical surface of the area where the endosseous dental implant abutment is located with Sirona Dental CAD/CAM fabricating devices. From the acquired data the position of the implant can be calculated. After an optical impression has been taken the scanbody is removed.

Sirona TiBase devices are compatible with following systems (Table 1):

**Table 1: Sirona TiBase Devices Compatibility**

Sirona TiBase	Compatible System		
	Manufacturer	System	Diameter
NBRS 3.5	Nobel Biocare	Replace® NP	3,5 mm
NBRS 4.3		Replace® RP	4.3 mm
NBRS 5.0		Replace® WP	5.0 mm
NBRS 6.0		Replace® 6.0	6.0 mm
NBB 3.4	Nobel Biocare	Brånemark®	3.4 mm
NBB 4.1		Brånemark®	4.1 mm
SSO 3.5	Straumann	Tissue level NN	3.5 mm
SSO 4.8		Tissue level RN	4.8 mm
SSO 6.5		Tissue level WN	6.5 mm
ATOS 3.5/4.0	Astra Tech	OsseoSpeed™	3.5 S / 4.0 S mm
ATOS 4.5/5.0		OsseoSpeed™	4.5 / 5.0 mm
FX 3.4	Friadent	Frialit® / Xive®	3.4 mm
FX 3.8		Frialit® / Xive®	3.8 mm
FX 4.5		Frialit® / Xive®	4.5 mm
FX 5.5		Frialit® / Xive®	5.5 mm
BO 3.4	Biomet 3i	Osseotite (Connec-tion type: Ex. Hex)	3.4 mm
BO 4.1		Osseotite (Connec-tion type: Ex. Hex)	4.1 mm
BO 5.0		Osseotite (Connec-tion type:	5.0 mm

Sirona TiBase	Compatible System		
	Manufacturer	System	Diameter
		Ex. Hex)	
ZTSV 3.5	Zimmer	Tapered Screw-Vent®	3.5 mm
ZTSV 4.5		Tapered Screw-Vent®	4.5 mm
ZTSV 5.7		Tapered Screw-Vent®	5.7 mm

### 5.1.2 Scientific Concept

The underlying scientific concept is the use of an already introduced technology of a titanium base abutment combined with individually CAD/CAM fabricated ceramic prosthetics.

### 5.1.3 Physical and Performance Characteristics

#### 5.1.3.1 Design

The TiBase devices have various diameters, are compatible with dedicated implant systems, and fit to compatible implants as provided in Table 2.

**Table 2: Implant Compatibility**

<b>TiBase</b>	<b>Implant- Manufacturer</b>	<b>Implant-System</b>	<b>510(k) Implant</b>
NBRS	Nobel Biocare	Replace	K020646
NBB	Nobel Biocare	Branemark	K022562
FX	Friadent	Xive	K013867
BO	Biomet 3i	Osseotite	K980549
ATOS	Astra Tech	OsseoSpeed	K091239
ZTSV	Zimmer	Tapered Screw-Vent	K061410
SSO	Straumann	SynOcta	K061176

#### **5.1.3.2 Material Used**

TiBase and abutment screw are made of Ti6Al4V.

#### **5.1.3.3 Physical Properties**

TiBase material composition and mechanical properties comply with ISO 5832-3:1996, Implants for surgery -- Metallic materials -- Part 3: Wrought titanium 6-aluminium 4-vanadium alloy

Sirona TiBase devices are compatible with systems listed in Table 1.

## **5.2 inCoris ZI meso**

### **5.2.1 Device Description**

The inCoris ZI meso offerings are blocks of various sizes from which individual dental mesostructures are grinded by milling machines (inLab MCXL, CEREC MCXL). The mesostructure is a part of a 2 part endosseous dental implant abutment which comprises a titanium base and a zirconium oxide mesostructure. The connection geometries are prefabricated.

## 5.2.2 Scientific Concept

The underlying scientific concept is the use of an already introduced technology of a titanium base abutment combined with individually CAD/CAM fabricated ceramic prosthetics.

## 5.2.3 Physical and Performance Characteristics

### 5.2.3.1 Design

The inCoris ZI meso are blocks of various sizes. The marketed ceramic is pre-sintered. One side of a block is mounted to a mandrel that will be inserted in the spindle's clamping chuck of the grinding machine. The connection geometry to titanium bases is prefabricated, i.e. already included in the shipped block. Connection geometries fit on Camlog (type K2244.xxxx) and Sirona (Tibase) titanium bases (Table 3 and Table 4).

**Table 3: Sirona inCoris ZI meso - TiBase Devices Compatibility**

Titanium Base			Ceramic Block		
TiBase (Sirona)	REF	Diameter	inCoris ZI meso (Sirona)	REF	Color inCoris ZI meso (Sirona)
NBRS 3.5	6282474	3,5 mm	inCoris ZI meso L	62 31 810	F0.5
				62 31 836	F2
NBRS 4.3	6282482	4.3 mm	inCoris ZI meso L	62 31 810	F0.5
				62 31 836	F2
NBRS 5.0	6282490	5.0 mm	inCoris ZI meso L	62 31 810	F0.5
				62 31 836	F2
NBRS 6.0	6282508	6.0 mm	inCoris ZI meso L	62 31 810	F0.5
				62 31 836	F2
NBB 3.4	6282516	3.4 mm	inCoris ZI meso L	62 31 810	F0.5
				62 31 836	F2
NBB 4.1	6282524	4.1 mm	inCoris ZI	62 31 810	F0.5

Titanium Base			Ceramic Block		
TiBase (Sirona)	REF	Diameter	inCoris ZI meso (Sirona)	REF	Color inCoris ZI meso (Sirona)
			meso L	62 31 836	F2
SSO 3.5	6284231	3.5 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
SSO 4.8	6284249	4.8 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
SSO 6.5	6284256	6.5 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
ATOS 3.5/4.0	6282532	3.5 S / 4.0 S mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
ATOS 4.5/5.0	6282540	4.5 / 5.0 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
FX 3.4	6282433	3.4 mm	inCoris ZI meso S	62 31 802 62 31 828	F0.5 F2
FX 3.8	6282441	3.8 mm	inCoris ZI meso S	62 31 802 62 31 828	F0.5 F2
FX 4.5	6282458	4.5 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
FX 5.5	6282466	5.5 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
BO 3.4	6282557	3.4 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
BO 4.1	6282565	4.1 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2

Titanium Base			Ceramic Block		
TiBase (Sirona)	REF	Diameter	inCoris ZI meso (Sirona)	REF	Color inCoris ZI meso (Sirona)
BO 5.0	6282573	5.0 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
ZTSV 3.5	6282581	3.5 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
ZTSV 4.5	6282599	4.5 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2
ZTSV 5.7	6282607	5.7 mm	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2

**Table 4: Sirona inCoris ZI meso - Camlog Devices Compatibility**

Titanium Base			Ceramic Block		
Camlog	REF	Dia- meter	inCoris ZI meso (Sirona)	REF	Color inCoris ZI meso (Sirona)
K2244.3348	K2244.3348	3.3	inCoris ZI meso S	62 31 802 62 31 828	F0.5 F2
K2244.3848	K2244.3848	3.8	inCoris ZI meso S	62 31 802 62 31 828	F0.5 F2
K2244.4348	K2244.4348	4.3	inCoris ZI meso S	62 31 802 62 31 828	F0.5 F2
K2244.5048	K2244.5048	5.0	inCoris ZI meso L	62 31 810 62 31 836	F0.5 F2

Titanium Base			Ceramic Block		
Camlog	REF	Dia- meter	inCoris ZI meso  (Sirona)	REF	Color inCoris ZI meso  (Sirona)
K2244.6048	K2244.6048	6.0	inCoris ZI meso L	62 31 810  62 31 836	F0.5  F2

### 5.2.3.2 Material Used

The inCoris ZI meso are pre-sintered zirconium oxide ceramic blocks. The metal block holder is made of aluminum. The material is composed of:

ZrO<sub>2</sub>+HfO<sub>2</sub>+Y<sub>2</sub>O<sub>3</sub> > 99.0%  
 Al<sub>2</sub>O<sub>3</sub> < 0.5%  
 Other oxides < 0.5%

### 5.2.3.3 Physical Properties

The final technical data of inCoris ZI meso are (after final sintering):

Density: 6.06 g/cm<sup>3</sup>  
 Coefficient of thermal expansion (CTE): 11.0\*10<sup>-6</sup> K<sup>-1</sup>  
 Flexural strength: > 900 MPa  
 Fracture toughness (KIC): 5.9 MPa·m<sup>1/2</sup>

## 5.3 Sirona Dental CAD/CAM Design and fabrication Devices

### 5.3.1 Device Description

The Sirona Dental CAD/CAM Design and fabricating devices for processing mesostructures perform and consist of

- Optical acquisition or recording of the topographical characteristics of dental impressions, or stone models using the devices Acquisition unit CEREC 3, CEREC AC and stationary scanning system inEos Blue

- 
- Design of mesostructures and processing the acquired or recorded data for this purposes using Sirona Dental CAD/CAM Software which runs on a CEREC 3, CEREC AC or PC. Design is performed by a dentist or dental technician
  - Milling of the mesostructure using CEREC MCXL or inLab MCXL milling machines from ceramic blocks intended for dental restorations and mesostructures

This Sirona Dental CAD/CAM Design and fabricating devices also process other dental restorations like crowns, bridge-frameworks, inlays, onlays and regulated under 21 CFR 872.3661, Optical Impression Systems for CAD/CAM, for such intended use.

### **5.3.2 Scientific Concept**

The underlying scientific concept is the use of CAD/CAM technology for the optical acquisition or recording of the topographical characteristics of dental impressions, or stone models, the design of individual mesostructures using recorded data (CAD), and eventually fabricating (milling) these designed mesostructures (CAM).

### **5.3.3 Physical and Performance Characteristics**

#### **5.3.3.1 Design**

Acquisition unit: the device consists of a camera for taking topographical characteristics of dental impressions optically. The recorded data are used for the design of individual mesostructures using specific CAD techniques for the dental field.

Fabricating devices: the devices mill the individual designed mesostructures from incoris ZI meso blocks. For this purpose the chucked block and the milling tools move according to the prescribed trajectories respectively the shape which is intended to be milled.

#### **5.3.3.2 Materials Used**

Not applicable.

### 5.3.3.3 Physical Properties

---

Not applicable.

## 6 Summary of the technological characteristics

### 6.1 TiBase

All proposed and predicates titanium bases and screws are made of Ti6Al4V, medical grade 5. Connection interfaces to the implants are identical for each individual diameter and connection type. Connection interfaces for dental restorations differ in that proposed devices have a notch in addition.

An extensive list is provided in Table 5.

### 6.2 inCoris ZI meso

Proposed and predicates device are to be bonded to titanium bases for supporting further dental restorations.

InCoris ZI meso material properties are practically identical to the properties available for the predicate device and are made of zirconium oxide and have practically the same composition. Because the composition of inCoris ZI meso blocs fulfills the standard ISO 13356:1997, "Implants for surgery, Ceramic materials based on yttria-stabilized tetragonal zirconia (Y-TZP)" they are classified as biocompatible by this standard..

The devices share the same thickness restrictions to be maintained.

Shape and bonding surface of connection interfaces to Camlog and Sirona Tibase are similar for Sirona inCoris ZI meso blocs and Camlog ceramic sleeves.

Same bonding material is used for proposed and predicate device.

An extensive list is provided in Table 6.

### 6.3 Sirona Dental CAD/CAM Design and fabrication Devices

Both, Sirona Dental CAD/CAM Design and fabricating devices and the LAVA system, uses their own devices for taking optical impressions to record topographical characteristics of teeth, dental impressions, or stone models for use in the computer aided design and fabricates dental restorative prosthetic devices used

---

in conjunction with endosseous dental implant abutments, i.e. it is an accessory to it. Both systems have the feature to transfer data of the optical impression to a remote milling machine via internet or exportation/importation of milling data.

An extensive list is provided in Table 7.

Table 5: Comparison of Sirona TiBase to Predicate Devices

Proposed Device TiBase	Predicate Device					Intended Use Differences	Abutment and Screw made of Ti6Al4V	Identical connection geometry to abutments	Connection geometry Type	Connection geometry Anti-rotational features	Screw geometry
	Manufacturer	System	Diameter	Titanium Base	Predicate Devices K-Number						
NBRS 3.5	Nobel Biocare	Replace® NP	3.5 mm	Nobel Biocare product catalog page 14. Product-No. 32376	K091756	no angulation correction to the implant axis	yes	yes	Internal 3 te-nons	yes	same
NBRS 4.3	Nobel Biocare	Replace® RP	4.3 mm	Product-No. 32377	K091756	no angulation correction to the implant axis	yes	yes	Internal 3 te-nons	yes	same
NBRS 5.0	Nobel	Replace® WP	5.0 mm	Product-No. 32378	K091756	no angulation	yes	yes	Internal	yes	same

October 20, 2010

Sirona Dental 510(k) Summary

Sirona Dental CAD/CAM System

Proposed Device TiBase	Predicate Device						Intended Use/ Differences	Abutment and Screw made of Ti6Al4V	Identical connection geometry to abutments	Con- nec- tion geo- metry Type	Con- nec- tion geo- metry Anti- rota- tional fea- tures	Screw geo- metry
	Manu- facturer	System	Diameter	Titanium Base	Predicate Devices K-Number							
	Biocare					correction to the implant axis			nal 3 te- nons			
NBR 6.0	Nobel Biocare	Replace® 6.0	6.0 mm	Product-No. 32375	K091756	no angulation correction to the implant axis	yes	yes	Inter- nal 3 te- nons	yes	same	
NBB 3.4	Nobel Biocare	Brånemark®	3.4 mm	Nobel Biocare product catalog page 14. Product-No. 32396	K091756	no angulation correction to the implant axis	yes	yes	Exter- nal Hexa- gonal	yes	same	
NBB 4.1	Nobel	Brånemark®	4.1 mm	Product-No. 32397	K091756	no angulation	yes	yes	Exter- nal	yes	same	

October 20, 2010

Sirona Dental 510(k) Summary

Sirona Dental CAD/CAM System

Proposed Device TiBase	Predicate Device				Intended Use/ Differences	Abutment and Screw made of Ti6Al 4V	Identical connection geometry to abutments	Connection geometry Type	Connection geometry	Screw geometry
	Manufacturer	System	Diameter	Titanium Base	Predicate Devices K-Number					
	Biocare					correction to the implant axis		nal Hexagonal		
SSO 3.5	Straumann	Tissue level NN	3.5 mm	Straumann product catalog page 51. Product-No. 048.505	K081005	no angulation correction to the implant axis	yes	external octagonal	yes	same
SSO 4.8	Straumann	Tissue level RN	4.8 mm	Product-No. 048.600 (p 55)	K081005	no angulation correction to the implant axis	yes	Internal Octagonal	yes	same
SSO 6.5	Straumann	Tissue level WN	6.5 mm	Product-No.	K081005	no angulation	yes	Inter-	yes	same

October 20, 2010

Sirona Dental 510(k) Summary

Sirona Dental CAD/CAM System

Proposed Device TiBase	Predicate Device				Intended Use/Differences	Abutment and Screw made of Ti6Al 4V	Identical connection geometry to abutments	Connection geometry Type	Connection geometry	Screw geometry
	Manufacturer	System	Diameter	Titanium Base						
				048.606 (p 62)	correction to the implant axis			nal Octagonal		
ATOS 3.5/4.0	Astra Tech	OsseoSpeed™	3.5 S / 4.0 S mm	Product-No. 24285	no angulation correction to the implant axis	yes	yes	Inter- nal Hexa- gonal	yes	same
ATOS 4.5/5.0	Astra Tech	OsseoSpeed™	4.5 / 5.0 mm	Product-No. 24235	no angulation correction to the implant axis	yes	yes	Inter- nal Hexa- gonal	yes	same
FX 3.4	Friadent	Frialit® / Xive®	3.4 mm	Friadent product catalog page 32.	no angulation correction to the implant	yes	yes	Inter- nal	yes	same

Sirona Dental 510(k) Summary

October 20, 2010

Sirona Dental CAD/CAM System

Proposed Device TiBase	Predicate Device				Intended Use/ Differences	Abutment and Screw made of Ti6Al 4V	Identical connection geometry to abutments	Connection geometry Type	Connection geometry Anti-rotational features	Screw geometry
	Manufacturer	System	Diameter	Titanium Base						
				Product-No. 46-2132	axis			Hexagonal		
FX 3.8	Friadent	Frialit® / Xive®	3.8 mm	Product-No. 46-2142	no angulation correction to the implant axis	yes	yes	Internal Hexagonal	yes	same
FX 4.5	Friadent	Frialit® / Xive®	4.5 mm	Product-No. 46-2152	no angulation correction to the implant axis	yes	yes	Internal Hexagonal	yes	same
FX 5.5	Friadent	Frialit® / Xive®	5.5 mm	Product-No. 46-2162	no angulation correction to the implant	yes	yes	Internal	yes	same

October 20, 2010

Sirona Dental 510(k) Summary

Sirona Dental CAD/CAM System

Proposed Device TiBase	Predicate Device				Intended Use/Differences	Abutment and Screw made of Ti6Al 4V	Identical connection geometry to abutments	Connection geometry Type	Connection geometry	Anti-rotational features	Screw geometry
	Manufacturer	System	Diameter	Titanium Base							
					axis				Hexagonal		
BO 3.4	Biomet 3i	Osseotite (Connection type: Ex. Hex)	3.4 mm	3i Biomet product catalog page 25. Product-No. MAP32G	no angulation correction to the implant axis	yes	yes	External Hexagonal	yes	same	same
BO 4.1	Biomet 3i	Osseotite (Connection type: Ex. Hex)	4.1 mm	Product-No. APP452G	no angulation correction to the implant axis	yes	yes	External Hexagonal	yes	same	same
BO 5.0	Biomet 3i	Osseotite (Connection type: Ex. Hex)	5.0 mm	Product-No. WPP552G	no angulation correction to the implant axis	yes	yes	External Hexagonal	yes	same	same

October 20, 2010

Sirona Dental 510(k) Summary

Sirona Dental CAD/CAM System

Proposed Device TiBase	Predicate Device					Intended Use/ Differences	Abutment and Screw made of Ti6Al4V	Identical connection geometry to abutments	Con- nec- tion geo- metry Type	Con- nec- tion geo- metry Anti- rota- tional fea- tures	Screw geo- metry
	Manu- facturer	System	Diameter	Titanium Base	Predicate Devices K-Number						
ZTSV 3.5	Zimmer	Tapered Screw- Vent®	3.5 mm	Zimmer product catalog page 9. Product-No. ZOA342S	K060880	no angulation correction to the implant axis	yes	yes	Inter- nal Hexa- gonal	yes	same
ZTSV 4.5	Zimmer	Tapered Screw- Vent®	4.5 mm	Product-No. ZOA442S	K060880	no angulation correction to the implant axis	yes	yes	Inter- nal Hexa- gonal	yes	same
ZTSV 5.7	Zimmer	Tapered Screw- Vent®	5.7 mm	Product-No. ZOA562S	K060880	no angulation correction to the implant axis	yes	yes	Inter- nal Hexa- gonal	yes	same

October 20, 2010

Sirona Dental 510(k) Summary

Sirona Dental CAD/CAM System

Proposed Device TiBase	Predicate Device			Intended Use/ Differences	Abutment and Screw made of Ti6Al 4V	Identical connection geometry to abutments	Connection geometry Type	Connection geometry	Screw geometry
Manufacturer	System	Diameter	Titanium Base	Predicate Devices K-Number					
					gonal				

October 20, 2010

Sirona Dental 510(k) Summary

Sirona Dental CAD/CAM System

**Table 6: Comparison of InCoris ZI meso to Predicate Device**

	<b>InCoris ZI meso</b>	<b>Camlog Ceramic sleeves K032448</b>
<b>Intended use</b>	inCoris ZI meso blocks are used in manufacturing individually designed mesostructures, which are glued to a fitting titanium base after milling and sintering.	ALTATEC BIOTECHNOLOGIES implants are indicated for single tooth replacement, as immediate abutments on long span to bridgework, as distal abutments on free-end edentulous areas to be restored with fixed bridgework, to support overdentures in totally or partially edentulous arches, and as abutments supporting a full arch fixed prosthesis in the totally edentulous mandible or maxilla.  The Ti-Ceramic abutment is used with the mentioned implant supported restorations where high aesthetics are desired.
<b>Application</b>	inCoris ZI mesostructures can only be used for the intended titanium bases or implants. Allocation of the connection size to the respective titanium base can be determined by the scanbody set of the respective implant system.  Please observe the indications and contraindications of the implant.	<ul style="list-style-type: none"> <li>• Bonded/cemented all-ceramic crown</li> <li>• Abutments with a 3.3 mm diameter for use only in the replacement of upper lateral incisors and lower incisors</li> <li>• Direct veneering with zirconium oxide dental ceramics</li> </ul>
<b>Contra-Indications</b>	<ul style="list-style-type: none"> <li>• Insufficient oral hygiene</li> <li>• Insufficient space available</li> <li>• Bruxism</li> <li>• For restorations with angulation correction to the implant axis</li> <li>• For individual tooth</li> </ul>	<ul style="list-style-type: none"> <li>• Primary abutment splinting</li> <li>• Restorations with an angulation correction of more than 20° to the implant axis</li> <li>• Single restorations with attached free end pontic</li> </ul>

	InCoris-ZI-meso	Camlog Ceramic sleeves K032448
	restorations with free end saddle <ul style="list-style-type: none"> <li>For restorations with a length to implant length ratio of more than 1:1.25</li> </ul>	<ul style="list-style-type: none"> <li>Restorations with a length relationship to the implant length greater than 1:1.25</li> </ul>
Technical Data		
Block-Material Composition	<ul style="list-style-type: none"> <li>ZrO<sub>2</sub>+HfO<sub>2</sub>+Y<sub>2</sub>O<sub>3</sub> &gt; 99.0%</li> <li>Al<sub>2</sub>O<sub>3</sub> &lt; 0.5%</li> <li>Other oxides &lt; 0.5%</li> </ul>	<ul style="list-style-type: none"> <li>ZrO<sub>2</sub>+HfO<sub>2</sub>+Y<sub>2</sub>O<sub>3</sub>: &gt; 99.0%</li> <li>Y<sub>2</sub>O<sub>3</sub>: 4.5 – 5.4%</li> <li>HfO<sub>2</sub>: &lt; 5%</li> <li>Al<sub>2</sub>O<sub>3</sub>: &lt; 0.5%</li> <li>Fe<sub>2</sub>O<sub>3</sub>: --</li> </ul>
Density (sintered)	6.06 g/cm <sup>3</sup>	> 6.00 g/cm <sup>3</sup>
Coefficient of thermal expansion (CTE)	11.0*10 <sup>-6</sup> K <sup>-1</sup> (20 °C - 500 °C)	No specification
Flexural strength	> 900MPa	≥ 800MPa
Fracture toughness (K <sub>IC</sub> )	5.8 MPa·m <sup>1/2</sup>	No specification
Grain Size	about 0.5 μm	< 0.6 μm
Anti-rotational feature	Notch	Notch
Bonding Material	Panavia F 2.0 (www.kuraray-dental.de)	Panavia F 2.0 (www.kuraray-dental.de)

**Table 7: Comparison of Sirona Dental CAD/CAM Design and fabricating Devices to Predicate Devices**

	<b>Sirona Dental CAD/CAM Design and fabricating Devices</b>	<b>Lava Software (K062493)</b>
Used for	<p>The Sirona Dental CAD/CAM-System is indicated for taking optical impressions to record the topographical characteristics of teeth, dental impressions, or stone models by computer aided design and fabrication in patients that require dental restorative prosthetic devices. The system also features the processing of mesostructures, a dental restorative prosthetic device used in conjunction with endosseous dental implant abutments.</p> <p>Sirona Dental CAD/CAM Design and fabricating devices feature the processing of mesostructures, a dental restorative prosthetic device used in conjunction with endosseous dental implant abutments, i.e. it is an accessory to it. Devices which feature the processing of mesostructures comprises CEREC3, CEREC AC, inEos, inEos Blue, CEREC MCXL and inLab MCXL</p>	<p>The Lava software is used with 3M ESPE's Lava system, an all-ceramic system for the CAD/CAM fabrication of dental restorations such as inlays, onlays, veneers, crowns and bridges. The software controls the measuring process, processing of the measurement data (3D-CAD tool), and export of the data to the milling machine. In addition, various patient and case information elements can be entered. Other functions are available for verification and service of the measuring system. The Lava software also facilitates the transfer of 3D data from a scanner to a remote milling machine via internet.</p>
Used with	Sirona Dental CAI/CAM Hardware	3M ESPE Lava System CAI/CAM Hardware
Controlling of recording process (CAI) (optical impression)	Yes	Yes
Processing the recorded data (data of optical	Yes	Yes

	<b>Sirona Dental CAD/CAM Design and fabricating Devices</b>	<b>Lava Software (K062493)</b>
impression) (CAD)		
Export of milling data to milling machine	Yes	Yes
Administration of patient data	Yes	Yes
Further functions	Calibration of CAI/CAM hardware	Verification and service of measuring system
Online capability	Option to upload/download the data from a web portal (Cerec Connect), to have CAI and CAM operating on two different locations connected via Internet	Facilitates the transfer of the 3D data from a scanner to a remote milling machine via Internet
Scan Implant Interface/surface	Yes (or with mounted scanbody)	Yes
Scan custom wax-up	Yes	Yes
Preparation of customized restoration ("meso-structure") to be mounted on the abutment	Yes	Yes
Bond of milled zirconia/ceramic customized meso-structure to metal abutment	Yes	Yes
Create of fitting crown to be mounted on top of meso-structure	Yes	Yes
Used with	Sirona Dental CAI/CAM Hardware	3M ESPE Lava System CAI/CAM Hardware
Used for	CAD creation of dental restorations including inlays,	CAD creation of dental restorations including inlays,

	<b>Sirona Dental CAD/CAM Design and fabricating Devices</b>	<b>Lava Software (K062493)</b>
	onlay, veneers, crowns, bridges and meso-structure to be mounted on top of abutments	onlay, veneers, crowns, bridges and meso-structure to be mounted on top of abutments
Controlling of measurement process (CAI)	Yes	Yes
Processing the measurement data (CAD)	Yes	Yes
Export to milling machine	Yes	Yes
Administration of patient data	Yes	Yes
Further functions	Calibration of CAI/CAM hardware	Verification and service of measuring system

### 7 Nonclinical Testing

According to FDA Guidance "Class II Special Controls Guidance Document: Root-form Endosseous Dental Implants and Endosseous Dental Implant Abutments", May 12, 2004, for no angled abutments fatigue testing is not required for straight abutments and has not been performed. A reverse-engineering analysis has demonstrated that the devices are identical with their predicates and compatible with their mating implants.

### 8 Clinical Testing

Clinical testing is not required and has not been performed.

### 9 Conclusion

Based on a comparison of intended use, indications, construction materials, principal of operations, features and technical data, the Sirona Dental CAD/CAM System which comprises of titanium bases TiBase, inCoris ZI meso blocks and Sirona Dental CAD/CAM Design and fabricating devices are safe and effective their

intended use and perform as well as and are substantially equivalent to their  
Predicate Devices.

---



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room -WO66-G609  
Silver Spring, MD 20993-0002

Ms. Fritz Kolle  
Quality Management/Regulatory Affairs  
Sirona Dental Systems GmbH  
Fabrikstrasse 31  
64625 Bensheim  
GERMANY

OCT 22 2010

Re: K100152

Trade/Device Name: Sirona Dental CAD/CAM System  
Regulation Number: 21 CFR 872.3630  
Regulation Name: Endosseous Dental Implant Abutment  
Regulatory Class: II  
Product Code: NHA  
Dated: October 6, 2010  
Received: October 8, 2010

Dear Ms. Kolle:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to

<http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Anthony D. Watson, B.S., M.S., M.B.A.  
Director

Division of Anesthesiology, General Hospital,  
Infection Control and Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(k) Number (if known): K100152

OCT 22 2010

Device Name: Sirona Dental CAD/CAM System

Indications for Use:

The Sirona Dental CAD/CAM System is intended for use in partially or fully edentulous mandibles and maxillae in support of single or multiple-unit cement retained restorations. The system consists of three major parts: TiBase, InCoris mesostructure, and CAD/CAM software. Specifically, the InCoris mesostructure and TiBase components make up a two-piece abutment which is used in conjunction with endosseous dental implants to restore the function and aesthetics in the oral cavity. The InCoris mesostructure may also be used in conjunction with the Camlog Titanium base CAD/CAM (types K2244.xxxx) (K083496) in the Camlog Implant System. The CAD/CAM software is intended to design and fabricate the InCoris mesostructure. The InCoris mesostructure and TiBase two-piece abutment is compatible with the following implants systems:

- *Nobel Biocare Replace (K020646)*
- *Nobel Biocare Branemark (K022562)*
- *Friadent Xive (K013867)*
- *Biomet 3i Osseotite (K980549)*
- *Astra Tech Osseospeed (K091239)*
- *Zimmer Tapered Screw-Vent (K061410)*
- *Straumann SynOcta (K061176)*

Prescription Use  X   
(Part 21 CFR 801 Subpart D)

OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NECESSARY)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Susan Purine

(Division Sign-Off)

Sirona Dental Systems 510(k)  
Additional Information K100152

October 20, 2009  
Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

510(k) Number: K100152