

1C101195  
JUL - 6 2010

## 510(k) Summary of Safety and Effectiveness

*This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.*

### Introduction

According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence.

### Submitter name, Address, and Contact

Lin-Zhi International, Inc.  
670 Almanor Avenue  
Sunnyvale, CA 94085  
Phone: (408) 732-3856  
Fax: (408) 732-3849  
e-mail: bclin@lin-zhi.com

Contact: Bernice Lin, Ph.D.  
V.P. of Operations

### Device Name and Classification

Classification Name: Enzyme Immunoassay, Opiates  
Class II, DJG (91 Toxicology),  
21 CFR 862.3650

6-Acetylmorphine Calibrators,  
Class II, DLJ (91 Toxicology),  
21 CFR 862.3200

6-Acetylmorphine Controls,  
Classification LAS (91 Toxicology),  
21 CFR 862.3280

Common Name: Homogeneous 6-Acetylmorphine Enzyme Immunoassay  
Proprietary Name: LZI 6-Acetylmorphine Enzyme Immunoassay,  
LZI 6-Acetylmorphine Drugs of Abuse (DAU) Calibrators  
LZI 6-Acetylmorphine Drugs of Abuse (DAU) Controls

## Legally Marketed Predicate Device(s)

The LZI 6-Acetylmorphine Enzyme Immunoassay (EIA) is substantially equivalent to the Microgenics CEDIA® DAU 6-Acetylmorphine Assay, Calibrators and Controls for Hitachi 717 Systems (K001178) manufactured by Microgenics Corp. The LZI 6-Acetylmorphine Enzyme Immunoassay is identical or similar to its predicate in terms of intended use, method principle, device components, and clinical performance.

## Device Description

The LZI 6-Acetylmorphine assay is a homogeneous enzyme immunoassay with ready-to-use liquid reagent. The assay is based on competition between drug in the sample and drug labeled with the enzyme glucose-6-phosphate dehydrogenase (G6PDH) for a fixed amount of antibody in the reagent. Enzyme activity decreases upon binding to the antibody, and the drug concentration in the sample is measured in terms of enzyme activity. In the absence of drug in the sample, 6-Acetylmorphine-labeled G6PDH conjugate is bound to antibody, and the enzyme activity is inhibited. On the other hand, when free drug is present in the sample, antibody would bind to free drug, the unbound 6-Acetylmorphine-labeled G6PDH then exhibits its maximal enzyme activity. Active enzyme converts nicotinamide adenine dinucleotide (NAD) to NADH, resulting in an absorbance change that can be measured spectrophotometrically at 340 nm.

## Intended Use

The LZI 6-Acetylmorphine Enzyme Immunoassay is intended for the qualitative and semi-quantitative determination of 6-Acetylmorphine in human urine, at a cutoff value of 10 ng/mL. The assay is designed for professional use with a number of automated clinical chemistry analyzers.

The LZI 6-Acetylmorphine Drugs of Abuse (DAU) Calibrators are for use as calibrators in the qualitative and semi-quantitative calibration of the LZI 6-Acetylmorphine Enzyme Immunoassay.

The 6-Acetylmorphine Drugs of Abuse (DAU) Controls are for use as assayed quality control materials to monitor the precision of the LZI 6-Acetylmorphine Enzyme Immunoassay.

***The assay provides only a preliminary analytical result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas or liquid chromatography/mass spectrometry (GC/MS or LC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be exercised with any drug of abuse test result, particularly when the preliminary test result is positive.***

## Performance Characteristics Summary: Hitachi 717 Analyzer

### Precision:

Precision: Semi-Quantitative, ng/mL

N=88 (ng/mL)	Within Run			Total Precision		
	Mean	SD	% CV	Mean	SD	% CV
0 ng/mL	0.2	0.19	N/A	0.2	0.3	169.4
2.5 ng/mL	2.8	0.32	11.6	2.8	0.3	12.6
5.0 ng/mL	5.1	0.28	5.5	5.1	0.4	7.2
7.5 ng/mL	7.4	0.37	5.0	7.4	0.4	5.1
10 ng/mL	9.7	0.34	3.5	9.7	0.4	4.4
12.5 ng/mL	12.2	0.4	3.1	12.2	0.4	3.6
15 ng/mL	14.6	0.47	3.3	14.6	0.5	3.7
17.5 ng/mL	17.7	0.44	2.5	17.7	0.5	2.9
20 ng/mL	19.9	0.51	2.6	19.9	0.6	2.9

### Semi-Quantitative Positive/Negative Results:

10 ng/mL Cutoff Result:		Within Run		Total Precision	
Sample Concentration	% of Cutoff	Number of Determination	Immunoassay Result	Number of Determination	Immunoassay Result
0 ng/mL	-100.0%	22	22 Negative	88	88 Negative
2.5 ng/mL	-75.0%	22	22 Negative	88	88 Negative
5.0 ng/mL	-50.0%	22	22 Negative	88	88 Negative
7.5 ng/mL	-25.0%	22	22 Negative	88	22 Pos/66 Neg
10 ng/mL	100.0%	22	6 Pos/16 Neg	88	88 Positive
12.5 ng/mL	+25.0%	22	22 Positive	88	88 Positive
15 ng/mL	+50.0%	22	22 Positive	88	88 Positive
17.5 ng/mL	+75.0%	22	22 Positive	88	88 Positive
20 ng/mL	+100.0%	22	22 Positive	88	88 Positive

### Precision: Qualitative, mA/min

N=88 (mA/min)	Within Run			Total Precision		
	Mean	SD	% CV	Mean	SD	% CV
0 ng/mL	358.6	2.2	0.6	358.6	2.4	0.7
2.5 ng/mL	380.0	2.5	0.7	380.0	2.8	0.7
5.0 ng/mL	399.4	2.2	0.6	399.4	2.4	0.6
7.5 ng/mL	419.5	2.4	0.6	419.5	2.7	0.6
10 ng/mL	437.4	2.2	0.5	437.4	2.6	0.6
12.5 ng/mL	454.8	2.0	0.4	454.8	3.0	0.7
15 ng/mL	470.4	3.8	0.8	470.4	4.6	1.0
17.5 ng/mL	489.2	2.3	0.5	489.2	3.0	0.6
20 ng/mL	500.1	2.2	0.4	500.1	2.8	0.6

**Qualitative Positive/Negative Results:**

10 ng/mL Cutoff Result:		Within Run		Total Precision	
Sample Concentration	% of Cutoff	Number of Determination	Immunoassay Result	Number of Determination	Immunoassay Result
0 ng/mL	-100.0%	22	22 Negative	88	88 Negative
2.5 ng/mL	-75.0%	22	22 Negative	88	88 Negative
5.0 ng/mL	-50.0%	22	22 Negative	88	88 Negative
7.5 ng/mL	-25.0%	22	22 Negative	88	88 Positive
10 ng/mL	100.0%	22	10 Pos/12 Neg	88	41 Pos/47 Neg
12.5 ng/mL	+25.0%	22	22 Positive	88	88 Positive
15 ng/mL	+50.0%	22	22 Positive	88	88 Positive
17.5 ng/mL	+75.0%	22	22 Positive	88	88 Positive
20 ng/mL	+100.0%	22	22 Positive	88	88 Positive

**Limit of Detection:**

The lowest concentration that can be differentiated from the negative urine with 95% confidence is determined as 2 ng/mL.

**Linearity:**

Hitachi 717 Instrument: 2 - 40 ng/mL

When comparing the result (y) and target (x) value, using the least squares regression technique, the regression equation and correlation are as follow:

$$y = 0.9632x + 0.2491, r^2 = 0.9961$$

**Method Comparison: Clinical Samples**

From a total of eighty (80) clinical unaltered samples:

10 ng/mL Cutoff

(100 % agreement with positive, 93 % agreement with negative samples)

**Endogenous Compound Interference and Specificity - Cross-Reactivity:**

No significant undesired cross reactants or endogenous substance interference was observed. See product insert for list of compounds tested.

**Summary:**

The information provided in this pre-market notification demonstrates that the LZI 6-Acetylmorphine Enzyme Immunoassay is substantially equivalent to the legally marketed predicate device for its general intended use. Substantial equivalence was demonstrated through comparison of intended use and physical properties to the commercially available predicate device as confirmed by liquid chromatography/mass spectrometry (LC/MS), an independent analytical method. The information supplied in this pre-market notification provides reasonable assurance that the LZI 6-Acetylmorphine Enzyme Immunoassay is safe and effective for its stated intended use.



DEPARTMENT OF HEALTH & HUMAN SERVICES

---

Lin-Zhi International, Inc.  
c/o Berince Lin, Ph.D.  
V.P. of Operations  
670 Almanor Ave  
Sunnyvale, CA 94085

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Mail Center -- WO66-0609  
Silver Spring, MD 20993-0002

Re: k101195

Trade/Device Name: 6-Acetylmorphine Enzyme Immunoassay and 6-Acetylmorphine  
Calibrators and Controls

Regulation Number: 21 CFR 862.3650

Regulation Name: Opiate test system

Regulatory Class: Class II

Product Code: DJG, DLJ and LAS

Dated: April 28, 2010

Received: April 29, 2010

**JUL 06 2010**

Dear Dr. Lin:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); and good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820).

Page 2

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (301) 796-5450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding postmarket surveillance, please contact CDRH's Office of Surveillance and Biometric's (OSB's) Division of Postmarket Surveillance at (301) 796-5760. For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or ( 301 ) 796-5680 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Courtney C. Harper, Ph.D.  
Director  
Division of Chemistry and Toxicology  
Office of *In Vitro* Diagnostic Device  
Evaluation and Safety  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K101195

Device Name: 6-Acetylmorphine Enzyme Immunoassay  
6-Acetylmorphine Calibrators and Controls

### Indications for Use:

The 6-Acetylmorphine Enzyme Immunoassay is intended for the qualitative and semi-quantitative determination of 6-Acetylmorphine in human urine, at a cutoff value of 10 ng/mL. The assay is designed for professional use with a number of automated clinical chemistry analyzers.

The 6-Acetylmorphine Drugs of Abuse (DAU) Calibrators are for use as calibrators in the qualitative and semi-quantitative calibration of the 6-Acetylmorphine Enzyme Immunoassay.

The 6-Acetylmorphine Drugs of Abuse (DAU) Controls are for use as assayed quality control materials to monitor the precision of the 6-Acetylmorphine Enzyme Immunoassay.

*The assay provides only a preliminary analytical result. A more specific alternative chemical method must be used in order to obtain a confirmed analytical result. Gas or liquid chromatography/mass spectrometry (GC/MS or LC/MS) is the preferred confirmatory method). Clinical consideration and professional judgment should be exercised with any drug of abuse test result, particularly when the preliminary test result is positive.*

Prescription Use  AND/OR  
(Part 21 CFR 801 Subpart D)

Over-The Counter Use   
(21 CFR 807 Subpart C)

---

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)



Division Sign-Off  
Office of In Vitro Diagnostic Device Evaluation and Safety  
510(k) K101195