

K101682

**510(k) Summary  
for the S 100 Pedicle Screw System**

In accordance with 21 CFR 807.92 of the Federal Code of Regulations  
the following 510(k) summary is submitted for the S 100 Pedicle Screw System

DEC - 7 2010

**Date Prepared:** November 8, 2010

**1. Submitter:**

Renovis Surgical Technologies  
1901 W. Lugonia Ave, Ste 340  
Redlands, CA 92374

**2. Contact Person:**

Anthony DeBenedictis

**2. Trade name:**

S 100 Pedicle Screw System

**Common Name:**

orthosis, pedicle screw system

**Classification Name:**

orthosis, spondylolisthesis spinal fixation/ MNH/ 888.3070

orthosis, spinal pedicle fixation/ MNI/ 888.3070

Class II

**3. Predicate or legally marketed devices which are substantially equivalent:**

ZODIAC™ Polyaxial Spinal Fixation - K042673 (Alphatec Spine Co)

Sequoia – K082032 (Abbott Spine)

Moss Miami - K000536 (DePuy)

Synergy VLS – K000236 (DePuy)

Rogozinski K884263 - (Smith & Nephew)

**4. Description of the device:**

As a posterior pedicle screw system designed for temporary stabilization of the anterior spine during the development of spinal fusion, the Renovis S 100 Pedicle Screw System is comprised of polyaxial pedicle screws, rods, and crosslinks. The S 100 System can be used for single or multiple level fixations

The screws are a top loading tulip design and are available in multiple diameters and lengths. Reduction screws are available for cases of spondylolisthesis where the short arms of the tulip of the standard screw are not long enough to engage the rod. The rods are available in straight and pre-lordosed (curved) configurations. The system also has variable and fixed crosslinks.

**Materials:**

The components are manufactured from titanium alloy (ASTM F136), CP titanium Grade 4 (ASTM F67) and CoCrMo (ASTM F1537).

**Function:**

The S 100 Pedicle Screw System is intended to provide immobilization and stabilization of spinal segments until fusion takes place.

**5. Substantial equivalence claimed to predicate devices:**

The S 100 Pedicle Screw System is substantially equivalent to the ZODIAC™ Polyaxial Spinal Fixation in terms of intended use, design, materials used, mechanical safety and performances. The table below compares the features and characteristics of the S 100 Pedicle Screw System to this predicate device. The S 100 Pedicle Screw is similar to the Sequoia Screw - K082032 (Abbott Spine) in terms of assembly and the Moss Miami - K000536 (DePuy), the Synergy VLS – open K000236 (DePuy) and the Rogozinski K884263 - (Smith & Nephew) in terms of strength.

sheet 1 of 2

Items \ Device Name	Renovis S 100 Pedicle Screw System	ZODIAC™ Polyaxial Spinal Fixation
Sponsor	Renovis Surgical Technologies	Alphatec Spine Co.
510(k) Number	--	K042673
Device Classification Name	orthosis, spinal pedicle fixation orthosis, spondylolisthesis spinal fixation spinal interlaminar fixation orthosis	orthosis, spinal pedicle fixation orthosis, spondylolisthesis spinal fixation spinal interlaminar fixation orthosis
Product Code	MNI, MNH, KWP	MNI, MNH, KWP
Indications for Use	Per MNI, MNH, KWP	Per MNI, MNH, KWP
Material	Ti-6Al-4V per ASTM F136 CP Ti Grade 4 per ASTM F67 CoCrMo per ASTM F1537	Ti-6Al-4V per ASTM F136 CP Ti per ASTM F67 CoCrMo per ASTM F1537
Top loading	Yes	Yes
Solid	Yes	Yes
Cannulated	Yes	Yes
Monoaxial	No	Yes
Polyaxial	Yes	Yes
Crosslink	Yes	Yes

**6. Intended Use:**

The Renovis S 100 Pedicle Screw System is intended to provide immobilization and stabilization of spinal segments in skeletally mature patients as an adjunct to fusion in the treatment of the following acute and chronic instabilities or deformities of thoracic, lumbar, and sacral spine: fracture, dislocation, failed previous fusion (pseudoarthrosis), spinal stenosis, degenerative spondylolisthesis with objective evidence of neurological impairment, spinal deformations such as scoliosis or kyphosis and loss of stability due to tumors.

When used as a pedicle screw system, the Renovis S 100 Pedicle Screw System is intended for the treatment of severe spondylolisthesis (Grade 3 and 4) of the L5-S1 vertebrae in skeletally mature patients receiving fusion by autogenous bone graft having implants attached to the lumbar and sacral spine (L3 to sacrum) with removal of the implants after the attainment of a solid fusion.

**7. Summary of Nonclinical Tests:**

The following tests were conducted:

- Testing per ASTM F1717
- Testing per ASTM F1798
- Component dissociate testing

The results of this testing was compared to a predicate system, with the results being equal or higher than the predicate system.

**8. Clinical Test Summary:**

No clinical studies were performed

**9. Conclusions Nonclinical and Clinical:**

The Renovis S 100 Pedicle Screw System is substantially equivalent to the predicate devices in terms of indications for use, design, material, function and method of assembly.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room W-O66-0609  
Silver Spring, MD 20993-0002

Renovis Surgical Technologies, LLC  
% The Orthomedix Group, Inc.  
Mr. J.D. Webb  
1001 Oakwood Boulevard  
Round Rock, Texas 78681

DEC - 7 2010

Re: K101682

Trade/Device Name: Renovis S 100 Pedicle Screw System  
Regulation Number: 21 CFR 888.3070  
Regulation Name: Pedicle screw spinal system  
Regulatory Class: Class II  
Product Code: MNI, MNH, KWP  
Dated: November 08, 2010  
Received: November 15, 2010

Dear Mr. Webb:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act

or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Mark N. Melkerson  
Director  
Division of Surgical, Orthopedic  
And Restorative Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K101682

Device Name: Renovis S 100 Pedicle Screw System

DEC - 7 2010

### Indications for Use:

The Renovis S 100 Pedicle Screw System is intended to provide immobilization and stabilization of spinal segments in skeletally mature patients as an adjunct to fusion in the treatment of the following acute and chronic instabilities or deformities of thoracic, lumbar, and sacral spine: fracture, dislocation, failed previous fusion (pseudoarthrosis), spinal stenosis, degenerative spondylolisthesis with objective evidence of neurological impairment, spinal deformations such as scoliosis or kyphosis and loss of stability due to tumors.

When used as a pedicle screw system, the Renovis S 100 Pedicle Screw System is intended for the treatment of severe spondylolisthesis (Grade 3 and 4) of the L5-S1 vertebrae in skeletally mature patients receiving fusion by autogenous bone graft having implants attached to the lumbar and sacral spine (L3 to sacrum) with removal of the implants after the attainment of a solid fusion.

Prescription Use X  
(Part 21 CFR 801 Subpart D)


AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
\_\_\_\_\_  
(Division Sign-Off)  
Division of Surgical, Orthopedic,  
and Restorative Devices

510(k) Number K101682

sheet 1 of 1