



K102159  
NOV 17 2010

**510(k) Summary**  
**21 CFR 807.92(a)**

**PowerPICC® SV Catheter**

---

<b>General Provisions</b>	Submitter Name:	Bard Access Systems, Inc.
	Address:	605 North 5600 West Salt Lake City, UT 84116
	Contact Person:	Jessica Agnello Regulatory Affairs Specialist
	Telephone Number:	(801) 522-5651
	Fax Number:	(801) 522-5425
	Date of Preparation:	November 4, 2010

---

<b>Subject Device</b>	Trade Name:	<b>PowerPICC® SV Catheter</b>
	Common Name:	Peripherally Inserted Central Catheter (PICC)
	Classification Name:	Percutaneous, implanted, long-term intravascular catheter LJS, 21 CFR § 880.5970 Class II

---

<b>Predicate Devices</b>	Trade Name:	<b>Poly Per-Q-Cath® PICC Catheter</b>
	Classification Name:	Percutaneous, implanted, long-term intravascular catheter, LJS, 21 CFR §880.5970 Class II
	Premarket Notification:	K034019, K051991
	Manufacturer:	Bard Access Systems, Inc.
	Trade Name:	<b>6F Triple Lumen Poly Per-Q-Cath® PICC Catheter</b>
	Classification Name:	Percutaneous, implanted, long-term intravascular catheter, LJS, 21 CFR §880.5970 Class II
	Premarket Notification:	K043502, K051991
	Manufacturer:	Bard Access Systems, Inc.
	Trade Name:	<b>4F Single Lumen PowerPICC® Catheter</b>
	Classification Name:	Percutaneous, implanted, long-term intravascular catheter, LJS, 21 CFR §880.5970 Class II
	Premarket Notification:	K070996
	Manufacturer:	Bard Access Systems, Inc
Trade Name:	<b>5F Single Lumen PowerPICC® Catheter</b>	
Classification Name:	Percutaneous, implanted, long-term intravascular catheter, LJS, 21 CFR §880.5970 Class II	
Premarket Notification:	K033389, K051991	
Manufacturer:	Bard Access Systems, Inc	

---

Trade Name: 6F Triple Lumen **PowerPICC**<sup>®</sup> Catheter  
 Classification Name: Percutaneous, implanted, long-term intravascular catheter, LJS, 21 CFR §880.5970 Class II  
 Premarket Notification: K053501, K051991  
 Manufacturer: Bard Access Systems, Inc

Trade Name: **PowerLine**<sup>®</sup> Catheters  
 Classification Name: Percutaneous, implanted, long-term intravascular catheter, LJS, 21 CFR §880.5970 Class II  
 Premarket Notification: K093927  
 Manufacturer: Bard Access Systems, Inc

Trade Name: **Morpheus**<sup>®</sup> CT PICC Catheter  
 Classification Name: Percutaneous, implanted, long-term intravascular catheter, LJS, 21 CFR §880.5970 Class II  
 Premarket Notification: K072196  
 Manufacturer: AngioDynamics, Inc

**Device Description**

The subject **PowerPICC**<sup>®</sup> SV catheters are members of Bard Access Systems' **PowerPICC**<sup>®</sup> series of power injectable catheters. The **PowerPICC**<sup>®</sup> SV catheters are available in 3F single lumen and 4F dual lumen configurations. The catheter extension leg, luer hub and junction are printed with description markings to facilitate proper use of the device. The **PowerPICC**<sup>®</sup> SV catheters are provided in sterile kit configurations.

**Intended Use**

The **PowerPICC**<sup>®</sup> SV catheters are intended for short or long term peripheral access to the central venous system for intravenous therapy and blood sampling.

**Indications For Use**

The **PowerPICC**<sup>®</sup> SV catheter is indicated for short or long-term peripheral access to the central venous system for intravenous therapy, blood sampling, power injection of contrast media, and allows for central venous pressure monitoring. For central venous pressure monitoring, it is recommended that a catheter lumen of 20 gauge or larger be used.

Catheter Size	Maximum Flow Rate
3 F Single Lumen	1 mL/sec
4 F Dual Lumen	2.5 mL/sec

**Technological Characteristics**

Technological characteristics of the subject **PowerPICC**<sup>®</sup> SV catheters are equivalent with respect to the basic catheter design and function to those of the predicate Poly **Per-Q-Cath**<sup>®</sup> Catheters. Technical characteristics based on the subject device's power injection capabilities are comparable to the predicate 4F SL **PowerPICC**<sup>®</sup> catheter and the **Morpheus**<sup>®</sup> CT PICC catheter. Differences do not raise any new questions regarding safety and effectiveness.

---

Verification and validation activities were designed and performed to demonstrate that the subject **PowerPICC**<sup>®</sup> SV catheters met predetermined performance specifications. Tests were performed on sterilized, finished devices. The following guidance documents and standards in conjunction with in-house protocols were used to determine appropriate methods for evaluating the performance of the device:

**Safety &  
Performance  
Tests**

- *Guidance on Premarket Notification [510(k)] Submission for Short-Term and Long-Term Intravascular Catheters, March 16, 1995*
- *ISO 10555-1:2009, Sterile, single-use intravascular catheters, Part 1. General requirements*
- *BS/EN/ISO 10555-3:1996/Cor 1:2002, Sterile, single-use intravascular catheters, Part 3. Central venous catheters*
- *AAMI/ANSI/ISO 10993-1:2009, Biological Evaluation of Medical Devices Part 1: Evaluation and Testing, and the FDA Modified ISO 10993 Test Profile*
- *AAMI/ANSI/ISO 10993-7:2008, Biological Evaluation of Medical Devices Part 7: Ethylene Oxide Sterilization Results*
- *AAMI/ANSI/ISO 11135:2007, Sterilization of Healthcare Products Ethylene oxide -- Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices*

The subject **PowerPICC**<sup>®</sup> SV catheters met all predetermined acceptance criteria derived from the above mentioned references. Design validation was conducted on the subject **PowerPICC**<sup>®</sup> SV configuration and yielded acceptable results.

Risk management, including a failure modes and effects analysis (FMEA), of the subject device was conducted in accordance with ISO 14971:2007, *Medical Devices – Risk Management for Medical Devices*. No new types of safety or efficacy questions were identified for the subject **PowerPICC**<sup>®</sup> SV catheters.

---

**Summary of  
Substantial  
Equivalence**

Based on the indications for use, technological characteristics, and safety and performance testing, the subject **PowerPICC**<sup>®</sup> SV catheters met the minimum requirements for its intended use/indications for use and is substantially equivalent in design, materials, sterilization, principles of operation and indications for use to current commercially available catheters/cited predicates.

---



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room -WO66-G609  
Silver Spring, MD 20993-0002

Ms. Jessica Agnello  
Regulatory Affairs Specialist  
C.R. Bard, Incorporated  
605 North 5600 West  
Salt Lake City, Utah 84116

NOV 17 2010

Re: K102159  
Trade/Device Name: PowerPICC<sup>®</sup> SV Catheter  
Regulation Number: 21 CFR 880.5970  
Regulation Name: Percutaneous, Implanted, Long-Term Intravascular Catheter  
Regulatory Class: II  
Product Code: LJS  
Dated: November 4, 2010  
Received: November 5, 2010

Dear Ms. Agnello:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.


If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to

<http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Anthony D. Watson, B.S., M.S., M.B.A.  
Director  
Division of Anesthesiology, general Hospital,  
Infection Control and dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**Indications for Use**

NOV 17 2010

510(k) Number (if known): K102159

Device Name: PowerPICC® SV Catheter

**Indications for Use:**

The **PowerPICC®** SV catheter is indicated for short or long-term peripheral access to the central venous system for intravenous therapy, blood sampling, power injection of contrast media, and allows for central venous pressure monitoring. For central venous pressure monitoring, it is recommended that a catheter lumen of 20 gauge or larger be used.

Catheter Size	Maximum Flow Rate
3 F Single Lumen	1 mL/sec
4 F Dual Lumen	2.5 mL/sec

Prescription Use    
 (Part 21 CFR §801 Subpart D)

AND/OR

Over-The-Counter Use    
 (21 CFR §801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

R. C. Chagnon 11/17/10  
(Division Sign-Off) Page 1 of 1  
Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

510(k) Number: K102159