

Reverse Medical Corporation
 RePort™ Guide Catheter System 510(k) Submission

K102418

4. 510(k) Summary according to 807.92

Pursuant to Section 12, Part (a)(i)(3A) of the Safe Medical Devices Act of 1990, Reverse Medical Corporation is providing the summary of Substantial Equivalence for the RePort™ Guide Catheter System.

4.1 Sponsor /Applicant Name and Address

Reverse Medical Corporation
 13900 Alton Parkway
 Suite 123
 Irvine, CA 92618

4.2 Sponsor Contact Information

Amy Eskina
 Vice President, Regulatory Quality & Clinical Affairs
 Phone: 949-215-0660 x234
 Fax: 949-215-0661
 email: aeskina@reversemed.com

4.3 Date of Preparation of 510(k) Summary

August 13, 2010

4.4 Device Trade or Proprietary Name

RePort™ Guide Catheter System

4.5 Device Common/Usual or Classification Name

Catheter, Percutaneous (Product Code: DQY)

4.6 Identification of the Legally Marketed Devices to which Equivalence is Being Claimed:

Name of Predicate Device	Name of Manufacturer (Town, State)	510(k) Number
Penumbra Neuron™ Intracranial Access System	Penumbra, Inc. Alameda, CA	K070970
ENVOY® Guiding Catheter	Codman & Shurtleff, Inc Raynham, MA	K093184
Chaperon™ Guiding Catheter System	Microvention, Inc. Aliso Viejo, CA	K082385

4.7 Device Description

The RePort™ Guide Catheter System is a two-catheter system comprised of the RePort™ Sheath and the RePort™ Dilator. The RePort™ Sheath has a distal region which expands when the RePort™ Dilator is advanced. The RePort™ Guide Catheter System can be used with hydrophilic guidewires up to .038 inches in diameter to access the desired anatomy.

The proximal ends of the RePort™ Sheath and the RePort™ Dilator have a luer fitting to allow attachment of accessories and infusion of liquids through the catheter. The catheter is offered in various sizes to accommodate physician preferences and patient anatomy. A split sheath is provided in the package to provide support and facilitate the introduction of the RePort™ Sheath's distal expandable region into the catheter introducer sheath. The System is provided sterile, non-pyrogenic, and is intended for single use only.

4.8 Intended Use

The Reverse Medical RePort™ Guide Catheter System is indicated for the introduction of interventional/diagnostic devices into the peripheral, coronary and neuro vasculature.

4.9 Comparison to Predicate Devices

	Penumbra Neuron™ Intracranial Access System	ENVOY® Guiding Catheter;	Chaperon™ Guiding Catheter System	RePort™ Guide Catheter
510(k) Number	K070970	K093184	K082385	TBD
Classification	Class II, DQY	Class II, DQY	Class II, DQY	Class II, DQY
Indication	The ... is indicated for the introduction of interventional devices into the peripheral, coronary, and neuro vasculature.	The ... is intended for use in the peripheral, coronary, and neurovasculature for the intravascular introduction of interventional/diagnostic devices.	The ...is intended for general intravascular use, including the neuro and peripheral vasculature. The... can be used to facilitate introduction of diagnostic or therapeutic devices. The... are not intended for use in coronary arteries.	The ... is indicated for the introduction of interventional / diagnostic devices into the peripheral, coronary, and neuro vasculature.
Materials				
- Shaft Materials	PTFE lined nylon/polyurethane with stainless steel braid support	PTFE lined nylon/polyurethane with stainless steel braid support	PTFE lined polymeric catheter with braid support	PTFE lined nylon/polyurethane with stainless steel and nitinol braid support
- Proximal End Configuration	Luer Hub	Luer Hub	Luer Hub	Luer Hub
- Radiographic markers / Radiopacity	<ul style="list-style-type: none"> • Delivery Catheter: Platinum reinforced distal shaft for radiopacity • Select Catheter: shaft is radiopaque throughout 	Radiopaque distal tip	Radiopaque distal tip	<ul style="list-style-type: none"> • Sheath: Platinum/tungsten bands embedded at junction of expandable section; gold marker at tip of expandable section; • Dilator: shaft is radiopaque throughout
- Packaging	Catheter in polyethylene hoop inside PET/PE/Tyvek pouch inside SBS carton	Catheter attached to packaging card inside PET/PE/Tyvek pouch inside SBS carton	Catheter attached to packaging card inside PET/PE/Tyvek pouch inside SBS carton	Catheter in polyethylene hoop attached to packaging card inside PET/PE/Tyvek pouch inside SBS carton
Sterilization	EtO	EtO	EtO	EtO

4.10 Summary of Non-clinical Data

4.10.1 Biocompatibility and Sterilization

The RePort™ Guide Catheter System is classified as Externally Communicating Devices, Circulating Blood, Limited Contact (≤ 24 hours). Results of the testing demonstrate that the blood contacting materials are biocompatible.

Blood contacting materials were tested in accordance with the tests recommended in the FDA General Program Memorandum #G95-1 (5/1/95): Use of International Standard ISO 10993-1 guidelines "Biological Evaluation of Medical Devices Part 1: Evaluation and Testing". The RePort™ Guide Catheter System successfully passed all of the following biocompatibility tests:

Test	Method
Cytotoxicity	L929 MEM Elution Test
Sensitization	Kligman Maximization
Intracutaneous Reactivity (Irritation)	Intracutaneous Injection Test
Systemic Toxicity (Acute)	ISO Acute Systemic Injection Test
Haemocompatibility	Complement Activation
	Hemolysis
	Inactivated Partial Thromboplastin Time Test
	<i>In vivo</i> thrombogenicity
Pyrogenicity	USP Material Mediated Rabbit Pyrogen Test

Sterilization conditions have been validated according to ANSI / AAMI / ISO 11135, *Sterilization of Health Care Products-Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices* to provide a Sterility Assurance Level (SAL) of 10⁻⁶.

4.10.2 Design Verification (Bench-Top Testing)

The physical, mechanical and performance testing of the RePort™ Guide Catheter System demonstrates that the product is substantially equivalent to the currently marketed predicate devices. Design Verification testing was conducted to evaluate the physical and mechanical properties of the RePort™ Guide Catheter System. All studies were conducted in accordance with Reverse Medical Design Control procedures. All testing was performed on units which were twice sterilized and met all inspection criteria. Tests on the the RePort™ Guide Catheter System included:

Verification and Test Summary	Bench Tests	Result
	Dimensional and Visual Inspection	Met established criteria
	Guidewire Compatibility	Met established criteria
	Torque Response	Met established criteria
	Torque Strength	Met established criteria
	System Deployment Cycle Test	Met established criteria
	Kink Resistance	Met established criteria
	Flexibility Test	Met established criteria
	Tensile Strength	Met established criteria
	Catheter Leak Test (Liquid Leakage)	Met established criteria
	Catheter Leak Test (Air Leakage)	Met established criteria
	Dynamic Pressure Test	Met established criteria
	Static Burst Test	Met established criteria
	Aspiration Test	Met established criteria
	Hub Gauging	Met established criteria
	Corrosion Resistance	Met established criteria
	USP Particulate Test	Met established criteria
	Navigation and Accessibility Capabilities <i>in vitro</i>	Met established criteria
	<i>In vivo</i> Tests	Result
	System Deliverability, Compatibility, Visibility and Aspiration Performance	Met established criteria
	Acute histopathology of treated vessels	Met established criteria
	Biocompatibility testing	Met established criteria

The physical, mechanical and performance testing of the subject RePort™ Guide Catheter System demonstrate that the product is safe and effective for its labeled indications and is Substantially Equivalent to the currently marketed predicate devices.

4.11 Substantial Equivalence

The performance of the RePort™ Guide Catheter System in this submission demonstrates that the product is substantially equivalent to the performance of the predicate devices. The equivalence was shown through comparison of component materials and specification, performance and biocompatibility testing and sterilization validation.

The RePort™ Guide Catheter System is substantially equivalent in intended use, design, technology/principles of operation, materials and performance to the predicate devices. Differences between the devices do not raise any significant issues of safety or effectiveness.



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room –WO66-G609
Silver Spring, MD 20993-0002

APR 15 2011

Reverse Medical Corporation
c/o Mr. Mark Job
Regulatory Technology Services LLC
1394 25th Street NW
Buffalo, MN 55313

Re: K102418
Trade/Device Name: RePort Guide Catheter
Regulation Number: 21 CFR 870.1250
Regulation Name: Percutaneous Catheter
Regulatory Class: Class II
Product Code: DQY
Dated: March 31, 2011
Received: April 1, 2011

Dear Mr. Job:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the [Federal Register](#).


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21

CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Bram D. Zuckerman, M.D.
Director
Division of Cardiovascular Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

3. Indications for Use

510(k) Number (if known): K102418

Device Name: Reverse Medical RePort™ Guide Catheter System

Indications for Use:

The Reverse Medical RePort™ Guide Catheter System is indicated for the introduction of interventional/diagnostic devices into the peripheral, coronary and neuro vasculature.


Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over the Counter Use _____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)


(Division Sign-Off)
Division of Cardiovascular Devices

510(k) Number K102418