

510(k) Summary

This summary of 510(k) safety and effectiveness is being submitted in accordance with the requirements of the Safe Medical Device Act of 1990 and 21 CFR 807.92.

The assigned 510(k) number is: K102908.

1) Applicant Name and Address

Applicant: r2 Diagnostics, Inc.
Address: 1801 Commerce Drive
South Bend, IN 46628
Contact Person: Marc D. Goldford
Phone #: 574-288-4377
Fax #: 574-288-2272
Email: marc@r2diagnostics.com
Date of Preparation: 15 November 2010

2) Device Name(s)

Trade Name: NoFact IX
Classification Name: Plasma, Factor Deficient (21CFR 864.7290, Product Code GJT)

3) Predicate Device(s)

Stago IX Deficient Plasma (K933432)

4) Device Description(s)

NoFact IX Deficient Plasma is a human plasma immunodepleted of Factor IX and intended for the quantitative determination of Factor IX activity in citrated plasma from patients suspected of FIX deficiency. FIX activity is based on the activated partial thromboplastin time. For in vitro diagnostic use.

5) Intended Use(s)

NoFact IX Deficient Plasma is a human plasma immunodepleted of Factor IX and intended for the quantitative determination of Factor IX activity in citrated plasma from patients suspected of FIX deficiency. FIX activity is based on the activated partial thromboplastin time. For in vitro diagnostic use.

6) Technological Characteristic Summary

a. Comparison of the submitted kit and the predicate kit is summarized in the table below:

Comparison of submitted device NoFact IX to predicate device STA Deficient IX		
Similarities		
Item	Submitted Device	Predicate Device
Intended use	NoFact IX Deficient Plasma is a human plasma immunodepleted of Factor IX and intended for the quantitative determination of Factor IX activity in citrated plasma from patients suspected of FIX deficiency. FIX activity is based on the activated partial thromboplastin time. For in vitro diagnostic use.	STA Deficient IX is an immunodepleted human plasma intended for use in tests for the determination of factor IX activity in plasma by analyzers of the STA line suitable with this reagent.
Constituent material	Citrated human plasma immunodepleted of Factor IX.	Citrated human plasma immunodepleted of Factor IX.
Measurement principle	Diluted patient sample is mixed with factor IX deficient plasma and then tested with an APTT. In these conditions the clotting time of the mixture is dependent on the concentration of FIX in the patient sample.	Diluted patient sample is mixed with factor IX deficient plasma and then tested with an APTT. In these conditions the clotting time of the mixture is dependent on the concentration of FIX in the patient sample.
Format	Lyophilized plasma	Lyophilized plasma
Analyte being tested	Factor IX activity	Factor IX activity
Differences		
Item	Submitted Device	Predicate Device
Reconstituted Stability	2-8C: 8 hr RT: 4 hr	2-8C: not listed RT: not listed On-board STA Compact: 4 hrs

b. NoFact IX Deficient Plasma was compared to the predicate STA IX Deficient plasma using the Stago PTT-A FIX assay on the STA Compact. Plasma samples were assayed in parallel with the Stago PTT-A assay using the STA IX deficient plasma ("Stago results"), and also with the same assay but where the NoFact IX deficient plasma was substituted for the STA IX deficient plasma ("NoFact results"). A total of two hundred and thirty-three plasma samples in three laboratories

were analyzed. The linear regression equation of this comparison was: $y = 0.858 * x + 5.73$, with a coefficient of determination of 0.915.

The regression statistics by site were:

	All Labs N=233	Site 1 N=100	Site 2 N=80	Site 3 N=53
Slope	0.858	0.785	0.923	0.817
Intercept	5.729	6.168	5.328	11.206
r^2	0.915	0.977	0.907	0.830
r	0.956	0.988	0.988	0.911

- c. The Stago PTT-A IX assay was evaluated for precision with three separate lots of NoFact IX Deficient Plasma according to the CLSI guideline EP5-A2 "Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-2nd Edition":

Plasma	Mean FIX activity, %	% CV, Within-run (S-r)	% CV, Lot-to-Lot (S-lot)	% CV, Within-Device (S-device)
System N (Normal Control Plasma) N=240	109%	5.1%	1.1%	7.2%
System P (Abnormal Control Plasma) N=240	43%	4.6%	2.9%	7.7%
Low FIX pooled patient plasma N=120	14%	5.7%	2.2%	7.8%

NoFact IX Deficient Plasma provides acceptable precision.

NoFact IX Deficient Plasma is substantially equivalent to Stago IX Deficient Plasma.



Food and Drug Administration
10903 New Hampshire Avenue
Silver Spring, MD 20993

r² Diagnostics, Inc.
c/o Mr. Marc D. Goldford
Director, Research and Development
1801 Commerce Dr.
South Bend, IN 46628

DEC 19 2011

Re: k102908
Trade/Device Name: NoFact IX Deficient Plasma
Regulation Number: 21 CFR 864.7290
Regulation Name: Factor Deficiency Test
Regulatory Class: Class II
Product Code: GJT
Dated: December 7, 2011
Received: December 8, 2011

Dear Mr. Goldford:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into class II (Special Controls), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); and good manufacturing practice

requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Parts 801 and 809), please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (301) 796-5450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,

for 

Maria M. Chan, Ph.D.
Director
Division of Immunology and Hematology Devices
Office of *In Vitro* Diagnostic Device Evaluation and Safety
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K102908

Device Name: NoFact IX Deficient Plasma

Indications For Use: NoFact IX Deficient Plasma is a human plasma immunodepleted of Factor IX and intended for the quantitative determination of Factor IX activity in citrated plasma from patients suspected of FIX deficiency. FIX activity is based on the activated partial thromboplastin time. For in vitro diagnostic use.

Prescription Use X
(21 CFR Part 801 Subpart D)

And/Or Over the Counter Use _____
(21 CFR Part 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE; CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Device Evaluation and Safety (OIVD)



Division Sign-Off
Office of In Vitro Diagnostic Device
Evaluation and Safety

510(k) K102908