



#### Performance Testing & Verification Testing

- Kink Testing
- Tensile Testing
- High Pressure Injection Testing
- Hub Aspiration Testing
- Embolic Agent Infusion Compatibility Testing
- Package Integrity Testing
- Corrosion Testing
- Diagnostic Agent Compatibility Testing
- Trackability Testing
- Dimensional Inspection
- Coating Integrity
- Antegrade Flow Testing
- Infusion Efficiency
- Acute System Toxicity

#### Biocompatibility Testing

- **Pyrogenic**- Test is conducted based on USP, General Chapter, <151> Pyrogen Test. The procedure is recommended in ISO 10993-11
- **intra-cutaneous irritation** based on ISO 10993-10: **Toxicity** of the Paladin Catheter testing was based on International Organization for Standardization 10993-11
- **Hemolysis** was tested according to procedures based on ASTM F756, Standard Practices for Assessment of Hemolytic Properties of Materials and ISO 10993-4
- **Sensitization** was tested based on the requirements of ISO 10993-10
- **Particulate** – USP Standards
- **Cytotoxic effects** were tested following the guidelines of ISO 10993-5:
- **Complement System** was performed

**Performance/Safety:** A risk/hazard analysis was conducted according to EN ISO 14971 Medical Devices- Application of Risk management to medical devices. Performance characteristic for this indication for use were determined. It was then justified that the performance of the Surefire™ Infusion Catheter System is substantially equivalent to the performance and safety of the Terumo Radifocus Glidecath. A battery of tests were performed according to protocols based on the requirements of the following standards and was shown to meet the acceptance criteria that were determined to be applicable to the safety and efficacy of the device:

- ISO 10555-1 Sterile, single use intravascular catheters Part 1 General requirements.
- Surefire™ Infusion Catheter System Section 5-1 Premarket Notification Traditional 510(k) Submission 15 February 2011
- ISO 10555-2 Sterile, single use intravascular catheters Part 2 Angiographic catheters.
- ISO 10993-1 Biological Evaluation of medical Devices Part 1: Evaluation and Testing, and the FDA Modified ISO 10993 Test Profile
- ISO 11135-1 Medical Devices –Validation and Routine Control of Ethylene Oxide Sterilization.

**Additional Safety  
Information:**

Manufacturing controls include visual, functional, dimensional and sterility tests. Blood contacting materials were tested in accordance with the tests recommended in the FDA General program Memorandum. Use of International Standard ISO 10993, "Biological Evaluation of Medical Devices Part -1 Evaluation and testing".

**Discussion of  
Animal Data:**

The animal study data is submitted in this 510(k) and is on file at Surefire Medical.

**Substantial  
Equivalence:**

The Surefire™ Infusion Catheter System is substantially equivalent in intended use, design, technology/principles of operation to the predicates. The Guide sheath is substantially equivalent to the Terumo Radifocus GLIDECATH, cleared under K090040. The Microcatheter is substantially equivalent to the EmboCath Plus Infusion Microcatheter by BioSphere Medical cleared under K062126. Differences between the devices do not raise any significant issues of safety or effectiveness.

**Submitter  
Information:**

Prepared by: Cheryl Hastings  
VP Clinical and Regulatory Affairs

Prepared for: Surefire Medical, Inc.  
8601 Turnpike Drive  
Suite 206  
Westminster, CO 80031

Date: June 16, 2011



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room –WO66-G609  
Silver Spring, MD 20993-0002

Surefire Medical, Inc.  
c/o Ms. Cheryl Hastings  
VP Clinical and Regulatory Affairs  
8601 Turnpike Dr., Suite 206  
Westminster, CO 80031

JUN 24 2011

Re: K110459

Trade/Device Name: Surefire™ Infusion Catheter System  
Regulation Number: 21 CFR 870.1200  
Regulation Name: Diagnostic intravascular catheter  
Regulatory Class: Class II  
Product Code: DQO  
Dated: June 17, 2011  
Received: June 20, 2011

Dear Ms. Hastings:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical


Page 2 – Ms. Cheryl Hastings

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucml15809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read 'Bram D. Zuckerman', with a long horizontal line extending to the right.

Bram D. Zuckerman, M.D.  
Director  
Division of Cardiovascular Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**Indication for Use**

510(k) Number (if known): K110459

Device Name: Surefire™ Infusion Catheter System

Indication for Use: The SUREFIRE™ INFUSION CATHETER SYSTEM is intended for use in angiographic procedures. It delivers radiopaque media and therapeutic agents to selected sites in the peripheral vascular system.

Prescription Use  X  AND/OR Over-The-counter Use \_\_\_\_\_  
(part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH,

  
\_\_\_\_\_  
(Division Sign-Off)  
Division of Cardiovascular Devices  
510(k) Number  K110459

**Indication for Use**

510(k) Number (if known): K110459

Device Name: Surefire™ Infusion Catheter System

Indication for Use: The SUREFIRE™ INFUSION CATHETER SYSTEM is intended for use in angiographic procedures. It delivers radiopaque media and therapeutic agents to selected sites in the peripheral vascular system.

Prescription Use  X  AND/OR Over-The-counter Use \_\_\_\_\_  
(part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH,



(Division Sign-Off)  
Division of Cardiovascular Devices  
510(k) Number  K110459