

**510(k) Summary for AccuTrac Laser Fiber**

MAR 28 2011

**A. Sponsor**

Boston Scientific Corporation  
Urology and Women's Health Division  
100 Boston Scientific Way  
Marlborough, MA 01756

**B. Contact**

Lauren Anderson  
Senior Specialist, Regulatory Affairs  
508-683-4707  
lauren.anderson@bsci.com

or

Nichole Riek  
Manager, Regulatory Affairs  
508-683-4175  
riekn@bsci.com

**C. Device Name**

Trade name: AccuTrac™  
Common/usual name: Laser Instrument, Surgical, Powered  
Classification Name: GEX – Laser surgical instrument for use in general and  
plastic surgery and in dermatology  
21 CFR 878.4810, Class II

**D. Predicate Device(s)**

Trade name: AccuMax™ (Straight Fire Holmium Laser Fiber)  
Common/usual name: Laser Instrument, Surgical, Powered  
Classification Name: GEX – Laser surgical instrument for use in general and  
plastic surgery and in dermatology  
21 CFR 878.4810, Class II  
Premarket Notification: Boston Scientific, K082928

And

Trade name: Flexiva™ (Modified Straight Fire Laser Fiber)  
Common/usual name: Laser Instrument, Surgical, Powered  
Classification Name: GEX – Laser surgical instrument for use in general and  
plastic surgery and in dermatology  
21 CFR 878.4810, Class II  
Premarket Notification: Boston Scientific, K100078

**E. Device Description**

The AccuTrac Laser Fibers are fiber optic laser energy delivery devices  
consisting of a SMA-905 connector, strain relief, and a silica core fiber jacketed

with ethylene tetrafluoroethylene (ETFE). The AccuTrac fibers are equipped with a polished and reinforced ball-shaped TracTip. These fibers may be used in a variety of laser based surgical cases as an integral part of laser systems.

For use with Ho:YAG laser systems with a standard SMA-905 connector that have been cleared for surgical use. Recommended Ho:YAG lasers are Dornier and New Star. Please refer to the laser system User Manual for complete information regarding applications, contraindications, precautions and warnings.

#### **F. Intended Use**

The AccuTrac Laser Fibers are intended for use in laser-based surgical applications, including, but not limited to endoscopic, laparoscopic and open surgical procedures involving vaporization, ablation and fragmentation of calculi (urinary and biliary) and surgical procedures involving vaporization, ablation, coagulation, hemostasis, excision, resection and incision of soft and cartilaginous tissue.

The AccuTrac Laser Fibers are designed for use with holmium (Ho:YAG) lasers with a standard SMA 905 connector that have been cleared for surgical use.

#### **G. Technological Characteristics**

The AccuTrac Laser Fiber has the same technological characteristics (i.e. SMA 905 connector, length of fiber optic cable, and strain relief) as the predicate devices. It is equipped with a polished and reinforced ball-shaped output tip, whereas the predicate devices have a polished flat tip.

#### **H. Substantial Equivalence**

A direct comparison of key characteristics demonstrates that the proposed laser fiber is substantially equivalent to the predicate devices in terms of intended use, technological characteristics, and performance characteristics. The AccuTrac Laser Fiber is as safe, as effective, and performs as well as the predicate devices.

#### **I. Performance Testing (Bench and User Evaluation)**

Boston Scientific has conducted performance testing with samples aged at T=0 and T=13 months accelerated aging in support of the distal ball tip design change. The following testing was completed to evaluate the effects of the design change:

- Tip Fracture Resistance
- Scope Testing
- Power Rating/Output Efficiency Testing
- Aiming Beam Testing

The results of the performance testing demonstrate equivalence of the AccuTrac to the predicate AccuMax and Flexiva laser fibers. The AccuTrac fibers are considered safe and effective for their intended use.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room W-066-0609  
Silver Spring, MD 20993-0002

Boston Scientific Corporation  
% Ms. Lauren B. Anderson  
Senior Regulatory Affairs Specialist  
100 Boston Scientific Way, M21  
Marlborough, Massachusetts 01752

MAR 28 2011

Re: K110686

Trade/Device Name: AccuTrac™ Laser Fiber

Regulation Number: 21 CFR 878.4810

Regulation Name: Laser surgical instrument for use in general and plastic surgery and in dermatology

Regulatory Class: II

Product Code: GEX

Dated: March 24, 2011

Received: March 25, 2011

Dear Ms. Anderson:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act


Page 2 - Ms. Lauren B. Anderson

or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink, appearing to read "Mark N. Melkerson" with a stylized flourish at the end.

Mark N. Melkerson  
Director  
Division of Surgical, Orthopedic  
and Restorative Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**Indications for Use Statement**

**510(k) Number** To be determined.

**Device Name** AccuTrac™ Laser Fiber

**Indications For Use**

The AccuTrac Laser Fibers are intended for use in laser-based surgical applications, including, but not limited to, endoscopic, laparoscopic and open surgical procedures involving vaporization, ablation and fragmentation of calculi (urinary and biliary) and surgical procedures involving vaporization, ablation, coagulation, hemostasis, excision, resection and incision of soft and cartilaginous tissue. The AccuTrac Laser Fibers are designed for use with Ho:YAG lasers with a standard SMA-905 connector that have been cleared for surgical use.

Prescription Use  X   
(21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

*Nigel Dyer for me*  
(Division Sign-Off)  
Division of Surgical, Orthopedic,  
and Restorative Devices

510(k) Number  K110686