

## BioPlex 2200 MMRV IgG 510(k) Summary

510(k) Number K 111072

Date Prepared: March 25, 2011

### ***Introduction***

Bio-Rad Laboratories hereby submits this Special 510(k) in accordance with the requirements of SMDA 1990 and 21 CFR 807.92. This summary of 510(k) safety and effectiveness information provides detail as a basis for a determination of substantial equivalence for the BioPlex 2200 MMRV IgG.

### ***Submitter name, address and contact***

Submitter	Contact Person
Bio-Rad Laboratories, Inc BioPlex Division 5500 E. Second Street Benicia, CA 94510	Juang Wang Regulatory Affairs Representative Phone: (510)741-4609 Fax: (510)741-4650

### ***Device name and Classification***

Product Trade Name	BioPlex 2200 MMRV IgG on the BioPlex 2200 Multi-Analyte Detection System
Common Name	Multi-Analyte Detection System – MMRV IgG
Classification name	Multiplex immunoassay for measles virus, mumps virus, rubella, and varicella-zoster virus
Device Class	Class II
Classification Panel	Microbiology
Registration Number	886.3510
Product Code	OPL

### ***Legally Marketed Predicate Device***

BioPlex 2200 MMRV IgG, k091616

### ***Intended Use/Indications For Use***

The BioPlex<sup>®</sup> 2200 MMRV IgG kit is a multiplex flow immunoassay intended for the qualitative detection of IgG antibodies to Measles, Mumps, Rubella and Varicella-zoster virus (VZV) in human serum and EDTA or heparinized plasma.

The BioPlex 2200 MMRV IgG kit is intended for use with the Bio-Rad BioPlex 2200 System.

This kit is intended as an aid in the determination of serological status to Measles, Mumps, Rubella, and VZV. This kit is not intended for use in screening blood or plasma donors.

The performance of this assay has not been established for use in neonatal, pediatrics and immunocompromised patients, or for use at point of care facilities.

### ***Device Description***

The BioPlex 2200 MMRV IgG kit uses multiplex flow immunoassay for simultaneous detection and identification of many antibodies in a single tube. Four (4) different populations of magnetic dyed beads are coated with antigens to identify the presence of IgG class antibodies against Measles, Mumps, Rubella and Varicella-zoster. The BioPlex 2200 System combines an aliquot of patient sample, sample diluent, and bead set reagent into a reaction vessel. The mixture is incubated at 37°C. After a wash cycle, anti-human IgG antibody, conjugated to phycoerythrin (PE), is added to the dyed beads and this mixture is incubated at 37°C. The excess conjugate is removed in another wash cycle, and the beads are re-suspended in wash buffer. The bead mixture then passes through the detector.

The identity of the dyed beads is determined by the fluorescence of the dyes, and the amount of antibody captured by the antigen is determined by the fluorescence of the attached PE. Raw data is calculated in relative fluorescence intensity (RFI). Three additional control beads, an Internal Standard Bead (ISB), a Serum Verification Bead (SVB), and a Reagent Blank Bead (RBB) are present in each reaction mixture to verify detector response, the addition of serum to the reaction vessel, and the absence of significant non-specific binding in serum.

The instrument is calibrated using a set of three (3) distinct calibrator vials, supplied separately by Bio-Rad Laboratories.

### ***Similarities and Differences***

#### Similarities

Feature	Modified Device
Intended Use/Indications For Use	No Change
Kit components	No Change
Technical Specifications	No Change
Fundamental Scientific Technology	No Change

Differences

The only difference of the BioPlex 2200 MMRV IgG is to modify QC testing from each reagent pack to once per day as stated in the Instructions For Use (IFU) of the BioPlex 2200 MMRV IgG reagent kit.

Feature	Modified	Predicate
Frequency of Reagent Pack QC Testing	QC once per day or per new reagent pack lot	QC once per pack and per day

***Summary of Design Control Activities***

A Failure Mode and Effect Analysis (FMEA) was used to facilitate, capture, and quantify potential impacts of Low Signal Pack (LSP) occurrence. The Risk Priority Number (RPN) is a quantitative measure of the combined effects of severity, occurrence, and detection of potential risks. Specific mitigations are recommended that may include changes to the design or formulation if the RPN score exceeds a chosen threshold.

The Design Control Activities include Risk Analysis method to identify the verification and validation activities required, test used and acceptance criteria.

Based on the conclusion of the risk management report, the modified QC procedure fulfills the requirements of the specifications of the design control process. Therefore, the performance of the modified QC test frequency is substantially equivalent to the current cleared kit.



Food and Drug Administration  
10903 New Hampshire Avenue  
Silver Spring, MD 20993

Bio-Rad Laboratories, Inc.  
c/o Mr. Juang Wang  
Regulatory Affairs Representative  
BioPlex Division  
5500 E. Second Street  
Benicia, CA 94510

AUG 23 2011

Re: K111072

Trade/Device Name: BioPlex 2200 MMRV IgG on the BioPlex 2200 Multi-Analyte  
Detection System  
Regulation Number: 21 CFR 866.3510  
Regulation Name: Rubella Virus Serological Reagents  
Regulatory Class: Class II  
Product Code: OPL, LJB, LJY, LFY  
Dated: August 8, 2011  
Received: August 16, 2011

Dear Mr. Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into class II (Special Controls), it may be subject to such additional controls. Existing major regulations affecting your device can be found in Title 21, Code of Federal Regulations (CFR), Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); and good manufacturing practice requirements as set forth in the


Page 2-Mr. Juang Wang

quality systems (QS) regulation (21 CFR Part 820). This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Parts 801 and 809), please contact the Office of *In Vitro* Diagnostic Device Evaluation and Safety at (301) 796-5450. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Sally A. Hojvat, M.Sc., Ph.D.  
Director  
Division of Microbiology Devices  
Office of In Vitro Diagnostic Devices  
Evaluation and Safety  
Center for Devices and Radiological Health

Enclosure

# Indication(s) for use

510(k) Number (if known): k111072

Device Name: BioPlex<sup>®</sup> 2200 MMRV IgG Kit on the BioPlex<sup>®</sup> 2200 Multi Analyte Detection System  
BioPlex<sup>®</sup> 2200 MMRV IgG Calibrator Set  
BioPlex<sup>®</sup> 2200 MMRV IgG Control Set

## Indications for Use:

### The BioPlex<sup>®</sup> 2200 MMRV IgG kit

The BioPlex<sup>®</sup> 2200 MMRV IgG kit is a multiplex flow immunoassay intended for the qualitative detection of IgG antibodies to Measles, Mumps, Rubella and Varicella-zoster virus (VZV) in human serum and EDTA or heparinized plasma.

The BioPlex 2200 MMRV IgG kit is intended for use with the Bio-Rad BioPlex 2200 System.

This kit is intended as an aid in the determination of serological status to Measles, Mumps, Rubella, and VZV. This kit is not intended for use in screening blood or plasma donors.

The performance of this assay has not been established for use in neonatal, pediatrics and immunocompromised patients, or for use at point of care facilities.

### BioPlex<sup>®</sup> 2200 MMRV IgG Calibrator Set

The BioPlex<sup>®</sup> 2200 MMRV IgG Calibrator Set is intended for the calibration of the BioPlex<sup>®</sup> 2200 MMRV IgG Reagent Pack.

### BioPlex<sup>®</sup> 2200 MMRV IgG Control Set

The BioPlex<sup>®</sup> 2200 MMRV IgG Control Set is intended for use as an assayed quality control to monitor the overall performance of the BioPlex 2200 instrument and BioPlex 2200 MMRV IgG Reagent Pack in the clinical laboratory. The performance of the BioPlex 2200 MMRV IgG Control Set has not been established with any other Measles, Mumps, Rubella or VZV IgG assays.

For Prescription Use Only X  
(Part 21 CFR 801 Subpart D)

AND/OR Over-the-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(Please do not write below this line-Continue on another page if needed)

---

Concurrence of CDRH, Office of In Vitro Diagnostic Devices (OIVD)

~~Division Sign-Off~~  
~~Office of In Vitro Diagnostic Device~~  
~~Evaluation and Safety~~

~~510(k) \_\_\_\_\_~~

*Tamara Feldman*  
Division Sign-Off

Office of In Vitro Diagnostic  
Device Evaluation and Safety

510(k) K111072