

K111808

SEP 21 2011



Cardinal Health

**510(K) SUMMARY OF SAFETY AND EFFECTIVENESS
for STERILE POLYISOPRENE POWDER-FREE BLUE SURGICAL GLOVES
WITH NEU-THERA® COATING**

(A summary of safety and effectiveness information in accordance with the requirements of 21 CFR 807.92)

Applicant: Cardinal Health
1430 Waukegan Road
McGaw Park, IL 60085

Establishment Registration

Number: 1423537

Regulatory Affairs

Contact: Tatyana Bogdan, RAC
Telephone: 847-887-2325
Fax: 847-887-2717
E-mail: tatyana.bogdan-curvin@cardinalhealth.com

Summary Prepared: June 14, 2011

Trade Name: Sterile Polyisoprene Powder-Free Blue Surgical Gloves with Neu-Thera Coating

Common Name: Surgeon's Gloves

Classification Name: Surgeon's Gloves

Classification Panel: General and Plastic Surgery

Regulation: 21 CFR 878.4460

Product Code(s): KGO

Legally marketed device(s)

to which equivalence is claimed: Esteem® Blue Sterile Polyisoprene Powder-Free Surgical Gloves with Neu-Thera Coating (containing Chitosan, Provitamin B5, Gluconolactone and Glycerol) cleared under 510(k) K042574 (product code KGO)

Reason for 510(k)

Submission: Modification of a legally marketed device

Device Description: The proposed device is a disposable device intended for over the counter use and is provided sterile. It is not made with natural rubber latex. Instead, it is formulated from polyisoprene, which is synthetic rubber latex. This sterile polyisoprene powder-free surgical glove is manufactured using exact same material used in the currently cleared device, Esteem Blue glove (K042574) that has been legally marketed by Cardinal Health for many years. The glove is coated with emollient coating containing Glycerol, Gluconolactone, D-Sorbitol and

Provitamin-B. The glove is manufactured using molds that feature anti-slip finish, independent thumb, and tapered mechanically locking cuffs to help reduce cuff roll down. It is offered powder-free and sterile.

Intended Use:

This powder-free surgeon's glove is a disposable device made of synthetic rubber intended to be worn by operating room personnel to protect a surgical wound from contamination.

Summary of the technological characteristics of the device compared to the predicate device		
Characteristic	Modified Device Sterile Polyisoprene Powder-Free Blue Surgical Gloves with Neu-Thera Coating	Original (Predicate) Esteem Blue Sterile Polyisoprene Powder-Free Surgical Glove with Neu-Thera Coating (K042574)
Material Composition	Synthetic Rubber Latex - Polyisoprene	Synthetic Rubber Latex - Polyisoprene
Design	Single Use Sterile Powder-free Hand Specific Independent Thumb Beaded Cuff Lubricated	Single Use Sterile Powder-free Hand Specific Independent Thumb Beaded Cuff Lubricated
Coating Contents	Provitamin B, Gluconolactone, D-Sorbitol and Glycerol	Chitosan, Provitamin B, Gluconolactone, D-Sorbitol and Glycerol
Intended Use/ Indications for Use	Powder-Free Surgeon's Glove	Powder-Free Surgeon's Glove
Dimensions & Physical Properties	Meets ASTM D3577	Meets ASTM D3577
Freedom from Holes	AQL meets 21CFR 800.20 & ASTM D3577 requirements	AQL meets 21CFR 800.20 & ASTM D3577 requirements
Powder Residual	Meets requirements of ≤ 2.0 mg/glove for Powder-Free designation per ASTM D3577	Meets requirements of ≤ 2.0 mg/glove for Powder-Free designation per ASTM D3577

PERFORMANCE DATA			
SUMMARY OF NON-CLINICAL TESTS CONDUCTED FOR DETERMINATION OF SUBSTANTIAL EQUIVALENCE*			
Performance Test Summary-New Device			
Characteristic	Standard/Test/FDA Guidance	Results Summary	
Biocompatibility: Primary Skin Irritation Guinea Pig Maximization Physical Characteristics: Dimensions Physical Properties Freedom from Holes Powder Residual	ISO 10993-10 ISO 10993-10 ASTM D3577 ASTM D3577 21 CFR 800.20 & ASTM D3577 ASTM D3577 tested using ASTM standard D6124	Gloves are non-irritating. Gloves do not display any potential for sensitization. Meet requirements Meet requirements for rubber surgical gloves Tested in accordance with ASTM D5151 with acceptable results Gloves meet powder level requirements for "Powder-Free" designation per ASTM D3577. Results generated values < 2mg of residual powder per glove.	
Comparative Performance Information Summary			
Characteristic	Requirement	New Device	Predicate Device
Biocompatibility:	ISO 10993-1	Meets requirements	Meets requirements
Primary Skin Irritation	ISO 10993-10	Pass	Pass
Guinea Pig Maximization	ISO 10993-10	Pass	Pass
Dimensions	ASTM D3577	Meets requirements	Meets requirements
Physical Properties	ASTM D3577	Meets requirements	Meets requirements
Freedom from Holes	21CFR800.20, ASTM D3577	Meets requirements	Meets requirements
Powder Residual	ASTM D3577	Meets requirements	Meets requirements
SUMMARY OF CLINICAL TESTS CONDUCTED FOR DETERMINATION OF SUBSTANTIAL EQUIVALENCE AND/OR OF CLINICAL INFORMATION			
Clinical data is not required.			
CONCLUSIONS DRAWN FROM NON-CLINICAL AND CLINICAL DATA			
Non-clinical data demonstrates that Sterile Polyisoprene Powder-Free Blue Surgical Gloves with Neu-Thera Coating meet the technological characteristics of ASTM D3577 standard, and are as safe, as effective, and performed as well as the legally marketed devices identified in this summary.			



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room -W066-G609
Silver Spring, MD 20993-0002

Ms. Tatyana Bogdan
Regulatory Affairs Manager
Cardinal Health-Medical Products and Services
1430 Waukegan Road
McGaw Park, Illinois 60085

SEP 21 2011

Re: K111868
Trade/Device Name: Esteem[®] Blue with Neu-Thera[®]
Regulation Number: 21 CFR 878.4460
Regulation Name: Surgeon's Glove
Regulatory Class: I
Product Code: KGO
Dated: July 26, 2011
Received: July 27, 2011

Dear Ms. Bogdan:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Anthony D. Watson, B.S., M.S., M.B.A.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure



Indications for Use

510(k) Number (if known): K111868

Device Name: Esteem® Blue with Neu-Thera®

Device description: Sterile Polyisoprene Powder-Free Blue Surgical Gloves with Neu-Thera Coating (containing Glycerol, Gluconolactone, D-Sorbitol and Provitamin-B)

Indications for Use: This powder-free surgeon's glove is a disposable device made of synthetic rubber intended to be worn by operating room personnel to protect a surgical wound from contamination.

Prescription Use _____ AND/OR Over-The-Counter Use X
(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Elizabeth F. Lawrence-Willis

(Division Sign-Off)

Page ___ of ___

Division of Anesthesiology, General Hospital
Infection Control, Dental Devices