

**510(k) Summary (as required by 21 CFR 807.92(c))**    **K112090**

**Manufacturer Name and Address**

Bahadir USA Corp.  
275 West Hoffman Avenue  
Lindenhurst, NY 11757  
Contact: Ismail Kilic  
Tel: 631-608-8522

JUN - 7 2012

**Submitter / Contact Person**

Bahadir USA Corp.  
275 West Hoffman Avenue  
Lindenhurst, NY 11757  
Contact: Ismail Kilic  
Tel: 631-608-8522  
Email: [info@bahadirusa.com](mailto:info@bahadirusa.com)

**Date Prepared**

Revised July 8, 2011  
Revised February 27, 2012  
Revised April 2, 2012  
Revised May 22, 2012

**Name of Device**

BAHADIR STERILIZATION TRAYS (INCLUDING WIDE BODY)

**Classification Name**

Sterilization Wraps, Trays, Containers  
Class 2 – KCT

**Predicate Device**

Bahadir – K102146

**Description of the Device  
Intended Use**

Intended for use in hospitals and healthcare facilities to contain other medical devices that are to be sterilized and to allow sterilization of the enclosed medical devices using steam sterilizers. Sterilized devices may be stored and transferred in the container.

The devices included in this submission are to be used with a pre-vacuum cycle of 4 minutes at 270 degrees F with a dry time of 20 minutes.

The devices (natural aluminum color containers and lids) subject to this submission are as follows:

Full Size Y111.10W  
Full Size Y111.13W  
Full Size Y111.15W  
Full Size Y111.20W  
Full Size Y111.26W  
Full Size Y110.62W  
Full Size Y110.68W  
Full Size Y111.62W  
Full Size Y111.68W

¾ Size Y211.10W  
¾ Size Y211.13W  
¾ Size Y211.15W

½ Size Y311.10W  
½ Size Y311.13W  
½ Size Y311.15W  
½ Size Y311.20W  
½ Size Y311.26W

**(Full Size Max Load Wt. 24.75 lbs  
including basket & instruments)**  
Lumen Size min 1mm diameter  
X 300mm length

**(¾ Max Load Wt. 12.5 lbs  
including basket & instruments)**  
Lumen Size min 1mm diameter  
X 200mm length

**(½ Size Max Load Wt. 9.25 lbs  
including basket & instruments)**  
Lumen size min 1mm diameter  
X 200mm length

## DESCRIPTION

The devices included in this submission are to be used with a pre-vacuum cycle of 4 minutes at 270 degrees F with a dry time of 20 minutes.

The devices subject to this submission are as follows:

Full Size Y111.10W	¾ Size Y211.10W	½ Size Y311.10W
Full Size Y111.13W	¾ Size Y211.13W	½ Size Y311.13W
Full Size Y111.15W	¾ Size Y211.15W	½ Size Y311.15W
Full Size Y111.20W		½ Size Y311.20W
Full Size Y111.26W		½ Size Y311.26W
Full Size Y110.62W		
Full Size Y110.68W		
Full Size Y111.62W		
Full Size Y111.68W		

The devices are composed of anodized aluminum (this does not include colored anodized lids).

### Comparison to Predicate Devices

The containers are the identical design, materials, sterilization cycle, and indications for use. The only difference is the dimensions (sizes).

Feature	Bahadir	Bahadir K102146
Intended Use	Intended for use in hospitals and healthcare facilities to contain other medical devices that are to be sterilized. Containers allow sterilization of the enclosed medical devices, including surfaces and lumens, using high vacuum steam sterilizers for 270 F for 4 minutes with 20 minutes (minimum) dry time.	Intended for use in hospitals and healthcare facilities to contain other medical devices that are to be sterilized. Containers allow sterilization of the enclosed medical devices, including surfaces and lumens, using high vacuum steam sterilizers for 270 F for 4 minutes with 20 minutes (minimum) dry time.
Material	Aluminum Alloy, stainless steel handles, silicone seal, paper filter	Aluminum Alloy, stainless steel handles, silicone seal, paper filter
Filter	Paper Filter	Paper Filter
Sterilization Method	Steam	Steam
Configurations / Dimensions	Includes: Full Size Y111.10W (580x280x100mm) Full Size Y111.13W (580x280x135mm) Full Size Y111.15W (580x280x150mm) Full Size Y111.20W (580x280x200mm) Full Size Y111.26W (580x280x260mm)	Includes 1/2 size, 3/4 size, and full size  <u>Full size</u> 580mm x 280mm x 200mm 580mm x 280mm x 150mm 580mm x 280mm x 135mm 580mm x 280mm x 100mm

	<p>Full Size Y110.62W (600x400x120mm)</p> <p>Full Size Y110.68W (600x400x180mm)</p> <p>Full Size Y111.62W (600x400x120mm)</p> <p>Full Size Y111.68W (600x400x180mm)</p> <p>¾ Size Y211.10W (465x280x100mm)</p> <p>¾ Size Y211.13W (465x280x135mm)</p> <p>¾ Size Y211.15W (465x280x150mm)</p> <p>½ Size Y311.10W (285x280x100mm)</p> <p>½ Size Y311.13W (285x280x135mm)</p> <p>½ Size Y311.15W (285x280x150mm)</p> <p>½ Size Y311.20W (285x280x200mm)</p> <p>½ Size Y311.26W (285x280x260mm)</p>	<p><u>¾ Size</u></p> <p>465mm x 280mm x 150mm</p> <p>465mm x 280mm x 135mm</p> <p>465mm x 280mm x 100mm</p> <p><u>½ size</u></p> <p>285mm x 280mm x 200mm</p> <p>285mm x 280mm x 150mm</p> <p>285mm x 280mm x 135mm</p> <p>285mm x 280mm x 100mm</p>
Perforation	The units included in this submission include perforated lids and both non-perforated and perforated bottoms.	The units included in this submission include perforated lids and non-perforated bottoms.
Air Permeance	Permits entry of sterilization agent and prevents microbial migration during storage.	Permits entry of sterilization agent and prevents microbial migration during storage.
Intended to be re-used	Yes	Yes
Sealed	Yes	Yes
Conformance to AAMI ST 77	Yes	Yes
Gasket	Silicone Based (gray color)	Silicone Based (gray color)
Handles	Wider handles with silicone holder	Wider handles with silicone holder
Cycle Parameters	270 F at 4 minutes 20 minutes drying	270 F at 4 minutes 20 minutes drying
Use with cannulized instruments	Yes	No

The differences are as follows:

- The only difference is the dimensions (wide body size).
- Includes non-perforated and perforated bottom
- Can be used with cannulized instruments

### **Non-Clinical Tests Performed**

The subject devices were subjected to sterility testing and performance testing.

### **Summary**

In conclusion we believe the devices are substantially equivalent to the predicate devices and do not introduce new issues of safety or effectiveness.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room -WO66-G609  
Silver Spring, MD 20993-0002

Mr. Ismail Kilic  
President  
Bahadir USA Corporation  
275 West Hoffman Avenue  
Lindenhurst, New York 11757

JUN - 7 2012

Re: K112090  
Trade/Device Name: Bahadir USA Sterilization Trays  
Regulation Number: 21 CFR 880.6850  
Regulation Name: Sterilization Wrap  
Regulatory Class: II  
Product Code: KCT  
Dated: May 29, 2012  
Received: May 31, 2012

Dear Mr. Kilic:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Anthony D. Watson, B.S., M.S., M.B.A.  
Director  
Division of Anesthesiology, General Hospital,  
Infection Control and Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): **112090**

Device Name: Bahadir USA Sterilization Trays

Indications For Use:

Intended for use in hospitals and healthcare facilities to contain other medical devices that are to be sterilized and to allow sterilization of the enclosed medical devices using steam sterilizers. Sterilized devices may be stored and transferred in the container.

The devices included in this submission are to be used with a pre-vacuum cycle of 4 minutes at 270 degrees F with a dry time of 20 minutes.

The devices (natural aluminum color containers and lids) subject to this submission are as follows:

Full Size Y111.10W  
Full Size Y111.13W  
Full Size Y111.15W  
Full Size Y111.20W  
Full Size Y111.26W  
Full Size Y110.62W  
Full Size Y110.68W  
Full Size Y111.62W  
Full Size Y111.68W

¾ Size Y211.10W  
¾ Size Y211.13W  
¾ Size Y211.15W

½ Size Y311.10W  
½ Size Y311.13W  
½ Size Y311.15W  
½ Size Y311.20W  
½ Size Y311.26W

(Full Size Max Load Wt. 24.75 lbs  
including basket & instruments)  
Lumen Size min 1mm diameter  
X 300mm length

(¾ Max Load Wt. 12.5 lbs  
including basket & instruments)  
Lumen Size min 1mm diameter  
X 200mm length

(½ Size Max Load Wt. 9.25 lbs  
including basket & instruments)  
Lumen size min 1mm diameter  
X 200mm length

Prescription Use \_\_\_\_\_  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use  X   
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
(Division Sign-Off)

Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

Page 1 of 1

510(k) Number:  K112090