

JAN - 5 2012

510(k) Summary of Safety and Effectiveness

Proprietary Name: Stryker Universal Neuro 3 System

Common Name: Neuro Plating System

Classification Name and Reference: Preformed alterable cranioplasty plate
21 CFR §882.5320
Burr hole cover
21 CFR §882.5250
Cranioplasty plate fastener
21 CFR §882.5360

Proposed Regulatory Class: Class II

Product Codes: GWO – Preformed alterable cranioplasty plate
GXR – Burr hole cover
HBW – Cranioplasty plate fastener

For Information contact: Rob Yamashita
Sr. Regulatory Compliance Analyst
Stryker Craniomaxillofacial
750 Trade Centre Way, Suite 200
Portage, MI 49002
Phone: (269) 389-4258 / Fax: (269) 329-2235
rob.yamashita@stryker.com

Date Prepared: September 1, 2011

Indications for Use / Intended Use

The Stryker Universal Neuro 3 System is intended for reconstruction, stabilization and/or rigid fixation of non load-bearing areas subsequent to craniotomy, craniectomy and cranial fractures in adults and adolescents (age 12 and higher).

Contraindications

The Stryker Universal Neuro 3 System is contraindicated for the following:

- Use of plates in non-reducible and unstable fractures
- Patients with active local infections
- Patients with metal allergies and foreign body sensitivity

- Potentially non-compliant patients with mental or neurological conditions who are unwilling or incapable of following postoperative care instructions
- Patients with limited blood supply to, or insufficient quality of, bone
- Use of products in cases where the fixation of the products could result in their peripheral edge coming into contact with the dura mater
- Screws coming in contact with the dura mater
- Use of implants adjacent to developing paranasal sinuses

Technological Characteristics

The Stryker Universal Neuro 3 System is designed for a wide selection of solutions for cranial fixation. It consists of an implant module (a storage module that contains various versions and shapes of plates and screws) for the respective anatomical and indicated areas.

The low profile plates of the Stryker Universal Neuro 3 System provide rigid fixation of cranial flaps with decreased palpability. There is a comprehensive selection of burr hole covers, straight plates, gap plates, 3D-plates, shunt plates, and box plates to provide many fixation options. The malleable plates can be easily contoured by hand without instruments. The pre-shaped skull-base plates provide covers for standard craniectomies, obviating the need to cut or trim mesh.

Performance Data

Materials used for the Stryker Universal Neuro 3 System are the same as the predicate devices. This includes all three product codes (GWO, GXR and HBW). The titanium materials used for manufacturing of the Stryker Universal Neuro 3 implants are rated to be biocompatible according to ISO 10993-1. Cytotoxicity testing was performed according to ISO 10993-1, 10993-5, 10993-12, and 10993-18. The corrosion resistance of all Neuro 3 screws, plates and meshes were demonstrated.

The bending stability of the Universal Neuro 3 plates (product code GWO) and burr hole covers (product code GXR) were tested by following ASTM F 382-99. The Lerch test was passed by all plates.

For the screws (product code HBW), the testing was performed via ASTM F 543 – Standard Specification and Test Methods for Metallic Medical Bone Screws, 2007. Torque, depth and angle were measured. The screws passed the automated insertion test.

Additionally, we tested the fixation stability of our screws with pull out safety testing and the retention force between the screw and the screwdriver blade also utilizing ASTM F 543. All acceptance criteria were met.

Substantial Equivalence

The Stryker Universal Neuro 3 System has been verified and validated according to Stryker procedures for product design and development. The validation proves the safety and effectiveness of the system. The information provided by Stryker in this 510(k) application was found to be substantially equivalent with these predicate devices:

- Stryker Universal Neuro 2 System (Stryker, K031659)
- Stryker Micro Dynamic Mesh (Stryker, K983528)
- Synthes Neuro Plate and Screw System (Synthes, K022012)
- KLS-Martin Micro Osteosynthesis System (KLS-Martin, K944561/K944565)

The Stryker Universal Neuro 3 System, as stated above, consists of devices with three product codes: GWO (plates), GXR (burr hole covers), and HBW (screws). It has the same material composition and operating principles as its predicates mentioned above. The intended use is similar to the predicate systems with the only difference being the inclusion of adolescent use. Further, there may be slight differences in dimensions and shapes between the Stryker Universal Neuro 3 System and the predicate devices; however, the information provided in this 510(k) proves substantial equivalence to the predicate devices.









Specifically speaking, our Universal Neuro 3 burr hole covers and screws are substantially equivalent to Stryker Universal Neuro 2 System (K031659) and

Synthes Neuro Plate and Screw System (K022012). Our Universal Neuro 3 plates are substantially equivalent to the same two predicates named above (K031659 and K022012) plus Stryker Micro Dynamic Mesh (K983528). Also, we included the KLS-Martin Micro Osteosynthesis System (K944561/K944565) as a predicate because of its indicated pediatric use.


Below are three tables that outline our Universal Neuro 3 substantial equivalence.

L

SUBSTANTIAL EQUIVALENCE TABLE: PLATES (Product Code GWO)

	Universal Neuro II	Universal Neuro III	Synthes Matrix Neuro	
	<p>53-34164 STRAIGHT PLATE, 16 HOLE</p> <p>53-34804 STRAIGHT PLATE, 8 HOLE</p> <p>53-34408 STRAIGHT PLATE, 4 HOLE W/BAR</p> <p>53-34212 2 HOLE, RIGID, 12mm BAR</p> <p>53-34212 2 HOLE, W/ITAB, 12mm BAR</p> <p>53-34216 2 HOLE, 16mm BAR</p> <p>53-34228 BOX PLATE, SMALL</p> <p>53-34230 BOX PLATE, LARGE, W/ITAB</p> <p>53-34240 BOX PLATE, LARGE</p> <p>53-34300 RECTANGLE PLATE</p> <p>53-34812 GAP PLATE, 8 HOLE, SMALL</p> <p>53-34822 GAP PLATE, 8 HOLE, LARGE</p> <p>53-34808 DOUBLE-Y PLATE, 6-HOLE, W/BAR</p> <p>53-34200 X PLATE, 4 HOLE</p> <p>53-00203 MESH PLATE, ROUND, SMALL</p> <p>53-00204 MESH PLATE, ROUND, MEDIUM</p> <p>53-00465 MESH PLATE, ROUND, LARGE</p> <p>53-00242 TRANS/SUBYRITHINE PLATE, SMALL</p> <p>53-00246 TRANS/SUBYRITHINE PLATE, LARGE</p> <p>53-00224 TEMPORAL PLATE, MEDIUM</p> <p>53-00282 SUBOCIPITAL PLATE, SMALL</p> <p>53-00485 SUBOCIPITAL PLATE, LARGE</p>	<p>53-05164 STRAIGHT PLATE, 16 HOLE</p> <p>53-05204 STRAIGHT PLATE, 8</p> <p>53-05408 STRAIGHT PLATE, 4 HOLE W/BAR</p> <p>53-06212 2 HOLE, RIGID, 12mm BAR</p> <p>53-06212 2 HOLE, W/ITAB, 12mm BAR</p> <p>53-05216 2 HOLE, 16mm BAR</p> <p>53-05228 BOX PLATE, SMALL</p> <p>53-05230 BOX PLATE, LARGE, W/ITAB</p> <p>53-05240 BOX PLATE, LARGE</p> <p>53-05812 GAP PLATE, 8 HOLE, SMALL</p> <p>53-05822 GAP PLATE, 8 HOLE, LARGE</p> <p>53-05808 DOUBLE-Y PLATE, 6-HOLE, W/BAR</p> <p>53-05300 X PLATE, 4 HOLE</p> <p>54-00344 Dynamic mesh 40x40x0.3mm</p> <p>54-00346 Dynamic mesh 60x60x0.3mm</p> <p>54-00345 Dynamic mesh 120x120x0.3mm</p> <p>54-00645 Dynamic mesh 40x40x0.6mm</p> <p>54-00646 Dynamic mesh 60x60x0.6mm</p> <p>54-00647 Dynamic mesh 120x120x0.6mm</p>	<p>04-503 072 Adaption Plate, 20 holes</p> <p>04-503 071 Adaption Plate, 7 holes</p> <p>04-503 070 Adaption Plate, 5 holes</p> <p>04-503 067 Straight Plate, 2 holes, 12mm</p> <p>04-503 063 Cranial Plate, straight, with centre space, 12 mm, 4 holes</p> <p>04-503 065 Frame Plate, square, 4 holes, 14x14 mm</p> <p>04-503 066 Frame Plate, square, 4 holes, 18x18 mm</p> <p>04-503 073 Frame Plate, rectangular, 4 holes, 10x16 mm</p> <p>04-503 074 Strut Plate, 2x3 holes, 14x24 mm</p> <p>04-503 075 Strut Plate, 2x4 holes, 14x34 mm</p> <p>04-503 068 Double-Y-Plate, 6 holes, length 18 mm</p> <p>04-503 069 Double-Y-Plate, 6 holes, length 21 mm</p> <p>04-503 064 X-Plate, 4 holes</p> <p>04-503 090 Circular, 30x0.4mm, malleable</p> <p>04-503 091 Circular, 70x0.4mm, malleable</p> <p>04-503 092 Circular, 100x0.4mm, malleable</p> <p>04-503 093 Circular, 30x0.4mm, rigid</p> <p>04-503 094 Circular, 70x0.4mm, rigid</p> <p>04-503 095 Circular, 100x0.4mm, rigid</p> <p>04-503 086 Mastoid Plate, small</p> <p>04-503 087 Mastoid Plate, medium</p> <p>04-503 088 Mastoid Plate, large</p> <p>04-503 087 Temporal Mesh Plate, contoured, 0.4mm</p> <p>04-503 086 Crescent-shaped, small, 0.4mm, malleable</p> <p>04-503 087 Crescent-shaped, large, 0.4mm, malleable</p> <p>04-503 088 Crescent-shaped, small, 0.4mm, rigid</p> <p>04-503 089 Crescent-shaped, large, 0.4mm, rigid</p>	<p>na</p>
Environment	<p>Styker Leiberger GmbH & Co. KG, Facility in Joseph-Lang-Str. 22, 78570 Mühlheim an der Donau, GER</p>	<p>Styker Leiberger GmbH & Co. KG, Facility in Joseph-Lang-Str. 22, 78570 Mühlheim an der Donau, GER</p>	<p>Synthes (USA), 1101 Synthes Avenue, Monument, CO 80132</p>	<p>For Universal Neuro II: Equipment manufacturer (Synthes Station) on the equivalent machines using the equivalent machine environment</p>
Indications for use	<p>see to new S10k number</p>	<p>031658/ 056328</p>	<p>K02202</p>	<p>na</p>
Application Area	<p>Neuro (Cranial)</p>	<p>Craniofacial</p>	<p>Neuro (Cranial), Mastoid, Maxilla & Chin</p>	<p>For Universal Neuro II: Equivalent application area compared to Universal Neuro III For Synthes: Enhanced indication area including the one of Universal Neuro III</p>
Material	<p>Commercially Pure Titanium, Grade II and IV</p>	<p>Commercially Pure Titanium, Grade II and IV</p>	<p>Commercially Pure Titanium</p>	<p>Equivalent material, therefore all plates are equal in regard to the mechanical and chemical properties of their material</p>
Size				<p>The shape of the Universal Neuro II subject devices is equivalent to the range offered by both predicate devices.</p>
Counterbore			<p>na</p>	<p>The shape of the counterbore is widened by means of the diameter at the lower opening. Due to the Universal Neuro II screw head diameter of 2.7mm the screw is not endangered to fall through the plate hole</p>
Thickness	<p>0.6mm (1 "rigid" plate) and 0.4mm (all others)</p>	<p>0.6mm (1 "rigid" plate) and 0.5mm (all others)</p>	<p>0.4mm</p>	<p>The thickness of the 0.4mm subject devices is equivalent to the Synthes predicate devices. The 0.6mm subject devices is equivalent to the 0.5mm Universal Neuro II device</p>
Surface treatment	<p>Type III Anodization</p>	<p>Type III Anodization</p>	<p>na</p>	<p>Equivalent surface treatment of Universal Neuro II and III devices</p>
Size				<p>The shape of the Universal Neuro III Skull Base subject devices is equivalent to the ones offered by the predicate devices from Synthes.</p> <p>The outer dimensions of the Universal Neuro III Skull Base subject devices is within the range offered by both predicate devices</p>
Internal Mesh Pattern				<p>The pattern of the subject devices as well as all predicate devices is equivalent characterized by fixation holes surrounded by connecting bars</p> <p>The pattern holes are of equivalent dimensions between the Universal Neuro II subject devices and the Universal Neuro III predicate devices. The hole connecting bars of the Universal Neuro III subject devices are wider than the ones of the Universal Neuro II predicate devices offering more stability</p>
Thickness	<p>1.3mm and 0.4mm</p>	<p>0.3mm and 0.6mm</p>	<p>0.4mm and 0.6mm</p>	<p>The thickness of the Universal Neuro III Skull Base subject devices is within the range offered by both predicate devices</p>
Surface treatment	<p>Type III Anodization</p>	<p>Type III Anodization</p>	<p>na</p>	<p>Equivalent surface treatment of Universal Neuro II and III devices</p>
Surgical Technique/Preparation	<p>see intended use</p>	<p>see intended use</p>	<p>see intended use</p>	<p>see intended use</p>

SUBSTANTIAL EQUIVALENCE TABLE: BURR HOLE COVERS (Product Code GXR)

Device	Subject	Predecessor	Similarities / Differences
	Universal Neuro II	Synthes Matrix Neuro	n/a
	53-34507: 7mm Burr Hole Cover 53-34510: 10mm Burr Hole Cover 53-34514: 14mm Burr Hole Cover 53-34520: 20mm Burr Hole Cover 53-34524: 24mm Burr Hole Cover 53-34620: 20mm Shunt	53-05507: 7mm Burr Hole Cover 53-05510: 10mm Burr Hole Cover 53-05514: 14mm Burr Hole Cover 53-05520: 20mm Burr Hole Cover 53-05524: 24mm Burr Hole Cover 04-503.021: 12mm Burr Hole Cover 04-503.022: 15mm Burr Hole Cover 04-503.023: 17mm Burr Hole Cover 04-503.024: 24mm Burr Hole Cover 04-503.026: 12mm Shunt 04-503.027: 15mm Shunt 04-503.028: 17mm Shunt 04-503.029: 24mm Shunt	n/a
Manufacturer	Stryker Leibinger GmbH & Co. KG, Facility in Joseph-Lang Str. 22, 78570 Muehlheim an der Donau, GER	Stryker Leibinger GmbH & Co. KG, Facility in Joseph-Lang Str. 22, 78570 Muehlheim an der Donau, GER	For Universal Neuro II (Equivalent manufacturer (Stryker/ Synthes) on the equivalent machines using the subject's design environment.
Environment			
510(k)			
Intended Use	see to new 510k number	K031659	n/a
Indications for use	The Stryker® Universal Neuro II implant system is intended for reconstruction, stabilization and/or rigid fixation of non-load-bearing bony areas subsequent to craniotomy, craniectomy and cranial fractures in adults and adolescents (age 12 and higher).	The Stryker Leibinger Universal Neuro System is a low-profile selective frame system intended for osteotomy, craniectomy, craniotomy and rigid fixation of non-load bearing areas, fractures and reconstruction in non-load bearing areas.	The subject device is intended to be used in cranial areas which are similar to the predicate. But it is limited to non-load bearing cranial indications only, whereas the Synthes predicate indications cover facial/orbital or mandibular sites as well. Additionally, the subject is intended to be used in adolescents and adults whereas the predicates have no limitation to
Application Area	Neuro (Cranial)	Craniofacial	For Universal Neuro II: Equivalent application area compared to Universal Neuro III For Synthes: Enhanced application area including the one of Universal Neuro III
Material	Commercially Pure Titanium	Commercially Pure Titanium	Equivalent material, therefore all plates are equal in regard to the mechanical and chemical properties of their material.
Design	53-34507: 7mm 53-34510: 10mm 53-34514: 14mm 53-34520: 20mm 53-34524: 24mm	53-05507: 7mm 53-05510: 10mm 53-05514: 14mm 53-05520: 20mm	The subject devices are of equivalent choice of sizes as available for the predicate devices. The number for countersinks for screw insertion of the subject devices are equivalent to the ones of the Synthes predicate devices.
Design	53-34614: 14mm 53-34620: 20mm	n/a	The subject devices are of equivalent choice of sizes as available for the predicate devices. The number for countersinks for screw insertion of the subject devices are equivalent to the ones of the Synthes predicate devices.
Countersink			The shape of the countersink is widened by means of the diameter at the lower opening. Due to the Universal Neuro III screw head diameter of 2.7mm the screw is not endangered to fail through the plate hole.
Thickness	0.4mm	0.5mm	
Surface treatment	Type III Anodization	Type III Anodization	Equivalent surface treatment of Universal Neuro II and III devices.
Surgical Techniques/Site Preparation	see intended use	see intended use	see intended use



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room –WO66-G609
Silver Spring, MD 20993-0002

Mr. Rob Yamashita
Senior Regulatory Affairs Representative
Stryker Craniomaxillofacial
750 Trade Centre Way, Suite 200
Portage, MI 49002

JAN - 5 2012

Re: K112557

Trade/Device Name: Stryker Universal Neuro 3 System
Regulation Number: 21 CFR 882.5320
Regulation Name: Performed alterable cranioplasty plate
Regulatory Class: Class II
Product Code: GWO, GXR, HBW
Dated: December 28, 2011
Received: December 29, 2011

Dear Mr. Yamashita:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Malvina B. Eydelman, M.D.

Director

Division of Ophthalmic, Neurological
and Ear, Nose and Throat Devices

Office of Device Evaluation

Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number: K112557

Device Name: Stryker Universal Neuro 3 System

Indications For Use:

The Stryker Universal Neuro 3 System is intended for reconstruction, stabilization and/or rigid fixation of non-load-bearing bony areas subsequent to craniotomy, craniectomy and cranial fractures in adults and adolescents (age 12 and higher).

Contraindications:

The Stryker Universal Neuro 3 System is contraindicated for the following:

- Use of plates in non-reducible and unstable fractures
- Patients with active local infections
- Patients with metal allergies and foreign body sensitivity
- Potentially non-compliant patients with mental or neurological conditions who are unwilling or incapable of following postoperative care instructions
- Patients with limited blood supply to or insufficient quality of bone
- Use of products in cases where the fixation of the products could result in their peripheral edge coming into contact with the dura mater
- Screws coming in contact with the dura mater
- Use of implants adjacent to developing paranasal sinuses

Prescription Use (Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____ (21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

J. Brown

(Division Sign-Off)
Division of Ophthalmic, Neurological and Ear,
Nose and Throat Devices

Page 1 of 1

510(k) Number K112557