

MAR 29 2012



TRADITIONAL 510(k) PREMARKET SUBMISSION  
IMPERIAL SURGICAL LTD. BLANKET and SOLUTION WARMING CABINETS

**510(k) Summary**  
**For**  
**Imperial Surgical Ltd. BLANKET and SOLUTION WARMING CABINETS**

**September 14, 2011**

Manufactured by: Imperial Surgical Ltd.  
Contact Name: Suzy Bairos, Vice President  
Address: 850 Halpern Ave.  
Dorval, Quebec  
H9P 1G6  
Tel.: 514-631-7988 Ext. 21  
Fax: 514-631-9083

IMPERIAL SURGICAL LTD.

850 Halpern Ave., Dorval, Quebec, H9P 1G6 Canada

Toll Free (800) 661-5432 • T.514.631.7988 • F.514.631.9083 • Web: [www.surgmed.com](http://www.surgmed.com)



TRADITIONAL 510(k) PREMARKET SUBMISSION  
 IMPERIAL SURGICAL LTD. BLANKET and SOLUTION WARMING CABINETS

1. **Device Name/Codes**

Device Proprietary Name: Imperial Surgical Blanket and Solution Warming Cabinets

Common/Usual Name: Warming Cabinets

Classification Name: Warmer, Thermal, Infusion Fluid

Classification Code: LGZ

2. **Predicate Device**

Steris Amsco Warming Cabinet (K092823, December 12, 2009)

Enthermics EC-7701 Fluid Warming Cabinet, K993797, January 20, 2000

3. **Description of Device**

The Imperial Surgical Ltd. Blanket and Solution Warming Cabinets are designed to store and warm blankets, hospital linens, irrigation fluids and/or injection fluids in accordance with the recommended warming temperatures and storage time guidelines provided by the manufacturers of such products

The subject device is available with either single compartment or double compartment cabinets, all designed to ensure maximum storage capability. The 7000 series are made of 18 gauge stainless steel with the doors and shelves made of 20 gauge stainless steel. The 8000 series are of a similar design to the 7000 series, however the stainless steel doors have an inset double pane tempered glass to allow interior cabinet viewing. The small 9000 series are made of powder coated steel, and the doors have an inset Plexiglas also to allow interior cabinet viewing.

4. **Intended Use**

The Imperial Surgical Ltd. Blanket and Solution Warming Cabinets are designed to store and warm blankets, hospital linens, irrigation fluids and/or injection fluids in accordance with the recommended warming temperatures and storage time guidelines provided by the manufacturers of such products

IMPERIAL SURGICAL LTD.

850 Haipern Ave., Dorval, Quebec, H9P 1G6 Canada

Toll Free (800) 661-5432 • T.514.631.7928 • F.514.631.9083 • Web: www.surgmed.com



TRADITIONAL 510(k) PREMARKET SUBMISSION  
 IMPERIAL SURGICAL LTD. BLANKET and SOLUTION WARMING CABINETS

**5. Description of Safety and Substantial Equivalence**

The Imperial Surgical Ltd. Blanket and Solution Warming Cabinets are substantially equivalent to the two predicate devices in all material respects. Please refer to the table below comparing the technological characteristics of the proposed Imperial Surgical Blanket and Solution Warming Cabinets to the predicate devices.

**Summary of the Proposed Device and Predicate Devices Technological Characteristics**

| Features                       | Imperial Surgical Ltd. Blanket and Solution Warming Cabinets  | Amsco Warming Cabinet (K092823)  | EC-7701 Fluid Warming Cabinet (K993797)   |
|--------------------------------|---|--|---|
| Intended Use                   | The Imperial Surgical Ltd. Blanket and Solution Warming Cabinets are designed to store and warm blankets, hospital linens, irrigation fluids and/or injection fluids in accordance with the recommended warming temperatures and storage time guidelines provided by the manufacturers of such products | Amsco Warming Cabinet is for heating flasks solutions, blankets and similar clinical articles. | The Enthermics Medical Systems EC-7701 Fluid Warming Cabinet is designed for safely store and warm irrigation fluids or injection fluids in accordance with the recommended warming temperatures and storage times stated in the fluid manufacturer's labeling. |
| Heating System                 | Convection Electric heating element with circulating fan  | Steam heat or Electric heater blanket and fan blower (convection heating)                      | Fully insulated electrothermal cable array (Convection heating) Usually radiant heating   |
| Unit Configuration             | Single / dual   | Single/Dual  | Single  |
| Unit Depth                     | 22" or 28" for table top and full size cabinets. 16" depth on count top model.  | 18" or 24"   | 34"   |
| Model                          | Wall or Counter   | Wall or Counter  | Wall  |
| Interior and Exterior Surfaces | Stainless Steel interior and exterior on table top and full size models. Counter top model made with baked on powder coated steel outer shell with Stainless Steel interior.  | Stainless Steel  | Stainless Steel   |
| Installation                   | Free standing or recessed   | Open-Mounted or Recessed   | Free-standing (mobile) or Recessed  |



**TRADITIONAL 510(k) PREMARKET SUBMISSION**  
**IMPERIAL SURGICAL LTD. BLANKET and SOLUTION WARMING CABINETS**

| Features                 | Imperial Surgical Ltd. Blanket and Solution Warming Cabinets   | Amsco Warming Cabinet (K092823)  | EC-7701 Fluid Warming Cabinet (K993797)  |
|--------------------------|--|--|--|
| Door                     | Stainless steel or double pane Glass set in stainless steel frame on table top and full size models. Counter top model has Plexiglas mounted on to anodized aluminum door frame.   | Stainless Steel  | Stainless Steel (Glass)  |
| Cabinet Storage Capacity | Dual compartment and full size cabinets equipped with 3 stainless steel perforated shelves<br>Table top model equipped with 1 stainless steel perforated shelf   | Dual Compartment Model-<br>Two shelves:<br>15 flasks -18" upper<br>45 flasks -18" lower<br>20 flasks -24" upper<br>60 flasks -24" lower<br><br>Single Compartment Model-<br>One Shelf<br>15 flasks -18"<br>20 Flasks - 24" | The cabinet is equipped with three (3) white epoxy-coated wire baskets, each with 24 liter capacity.                                 |
| Cabinet Volume           | Dual compartment capacity 18.1 cu ft overall (4.3 cu ft upper - 13.8 cu ft lower)<br>Single compartment full size capacity 19.8 cu t<br>Single compartment shallow depth design capacity 14.6 cu ft<br>Table top model capacity 7.5 cu ft<br>Table top shallow design capacity 5.2 cu ft<br>Counter top version has capacity of 2.0 cu ft (approx) | 18" upper chamber - 3.1 cu ft<br>24" upper chamber - 4.2 cu ft<br>18" lower chamber - 8.9 cu ft<br>24" lower chamber - 12 cu ft  | Single compartment Apprx 7.8 cu ft internal capacity   |
| Controls                 | Electronic temp controller with LED display.<br>Illuminated power switch/breaker (red)<br>Amber neon indicator for Element on<br>Red neon indicator for trouble  | Thermostat/power switch/fuse with indicating light/color coded temperature selector.   | Electronic control consists of 4 digit LED display, on/off keys, integrated lock feature and a series of prompt sequence indicators. |
| Software                 | N/A  | N/A  | Not applicable   |



TRADITIONAL 510(k) PREMARKET SUBMISSION  
 IMPERIAL SURGICAL LTD. BLANKET and SOLUTION WARMING CABINETS

| Features                     | Imperial Surgical Ltd.<br>Blanket and Solution<br>Warming Cabinets   | Amsco Warming Cabinet<br>(K092823)  | EC-7701 Fluid Warming<br>Cabinet (K993797)   |
|------------------------------|--|---|--|
| Temperature Selection Range  | Can be factory programmed for Blankets:<br>30°C to 71°C (86°F to 160°F)<br>For Irrigation Fluids:<br>30°C to 66°C (86°F to 150°F)<br>For Injection Fluids:<br>30°C to 40°C (86°F to 104°F)<br>Must be specified at time of ordering. | 95°F to 150°F (35°C to 65°C)  | 90°F to (32°C) to 150°F to (66°C)  |
| Temperature Lock             | Settings are factory set and locked. The device can be unlocked in the field and user can change settings if they wish to.   | N/A   | The device allows the user to "lock" the mode (IRR or INJ) and temperature setting controls.   |
| Door Lock                    | Not available  | Available by SSQ (Special Sales Quote) only.  | Available as an option   |
| Over Temperature Alarm Point | Visual alarm if chamber temperature exceeds 6°C above set temperature. Internal sensor shuts off heating elements when over temperature occurs   | Visual alarm if unit has a chamber temperature greater than 12°F above set temperature. In the event of an over temp condition, sensors automatically turns off the heater (s). | Visual and audible alarm if unit has a chamber temperature greater than 10°F (5.5°C) above set temperature. In the event of an over temp condition, the heating system shuts down. |
| Voltage Requirements         | 120VAC 60Hz  | Electric Model: 110/120 Vac, 220/240 Vac nominal, 50/60 HZ. Steam Model: 120 Vac single phase.  | 125 Vac, 60 Hz, 1 ph   |

IMPERIAL SURGICAL LTD.

350 Halpern Ave., Dorval, Quebec, H9P 1G6 Canada

Toll Free (800) 661-5432 • T.514.631.7928 • F.514.631.9023 • Web: www.surgimed.com



TRADITIONAL 510(k) PREMARKET SUBMISSION  
IMPERIAL SURGICAL LTD. BLANKET and SOLUTION WARMING CABINETS

**6. Substantial Equivalent Assessment of non-clinical performance data.**

The Imperial Surgical Ltd. Blanket and Solution Warming Cabinets have been tested to confirm they are safe and effective for their intended use. The following safety and functional tests have been performed:

Electrical safety is demonstrated as each unit is field certified by CSA International using CSA Standard(s) C22.2 No(s). 0,0.4 as a guide and SPE-1000, UL1262, US Field Evaluation.

Performance testing has been conducted to evaluate the following:

- accuracy of the temperature controller
- even distribution of heat throughout the cabinet with limited temperature variance, and with multiple temperature readings within the accepted temperature variance of +/- 1.5°C.
- how quickly the cabinet will reach the desired temperature setting from a cold start.
- the recovery time required when the cabinet door is opened repeatedly for the temperature to return to the desired setting.
- time required for the cabinet to achieve the desired temperature setting loaded to 3/4 of its capacity (typical load expected during operation).

All tests met the expected specifications for the device. It is concluded that the electrical safety and performance testing conducted support the safety and effectiveness of the Imperial Surgical Ltd. Blanket and Solution Warming Cabinets.

**7. Conclusion**

Based on the information provided in this 510(k) premarket notification, the Imperial Surgical Ltd. Blanket and Solution Warming Cabinets are substantially equivalent in terms of safety and effectiveness to the predicate devices identified above.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Room - WO66-G609  
Silver Spring, MD 20993-0002

Imperial Surgical Limited  
C/O Ms. Kathryn Ronalds  
Associate Director, Regulatory Affairs (Devices)  
CanReg Incorporation  
4 Innovation Drive  
Dundas  
CANADA L9H 7P3

MAR 29 2012

Re: K112702  
Trade/Device Name: Blanket and Solution Warming Cabinets  
Regulation Number: Unclassified  
Regulation Name: None  
Regulatory Class: Unclassified  
Product Code: LGZ  
Dated: March 21, 2012  
Received: March 22, 2012

Dear Ms. Ronalds:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



 Anthony D. Watson, B.S., M.S., M.B.A.  
Director  
Division of Anesthesiology, General Hospital,  
Infection Control and Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

### INDICATIONS FOR USE

510(k) Number (if known):

Device Name: Blanket and Solution Warming Cabinets

Indications for Use:

The Imperial Surgical Ltd. Blanket and Solution Warming Cabinets are designed to store and warm blankets, hospital linens, irrigation fluids and/or injection fluids in accordance with the recommended warming temperatures and storage time guidelines provided by the manufacturers of such products.

Prescription Use \_\_\_\_\_ AND/OR Over-The-Counter Use ✓  
(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

RHC Chyma 3/28/12  
(Division Sign-Off)  
Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

510(k) Number: K112702

Page 1 of 1