

MAR 22 2012

Genii, Inc.
Traditional 510(k) Premarket Submission
gi 4000 Electrosurgical Generator

Section 005 – 510(k) Summary**Submission Correspondent:**

Name: Emergo Group, Inc.
www.emergogroup.com
Address: 611 West Fifth Street
Third Floor
Austin, TX 78701
Phone: (512) 327-9997
Fax: (512) 327-9998
Contact: Julie Powell
Email: jpowell@emergogroup.com

Submission Sponsor:

Genii, Inc.
2145 Woodlane Drive Suite 101-W
St. Paul, MN 55125
Tel. 651-501-4810
Fax. 651-501-4819

Date Prepared: 2 November 2011; Revised 23 February 2012; Revised 14 March 2012

Trade Name: Genii *gi 4000* Electrosurgical Generator

Regulation Number: 878.4400

Classification Name: Electrosurgical cutting and coagulation device and accessories.

Product Code: GEI

Classification Panel: General and Plastic Surgery

Regulatory Class: Class II

Device Description:

The *gi 4000* has seven monopolar outputs, a monopolar argon assisted coagulation output, and a bipolar output. Users make selections from a touch screen that displays stored default power setting start points, convenience and safety information, as well as indications as to current selection mode and return electrode monitoring status.

The *gi 4000* has an integrated peristaltic lavage pump designed to provide lavage washing in the gut, either by accessories connected directly to endoscopes which have dedicated washing channels, or by accessories which introduce lavage fluid into the biopsy port of an endoscope or via dedicated washing channels incorporated into a therapeutic accessory such as a bipolar

Genii, Inc.
Traditional 510(k) Premarket Submission
gi 4000 Electrosurgical Generator

endostasis probe. The system is operated by a dual pedal footswitch which controls the lavage pump and power delivery.

The argon output is designed for an argon assisted coagulation method. This method uses 99.999% pure argon gas to provide ionized RF arcs to the targeted tissue. The monopolar outputs are designed for cutting and coagulation methods. The unit is designed with a standard monopolar active cord receptacle (foot switch activated).

Endoscopic accessories connect to the monopolar receptacle via a widely available detachable active cord accessory. Both the monopolar and argon outputs require the use of a patient return electrode (grounding pad) to return the RF energy from the patient back to the generator. The unit is designed with a standard GI bipolar single jack receptacle. No patient return electrode (grounding pad) is needed for bipolar methods. The output power of the argon, monopolar and bipolar methods are independently controlled through a touch panel LCD display located on the front panel of the generator.

The argon gas for the argon output is exclusively supplied by the disposable Genii Argon Gas Canister. Each canister contains 49 liters of ultra-high purity (99.999% pure) argon gas (certified). When using the Argon Coagulation Method, the ArConnect is intended for single-use only and will be supplied as a non-sterile, disposable item. The ArConnect has been designed to allow the delivery of monopolar high frequency electrosurgical energy and argon gas to a single use argon coagulation probe sold separately. The ArConnect is supplied individually packaged in disposable packaging and is not to be reprocessed or re-used.

Indications for Use:

The *gi 4000* Electrosurgery Generator is intended to deliver electrosurgical outputs to perform cutting and/or coagulation and coagulation with or without delivery of argon gas in flexible endoscopic applications. The integrated lavage pump is intended to wash tissue or mucosa within the gut when used in conjunction with washing catheters or bipolar endoscopic probes, endoscope water jet channels, or endoscope working channels. The Genii ArConnect™ Conducting Adapter is intended for use with the *gi 4000* during argon assisted coagulation.

Predicate Devices:

1. Olympus AFU-100 Peristaltic Pump Unit (PPU); K073207
2. Erbe APC 300 Argon Plasma Coagulator and Accessories; K963189
3. Erbe ICC 200 Electrosurgical Workstation for Flexible Endoscopy; K933157
4. Erbe APC Connector Hose and Probes; K013348

Summary of Technical Comparisons:

The following table compares *gi 4000* Electrosurgery Generator to the predicate devices with respect to intended use, target population, technological characteristics and principles of operation, providing more detailed information regarding the basis for the determination of substantial equivalence.

Genii, Inc.
 Traditional 510(k) Premarket Submission
 gi 4000 Electrosurgical Generator

Manufacturer	Genii	Olympus		
Trade Name	<i>Gi 4000</i>	Olympus AFU-100	Erbe APC 300	Erbe ICC 200
510(k) Number	NA	K073207	K963189	K933157
Product Code	GEI	FEQ (with GEI)	GEI	GEI
Regulation Number	878.4400	876.1500 (with 878.4400)	878.4400	878.4400
Regulation Name	Electrosurgical cutting and coagulation device and accessories	Pump, air non manual, for endoscope (with Electrosurgical cutting and coagulation device and accessories)	Electrosurgical cutting and coagulation device and accessories	Electrosurgical cutting and coagulation device and accessories
Intended Use	The gi 4000 Electrosurgery Generator is intended to deliver electrosurgical outputs to perform cutting and/or coagulation and coagulation with or without delivery of argon gas in flexible endoscopic applications. The integrated lavage pump is intended to wash tissue or mucosa within the gut when used in conjunction with washing catheters or bipolar endoscopic probes, endoscope water jet channels, or endoscope working channels.	The Olympus AFU-100 peristaltic flushing pump is intended for irrigation of instruments or irrigation/flushing/cleansing of tissue surfaces and wounds supporting endoscopic diagnosis or therapy.	The Erbe APC 300 Argon Plasma Coagulator is an argon gas delivery system that is designed for coagulation when used in combination with Erbe ICC electrosurgery Generators and Erbe APC Probes.	The ICC 200 is intended to deliver high frequency electrical current for the cutting/or coagulation of tissue.
Pump Type	Peristaltic	Peristaltic	n/a	n/a
Flow Rate	0-735 ml/minute(tubing dependent)200 ml max through leur valves	0-600ml/minute (tubing dependent) 200 ml max through leur valves	n/a	n/a
Dimensions (Width x Height x Depth)	16" x 6.5" x 16" (19.5 X 7 X 19 max overall including pump head, water bottle holder and feet)	295 X 430 X115 mm (note millimeters not inches)	APC 16" x 4" x 14.5" (needs cart for gas tanks)	11" x 6" x 14.5"
Max Voltage Peak	4200 Vp.	n/a	n/a	4000 Vp
Frequency	315-550 KHz.	n/a	n/a	330-1,000 KHz (1 MHz).
Output Power Limit	60-200 watts.	n/a	n/a	120-200 watts
Gas Flow Rate	0.1-2.0 mL/minute.	n/a	0.1-9.0 mL/minute.	n/a
Gas Purity	99.999%.	n/a	99.998%.	n/a
Tank Size	49 liters.	n/a	2 X 934 liters.	n/a
Electrical Safety/other standards	IEC 60601-1 IEC 60601-1-2, IEC 60601-2-2, CISPR 11, IEC 62304, ISO 14971	IEC 60601-1, UL 60601-1, IEC 60601-1-2, IEC 60601-1-4	ANSI/AAMI HF-18/1993, IEC 60601-1-2	Not known
Wave Form • Mono Soft Coag • Mono Coag • Mono Blend Cut	Continuous Sine Modulated Alternating 50%DC	n/a	n/a	Continuous Sine Modulated n/a

Genii, Inc.
 Traditional 510(k) Premarket Submission
 gi 4000 Electrosurgical Generator

Manufacturer	Genii	Olympus		
Trade Name	Gi 4000	Olympus AFU-100	Erbe APC 300	Erbe ICC 200
<ul style="list-style-type: none"> •Mono Cut •Bipolar •Argon 	Continuous sine Continuous sine Modulated			Continuous sine Continuous sine Modulated
Max Voltage Peak <ul style="list-style-type: none"> •Mono Soft Coag •Mono Coag •Mono Blend Cut /PBC •Mono Pulse Cut •Mono Cut •Bipolar •Argon 	300 2200 1625/1050 1000 780 104 4200	n/a	n/a	190 2300-2500 n/a (EC) 380-600/190 (EC) 600 190 4000
Frequency (KHz) <ul style="list-style-type: none"> •Mono Soft Coag •Mono Coag •Mono Blend Cut/PBC •Mono Pulse Cut •Mono Cut •Bipolar •Argon 	350KHz 460-500KHz 350/350 350 350 350 460-500	n/a	n/a	330KHz 1,000KHz n/a 330 330 350 1,000KHz
Max Power <ul style="list-style-type: none"> •Mono Soft Coag •Mono Coag •Mono Blend Cut/ PBC •Mono Pulse Cut •Mono Cut •Bipolar •Argon 	120 Watts @ 100 Ω 120 Watts @ 500 Ω 120 Watts @ 500 Ω/ 120 Watts @ 500 Ω 120 Watts @ 500 Ω 200 Watts @ 500 Ω 60 Watts @ 400 Ω 120 Watts @ 500 Ω	n/a	n/a	120 Watts @ 125 Ω 120 Watts @ 350 Ω n/a 200 Watts @ 500 Ω 200 Watts @ 500 Ω 120 Watts @ 125 Ω 99 Watts @ 350 Ω
Duty Cycle <ul style="list-style-type: none"> •Mono Soft Coag •Mono Coag •Mono Cut •Bipolar 	100% 6% 100% 100%	n/a	n/a	100% 4% - 6% 100% 100%
Crest Factor <ul style="list-style-type: none"> •Mono Soft Coag •Mono Cut •Bipolar 	1.4 1.4 1.4	n/a	n/a	1.4 1.4 1.4
Power Curve <ul style="list-style-type: none"> •Mono Soft Coag •Mono Coag •Bipolar 	Narrow Wide Narrow	n/a	n/a	Narrow Wide Narrow
Pulse Speed <ul style="list-style-type: none"> •Mono Pulse Blend Cut •Mono Pulse Cut 	700ms ON/50ms OFF 700ms ON/50ms OFF	n/a	n/a	750ms ON/50ms OFF 750ms ON/50ms OFF

Genii, Inc.
 Traditional 510(k) Premarket Submission
 gi 4000 Electrosurgical Generator

The following table compares Genii ArConnect accessory to the predicate accessory device.

Manufacturer	Genii	
Trade Name	Genii ArConnect	Erbe APC Connector Hose
510(k) Number	NA	K013348
Product Code	GEI	GEI
Regulation Number	878.4400	878.4400
Regulation Name	Electrosurgical cutting and coagulation device and accessories	Electrosurgical cutting and coagulation device and accessories
Intended Use	The Genii ArConnect conducting adapter is intended for use with the Genii Gi4000 during argon assisted coagulation.	The APC Connector Hose and Probes are intended for use in argon plasma coagulation. The device is used to treat many conditions in endoscopy for various surgical procedures.
Material	PVC tubing Shore A75+3Clear, Cable 24AWG Terminal: brass with nickel and Gold plating; Insert: ABS PA707 with EX-33084 white; Overmold HDPE LH523, nature; PE white, #PEM864; Gas Tube: ABS PA707, with EX-33084 White; O ring: AS 568-008, Katon P/N: OR-NB50NL-008-BK, ShoreA: 50, BUNA-U; Heat shrink sleeve band: 82*98.5*	Polyamide, PEEK, Thermoflex, silicone.
Length	5 feet.	6.5 feet.
Back Flow	Integrated in the generator gas flow design	A separate back flow filter/valve is sold separately to be used with the hose.
How Supplied	Single use; disposable.	Re-useable.

Non-clinical Data – Bench Testing:

As part of demonstrating safety and effectiveness of the Genii *gi 4000* Electrosurgical Generator and in showing substantial equivalence to the predicate devices that are subject to this 510(k) submission, the following testing has been performed:

In vitro testing
 Electrical Safety
 EMC
 Biocompatibility
 Software Testing

In Vitro Testing: An *in vitro* tissue study was conducted to compare the lesions created by the subject *gi 4000* electrosurgical generator and the predicate Erbe ICC 200E with APC 300 system when using equivalent outputs. This *in vitro* study demonstrated substantial similarity between lesions created by the subject *gi 4000* device and the predicate Erbe ICC 200E with APC 300 system. The data demonstrate that the performance of the subject device is equivalent to the predicate devices.

Electrical Safety and EMC: The *gi 4000* Electrosurgical Generator has been tested and met the requirements of the Electrical Safety and EMC standards used.

Genii, Inc.
Traditional 510(k) Premarket Submission
gi 4000 Electrosurgical Generator

Biocompatibility Testing: Cytotoxicity, Irritation and Sensitization testing have been completed on the ArConnect Conducting Cable. The test results demonstrate that the ArConnect Conducting Cable is non-toxic.

Software Verification and Validation: Software verification and validation testing has been performed on the *gi 4000* Electrosurgical Generator.

Safety and Effectiveness:

By definition, a device is substantially equivalent to a predicate device when the device has the same intended use and the same technological characteristics as the previously cleared predicate device, or has the same intended use and different technological characteristics, and it can be demonstrated that the device is substantially equivalent to the predicate device, and that the new device does not raise different questions regarding its safety and effectiveness as compared to the predicate device.

It has been shown in this 510(k) submission that the differences between the Genii *gi 4000* Electrosurgical Generator and the predicate devices do not raise any questions regarding its safety and effectiveness. The Genii *gi 4000* Electrosurgical Generator, as designed and manufactured, therefore is determined to be substantially equivalent to the referenced predicate devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room - WO66-G609
Silver Spring, MD 20993-0002

Genii, Incorporated
% Emergo Group, Incorporated
Ms. Julie Powell
611 West 5th Street, Third Floor
Austin, Texas 78701

MAR 22 2012

Re: K113265

Trade/Device Name: gi 4000 Electrosurgical Generator
Regulation Number: 21 CFR 878.4400
Regulation Name: Electrosurgical cutting and coagulation device and accessories
Regulatory Class: Class II
Product Code: GEI
Dated: March 14, 2012
Received: March 15, 2012

Dear Ms. Powell:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Ms. Julie Powell

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Mark N. Melkerson
Director
Division of Surgical, Orthopedic
and Restorative Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

