

510(k) Summary for GC Aadva™ Abutments – Zirconia – Bo Ac Bh Sv*510(k) Summary***GC Aadva™ Abutments – Zirconia – Bo Ac Bh Sv**

Date Prepared: October 18, 2012

Submitter/Contact Person	H. Carl Jenkins The Wood Burditt Group FDA Regulatory Counseling 10 E. Scranton Avenue, Suite 201 Lake Bluff, IL 60044 (ph) (847) 234-7500 x 205 (fax) (847) 574-0728 (email) hcjenkins@woodburditt.com
Applicant	GC America, Inc. 3737 W. 127th Street Alsip, IL 60803 800.323.3386 x4042 708.897.4042 708.897.4031 (fax)
Manufacturer	GC CORPORATION. 76-1 HASUNUMA-CHO, ITABASHI-KU TOKYO 174-8585 JAPAN
Device Name	GC Aadva™ Abutments – Zirconia – Bo Ac Bh Sv
Common Name	Endosseous Dental Implant Abutment
Classification	Class II Procode NHA Regulation: 21 CFR 872.3630

Device Description

The "GC Aadvant[™] Abutments - Zirconia - Bo Ac Bh Sv" are abutments which are placed into a dental implant to provide support for a prosthetic restoration. The subject abutments are indicated for use with the following implant fixtures:

Abutment - Implant Compatibility Table

Types	Material	Antirotational Features	Implant - Engaged Fixture (Manufacturer)
IN-Bo-3.3mm #154	Zirconia : Y-TZP	Internal Square	Φ3.3 Bone Level Implant NC (Institut straumann AG) K062129
IN-Bo-4.1mm #155	Zirconia : Y-TZP	Internal Square	Φ4.1, Φ4.8 Bone Level Implant RC (Institut straumann AG) K062129
IN-Ac-3.5mm #156	Zirconia : Y-TZP	Internal Hexagon	Φ3.5 NobelActive NP (Nobel Biocare AB) K071370
IN-Ac-4.3mm #157	Zirconia : Y-TZP	Internal Hexagon	Φ4.3, Φ5.0 NobelActive RP (Nobel Biocare AB) K071370
IN-Bh-3.8mm #158	Zirconia : Y-TZP	Internal Hexagon	Φ3.8 Tapered Internal Implant (BioHorizons) K071638
IN-Bh-4.6mm #159	Zirconia : Y-TZP	Internal Hexagon	Φ4.6 Tapered Internal Implant (BioHorizons) K071638
IN-Sv-3.5mm #151	Zirconia : Y-TZP	Internal Hexagon	Φ3.7, Φ4.1 Tapered Screw-Vent (Zimmer Dental Inc.) K013227
IN Sv-4.5mm #152	Zirconia : Y-TZP	Internal Hexagon	Φ4.7 Tapered Screw-Vent (Zimmer Dental Inc.) K013227

The "GC Aadva™ Abutments – Zirconia – Bo Ac Bh Sv" abutment components comply with ISO 13356:2008 – Implants for Surgery – Ceramic materials based on yttria-stabilized tetragonal zirconia (Y-TZP). The abutments are mounted into the implant with a screw made of Titanium grade Ti-6Al-4V, which meets the requirements of ISO : 5832-3. Each of these abutments have the following design limitations:

<u>Types</u>	<u>Collar</u>	<u>Post</u>	<u>Angulation</u>	<u>Length</u>	<u>Width</u>	<u>Head space</u>
<u>IN-Bo-3.3mm #154</u>	<u>Max.6.0mm</u>	<u>Max.10.0mm</u>	<u>Max.20°</u>	<u>Max.12.5mm</u>	<u>Max. Ø 11.6mm</u>	<u>Max.11.8mm</u>
	<u>Min.2.5mm</u>	<u>Min.4.0mm</u>	<u>Min.0°</u>	<u>Min.6.5mm</u>	<u>Min. Ø 4.1mm</u>	<u>Min.1.0mm</u>
<u>IN-Bo-4.1mm #155</u>	<u>Max.6.0mm</u>	<u>Max.10.0mm</u>	<u>Max.20°</u>	<u>Max.12.5mm</u>	<u>Max. Ø 11.6mm</u>	<u>Max.11.8mm</u>
	<u>Min.2.5mm</u>	<u>Min.4.0mm</u>	<u>Min.0°</u>	<u>Min.6.5mm</u>	<u>Min. Ø 4.9mm</u>	<u>Min.1.0mm</u>
<u>I-Ac-3.5mm #156</u>	<u>Max.6.0mm</u>	<u>Max.10.0mm</u>	<u>Max.20°</u>	<u>Max.12.5mm</u>	<u>Max. Ø 11.6mm</u>	<u>Max.11.8mm</u>
	<u>Min.2.5mm</u>	<u>Min.4.0mm</u>	<u>Min.0°</u>	<u>Min.6.5mm</u>	<u>Min. Ø 4.1mm</u>	<u>Min.1.0mm</u>
<u>IN-Ac-4.3mm #157</u>	<u>Max.6.0mm</u>	<u>Max.10.0mm</u>	<u>Max.20°</u>	<u>Max.12.5mm</u>	<u>Max. Ø 11.6mm</u>	<u>Max.11.8mm</u>
	<u>Min.2.5mm</u>	<u>Min.4.0mm</u>	<u>Min.0°</u>	<u>Min.6.5mm</u>	<u>Min. Ø 4.9mm</u>	<u>Min.1.0mm</u>
<u>IN-Bh-3.8mm #158</u>	<u>Max.6.0mm</u>	<u>Max.10.0mm</u>	<u>Max.20°</u>	<u>Max.12.5mm</u>	<u>Max. Ø 11.6mm</u>	<u>Max.11.8mm</u>
	<u>Min.2.5mm</u>	<u>Min.4.0mm</u>	<u>Min.0°</u>	<u>Min.6.5mm</u>	<u>Min. Ø 4.1mm</u>	<u>Min.1.0mm</u>
<u>IN-Bh-4.6mm #159</u>	<u>Max.6.0mm</u>	<u>Max.10.0mm</u>	<u>Max.20°</u>	<u>Max.12.5mm</u>	<u>Max. Ø 11.6mm</u>	<u>Max.11.8mm</u>
	<u>Min.2.5mm</u>	<u>Min.4.0mm</u>	<u>Min.0°</u>	<u>Min.6.5mm</u>	<u>Min. Ø 4.9mm</u>	<u>Min.1.0mm</u>
<u>IN-Sv-3.5mm #151</u>	<u>Max.6.0mm</u>	<u>Max.10.0mm</u>	<u>Max.20°</u>	<u>Max.12.5mm</u>	<u>Max. Ø 11.6mm</u>	<u>Max.11.8mm</u>
	<u>Min.2.5mm</u>	<u>Min.4.0mm</u>	<u>Min.0°</u>	<u>Min.6.5mm</u>	<u>Min. Ø 4.1mm</u>	<u>Min.1.0mm</u>
<u>IN-Sv-4.5mm #152</u>	<u>Max.6.0mm</u>	<u>Max.10.0mm</u>	<u>Max.20°</u>	<u>Max.12.5mm</u>	<u>Max. Ø 11.6mm</u>	<u>Max.11.8mm</u>
	<u>Min.2.5mm</u>	<u>Min.4.0mm</u>	<u>Min.0°</u>	<u>Min.6.5mm</u>	<u>Min. Ø 4.9mm</u>	<u>Min.1.0mm</u>

Intended Use:

The "GC Aadva™ Abutments - Zirconia - Bo Ac Bh Sv" is intended for use with an endosseous dental implant to support a prosthetic device in a partially or completely edentulous patient. It is intended for use to support a tooth prosthesis, in the mandible or maxilla. The abutment screw is intended to secure the abutment to the endosseous dental implant.

Indications for use:

Indications for Use: GC Aadva™ Abutments - Zirconia - Bo Ac Bh Sv are dental implant abutments for use with partially or fully edentulous patients to restore chewing function by attachment to a dental implant fixture placed in the maxilla or mandible. Each abutment is accompanied by a screw in order to engage corresponding dental implant fixture.

Abutment - Implant Compatibility Table

Types	Material	Antirotational Features	Implant - Engaged Fixture (Manufacturer)
IN-Bo-3.3mm #154	Zirconia : Y-TZP	Internal Square	Φ3.3 Bone Level Implant NC (Institut straumann AG) K062129
IN-Bo-4.1mm #155	Zirconia : Y-TZP	Internal Square	Φ4.1, Φ4.8 Bone Level Implant RC (Institut straumann AG) K062129
IN-Ac-3.5mm #156	Zirconia : Y-TZP	Internal Hexagon	Φ3.5 NobelActive NP (Nobel Biocare AB) K071370

IN-Ac-4.3mm #157	Zirconia : Y-TZP	Internal Hexagon	Φ4.3, Φ5.0 NobelActive RP (Nobel Biocare AB) K071370
IN-Bh-3.8mm #158	Zirconia : Y-TZP	Internal Hexagon	Φ3.8 Tapered Internal Implant (BioHorizons) K071638
IN-Bh-4.6mm #159	Zirconia : Y-TZP	Internal Hexagon	Φ4.6 Tapered Internal Implant (BioHorizons) K071638
IN-Sv-3.5mm #151	Zirconia : Y-TZP	Internal Hexagon	Φ3.7, Φ4.1 Tapered Screw-Vent (Zimmer Dental Inc.) K013227
IN Sv-4.5mm #152	Zirconia : Y-TZP	Internal Hexagon	Φ4.7 Tapered Screw-Vent (Zimmer Dental Inc.) K013227

Prescription Use Only.

Substantial Equivalence:

The applicant device is substantially equivalent to the predicate devices in its intended use, indications for use, and design, as described below:

Predicate devices list

Product	Manufacturer	K Number
Atlantis Straumann Bone Level Abutment	Astra Tech Inc. 25 First Street Cambridge, Massachusetts	K083871

	02141 USA	
NobelActive Zirconia Abutment	Nobel Biocare USA, LLC. 22715 Savi Ranch Parkway, Yorba Londa CA 92887 USA	K072129
GC Aadvia Abutment	GC Corporation 76-1 Hasunuma-Cho, Itabashi-Ku Tokyo, Japan	K072100
Atlantis Abutment in Zirconia for BioHorizons Implant	Atlantis Components Inc. 25 First Street Cambridge, Massachusetts 02141 USA	K073540
Atlantis Abutment for BioHorizons Implant	Atlantis Components Inc. 25 First Street Cambridge, Massachusetts 02141 USA	K073258
ATLANTIS ABUTMENT FOR ZIMMER INTERFACE	Astra Tech Inc. 25 First Street Cambridge, Massachusetts 02141 USA	K053373

Materials used in blocks and screws

Component	Material	Predicate devices using same materials specified here
Block	Zirconia : Y-TZP ISO : 13356	Atlantis Straumann Bone Level Abutment (K083871)
	ZrO ₂ +HfO ₂ +Y ₂ O ₃ > 99.0 Y ₂ O ₃ : 4.5~5.4 Hf ₂ O ₃ <5	NobelActive Zirconia Abutment (K072129) GC AADVA ABUTMENT (K072100)

	Al ₂ O ₃ <0.5 Other oxides<0.5	Atlantis abutment in zirconia for 3i certain interface (K063734) Atlantis Abutment in Zirconia for BioHorizons Implant (K073540) ATLANTIS ABUTMENT FOR ZIMMER INTERFACE (K053373)
Screw	Ti-6Al-4V ISO : 5832-3 Ti : Bal. Al : 5.5~6.75 V : 3.5~4.5 Fe<0.3 O<0.2 C<0.08 N<0.05 H<0.015	Atlantis Straumann Bone Level Abutment (K083871) NobelActive Internal Connection (K071370) GC AADVA ABUTMENT (K072100) Atlantis Abutment in Zirconia for BioHorizons Implant (K073540) Atlantis Abutment for BioHorizons Implant (K073258) ATLANTIS ABUTMENT FOR ZIMMER INTERFACE (K053373)

Technological Characteristic	Applicant Device	Predicate Devices
Material	Zirconia : Y-TZP ISO : 13356	Zirconia : Y-TZP ISO : 13356
Performance Characteristics	Allows the prosthesis to be retained to the abutment; abutment screw is intended to	Allows the prosthesis to be retained to the abutment; abutment screw is intended to

	secure the abutment to the endosseous dental implant.	secure the abutment to the endosseous dental implant.
Intended Use	Intended for use with an endosseous dental implant to support a prosthetic device in a partially or fully edentulous patient to restore chewing function. Intended to be attached to a dental implant fixture placed in the maxilla or mandible.	Intended for use with an endosseous dental implant to support a prosthetic device in a partially or fully edentulous patient to restore chewing function. Intended to be attached to a dental implant fixture placed in the maxilla or mandible.

Summary of Non-Clinical Performance Testing:

Static and fatigue testing was conducted on the “worst case scenario” implant-abutment combination assemblies in accordance with FDA’s Guidance Document for Dental Implants and ISO 14801. Test results demonstrated that the “GC Aadva™ Abutments – Zirconia – Bo Ac Bh Sv” are compatible with the referenced implant fixtures and the implant-abutment assemblies support adequate static and fatigue test loads. Performance testing demonstrates that the device performs as intended and is as safe and effective as the cited predicates.

Substantial Equivalence Conclusion Statement:

Based on noted similarities and comparative traits in the indications for use, manufacturing materials, design and performance characteristics, and the fact that the applicant device and the predicate devices have demonstrated acceptable static and fatigue test loads during performance testing conducted in accordance with FDA’s Guidance Document for Dental Implants and ISO 14801, the “GC Aadva™ Abutments – Zirconia – Bo Ac Bh Sv” is substantially equivalent to the following predicate devices:

Product	Manufacturer	K Number
Atlantis Straumann Bone Level Abutment	Astra Tech Inc. 25 First Street Cambridge, Massachusetts 02141 USA	K083871
NobelActive Zirconia Abutment	Nobel Biocare USA, LLC. 22715 Savi Ranch Parkway, Yorba Londa CA 92887 USA	K072129
GC Aadva Abutment	GC Corporation 76-1 Hasunuma-Cho, Itabashi-Ku Tokyo, Japan	K072100

K113816

Atlantis Abutment in Zirconia for BioHorizons Implant	Atlantis Components Inc. 25 First Street Cambridge, Massachusetts 02141 USA	K073540
Atlantis Abutment for BioHorizons Implant	Atlantis Components Inc. 25 First Street Cambridge, Massachusetts 02141 USA	K073258
ATLANTIS ABUTMENT FOR ZIMMER INTERFACE	Astra Tech Inc. 25 First Street Cambridge, Massachusetts 02141 USA	K053373



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center – WO66-G609
Silver Spring, MD 20993-0002

January 30, 2013

GC America, Incorporated
C/O Mr. Carl H. Jenkins
The Wood Burditt Group
10 East Scranton Avenue, Suite 201
LAKE BLUFF IL 60044

Re: K113816

Trade/Device Name: GC Aadva™ Abutments – Zirconia – Bo Ac Bh Sv
Regulation Number: 21 CFR 872.3630
Regulation Name: Endosseous Dental Implant Abutment
Regulatory Class: II
Product Code: NHA
Dated: January 17, 2013
Received: January 25, 2013

Dear Mr. Jenkins:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.


If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

N
hor


Anthony D. Watson, B.S., M.S., M.B.A.
Director
Division of Anesthesiology, General Hospital,
Respiratory, Infection Control and
Dental Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

K 113816

Section 5 – Indications for Use Statement

Indications for Use

510(k) Number (if known): K113816

Device Name: GC Aadva™ Abutments – Zirconia – Bo Ac Bh Sv

Indications for Use: GC Aadva™ Abutments – Zirconia – Bo Ac Bh Sv are dental implant abutments for use with partially or fully edentulous patients to restore chewing function by attachment to a dental implant fixture placed in the maxilla or mandible. Each abutment is accompanied by a screw in order to engage corresponding dental implant fixture.

(Cont'd on next page)

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF
NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Susan Runner DDS, MA 2013.01.30
08:08:47 -05'00'

(Division Sign-Off)
Division of Anesthesiology, General Hospital
Infection Control, Dental Devices

510(k) Number: K113816

Indications for Use

510(k) Number (if known): K113816

Device Name: GC Aadva™ Abutments – Zirconia – Bo Ac Bh Sv

(Cont'd from previous page)

Abutment – Implant Compatibility Table

Types	Material	Antirotational Features	Implant - Engaged Fixture (Manufacturer)
IN-Bo-3.3mm #154	Zirconia : Y-TZP	Internal Square	Φ3.3 Bone Level Implant NC (Institut straumann AG) K062129
IN-Bo-4.1mm #155	Zirconia : Y-TZP	Internal Square	Φ4.1, Φ4.8 Bone Level Implant RC (Institut straumann AG) K062129
IN-Ac-3.5mm #156	Zirconia : Y-TZP	Internal Hexagon	Φ3.5 NobelActive NP (Nobel Biocare AB) K071370
IN-Ac-4.3mm #157	Zirconia : Y-TZP	Internal Hexagon	Φ4.3, Φ5.0 NobelActive RP (Nobel Biocare AB) K071370
IN-Bh-3.8mm #158	Zirconia : Y-TZP	Internal Hexagon	Φ3.8 Tapered Internal Implant (BioHorizons) K071638
IN-Bh-4.6mm #159	Zirconia : Y-TZP	Internal Hexagon	Φ4.6 Tapered Internal Implant (BioHorizons) K071638
IN-Sv-3.5mm #151	Zirconia : Y-TZP	Internal Hexagon	Φ3.7, Φ4.1 Tapered Screw-Vent (Zimmer Dental Inc.) K013227
IN Sv-4.5mm #152	Zirconia : Y-TZP	Internal Hexagon	Φ4.7 Tapered Screw-Vent (Zimmer Dental Inc.) K013227