

K121329

**510(k) Summary of Safety and Effectiveness**

AUG 24 2012

This summary of safety and effectiveness information is submitted in accordance with 21 CFR §807.92.

**1) Submitter's name, address, telephone number, contact person**

Submitted by:  
SuperSonic Imagine, S.A.  
Les Jardins de la Duranne – Bât. E & F  
510, rue René Descartes  
13857 Aix-en-Provence Cedex  
France  
Telephone: 011 33 442 99 24 24

Distributed by:  
SuperSonic Imagine, Inc.  
11714 North Creek Parkway N  
Suite 150  
Bothell, WA 98011  
North America  
Telephone: +1(425) 686 6380

**Corresponding Official:**

Jacques Souquet  
Chief Executive Officer  
Telephone: 011 33 442 99 24 35

Date: 2012/08/21

**2) Name of the device, including the trade or proprietary name if applicable, the common or usual name, and the classification name, if known:**

Common/Usual Name: Diagnostic Ultrasound System with Accessories  
Proprietary Name: Aixplorer®

Classification:

Regulatory Class: II

Classification Name:	21 CFR Section	Product Code
Ultrasonic Pulsed Doppler Imaging System	892.1550	90-IYN
Ultrasonic Pulsed Echo Imaging System	892.1560	90-IYO
Diagnostic Ultrasound Transducer	892.1570	90-ITX

**3) Substantially Equivalent/Predicate Devices**

AIXPLORER® Ultrasound Imaging System (K102041), cleared on 10/13/2010  
Siemens Acuson S2000TM Diagnostic Ultrasound System (K072786), cleared on 11/13/2007  
Philips iU22 Ultrasound System (K093563), cleared on 02/01/2010

#### 4) Description of Device

The SuperSonic Imagine AIXPLORER® system is a cart based ultrasound imaging system used to perform non-invasive diagnostic general purpose ultrasound imaging studies. The system contains a scan converter and can be coupled to a variety of linear, curved, micro-convex, and motorized linear array transducers to produce images, which are displayed on a LCD monitor. An adjustable control panel with integrated touch screen allows the user to perform an ultrasound exam quickly and efficiently in accordance with ALARA principles. The system also allows the user to perform measurements, capture images to digital memory or to an external device (such as a printer), and review diagnostic studies in the form of a report. The system functions in a manner identical to the predicate devices and transducers for the imaging modes: B-Mode, Color Flow, Pulsed Wave Doppler, Harmonic Imaging, Amplitude Doppler, 3D imaging and for ShearWave™ elastography.

#### 5) Intended Use

The SuperSonic Imagine AIXPLORER® ultrasound system and transducer are intended for general purpose pulse echo ultrasound imaging and Doppler fluid flow analysis of the human body.

#### 6) Indication for Use

The SuperSonic Imagine AIXPLORER® ultrasound system is indicated for use in the following applications: Abdominal, Small Organs, Musculoskeletal, Superficial Musculoskeletal, Vascular, Peripheral Vascular, GYN, Pelvic, Pediatric, Urology, Trans-rectal, Trans-vaginal and Neonatal Cephalic.

The system also provides the ability to measure anatomical structures (Abdominal, Small Organs, Musculoskeletal, Superficial Musculoskeletal, Peripheral Vascular, GYN, Pelvic, Pediatric, Urology, Trans-rectal, Trans-vaginal, Neonatal Cephalic).

#### 7) Safety Considerations

As a Track 3 ultrasound device, the SuperSonic Imagine AIXPLORER® ultrasound system is designed to comply with the "Standard For Real Time Display Of Thermal And Mechanical Acoustic Output Indices On Diagnostic Ultrasound Equipment" AIUM/NEMA 2004a published by the National Electrical Manufacturers Association as UD -3. With respect to limits on acoustic outputs, the SuperSonic Imagine AIXPLORER® ultrasound system complies with the FDA guideline limits set in the September 9, 2008, 510(k) diagnostic ultrasound guidance.

With regard to general safety, the SuperSonic Imagine AIXPLORER® ultrasound system scanner is designed to comply with IEC 60101 -1 (2005) Medical Electrical Equipment, Part 1: General Requirements for Basic Safety and Essential Performance, and IEC 60601 - 2-37 (2007): Particular Requirements For The Safety Of Ultrasonic Medical Diagnostic And Monitoring Equipment.

The device's acoustic output limits are:

Mechanical Index	1.9 (Maximum)
TIS/TIB	0.1 – 4.0 (Range)
ISPTA (d)	720 mW/cm <sup>2</sup>
ISPPA (d)	0 – 700 W/cm <sup>2</sup>

The limits are the same as predicate Track 3 devices. These considerations apply to all modes the system offers.

#### 8) Comparison to Predicate Devices

The SuperSonic Imagine AIXPLORER® system and transducers are substantially equivalent to the predicate devices with regard to intended use, imaging capabilities, safety and effectiveness.

- The systems are all intended for diagnostic ultrasound imaging and fluid flow analysis.
- The systems have the same clinical indications for use.

- The systems have the same B-Mode (grayscale imaging) and Doppler capabilities.
- The systems have similar capability in terms of harmonic imaging, spatial compound imaging, elastography imaging and other image post-processing features to improve the image quality and aid in clinical evaluation and diagnosis.
- The transducers are similar in materials, manufacture and clinical capability.
- The systems are manufactured with materials which have been evaluated and found to be safe for the intended use of the device.
- The systems have acoustic power levels which are below the applicable FDA limits.
- The systems have similar capability in terms of performing measurements, capturing digital images, reviewing and reporting studies.
- The systems have been found to be manufactured in compliance with approved electrical and physical safety standards.

**9) Nonclinical Performance Data**

Non-clinical testing was conducted per the following standards to support a determination of substantial equivalence to the predicate devices.

Reference Standard	Tests Performed
IEC 60601-1 3 <sup>rd</sup> Edition	All applicable electrical, basic safety and essential performance tests.
UL 60601-1 1 <sup>st</sup> Edition	All applicable electrical, basic safety and essential performance tests specific to the U.S.A.
IEC 60601-1-1 2 <sup>nd</sup> Edition	All applicable tests pertaining to Medical Electrical Systems.
IEC 60601-1-2 3 <sup>rd</sup> Edition	All applicable testing pertaining to electromagnetic compatibility.
IEC 60601-2-37 2 <sup>nd</sup> Edition	All applicable testing pertaining to the particular requirements for the safety of ultrasonic medical diagnostic and monitoring equipment.
NEMA UD 2 (Rev. 3)	All tests applicable in order to demonstrate compliance with the "Acoustic Output Measurement Standard for Diagnostic Ultrasound Equipment".
NEMA UD 3 (Rev. 2)	All tests applicable in order to demonstrate compliance with the "Standard For Real Time Display Of Thermal And Mechanical Acoustic Output Indices On Diagnostic Ultrasound Equipment".
ISO 10993-1	Applicable biocompatibility tests per FDA 510(k) Memorandum - #G95-1 – per the appropriate device category.

In addition to the referenced standards testing, SuperSonic Imagine conducted the following performance tests with respect to the Panoramic and High PRF features:

Panoramic Measurement Verification Tests
High PRF Measurement Verification Tests

The above testing confirmed that the Aixplorer System performs according to the stated intended use. All data fell within pre-determined product specifications and external standard requirements. Results of non-clinical testing confirmed the substantial equivalence of the Aixplorer System to the predicate device(s).

## **10) Conclusion**

The documentation provided demonstrates that:

- 1) The system and transducers are substantially equivalent to the predicate devices.
- 2) There are no new questions of safety and effectiveness concerning the SuperSonic Imagine AIXPLORER® ultrasound system and transducers.
- 3) The ultrasound device has been scientifically evaluated and has been demonstrated to be at least as safe and effective as the predicate devices cited in item 3.  
The system's acoustic power levels are below the applicable FDA limits.



10903 New Hampshire Avenue  
Silver Spring, MD 20993

AUG 24 2012

Jacques Souquet, Ph.D.  
Chief Executive Officer  
SUPERSONIC imagine  
Les Jardins de la Duranne  
510, rue Rene Descartes – Bat. E et F  
13 857 AIX-EN-PROVENCE CEDEX  
FRANCE

Re: K121329

Trade/Device Name: AIXPLORER® Ultrasound System  
Regulation Number: 21 CFR 892.1550  
Regulation Name: Ultrasonic pulsed doppler imaging system  
Regulatory Class: II  
Product Code: IYN, IYO, and ITX  
Dated: July 2, 2012  
Received: July 2, 2012

Dear Dr. Souquet:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

This determination of substantial equivalence applies to the following transducers intended for use with the AIXPLORER® Ultrasound System, as described in your premarket notification:

Transducer Model Number

SL15-4 (1D Linear Array)  
SC6-1 (Curved Array)  
SE12-3 (Endocavitary)

SLV16-5 (Motorized Linear)  
SL10-2 (Linear)  
SMC12-3 (Micro-Curved)

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus permits your device to proceed to market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

If you have any questions regarding the content of this letter, please contact Shahram Vaezy at (301) 796-6242.

Sincerely Yours



Janine M. Morris  
Acting Director  
Division of Radiological Devices  
Office of In Vitro Diagnostic Device  
Evaluation and Safety  
Center for Devices and Radiological Health

Enclosure(s)

## Diagnostic Ultrasound Indications for Use

510(k) number (if known):

Device Name: AIXPLORER® Ultrasound System

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Indications for Use:

The SuperSonic Imagine AIXPLORER® ultrasound system is indicated for use in the following applications: Abdominal, Small Organs, Musculoskeletal, Superficial Musculoskeletal, Vascular, Peripheral Vascular, GYN, Pelvic, Pediatric, Urology, Trans-rectal, Trans-vaginal and Neonatal Cephalic.

The system also provides the ability to measure anatomical structures (Abdominal, Small Organs, Musculoskeletal, Superficial Musculoskeletal, Peripheral Vascular, GYN, Pelvic, Pediatric, Urology, Trans-rectal, Trans-vaginal, Neonatal Cephalic).


Prescription Use XX OR Over-The-Counter Use \_\_\_\_\_

(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page 1 of 11

  
\_\_\_\_\_  
(Division Sign-Off)  
Division of Radiological Devices  
A121329  
5/20/11

# Diagnostic Ultrasound Indications for Use

510(k) number (if known):

Device Name: AIXPLORER® Ultrasound System

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation						
General (Track 1 Only)	Specific (Tracks 1 & 3)	B	M	PWD	CWD	Color Doppler	Combined (Specify)	Other* (Specify)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal							
	Abdominal (including urology)	P		P		P	P, 1, 2, 3, N, 4	P, 5, 6, 7, 8, 10, N, 9
	Intra-operative (Specify)							
	Intra-operative (Neuro)							
	Laparoscopic							
	Pediatric	P		P		P	P, 1, 2, 3, N, 4	P, 5, 6, 7, 8, 10, N, 9
	Small Organ (Breast, Thyroid, Testicle, Prostate, penis, etc...)	P		P		P	P, 1, 2, 3, N, 4	P, 5, 6, 7, 8, 10, N, 9
	Neonatal Cephalic	N		N		N	N, 1, 2, 3, 4	N, 5, 6, 7, 9
	Adult Cephalic							
	Trans-rectal	P		P		P	P, 1, 2, 3, N, 4	P, 5, 6, 7, 8
	Trans-vaginal	P		P		P	P, 1, 2, 3, N, 4	P, 5, 6, 7, 8
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skeletal (Conventional)	P		P		P	P, 1, 2, 3, N, 4	P, 5, 6, 7, 8, 10, N, 9
	Musculo-skeletal (Superficial)	P		P		P	P, 1, 2, 3, N, 4	P, 5, 6, 7, 8, 10, N, 9
	Intravascular							
GYN	P		P		P	P, 1, 2, 3, N, 4	P, 5, 6, 7, 8, 10	
Pelvic	P		P		P	P, 1, 2, 3, N, 4	P, 5, 6, 7, 8, 10	
Other (Specify)								
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Intravascular (Cardiac)							
	Trans-esoph. (Cardiac)							
	Intra-cardiac							
Other (Specify)								
Peripheral	Peripheral vessel	P		P		P	P, 1, 2, 3, N, 4	P, 5, 6, 7, 8, 10, N, 9
Vessel	Other (Specify)	P		P		P	P, 1, 2, 3, N, 4	P, 5, 6, 7, 8, N, 9

N = new indication; P = previously cleared by FDA (K102041)

**Additional Comments:**

- 1: Combined modes include: B+ Color Flow
- 2: Combined modes include: B+ ShearWave™ Elastography
- 3: Combined modes include: B+ Pulsed Wave
- 4: Combined modes include: B+ Pulsed Wave + Color Flow
- 5: Harmonic Imaging
- 6: Spatial Compounding
- 7: ShearWave™ Elastography
- 8: Imaging Guidance for Biopsies
- 9: Panoramic Imaging
- 10: 3D Imaging

Prescription Use   XX   OR Over-The-Counter Use       

(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

[Signature]

(Division Sign-Off)  
Division of Radiological Devices

K121329

# Diagnostic Ultrasound Indications for Use

510(k) Number (if known):

Device Name: SL15-4 transducer (1D Linear Array Transducer)

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation							
General (Track 1 Only)	Specific (Tracks 1 & 3)	B	M	PWD	CWD	Color Doppler	Combined (Specify)	Other* (Specify)	
Ophthalmic	Ophthalmic								
Fetal Imaging & Other	Fetal								
	Abdominal	P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8, N, 9	
	Intra-operative (Specify)								
	Intra-operative (Neuro)								
	Laparoscopic								
	Pediatric	P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8, N, 9	
	Small Organ (for example Breast, Thyroid, Testicle, Prostate, Penis)	P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8, N, 9	
	Neonatal Cephalic	N		N		N	N, 1, 2, 3, 4	N, 5, 6, 7, 9	
	Adult Cephalic								
	Trans-rectal								
	Trans-vaginal								
	Trans-urethral								
	Trans-esoph. (non-Card.)								
	Musculo-skeletal (Conventional)	P		P			P	P 1, 2, 3, N, 4	P 5, 6, 7, 8, N, 9
	Musculo-skeletal (Superficial)	P		P			P	P 1, 2, 3, N, 4	P 5, 6, 7, 8, N, 9
	Intravascular								
	GYN								
	Pelvic								
Other (Specify)									
Cardiac	Cardiac Adult								
	Cardiac Pediatric								
	Intravascular (Cardiac)								
	Trans-esoph. (Cardiac)								
	Intra-cardiac								
	Other (Specify)								
Peripheral	Peripheral vessel	P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8, N, 9	
Vessel	Other (Specify)								

N = new indication; P = previously cleared by FDA (K102041)

Additional Comments:


- 1: Combined modes include: B+ Color Flow
- 2: Combined modes include: B+ ShearWave™ Elastography
- 3: Combined modes include: B+ Pulsed Wave
- 4: Combined modes include: B+ Pulsed Wave + Color Flow
- 5: Harmonic Imaging
- 6: Spatial Compounding
- 7: ShearWave™ Elastography
- 8: Imaging Guidance for Biopsies
- 9: Panoramic Imaging
- 10: 3D Imaging

Prescription Use XX OR Over-The-Counter Use \_\_\_\_\_

(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
 (Division Sign-Off)  
 Division of Radiological Devices  
 510k 15121329

# Diagnostic Ultrasound Indications for Use

510(k) Number (if known):

Device Name: SC6-1 transducer (curved array transducer)

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation						
General (Track 1 Only)	Specific (Tracks 1 & 3)	B	M	PWD	CWD	Color Doppler	Combined (Specify)	Other* (Specify)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal							
	Abdominal (including urology)	P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8, N, 9
	Intra-operative (Specify)							
	Intra-operative (Neuro)							
	Laparoscopic							
	Pediatric	P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8, N, 9
	Small Organ (Breast, Thyroid, Testicle, Prostate, penis, etc...)	P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skeletal (Conventional)	P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8, N, 9
	Musculo-skeletal (Superficial)	P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8, N, 9
Intravascular								
GYN	P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8	
Pelvic	P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8	
Other (Specify)								
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Intravascular (Cardiac)							
	Trans-esoph. (Cardiac)							
	Intra-cardiac							
Other (Specify)								
Peripheral	Peripheral vessel	P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8, N, 9
Vessel	Other (Specify)	P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8, N, 9

N = new indication; P = previously cleared by FDA (K102041)

**Additional Comments:**

- 1: Combined modes include: B+ Color Flow
- 2: Combined modes include: B+ ShearWave™ Elastography
- 3: Combined modes include: B+ Pulsed Wave
- 4: Combined modes include: B+ Pulsed Wave + Color Flow
- 5: Harmonic Imaging
- 6: Spatial Compounding
- 7: ShearWave™ Elastography
- 8: Imaging Guidance for Biopsies
- 9: Panoramic Imaging
- 10: 3D Imaging

Prescription Use   XX   OR Over-The-Counter Use       

(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

[Signature]

(Division Sign-Off)  
Division of Radiological Devices  
DVR

510k B121329

# Diagnostic Ultrasound Indications for Use

510(k) Number (if known):

Device Name: SE12-3 transducer (endocavitary transducer)

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation							
General (Track 1 Only)	Specific (Tracks 1 & 3)	B	M	PWD	CWD	Color Doppler	Combined (Specify)	Other* (Specify)	
Ophthalmic	Ophthalmic								
Fetal Imaging & Other	Fetal								
	Abdominal								
	Intra-operative (Specify)								
	Intra-operative (Neuro)								
	Laparoscopic								
	Pediatric								
	Small Organ (for example Breast, Thyroid, Testicle, Prostate, penis, etc...)	P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8	
	Neonatal Cephalic								
	Adult Cephalic								
	Trans-rectal	P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8	
	Trans-vaginal	P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8	
	Trans-urethral								
	Trans-esoph. (non-Card.)								
	Musculo-skeletal (Conventional)								
	Muscuio-skeletal (Superficial)								
	Intravascular								
	GYN		P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8
	Pelvic		P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8
	Other (Specify)								
	Cardiac	Cardiac Adult							
Cardiac Pediatric									
Intravascular (Cardiac)									
Trans-esoph. (Cardiac)									
Intra-cardiac									
Other (Specify)									
Peripheral Vessel	Peripheral vessel								
Vessel	Other (Specify)	P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8	

N = new indication; P = previously cleared by FDA (K102041)

**Additional Comments:**


- 1: Combined modes include: B+ Color Flow
- 2: Combined modes include: B+ ShearWave™ Elastography
- 3: Combined modes include: B+ Pulsed Wave
- 4: Combined modes include: B+ Pulsed Wave + Color Flow
- 5: Harmonic Imaging
- 6: Spatial Compounding
- 7: ShearWave™ Elastography
- 8: Imaging Guidance for Biopsies
- 9: Panoramic Imaging
- 10: 3D Imaging

Prescription Use XX OR Over-The-Counter Use \_\_\_\_\_

(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
 (Division Sign-Off)  
 Division of Radiological Devices  
 ODP  
 510k K121389

# Diagnostic Ultrasound Indications for Use

510(k) Number (if known):

Device Name: SLV16-5 transducer (motorized linear transducer)

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation							
General (Track 1 Only)	Specific (Tracks 1 & 3)	B	M	PWD	CWD	Color Doppler	Combined (Specify)	Other* (Specify)	
Ophthalmic	Ophthalmic								
Fetal Imaging & Other	Fetal								
	Abdominal	P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8, 10, N, 9	
	Intra-operative (Specify)								
	Intra-operative (Neuro)								
	Laparoscopic								
	Pediatric	P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8, 10, N, 9	
	Small Organ (for example Breast, Thyroid, Testicle, Prostate, penis, etc...)	P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8, 10, N, 9	
	Neonatal Cephalic								
	Adult Cephalic								
	Trans-rectal								
	Trans-vaginal								
	Trans-urethral								
	Trans-esoph. (non-Card.)								
	Musculo-skeletal (Conventional)	P		P			P	P 1, 2, 3, N, 4	P 5, 6, 7, 8, 10, N, 9
	Musculo-skeletal (Superficial)	P		P			P	P 1, 2, 3, N, 4	P 5, 6, 7, 8, 10, N, 9
Intravascular									
GYN									
Pelvic									
Other (Specify)									
Cardiac	Cardiac Adult								
	Cardiac Pediatric								
	Intravascular (Cardiac)								
	Trans-esoph. (Cardiac)								
	Intra-cardiac								
	Other (Specify)								
Peripheral	Peripheral vessel	P		P		P	P 1, 2, 3, N, 4	P 5, 6, 7, 8, 10, N, 9	
Vessel	Other (Specify)								

N = new indication; P = previously cleared by FDA (K102041)

Additional Comments:

- 1: Combined modes include: B+ Color Flow
- 2: Combined modes include: B+ ShearWave™ Elastography
- 3: Combined modes include: B+ Pulsed Wave
- 4: Combined modes include: B+ Pulsed Wave + Color Flow
- 5: Harmonic Imaging
- 6: Spatial Compounding
- 7: ShearWave™ Elastography
- 8: Imaging Guidance for Biopsies
- 9: Panoramic Imaging
- 10: 3D Imaging

Prescription Use   XX   OR Over-The-Counter Use           

(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
 (Division Sign-Off)  
 Division of Radiological Devices  
 510k K121329

# Diagnostic Ultrasound Indications for Use

510(k) Number (if known):

Device Name: SL10-2 transducer (linear transducer)

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation						
General (Track 1 Only)	Specific (Tracks 1 & 3)	B	M	PWD	CWD	Color Doppler	Combined (Specify)	Other* (Specify)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal							
	Abdominal	N		N		N	N, 1, 2, 3, 4	N, 5, 6, 7, 8, 9
	Intra-operative (Specify)							
	Intra-operative (Neuro)							
	Laparoscopic							
	Pediatric	N		N		N	N, 1, 2, 3, 4	N, 5, 6, 7, 8, 9
	Small Organ (for example Breast, Thyroid, Testicle, Prostate, penis, etc...)	N		N		N	N, 1, 2, 3, 4	N, 5, 6, 7, 8, 9
	Neonatal Cephalic	N		N		N	N, 1, 2, 3, 4	N, 5, 6, 7, 9
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skeletal (Conventional)	N		N		N	N, 1, 2, 3, 4	N, 5, 6, 7, 8, 9
	Musculo-skeletal (Superficial)	N		N		N	N, 1, 2, 3, 4	N, 5, 6, 7, 8, 9
	Intravascular							
	GYN							
Pelvic								
Other (Specify)								
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Intravascular (Cardiac)							
	Trans-esoph. (Cardiac)							
	Intra-cardiac							
Other (Specify)								
Peripheral	Peripheral vessel	N		N		N	N, 1, 2, 3, 4	N, 5, 6, 7, 8, 9
Vessel	Other (Specify)	N		N		N	N, 1, 2, 3, 4	N, 5, 6, 7, 8, 9

N = new indication; P = previously cleared by FDA (K102041)

Additional Comments:


- 1: Combined modes include: B+ Color Flow
- 2: Combined modes include: B+ ShearWave™ Elastography
- 3: Combined modes include: B+ Pulsed Wave
- 4: Combined modes include: B+ Pulsed Wave + Color Flow
- 5: Harmonic Imaging
- 6: Spatial Compounding
- 7: ShearWave™ Elastography
- 8: Imaging Guidance for Biopsies
- 9: Panoramic Imaging
- 10: 3D Imaging

Prescription Use XX OR Over-The-Counter Use \_\_\_\_\_

(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
 (Division Sign-Off)  
 Division of Radiological Devices  
 510k 5121329

# Diagnostic Ultrasound Indications for Use

510(k) Number (if known):

Device Name: SMC12-3 transducer (micro-curved transducer)

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

Clinical Application		Mode of Operation						
General (Track 1 Only)	Specific (Tracks 1 & 3)	B	M	PWD	CWD	Color Doppler	Combined (Specify)	Other* (Specify)
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal							
	Abdominal	N		N		N	N, 1, 2, 3, 4	N, 5, 6, 7, 8, 9
	Intra-operative (Specify)							
	Intra-operative (Neuro)							
	Laparoscopic							
	Pediatric	N		N		N	N, 1, 2, 3, 4	N, 5, 6, 7, 8, 9
	Small Organ (for example Breast, Thyroid, Testicle, Prostate, penis, etc...)	N		N		N	N, 1, 2, 3, 4	N, 5, 6, 7, 8, 9
	Neonatal Cephalic	N		N		N	N, 1, 2, 3, 4	N, 5, 6, 7, 9
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Trans-esoph. (non-Card.)							
	Musculo-skeletal (Conventional)	N		N		N	N, 1, 2, 3, 4	N, 5, 6, 7, 8, 9
	Musculo-skeletal (Superficial)	N		N		N	N, 1, 2, 3, 4	N, 5, 6, 7, 8, 9
	Intravascular							
GYN								
Pelvic								
Other (Specify)								
Cardiac	Cardiac Adult							
	Cardiac Pediatric							
	Intravascular (Cardiac)							
	Trans-esoph. (Cardiac)							
	Intra-cardiac							
Other (Specify)								
Peripheral	Peripheral vessel	N		N		N	N, 1, 2, 3, 4	N, 5, 6, 7, 8, 9
Vessel	Other (Specify)	N		N		N	N, 1, 2, 3, 4	N, 5, 6, 7, 8, 9

N = new indication; P = previously cleared by FDA (K102041)

Additional Comments:


- 1: Combined modes include: B+ Color Flow
- 2: Combined modes include: B+ ShearWave™ Elastography
- 3: Combined modes include: B+ Pulsed Wave
- 4: Combined modes include: B+ Pulsed Wave + Color Flow
- 5: Harmonic Imaging
- 6: Spatial Compounding
- 7: ShearWave™ Elastography
- 8: Imaging Guidance for Biopsies
- 9: Panoramic Imaging
- 10: 3D Imaging

Prescription Use XX OR Over-The-Counter Use \_\_\_\_\_

(Part 21 CFR 801 Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

  
 (Division Sign-Off)  
 Division of Radiological Devices  
 510k B121329