

510(K) Summary

SEP 7 2012

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

The assigned 510(k) number is: K121353

1. Submitter's Identifications:

Well Life Healthcare Limited
1F., No. 16, Lane 454, Jungjeng Road, Yunghe City, Taipei County 234, Taiwan, ROC

Contact: Jenny Hsieh
Telephone: + 886 2 2928 2112

Date of Summary Preparation: August 17, 2012

2. Name of the Device:

Trade/Device Name: Mini Patch, model : WL-2301A & WL-2301B.
Regulation Number: 21 CFR 882.5890
Regulation Name: Stimulator, Nerve, Transcutaneous For Pain Relief.
Regulatory Class: II -
Product Code: NUH

3. Information of the 510(k) Cleared Device (Predicate Device): Well-Life OTC TENS for Arm & Leg Pain Relief. (K091757)

4. Device Description:

The Mini Patch, models WL-2301A & WL-2301B are transcutaneous electrical nerve stimulator used for pain relief by applying an electrical current to electrodes, which are attached on the skin. The output and waveform characteristic is fixed for every model, only the intensity is adjustable.

The Mini Patch, models WL-2301A & WL-2301AB, consist mainly of two parts: the stimulus generator, electrode. The stimulus generator generates the output current specified as the input of controller. The output port transmits the output current to the electrode, which is attached to the patient's skin so as to transmit this stimulus current to the patient for pain relief.

The stimulation mode for Mini Patch is burst mode with fixed pulse width, pulse rate, frequency, and fixed timer, only amplitude is adjustable. This operation way is considered the simplification from the burst mode of comparison clear model, WL-2302(K020020). Every model of Mini-TENS has its individual stimulation operation cycle.

5. Intended Use:

The Mini Patch, model: WL-2301A & WL-2301B Transcutaneous Electrical Nerve Stimulators are intended for temporary relief of pain associated with sore and aching muscles in the upper and lower extremities (arm and/or leg) due to strain from exercise or normal household and work activities.

6. Discussion of Non-Clinical Tests Performed Determination of Substantial Equivalence are as follows:

Compliance to applicable voluntary standards includes ANSI/AAMI, NS4-1985, as well as IEC 60601-1, and IEC 60601-1-2 requirement.

In addition to the compliance of voluntary standards, the software verification has been carried out according to the FDA software guidance.

7. Comparison of Significant device features

Comparison Features	Models to be Compared	
	OTC TENS for Arm & Leg	Mini Patch
Model	WL-2407	WL-2301A & WL-2301B
510(K) No.	K091757	K121353
Prescription or OTC	OTC	Same
Indication for use	intended for temporary relief of pain associated with sore and aching muscles in the upper and lower extremities (arm and/or leg) due to strain from exercise or normal household and work activities	Same
FDA product code	NUH	Same
Electrode Used	K082065, Snap Type self adhesive electrode(5 X 5 cm)	K082065, Snap Type self adhesive electrode (2.925 x 5.6 cm)
Manufacturer	Well-Life	Same

1K121353

8. Comparison of Unit Characteristics & Output Specification

Parameter		Predicate Device	New Device	
510(K) Number		K091757	K121353	K121353
Device Name and Model		WL-2407	WL-2301A	WL-2301B
Manufacturer		Well-Life	Same	Same
Power Source(s)		1.5Vx3 (AAA Size)	3Vx1	3Vx1
- Method of Line current Isolation		Type BF	Same	Same
- Patient Leakage Current		---	---	---
- Normal condition (uA)		Under 0.1	Same	Same
- single Fault condition (uA)		Under 0.5	Same	Same
Average DC current through electrodes when device is on but no pulses are being applied (uA)		Not applicable	Same	Same
Number of Output Modes		8	3	3
Number of Output Channels:	Synchronous or Alternating?	Synchronous	Same	Same
	Method of Channel Isolation	Output Coil	Same	Same
Regulated Current or Regulated Voltage?		Voltage	Same	Same
Software/Firmware/Microprocessor control?		Yes	Same	Same
Automatic Overload Trip?		No	Same	Same
Automatic No-Load Trip?		Yes	Same	Same
Automatic Shut Off?		Yes	Same	Same
User Override control?		No	Same	Same
Indicator Display:	On/Off Status?	Yes	Same	Same
	Low Battery?	Yes	Same	Same
	Voltage/Current Level?	Yes	Same	Same
Timer Range (Minutes)		5-60	20	20
Compliance with Voluntary Standards?		ANSI/AAMI NS4, IEC 60601-1 & IEC 60601-1-2.	Same	Same
Compliance with 21 CFR 898?		Yes	Same	Same
Weight (g)		43	20	20
Dimensions (mm.) [W x H x D]		90x52.5x19.4	59x30x12.15	60.2x31.9x14.4
Housing Materials and construction		ABS	Same	Same

Parameter		Predicate Device	New Device	
Mode or Program Name		WL-2407	WL-2301A	WL-2301B
Waveform (e.g., pulsed monophasic, biphasic)		Biphasic	Same	Same
Shape (e.g., rectangular, spike, rectified sinusoidal)		Rectangular (with deformation on the top of each pulse)	Same	Same
Maximum Output Voltage (volts) (+/- 20 %)	40V	40V	30V	30V
	@500Ω	@500Ω	@500Ω	@500Ω
	62V	62V	50V	50V
	@2kΩ	@2kΩ	@2kΩ	@2kΩ
Maximum Output Current (specify units) (+/- 20 %)	104V	104V	80V	80V
	@10kΩ	@10kΩ	@10kΩ	@10kΩ
	80mA	80mA	60mA	60mA
	@500Ω	@500Ω	@500Ω	@500Ω
	31mA	31mA	25mA	25mA
	@2kΩ	@2kΩ	@2kΩ	@2kΩ
	10.4mA	10.4mA	8.0mA	8.0mA
	@10kΩ	@10kΩ	@10kΩ	@10kΩ
Duration of primary (despolarizing) phase (usec)		260 Max	250 Max	250 Max
Pulse Duration (usec)		500.26 ms/Max.	33.58 ms/Max.	33.58 ms/Max.
Frequency (Hz) [or Rate (pps)]		60 Max	60 Max	60 Max
For multiphasic waveforms only:	Symmetrical phases?	No	Same	Same
	Phase Duration (include units), (Stage range, if applicable), (both phases, if asymmetrical)	Not applicable	Same	Same

9. Significant output characteristics comparison table:

Comparison feature		Model		Remarks
		WL-2407	WL-2301/2301B	
Net charge		0	0	Note 1
Max. phase charge		20.8 uc	15uc	
Max. current Density		0.050 mA/cm ²	0.055 mA/cm ²	
Max. Average current (RMSA)	500 Ω	38.43mA	29.394 mA	
	2K Ω	14.892mA	12.248 mA	
	10 Ω	4.996mA	3.9192 mA	
Max. Power Density		0.0020Watts/ cm ²	0.0017Watts/ cm ²	
Burst Mode		No	6Sec/On & 1Sec/Off	

Note: Net Charge is Zero for the symmetrical stimulation biphasic wave

10. Conclusions

The Well-Life Mini Patch, model: WL-2301A & WL-2301B Transcutaneous Electrical Nerve Stimulators has the same intended use and the similar technological characteristics as the cleared device- Well-Life OTC TENS for Arm & Leg Pain Relief. (K091757) . Moreover, verification and validation tests contained in this submission demonstrate that the difference in the submitted demonstrate that the difference in the submitted models could maintain the same safety and effectiveness as that of cleared device.

In the other words, those engineering difference do not: (1) affect the intended use or (2) alter the fundamental scientific technology of the device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Room -WO66-G609
Silver Spring, MD 20993-0002

Well Life Healthcare, Ltd
% Ms. Jenny Hsieh
1Fl, No.16, Lane 454, Jungjeng Rd.
Yunghe City, Taipei County
Taiwan, R.O.C.

SEP 7 2012

Re: K121353
Trade/Device Name: Mini Patch, Model WL-2301A & WL-2301B
Regulation Number: 21 CFR 882.5890
Regulation Name: Transcutaneous Electrical Nerve Stimulator for Pain Relief
Regulatory Class: II
Product Code: NUH
Dated: August 17, 2012
Received: August 21, 2012

Dear Ms. Hsieh:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

for

Kesia Alexander

Malvina B. Eydelman, M.D.
Director
Division of Ophthalmic, Neurological,
and Ear, Nose and Throat Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications For Use

510(k) Number (if known): K121353

Device Name: Mini Patch, model : WL-2301A & WL-2301B

Indications For Use:

The Mini Patch, model: WL-2301A & WL-2301B Transcutaneous Electrical Nerve Stimulators are intended for temporary relief of pain associated with sore and aching muscles in the upper and lower extremities (arm and/or leg) due to strain from exercise or normal household and work activities


Prescription Use _____
(Part 21 CFR 801 Subpart D)

OR

Over-The-Counter Use √
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)

Division of Ophthalmic, Neurological and Ear,
Nose and Throat Devices

510(k) Number _____

K121353