

k122467

3M ESPE  
Dental Products

3M Center  
St. Paul, MN 55144-1000  
651 733 1110

NOV 20 2012



### 510(k) Summary

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of 21 CFR 807.92.

**510(k) Submitter**..... 3M Company  
3M ESPE Dental Products  
3M Center, 275-2W-08  
2510 Conway Avenue  
St. Paul, MN 55144-1000 USA

**Contact person**..... Ginger Cantor, RAC  
Regulatory Affairs Advanced Specialist  
Phone: (651) 733-1317  
Fax: (651) 737-9665  
[gcantor@mmm.com](mailto:gcantor@mmm.com)

**Date Summary was Prepared**.....November 16, 2012

**Trade Name**.....3M True Definition Scanner (Model G5)

**Common Name(s)**.....Intra-Oral Scanner

**Recommended Classification**.....Optical Impression Systems for  
CAD/CAM, Class II device,  
(21 CFR § 872.3661, Product Code: NOF)

**Predicate Devices:**

3M™ ESPE™ Lava Chairside Oral Scanner (C.O.S.), K073199

3M™ Unitek™ Lava Chairside Oral Scanner (C.O.S.), K081961

**Description of Device:**

The 3M True Definition Scanner is a digital impression generating system consisting of a computer system on a mobile cart, a lightweight scanning wand, and embedded software (including firmware). Accessory items include a contrast powder to be applied to the patient's teeth and/or oral anatomy and a battery powered powder sprayer.

The computer system consists of a commercial off-the-shelf personal computer (PC) and a touch screen monitor.

The scanning wand is a hand held optical device that captures high-resolution video images, in real time, as the patient is being scanned. The wand contains an optical system comprised of low intensity LED's, a light sensor, a lens and supporting electronics. The wand is connected to the cart via a high speed data transfer cable. The wand is designed to be easily maneuvered inside the patient's mouth and captures video imagery at 20 frames per second. Those images are converted to 3D data sets and displayed in real time.

The software contains high-speed image processing algorithms, real-time modeling, case management, and archival functionality.

3M True Definition Scanner facilitates a digital workflow which reduces or eliminates many steps traditionally required by the dentist and lab, including tray selection, plaster pouring, material dispensing, base & pin, material setting, die cutting, trimming, articulation, packaging and shipping.

**Indications for Use:**

3M True Definition Scanner is an optical impression system used to record the topographical characteristics of the dentition and/or full arch and preparation areas (including features such as implant scan locator fixtures, braces, brackets, etc.). In addition it can record the topographical characteristics of the oral anatomy (such as soft tissue, gingivae and palate).

The three dimensional (3D) model generated from the scan may be further used for the design and manufacturing of dental restorations including implant supported prostheses and partial frameworks, and can be used to design and manufacture physical models of the teeth.

It may also be used in conjunction with production of orthodontic appliances, retainers and accessories.

**Electrical Safety/Electromagnetic Compatibility (EMC)**

3M ESPE provided data in the 510(k) to demonstrate the device meets the requirements of ANSI/AAMI ES60601-1: 2005 (i.e., IEC 60601-1: 2005), *Medical electrical equipment – Part 1: General requirements for basic safety and essential performance*, IEC60601-1:1998 + A1:1992 + A1:1995, *Medical electrical equipment – Part 1: General requirements for safety*, and IEC 60601-1-2: 2007, *Medical Electrical Equipment –Part 1: General Requirements for Safety: Collateral Standard: Electromagnetic Compatibility Requirements and Tests*.

In addition, information was provided in the 510(k) to demonstrate the 3M True Definition Scanner meets the requirements of IEC62471:2006, *Photobiological safety of lamps and lamp systems*.

3M concludes that the 3M True Definition Scanner is safe for its intended use/indications for use.

**Biocompatibility**

3M provided data in this 510(k) to demonstrate the 3M True Definition Scanner meets the requirements of relevant parts of ISO10993-1(2009), *Biological Evaluation of Medical Devices—Part 1: Evaluation and Testing within a Risk Management Process* and that the device meets the requirements of ISO 7405: 2008, *Dentistry -- Evaluation of biocompatibility of medical devices used in dentistry*.

**Non-Clinical Testing**

3M successfully performed software verification and validation activities including testing of the device's accuracy, case processing, data management, fault management, network access and remote service, patient, doctor and prescription management, scanning operations and workflows to laboratories and mills. Results of accuracy testing are shown in the Substantial Equivalence summary below.

**Substantial Equivalence:**

3M ESPE's evaluation of the substantial equivalence of the True Definition Scanner to the Lava C.O.S. predicates was based on a comparison of device classification, intended use and indications for use, contraindications, labeling, technological attributes, accuracy, and accessories. 3M also compared risks and End-of Useful Life for both devices.

<b>Comparison of Intended Use/Indications for Use/Contraindications</b>			
<b>Feature</b>	<b>3M True Definition Scanner (FDA 510(k)# K122467)</b>	<b>3M ESPE Lava C.O.S. (FDA 510(k)# K073199)</b>	<b>3M Unitek Lava C.O.S. (FDA 510(k)# K081961)</b>
Procode	NOF	NOF	NOF
Classification	Optical Impression System for CAD/CAM 21 CFR §872.3661	Optical Impression System for CAD/CAM 21 CFR §872.3661	Optical Impression System for CAD/CAM 21 CFR §872.3661
Indications for Use	<p>3M True Definition Scanner is an optical impression system used to record the topographical characteristics of the dentition and/or full arch and preparation areas (including items such as implant scan locator fixtures, braces, brackets, etc.). In addition, it can record the topological characteristics of the oral anatomy (such as soft tissue, gingivae, and palate).</p> <p>The three dimensional (3D) model generated from the scan may be further used for the design and manufacturing of dental restorations including implant supported prostheses and partial frameworks, and can be used to design and manufacture physical models of the teeth. It may also be used in conjunction with the production of orthodontic appliances, retainers and accessories.</p>	<p>The 3M ESPE Lava Chairside Oral Scanner (C.O.S.) is an optical impression system (CAD/CAM) used to record the topographical characteristics of teeth.</p>	<p>The 3M Unitek Lava Chairside Oral Scanner (C.O.S.) is an optical impression system (CAD/CAM) used to record the topographical characteristics of teeth.</p> <p>Data generated from the 3M Unitek Lava Chairside Oral Scanner may be used in conjunction with the production of orthodontic appliances, retainers and accessories.</p>
Contra-indications	None	None	None

<b>Technological Comparison</b>			
<b>Item</b>	<b>3M True Definition Scanner (K122467)</b>	<b>Lava C.O.S. (510(k)# K073199)</b>	<b>Lava C.O.S. (510(k)# K081961)</b>
Operating System	Fedora 15 (Linux)	Custom Gentoo Linux	Custom Gentoo Linux
Computer Hardware	Intel XEON processor – closed system	AMD based processor – closed system	AMD based processor – closed system
Internet Communications	Wireless 802.11 b/g/n	Wireless 802.11 b/g/n	Wireless 802.11 b/g/n
Scanning Wand Computer Interface	USB 3.0 Interface	Proprietary serial interface to custom PCI adapter	Proprietary serial interface to custom PCI adapter
Scanning Wand Imaging Sensor(s)	Single	Three	Three
Scanning Wand LEDs (# and wavelength)	6 LEDs (465 nm wavelength)	192 LEDs (465 nm wavelength)	192 LEDs (465 nm wavelength)
Isolation Transformer	Powertronix Isolation Station™	Custom built	Custom built

<b>Accuracy Comparison</b>		
<b>Parameter</b>	<b>3M True Definition Scanner (K122467)</b>	<b>Lava COS (K073199 and K081961)</b>
<b>Initial Average Error (%)</b>	0.13%	0.57%
<b>Initial Max Error (%)</b>	0.22%	0.74%

<b>Accessories Comparison</b>				
<b>Item</b>	<b>3M True Definition Scanner (K122467)</b>	<b>3M ESPE Lava C.O.S. (510(k)# K073199)</b>	<b>3M Unitek Lava C.O.S. (510(k) #081961)</b>	<b>Comments</b>
Contrast patterning Powder	Composition: <ul style="list-style-type: none"> <li>• titanium dioxide</li> <li>• amorphous silica</li> <li>• aluminium hydroxide</li> <li>• synthetic amorphous silica</li> </ul>	Composition: <ul style="list-style-type: none"> <li>• titanium dioxide</li> <li>• amorphous silica</li> <li>• aluminium hydroxide</li> <li>• synthetic amorphous silica</li> </ul>	Composition: <ul style="list-style-type: none"> <li>• titanium dioxide</li> <li>• amorphous silica</li> <li>• aluminium hydroxide</li> <li>• synthetic amorphous silica</li> </ul>	No change- the same powder described in K073199 was used in K081961 and will be used for the 3M True Definition Scanner
Powder delivery device	Battery (9V) powered sprayer	Battery (9V) powered sprayer	Battery (9V) powered sprayer	No change – the same sprayer described in K073199 was used in K081961 and will be used for the 3M True Definition Scanner

3M ESPE concludes that 3M True Definition Scanner is substantially equivalent to the named predicate devices (3M ESPE Lava C.O.S. (K073199) and 3M Unitek Lava C.O.S (K081961)). Any differences between the 3M True Definition Scanner and the predicate(s) do not raise new questions of safety or efficacy.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-002

November 20, 2012

3M Company  
C/O Mr. Mark Job  
Regulatory Technology Services, Limited Liability Company  
1394 25<sup>TH</sup> Street, North West  
Buffalo, Minnesota 55313

Re: K122467

Trade/Device Name: 3M True Definition Scanner (Model G5)  
Regulation Number: 21 CFR 872.3661  
Regulation Name: Optical Impression Systems for CAD/CAM  
Regulatory Class: II  
Product Code: NOF  
Dated: November 2, 2012  
Received: November 5, 2012

Dear Mr. Job:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.



## Indications for Use

**510(k) Number (if known):** K122467

**Device Name:** 3M True Definition Scanner (Model G5)

### Indications for Use:

3M True Definition Scanner is an optical impression system used to record the topographical characteristics of the dentition and/or full arch and preparation areas (including features such as implant scan locator fixtures, braces, brackets, etc.). In addition it can record the topographical characteristics of the oral anatomy (such as soft tissue, gingivae and palate).

The three dimensional (3D) model generated from the scan may be further used for the design and manufacturing of dental restorations including implant supported prostheses and partial frameworks, and can be used to design and manufacture physical models of the teeth.

It may also be used in conjunction with production of orthodontic appliances, retainers and accessories.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE OF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Susan Runner DDS, MA 2012.11.20 12:13:34  
05'00'

**(Division Sign-Off)**  
**Division of Anesthesiology, General Hospital**  
**Infection Control, Dental Devices**

**510(k) Number:**   K122467