

## 510(k) Summary

**Submission Date:** 26 September 2012

**Submitter:** Bracco Diagnostics, Inc.  
107 College Road  
Princeton, NJ 08540 USA

AUG 06 2013

**Submitter and Official Contact:** Bracco Diagnostics Inc.  
532 Broadhollow Road, Suite 126  
Melville, NY 11747

Ms. Tracey Alexander  
Director Regulatory Affairs  
Phone: +1 (609) 524-2865  
Fax: +1 (631) 847-3904  
Email: [tracey.alexander@diag.bracco.com](mailto:tracey.alexander@diag.bracco.com)

**Trade Name:** Bracco Diagnostics, Inc. CO<sub>2</sub> ENDOSCOPIC INSUFFLATOR IN LINE TUBING SET

**Common Name:** Tubing Set

**Classification Name:** Endoscope and accessories

**Classification Regulation:** 21 CFR §876.1500

**Product Code:** FCX, FDS

<b>Substantially Equivalent Devices:</b>	<b>BDI Model</b>	<b>Predicate 510(k) Number</b>	<b>Predicate Manufacturer and Model</b>
	Bracco Diagnostics, Inc. CO <sub>2</sub> ENDOSCOPIC INSUFFLATOR IN LINE TUBING SET	K053008	E-Z-EM, Inc. (now BDI) CO <sub>2</sub> Efficient Endoscopic Insufflator
		K954451	Olympus America, Inc. EVIS 140 System (MH-970 Water Container Accessory)

**Device Description:** The Bracco Diagnostics Inc. (BDI) CO<sub>2</sub> ENDOSCOPIC INSUFFLATOR IN LINE TUBING SET (ILTS) is intended for use with BDI CO<sub>2</sub> Endoscopic Insufflators and OEM endoscope systems, in this specific case the Olympus® 140, 160, 180, and 190 series endoscopes.

The ILTS is supplied sterilized using Gamma radiation to a Sterility Assurance Level (SAL) of 10<sup>-6</sup>.

The ILTS is individually packaged At the time of this submission, accelerated aging tests confirmed a one (1) year shelf life.

**Intended Use:** The CO<sub>2</sub> ENDOSCOPIC INSUFFLATOR IN LINE TUBING SET is intended to connect a CO<sub>2</sub> source (insufflator), and a sterile water source (water bottle), to an endoscope to supply CO<sub>2</sub> during gastrointestinal endoscopic procedures.

**Technology Comparison:** The ILTS employs the same technological characteristics as the predicate devices.

**Characteristic Comparisons:**

	<i>E-Z-EM, Inc. (Bracco Diagnostics, Inc.) CO<sub>2</sub> ENDOSCOPIC INSUFFLATOR Tubing Set (K053008)</i>	<i>Olympus MH-970 Water Container (Accessory within the Olympus EVIS 140 Video Endoscopy System) (K954451)</i>	<i>Bracco Diagnostics, Inc. CO<sub>2</sub> ENDOSCOPIC INSUFFLATOR IN LINE TUBING SET</i>
<b>Intended Use</b>	The CO <sub>2</sub> ENDOSCOPIC INSUFFLATOR IN LINE TUBING SET is intended to be used with a CO <sub>2</sub> source and insufflator, along with a sterile water source, to supply CO <sub>2</sub> and sterile water to an endoscope during gastrointestinal endoscopic procedures.	The Olympus MH-970 Water Container is intended to be used with a CO <sub>2</sub> source and insufflator, along with a sterile water source, to supply CO <sub>2</sub> and sterile water to an endoscope during endoscopic procedures.	The CO <sub>2</sub> ENDOSCOPIC INSUFFLATOR IN LINE TUBING SET is intended to connect a CO <sub>2</sub> source (insufflator), and a sterile water source (water bottle), to an endoscope to supply CO <sub>2</sub> during gastrointestinal endoscopic procedures.
<b>Compatibility</b>	Connects to OEM endoscope systems that are equipped with an external CO <sub>2</sub> luer connector.	Connects to Olympus® 140, series endoscopes.	Same for both predicates
<b>Sterility Method</b>	Non-sterile.	Glutaraldehyde; Autoclaveable (Provided non-sterile)	Radiation (Provided sterile)
<b>Sterility Assurance Level</b>	Not applicable	Unknown	10 <sup>-6</sup>

<i>Disposable or Reusable</i>	Disposable	Reusable	Disposable
<i>Packaging</i>	Poly bag	Unknown	10" x 14" Tyvek®, peelable pouch or similar
<i>Shelf-Life</i>	Six (6) years	Unknown	One (1) year initially

**Performance Testing:****Sterilization**

The ILTS is gamma radiation sterilized and was validated to a sterility assurance level of  $10^{-6}$  in accordance with the following standards:

- *ISO 11137-1: 2006, Sterilization of health care products – Requirements for validation and routine control – Radiation sterilization; and*
- *ISO 11137-2:2006 (currently under revision), Sterilization of health care products -- Radiation -- Part 2: Establishing the sterilization dose.*

Verification results indicate that the ILTS complies with the standards.

**Shelf-Life**

The ILTS is sterilized and its packaging was validated in accordance with the following standards:

- *ISO 11607-1: 2006 Packaging for terminally sterilized medical devices – Part 1: requirements for materials, sterile barrier systems and packaging systems; and*
- *ISO 11607-2: 2006 Packaging for terminally sterilized medical devices – Part 2: Validation requirements for forming, sealing and assembly processes.*

Verification results indicate that the ILTS complies with the standards.

**Biocompatibility**

ILTS patient contact materials were verified in accordance with the following standard:

- *ISO 10993-1: 2009, Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process.*

Verification results indicated that the ILTS patient contact materials comply with the standard.

*Performance Testing* – *Bench* The ILTS was tested for performance in accordance with its predetermined specifications. Testing included:

- Device compatibility;
- Usability;
- Environmental;
- Delivery flow rate;
- Delivery pressure;
- Leakage;
- Mechanical integrity;
- Shipping and transportation; and
- Labeling.

Test results indicate that the ILTS complies with its predetermined specification.

*Conclusion*

Verification and validation activities were conducted to establish the performance and safety characteristics of the ILTS. The results of these activities demonstrate that the ILTS is safe and effective when used in accordance with its intended use and labeling.

Therefore, the ILTS is considered substantially equivalent to the predicate devices.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center -- WO66-G609  
Silver Spring, MD 20993-0002

August 6, 2013

Bracco Diagnostics, Inc.  
% Tracey Alexander  
Director Regulatory Affairs  
532 Broadhollow Road, Suite 126  
Melville, NY 11747

Re: K123047  
Trade/Device Name: Bracco Diagnostics, Inc. CO<sub>2</sub> ENDOSCOPIC INSUFFLATOR IN  
LINE TUBING SET  
Regulation Number: 21 CFR 876.1500  
Regulation Name: Endoscope and accessories  
Regulatory Class: II  
Product Code: FCX, FDS  
Dated: July 12, 2013  
Received: July 26, 2013

Dear Tracey Alexander,

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA).

You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**Herbert P. Lerner -S**

for

Benjamin R. Fisher, Ph.D.  
Director  
Division of Reproductive, Gastro-Renal,  
and Urological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K 123047

---

Device Name: Bracco Diagnostics, Inc. CO<sub>2</sub> ENDOSCOPIC INSUFFLATOR  
IN LINE TUBING SET

Indications for Use: The CO<sub>2</sub> ENDOSCOPIC INSUFFLATOR IN LINE TUBING  
SET is intended to connect a CO<sub>2</sub> source (insufflator), and a  
sterile water source (water bottle), to an endoscope to supply  
CO<sub>2</sub> during gastrointestinal endoscopic procedures.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF  
NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

**Herbert P. Lerner -S**

---

(Division Sign-Off)

Division of Reproductive, Gastro-Renal, and  
Urological Devices

510(k) Number           K123047