

K123131

FEB 08 2013

## 510(k) Summary

510(k) Summary

ST AIA-PACK 25-OH Vitamin D

FEB 08 2013

**Date:** October 3, 2012  
Amended: January 7, 2013

**Submitter:** Tosoh Bioscience, Inc  
3600 Gantz Road  
Grove City, OH 43123

**Contact Person:** Robert L. Wick  
Regulatory Specialist  
6000 Shoreline Ct., Ste. 101  
South San Francisco, CA 94080  
Phone: 650-636-8117  
Fax: 650-636-8121  
Email: [Robert.Wick@Tosoh.com](mailto:Robert.Wick@Tosoh.com)

**Device Name:** ST AIA-PACK 25-OH Vitamin D  
**Classification:** Class II  
MRG  
Clinical Chemistry  
21 CFR 862.1825

**Device Name:** ST AIA-PACK 25-OH Vitamin D Calibrator Set  
**Classification:** Class II  
JIT  
Clinical Chemistry  
21 CFR 862.1150

**Device Name:** AIA-PACK 25-OH Vitamin D Control Set  
**Classification:** Class II  
JJX  
Clinical Chemistry  
21 CFR 862.1660

**Predicate Device:** k983617  
DiaSorin Inc.  
25-Hydroxyvitamin D <sup>125</sup>I RIA Kit

## 510(k) Summary

### ST AIA-PACK 25-OH Vitamin D

According to the requirements of 21 CFR 807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence.

#### Device Description:

The ST AIA-PACK 25-OH Vitamin D is a one-step delayed competitive enzyme immunoassay which, after sample pretreatment, is performed entirely in the ST AIA-PACK 25-OH Vitamin D test cup. Sample pretreatment reagents (containing sodium hydroxide) disassociate 25-OH vitamin D from its binding proteins in the test sample. 25-OH vitamin D present in the pretreated sample is bound to 25-OH vitamin D-specific monoclonal antibody immobilized on magnetic beads. After that, the enzyme-labeled 25-OH vitamin D is added to the reaction mixture. The enzyme-labeled 25-OH vitamin D competes with 25-OH vitamin D for binding to the antibody on magnetic beads in the reaction mixture.

After the second incubation, the magnetic beads are washed to remove the unbound enzyme-labeled 25-OH vitamin D and are then incubated with a fluorogenic substrate, 4-methylumbelliferyl phosphate (4MUP). The amount of enzyme-labeled 25-OH vitamin D that binds to the beads is inversely proportional to the 25-OH vitamin D concentration in the test sample. A standard curve is constructed, and unknown 25-OH vitamin D concentrations are calculated using this curve.

The following products are required to use the ST AIA-PACK 25-OH Vitamin D P/N 025234 :

ST AIA-PACK 25-OH Vitamin D Calibrator Set	025334
ST AIA-PACK 25-OH Vitamin D Pretreatment Set	025734
AIA-PACK 25-OH Vitamin D Control Set	025434

#### Device Intended Use:

ST AIA-PACK 25-OH Vitamin D is designed for in vitro diagnostic use only for the quantitative measurement of total 25-hydroxyvitamin D (25-OH vitamin D) in human serum, Na heparinized or EDTA plasma on TOSOH AIA System Analyzers. The Tosoh ST AIA-PACK 25-OH Vitamin D is intended as an aid in the determination of Vitamin D sufficiency.

#### Calibrators:

The ST AIA-PACK 25-OH Vitamin D Calibrator Set is intended for In Vitro Diagnostic Use Only for the calibration of the ST AIA-PACK 25-OH Vitamin D assay.

#### Controls:

The AIA-PACK 25-OH Vitamin D Control Set is intended for IN VITRO DIAGNOSTIC USE ONLY for performing quality control procedures with the ST AIA-PACK 25-OH Vitamin D assay.

**Substantial Equivalence:**

**Comparison between the Tosoh ST AIA-PACK 25-OH Vitamin D and DiaSorin 25-Hydroxyvitamin D <sup>125</sup>I RIA Kit**

**Similarities**

<b>Parameter</b>	<b>ST AIA-PACK 25-OH Vitamin D</b>	<b>DiaSorin 25-Hydroxyvitamin D <sup>125</sup>I RIA Kit</b>
Intended use	ST AIA-PACK 25-OH Vitamin D is designed for in vitro diagnostic use only for the quantitative measurement of total 25-hydroxyvitamin D (25-OH vitamin D) in human serum, Na heparinized or EDTA plasma on TOSOH AIA System Analyzers. The Tosoh ST AIA-PACK 25-OH Vitamin D is intended as an aid in the determination of Vitamin D sufficiency.	This kit is intended for the quantitative determination of the 25-hydroxy vitamin D (25-OH-D) and other hydroxylated vitamin D metabolites in human serum or plasma to be used in the assessment of vitamin D sufficiency. Assay results should be used in conjunction with other clinical and laboratory data to assist the clinician in making individual patient management decisions in an adult population.
Specimen type	EDTA plasma, Na heparinized plasma or serum	Serum or plasma

**Calibrator Set**

<b>Parameter</b>	<b>ST AIA-PACK 25-OH Vitamin D Calibrator Set</b>	<b>25-OH-D3 Calibrators for Diasorin <sup>125</sup>I RIA Kit</b>
Intended use	The ST AIA-PACK 25-OH Vitamin D Calibrator Set is intended for IN VITRO DIAGNOSTIC USE ONLY for the calibration of the ST AIA-PACK 25-OH Vitamin D assay.	The 25-OH-D3 Calibrators is used to calibrate the Diasorin <sup>125</sup> I RIA Kit.
Matrix	Human serum	Human serum
Number of Calibrators	6	6

**Control Set**

<b>Parameter</b>	<b>AIA-PACK 25-OH Vitamin D Control Set</b>	<b>25-OH-D Controls for Diasorin <sup>125</sup>I RIA Kit</b>
Intended use	The AIA-PACK 25-OH Vitamin D Control Set is intended for IN VITRO DIAGNOSTIC USE ONLY for performing quality control procedures with the ST AIA-PACK 25-OH Vitamin D assay.	The 25-OH-D Controls is used for performing quality control with the Diasorin <sup>125</sup> I RIA Kit .
Matrix	Human serum	Human serum
Number of levels	2	2

**Differences**

Parameter	ST AIA-PACK 25-OH Vitamin D	DiaSorin 25-Hydroxyvitamin D <sup>125</sup> I RIA Kit
Assay Technology	Fluorescence Immunosassay	Radioimmunoassay
Incubation Time	10 minute cycle	110 Minute Cycle
Reference range	10.8 to 54.75 ng/mL	9.0 – 37.6 ng/mL
Limit of detection	2.6 ng/mL	1.5 ng/mL
Assay Range	4.0 ng/mL – 120.0 ng/mL	4.8 ng/mL – 100.0 ng/mL

**Calibrator Set**

Parameter	ST AIA-PACK 25-OH Vitamin D Calibrator Set	25-OH-D3 Calibrators for Diasorin <sup>125</sup> I RIA Kit
Format	Lyophilized	Ready to use liquid
Calibrator Range	0-165 ng/mL	0-100 ng/mL

**Control Set**

Parameter	AIA-PACK 25-OH Vitamin D Control Set	25-OH-D Controls for Diasorin <sup>125</sup> I RIA Kit
Format	Lyophilized	Ready to use liquid
Controls Levels	Two levels at approximately 20 and 80 ng/mL of 25-OH vitamin D	Two levels at the low normal (9.0 ng/mL) and high normal (37.6 ng/mL)

**PERFORMANCE CHARACTERISTICS**

**Precision**

The precision for ST AIA-PACK 25-OH Vitamin D was determined based on guidance from CLSI Protocol EP5-A2.

Within- run precision was assessed by assaying three levels of pooled and spiked Na heparinized plasma, EDTA plasma and serum specimens with each lot number and instrument. Estimates of total and within-run precision were obtained from measurements of 2 replicates in a single run, 2 times a day for 20 non-consecutive days for each specimen/reagent/ instrument combination. This equaled a total of 40 runs and 80 determinants for each data set. There were three (3) data sets at each level for each specimen type.

**Intra-assay (within run) Precision**

<b>Sample</b>	<b>Mean (ng/mL)</b>	<b>Standard Deviation (ng/mL)</b>	<b>Coefficient of Variation (%)</b>
Serum A1	16.0	1.2	7.3
Serum B1	36.0	1.0	2.7
Serum C1	94.1	1.5	1.6
Serum A2	21.8	0.7	3.3
Serum B2	44.4	1.6	3.6
Serum C2	98.0	1.9	2.0
Serum A3	21.8	0.9	4.3
Serum B3	46.5	1.8	3.9
Serum C3	103.7	2.1	2.0
Na Hep Plasma A1	15.7	0.8	4.8
Na Hep Plasma B1	32.9	0.8	2.4
Na Hep Plasma C1	75.4	1.2	1.6
Na Hep Plasma A2	21.5	1.0	4.5
Na Hep Plasma B2	47.0	1.2	2.6
Na Hep Plasma C2	102.8	1.9	1.9
Na Hep Plasma A3	21.8	1.1	5.2
Na Hep Plasma B3	50.0	1.6	3.1
Na Hep Plasma C3	108.9	1.8	1.7
EDTA Plasma A1	13.7	0.7	4.9
EDTA Plasma B1	34.7	0.9	2.7
EDTA Plasma C1	71.1	1.2	1.7
EDTA Plasma A2	24.4	0.8	3.3
EDTA Plasma B2	47.4	0.9	2.0
EDTA Plasma C2	105.2	1.5	1.4
EDTA Plasma A3	25.6	1.1	4.2
EDTA Plasma B3	50.3	1.0	1.9
EDTA Plasma C3	110.9	1.3	1.2

**Total Precision**

Sample	Mean (ng/mL)	Standard Deviation (ng/mL)	Coefficient of Variation (%)
Serum A1	16.0	1.2	7.4
Serum B1	36.0	1.3	3.7
Serum C1	94.1	2.2	2.3
Serum A2	21.8	0.8	3.8
Serum B2	44.4	1.5	3.3
Serum C2	98.0	2.0	2.0
Serum A3	21.8	1.3	5.8
Serum B3	46.5	1.8	3.9
Serum C3	103.7	2.4	2.3
Na Hep Plasma A1	15.7	1.0	6.4
Na Hep Plasma B1	32.9	1.4	4.1
Na Hep Plasma C1	75.4	2.2	2.9
Na Hep Plasma A2	21.5	1.0	4.4
Na Hep Plasma B2	47.0	1.5	3.3
Na Hep Plasma C2	102.8	2.7	2.7
Na Hep Plasma A3	21.8	1.1	5.2
Na Hep Plasma B3	50.0	2.0	4.0
Na Hep Plasma C3	108.9	3.3	3.0
EDTA Plasma A1	13.7	0.9	6.6
EDTA Plasma B1	34.7	1.3	3.9
EDTA Plasma C1	71.1	1.9	2.7
EDTA Plasma A2	24.4	1.0	4.1
EDTA Plasma B2	47.4	1.4	2.9
EDTA Plasma C2	105.2	2.3	2.2
EDTA Plasma A3	25.6	1.2	4.7
EDTA Plasma B3	50.3	1.6	3.2
EDTA Plasma C3	110.9	2.8	2.5

**Controls Precision****Intra-assay (within run) Precision**

	Mean (ng/mL)	Pooled SD (ng/mL)	CV (%)
Level 1	20.0	0.79	3.9
Level 2	79.1	1.63	2.1

Total Precision

	Mean (ng/mL)	Pooled SD (ng/mL)	CV (%)
Level 1	20.0	1.16	5.8
Level 2	79.1	2.22	2.8

**Linearity:**

The linearity for ST AIA-PACK 25-OH Vitamin D was determined, based on guidance from CLSI protocol EP6-A. The linearity was measured on the AIA-2000 instrument and has been demonstrated to be linear from 4 ng/mL to 150 ng/mL.

**Correlation**

The methods comparison study was determined based on guidance from (EP9-A2).

a. Method Comparison

A total of 156 unaltered serum specimens were assayed in singleton utilizing the ST AIA-PACK 25-OH Vitamin D assay on the AIA-2000 analyzer and the 25-Hydroxyvitamin D <sup>125</sup>I RIA Kit (Diasorin). The regression analysis for the correlation between the alternate method (x) and the ST AIA-PACK 25-OH Vitamin D (y) is as follows:

	Deming	Regular
Slope:	0.934 (0.884 to 0.984)	0.885 (0.836 to 0.935)
Intercept:	2.53 (0.75 to 4.31)	3.97 (2.21 to 5.73)
Corr Coef (R):	0.944	
Result Ranges:	Diasorin 4.4 to 103.0 ng/mL Tosoh 6.2 to 103.2 ng/mL	

\*95% Confidence Intervals are shown in parentheses

b. Matrix Comparison

The correlation between serum (x) and Na heparinized plasma (y) on the ST AIA-PACK 25-OH Vitamin D was carried out using 115 unaltered specimens.

	Deming	Regular
Slope:	0.993 (0.974 to 1.012)	0.988 (0.978 to 1.007)
Intercept:	-0.390 (-1.13 to 0.36)	-0.20 (-0.94 to 0.55)
Corr Coef (R):	0.995	
Result Ranges:	Na Heparin 11.3 to 112.8 ng/mL Serum 11.0 – 111.1 ng/mL	

\*95% Confidence Intervals are shown in parentheses

The correlation between serum (x) and EDTA Plasma (y) on the ST AIA-PACK 25-OH Vitamin D was carried out using 115 unaltered specimens.

	Deming	Regular
Slope:	1.041 (1.020 to 1.062)	1.034 (1.013 to 1.056)
Intercept:	0.091 (0.09 to 1.73)	1.14 (0.032 to 1.96)
Corr Coef (R):	0.994	
Result Ranges:	EDTA 11.6 to 111.9 ng/mL Serum 11.0 – 111.1 ng/mL	

\*95% Confidence Intervals are shown in parentheses

### Cross Reactivity

The following cross reactants were tested on the ST AIA-PACK 25-OH Vitamin D assay by spiking known quantities of the cross reactant in serum specimen.

Cross-reactants	Cross-reactivity* [mol%]
25-OH Vitamin D2 (30 ng/mL)	101.1
25-OH Vitamin D3 (158.2 ng/mL)	99.2
3-epi 25-OH Vitamin D2 (30 ng/mL)	131.3
3-epi 25-OH Vitamin D3 (30 ng/mL)	107.7
24,25-(OH) <sub>2</sub> Vitamin D2 (100 ng/mL)	5.2
24,25-(OH) <sub>2</sub> Vitamin D3 (100 ng/mL)	18.1
Vitamin D2 (1000 ng/mL)	0.5
Vitamin D3 (1000 ng/mL)	<0.1
1, 25-(OH) <sub>2</sub> Vitamin D2 (1.0 ng/mL)	<0.1
1, 25-(OH) <sub>2</sub> Vitamin D3 (1.0 ng/mL)	<0.1
Paricalcitol (2.0 ng/mL)	<0.1

Tosoh Bioscience, Inc.

**Conclusion:**

The ST AIA-PACK 25-OH Vitamin D recovers  $\geq 100\%$  for both 25-OH Vitamin D2 and 3-epi- 25-OH Vitamin D3 at 30 ng/mL.

The ST AIA-PACK 25-OH Vitamin D recovers approximately 100 % for 25-OH Vitamin D3 at 158.2 ng/mL.

The ST AIA-PACK 25-OH Vitamin D had a cross reactivity of 5.2% at 100 ng/mL of 24,25-Dihydroxyvitamin-D2.

The ST AIA-PACK 25-OH Vitamin D did not significantly cross react with 24, 25-(OH)<sub>2</sub> Vitamin D3 at 100 ng/mL.

The ST AIA-PACK 25-OH Vitamin D did not significantly cross react with Vitamin D3 and Vitamin D2 at 1000 ng/mL.

The ST AIA-PACK 25-OH Vitamin D did not significantly cross react with 1, 25-(OH)<sub>2</sub> Vitamin D3, 1, 25-(OH)<sub>2</sub> Vitamin D2 and Paricalcitol at 2.0 ng/mL.

**Reference Ranges**

Number of Samples (n) 233  
Reference Interval 10.8 - 54.75 ng/mL

The central 95% of the reference range was used to determine the reference interval.

**Interference**

Interference is defined, for the purposes of this study, with recovery outside of 100 +/- 10% of the known concentration of the specimen after the following substances are added to human specimens. Three studies were conducted using EDTA plasma, Na heparinized plasma and serum.

- Added hemoglobin (up to 11.4 mg/dL) does not interfere with the assay.
- Conjugated bilirubin (up to 18.4 mg/dL) and free bilirubin (up to 16.8 mg/dL) does not interfere with the assay.
- Lipemia, as indicated by triglyceride concentration (up to 357 mg/dL), does not interfere with the assay.

Tosoh Bioscience, Inc.

- Protein, as indicated by human albumin concentration (up to 8.5 g/dL), does not interfere with the assay.
- Ascorbic acid (up to 20 mg/dL) does not interfere with the assay.
- EDTA•2K (up to 10 mg/mL) does not interfere with the assay.
- Sodium Heparin (up to 100 U/mL) does not interfere with the assay.
- Rheumatoid factor (up to 500 IU/mL) does not interfere with the assay.

#### Limit of Detection and Limit of Quantification

The LoB and LoQ of the ST AIA-PACK 25-OH Vitamin D was determined based on CLSI guideline EP17-A.

The LoB was determined by 60 measurements of 3 different blank specimens. The LoB was the value at the 95th percentile. In this case the LoB was determined to be 1.6 ng/mL.

The LoD was determined by 10 measurements of 6 low level samples. The sample range was chosen to be between LoB and 4xLoB and 3 lots of reagents were utilized. The LoD was determined to be 2.6 ng/mL.

The LoQ was determined by 10 measurements of 18 low level samples and the CV% for each sample was calculated. The samples were selected to be in the range of LoB to 4xLoB. A precision profile was plotted for concentration vs. the CV%. The functional sensitivity at 20% CV was calculated and determined as the LoQ. The LoQ was determined to be 2.9 ng/mL.

The reportable range for the assay is 4.0 – 120 ng/mL.

#### Standards:

Number	FDA Recognition Number	Revision Date	Title
EP5-A2	7-110	10/31/2005	Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition
EP6-A	7-193	03/18/2009	Evaluation of the Linearity of Quantitative Measurement Procedures: a Statistical Approach: Approved Guideline
C28-A3	7-202	09/08/2009	Defining, Establishing, and verifying Reference Intervals in the Clinical Laboratory; Approved Guideline – Third Edition
EP17-A	7-194	03/18/2009	Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline
EP9-A2	7-92	03/08/2004	Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Second Edition
EP21-A	7-174	03/18/2009	Estimation of Total Analytical Error for

Tosoh Bioscience, Inc.

			Clinical Laboratory Methods; Approved Guideline
--	--	--	--

Tosoh Bioscience, Inc.

**Conclusion:**

The Tosoh Bioscience, Inc. ST AIA-PACK 25-OH Vitamin D is substantially equivalent to the DiaSorin Inc. (k)983617-Hydroxyvitamin D <sup>125</sup>I RIA Kit for in vitro diagnostic use only for the quantitative measurement of total 25-hydroxyvitamin D (25-OH vitamin D) in human serum, Na heparinized or EDTA plasma on TOSOH AIA System Analyzers. The Tosoh ST AIA-PACK 25-OH Vitamin D is intended as an aid in the determination of Vitamin D sufficiency



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

February 8, 2013

Tosoh Bioscience, Inc  
c/o Robert L. Wick  
6000 Shoreline Court, Suite 101  
South San Francisco, CA 94080

Re: k123131

Trade/Device Name: ST AIA-PACK 25-OH Vitamin D,  
ST AIA-PACK 25-OH Vitamin D Calibrator Set,  
AIA-PACK 25-OH Vitamin D Control Set

Regulation Number: 21 CFR 862.1825

Regulation Name: Vitamin D test system

Regulatory Class: II

Product Code: MRG, JIT, JJX

Dated: January 08, 2013

Received: January 09, 2013

Dear Mr. Wick:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2—Mr. Wick

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

 Carol G. Benson for

Courtney H. Lias, Ph.D.  
Director  
Division of Chemistry and Toxicology Devices  
Office of In Vitro Diagnostics  
and Radiological Health  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): k123131

Device Name:

ST AIA-PACK 25-OH Vitamin D

ST AIA-PACK 25-OH Vitamin D Calibrator Set

AIA-PACK 25-OH Vitamin D Control Set

Indications For Use:

ST AIA-PACK 25-OH Vitamin D is designed for in vitro diagnostic use only for the quantitative measurement of total 25-hydroxyvitamin D (25-OH vitamin D) in human serum, Na-heparinized plasma or EDTA plasma on TOSOH AIA System Analyzers. The Tosoh ST AIA-PACK 25-OH Vitamin D is intended as an aid in the determination of Vitamin D sufficiency.

The ST AIA-PACK 25-OH Vitamin D Calibrator Set is intended for in vitro diagnostic use only for the calibration of the ST AIA-PACK 25-OH Vitamin D assay.

The AIA-PACK 25-OH Vitamin D Control Set is intended for in vitro diagnostic use only for performing quality control procedures with the ST AIA-PACK 25-OH Vitamin D assay.

Prescription Use  X   
(21 CFR Part 801 Subpart D)

And/Or

Over the Counter Use    
(21 CFR Part 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE; CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Diagnostic Device Evaluation and Safety (OIVD)

Yung W. Chan -S

Division Sign-Off  
Office of In Vitro Diagnostic Device  
Evaluation and Safety

510(k)  k123131