

TRANSTEK

JAN 09 2013

Section 5 - 510(k) Summary

Date of Summary Preparation: 11/30/2012

1. Submitter's Identifications

Submitter's Name: ZHONGSHAN TRANSTEK ELECTRONICS CO., LTD.
Address: Jin'an Road, Minzhong, Zhongshan City, Guangdong, China
Contact Person: Lisa Li
Contact Email Address: lisha1@transtek.cn
Telephone: 086-760-88282982 ext. 876
Fax: 086-760-85339231

2. Correspondent's Identifications

Correspondent's Name: A03 Lab of BTS
Address: No.1 Fanghua Street, Hi-tech Zone, Chengdu City, Sichuan, China
Contact Person: Leo Wang
Contact Email Address: leo.w@hibts.com
Telephone: 086-28-86083300
Fax: 086-20-80727399

3. Name of the Device

Device Classification Name: Analyzer, Body Composition (Impedance Plethysmograph)
Product Name: Glass Body Fat Analyzer
Trade/Proprietary Name: TRANSTEK Glass Body Fat Analyzer
Model: LS206-E
Classification Panel: Cardiovascular
Common/Usual Name: Body Composition Analyzer/Scales
Product Code: MNW
Device Classification: Class II
Contraindications: Do not use the Analyzer if you have a pacemaker or other internal medical device.

4. The Predicate Devices

TRANSTEK, Glass Body Fat Analyzer, Model GBF-950-D, K112932

5. Device Description

Transtek Glass Body Fat Analyzer uses BIA (Bioelectrical Impedance Analysis) technology which passes an electrical current through the body to estimate body fat mass, total body water, muscle mass and bone mass. The electrical current is low and may not be felt. The current passes freely through the

TRANSTEK

fluids contained in muscle tissue, but encounters difficulty/resistance when it passes through fat tissue. This resistance of the fat tissue to the current is termed 'Bioelectrical Impedance', and is accurately measured by Glass Body Fat Analyzer LS206-E.

This method simultaneously calculates your personal weight, body fat, total body water, muscle mass and bone mass, giving you a more accurate reading of your overall health and fitness.

This scale stores the personal data of up to 4 users. As well as being an analyzer, this device can be used as a conventional weight scale.

Transtek Glass Body Fat Analyzer LS206-E embeds a Wireless network connections module that allows it to connect to nearby receiving end (such as specific equipment that named Bridge) which is connected to the Internet. Once measurement is over, the LCD of device displays results. And the device will start to send out data. The Bridge receive / storage, and transmission data to Internet server. Thus users can receive, and display/storage, measurement results from LS206-E unit through their end devices (e.g. PC, cellular, tablet) that connected Internet.

6. Intended Use of Device

The TRANSTEK Glass Body Fat Analyzer measure weight and uses bioelectrical impedance analysis (BIA) technology to estimate body fat, total body water percentage, muscle mass, and bone mass in generally healthy adults 18 years of age or older. It is intended for use in the home/domestic setting only.

7. Design Control Activities and Performance Tests Summary

Design control activities for this modification were performed and bench tests have been done. Those performance tests, risk management, and design verification tests provide demonstration that the difference does not raise any new questions of safety and effectiveness.

LS206-E conforms to the following standards:

ISO14971, Risk management to medical devices

IEC60601-1, Electrical safety; IEC60601-1-2, Electromagnetic compatibility

FCC Part 15, EMI tests of FCC Radiation & RF rules and regulations

Explanation: The new wireless function does not affect body analyzer measurement function. Therefore we have not done the Clinical test.

8. Summary of Substantial Equivalence

8.1 Differences between proposed device and the predicate device

The only significant function difference between the two devices is that LS206-E add-on a wireless data communication, what user option, which can transmit measurement results to those end devices which connected Internet.

More modification details are described in this submission.

TRANSTEK

8.2 Discussion

The Transtek Glass Body Fat Analyzer LS206-E has identical indication for use, fundamental scientific technology, energy type, dimensional specifications, and similar performance specifications, software/firmware, functions, labeling to the predicate device.

The only function difference between LS206-E and the predicate device is that the modified device provides user an optional wireless data transmission. It is an add-on function that is entirely independent from the body analyzer function, which does not rely on the wireless connection to carry out a bioelectrical impedance analysis and display its results. Thus the wireless data transmission function does not affect the safety and effectiveness of the body analyzer function.

9. Conclusions

The Transtek Glass Body Fat Analyzer LS206-E is substantially equivalent to the predicate device GBF-950-D by having the identical indication for use, identical technologies, and an add-on function which does not impact the safety and effectiveness of the device.

--- End of this section ---



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

January 9, 2013

ZHONGSHAN TRANSTEK ELECTRONICS CO., LTD.
% Mr. Leo Wang
Senior Consultant
A03 Lab of BTS
No. 1 Fanghua Street, Hi-tech District
CHENGDU SICHUAN
CHINA 610041

Re: K123781
Trade/Device Name: Transtek Glass Body Fat Analyzer
Models: LS206-E
Regulation Number: 21 CFR§ 870.2770
Regulation Name: Impedance plethysmograph
Regulatory Class: II
Product Code: MNW
Dated: December 10, 2012
Received: December 10, 2012

Dear Mr. Wang:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies.

You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Herbert R. Lerner

for
Benjamin R. Fisher, Ph.D.
Director
Division of Reproductive, Gastro-Renal,
and Urological Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

TRANSTEK

Section 4 - Indications for Use

510(k) Number (if known): K123781

Device Name:

Transtek Glass Body Fat Analyzer
Models: LS206-E

Indications for Use:

The Transtek Glass Body Fat Analyzer measure weight and uses bioelectrical impedance analysis (BIA) technology to estimate body fat, total body water percentage, bone mass, and muscle mass in generally healthy adults 18 years of age or older.

It is intended for use in the home/domestic setting only.

Prescription Use _____
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use X
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Herbert R. Lerner

(Division Sign-Off)
Division of Reproductive, Gastro-Renal, and
Urological Devices
510(k) Number K123781