

K130014



**SAKAE CORPORATION**

239 Onishi, Fujioka-shi, Gunma, 370-1401 Japan  
Phone: +81-274-52-3126, Fax: +81-274-52-4240

**510(k) Summary**

**APR 04 2014**

***1 Applicant***

SAKAE CORPORATION  
239 Onishi, Fujioka, Gunma, 370-1401 JAPAN  
Phone: (81)-274-52-3126 Fax: (81)-274-52-4240

***2 Contact Information:***

Erika B. Ammirati, R. A. C., MT (ASCP)  
Regulatory Consultant to SAKAE CORPORATION  
Phone: (650) 949-2768 E-mail: erikaba@pacbell.net

***3 Device Trade Name:***

A1c GEAR System

***4 Device Common Name:***

Glycated hemoglobin assay and discrete photometric chemistry analyzer for clinical use

***5 Manufacturer Address***

SAKAE CORPORATION  
239 Onishi, Fujioka, Gunma, 370-1401 JAPAN  
Phone: (81)-274-52-3126 Fax: (81)-274-52-4240

***6 Device Classification:***

SAKAE CORPORATION A1c GEAR System (new) is a Class II device and reagent, and is classified by FDA under 21 CFR 864.7470 Glycosylated hemoglobin assay and the FDA Product Code is LCP.

A1c GEAR (the instrument) is a Class I device and is classified by FDA under 21 CFR 862.2160, Discrete photometric chemistry analyzer for clinical use, and the FDA Product Code is JJE.

***7 Device Description:***

The A1c GEAR instrument is a fully automated desktop electric spectrophotometer that measures %HbA1c in human whole blood using a dedicated reagent (MEDIDAS HbA1c). The system illuminates a 660 nm LED (Light Emitting Diode) through the test material and quantitatively measures the percent of hemoglobin A1c in the total hemoglobin (%HbA1c) by means of light absorbance changes and a non-linear calibration curve. The system includes the Hemoglobin A1c Analyzer (A1c GEAR), thermal printer, barcode reader, power cable, and fan filter.



MEDIDAS HbA1c is composed of a test cartridge, capillary, pipette tip and master calibration card. The cartridge is pre-filled with reagent; latex (reagent R1), antibody (reagent R2), and sample dilute solution.

**8 Indications for Use**

The A1c GEAR System is intended for in vitro diagnostic use only for the quantitative measurement of the percent hemoglobin A1c (%HbA1c) from finger-stick blood or venous whole blood collected in either EDTA or sodium fluoride (NaF) for clinical laboratory and point of care use. The measurement of HbA1c is recommended to monitor long-term glycemic control of persons previously diagnosed with diabetes mellitus. This test is not for screening or diagnosis of diabetes.

**9 Limitations**

- This test should not be used in monitoring daily glucose control.
- Should not be used to replace daily home testing of urine and blood glucose levels.
- Should not be used for analyzing samples from patients with conditions causing shortened red blood cell survival, such as hemolytic diseases, pregnancy and significant acute or chronic blood loss.

**10 Expected Values and Reference (non-diabetic) Level**

The American Diabetes Association (ADA) expected value range is 4.0-6.0% HbA1c for people without diabetes.

The American Diabetes Association's (ADA) most recent Clinical Practice recommendation for diabetes specified a treatment goal of less than 7% and suggests additional action when HbA1c is above 8%

HbA1c Value	Glycemic Goal
<8% HbA1c	Less stringent
<7% HbA1c	General (Non-Pregnant Adults)
<6.5% HbA1c	More stringent

American Diabetes Association Standards of Medical Care in Diabetes 2012, 35 (Supplement1), S11-S63

**11 Predicate Device**

DCA Vantage, a test system for hemoglobin A1c by Siemens Medical Solutions Diagnostics, K071466. The predicate device has the same intended use, uses substantially the same assay methodology, and is substantially equivalent to the A1c GEAR System.

**Device Comparison Chart**

<b>Manufacturer</b>	<b>SAKAE CORPORATION</b>	<b>Siemens Medical Solutions Diagnostics</b>
<b>Trade Name</b>	<b>A1c GEAR System</b>	<b>DCA Vantage</b>
<b>510(k) Number</b>	K130014	K071466
<b>Product Code</b>	LCP	LCP
<b>Regulation Number</b>	864.7470	864.7470
<b>Indications for use:</b>	Quantitative measurement of percent hemoglobin A1c in human whole blood	Quantitative measurement of percent hemoglobin A1c in human whole blood
<b>Methodology</b>	Immuno-turbidimetric	Immuno-turbidimetric (inhibition)
<b>Sample</b>	Finger-stick blood or venous whole blood collected in K2 EDTA or sodium fluoride	Finger-stick blood or venous whole blood collected in EDTA, heparin, fluoride/oxalate, and citrate
<b>Visual Display</b>	LCD	LCD
<b>Hemolysate preparation</b>	Automatic	Automatic
<b>Detection Method</b>	Transmission	Transmission
<b>Calibration</b>	User; calibration card	User; calibration card
<b>Recommended testing environment</b>	Professional use; point of care	Professional use; point of care
<b>Throughput</b>	6-7 minutes per sample	6-7 minutes per sample
<b>Analytical Range</b>	4.3-12.5%	2.5-14.0%
<b>Reagent Storage</b>	2-8 degrees Celsius (36-46 degrees Fahrenheit)	2-8 degrees Celsius (36-46 degrees Fahrenheit)
<b>Accuracy (Comparison)</b>	Versus HPLC method $Y=1.03x-0.33$ , $R=0.99$ $N=158$	Versus DCCT reference method (HPLC) $Y=1.02x-0.00$ , $R=0.98$ $N=100$
<b>NGSP Certification Status</b>	Certified	Certified
<b>Complies with IEC 60601-1</b>	Yes	Yes
<b>Complies with IEC 60601-1-2</b>	Yes	Yes

**12 Performance Data****12.1 Linearity**

Linearity of the A1c GEAR System was verified with the use of two whole blood samples collected into EDTA tubes. Sample low (L: 4.0% HbA1c, result from ion-exchange HPLC) and high (H: 15.0% HbA1c, result from ion-exchange HPLC) were mixed in different proportions to obtain a series of 11 samples. All samples were measured in triplicate. Recovery rate was used as an indicator for the degree of the deviation of expected values. The linear regression analysis was performed.

**Table 1 Linearity of the A1c GEAR System. Y: observed value, X: expected value,  $r^2$ : squared coefficient of correlation, recovery (%) = observed value / expected value x 100.**

Range (% HbA1c)	Regression line	$r^2$	Recovery (%)
4.0 - 13.1	$y = 0.98x + 0.19$	1.00	98 - 103

**12.2 Method Comparisons (venous to venous sampling, in-house)**

Method comparison studies were performed with three comparison methods; two different ion-exchange HPLC methods and one point of care (POC, DCA Vantage) method were each compared to the A1c GEAR System. Venous whole blood collected into EDTA tubes were prepared from donors and analyzed.

**Table 2 Linear regression analysis data of method comparison. Y: A1c GEAR, X: comparison method, N: number of samples,  $r^2$ : squared coefficient of correlation.**

Comparison Method	N	HbA1c (%)	Regression Line	$r^2$
HPLC 1	158	4.6-10.6	$y = 1.03x - 0.33$	0.98
HPLC 2	40	4.2-9.8	$y = 0.99x + 0.31$	0.98
Another POC analyzer (DCA Vantage)	60	4.7-11.7	$y = 0.95x - 0.12$	0.99

**12.3 Matrix Comparison**

Matrix comparison studies were performed to evaluate the effect of the sample matrix. A finger-stick sample and venous whole blood samples with anticoagulants EDTA or sodium fluoride (NaF), were collected from each donor and analyzed with the A1c GEAR System.

**Table 3 Linear regression analysis data of matrix comparison. N: number of samples,  $r^2$ : squared coefficient of correlation.**

Matrix	N	HbA1c (%)	Regression line	$r^2$
finger (y) vs. EDTA-venous (x)	78	4.3-9.0	$y = 0.96x + 0.15$	0.99
NaF-venous (y) vs. finger (x)	46	4.8-8.8	$y = 1.04x - 0.06$	0.99
NaF-venous (y) vs. EDTA-venous (x)	81	5.3-10.9	$y = 1.01x + 0.01$	0.99

**12.4 Precision**

Precision studies were performed at both internal and external sites. The studies followed CLSI (Clinical and Laboratory Standards Institute) Guideline EP5-A2.



Within-run (repeatability), between-day, and total precision were determined for two control materials (control L and H) and three EDTA whole blood samples at the internal site, and with two control materials (control L and H) and two EDTA whole blood samples at the external site. The samples were analyzed for 20 days, in duplicate, twice a day (n = 80).

**Table 4 Results from the internal site. %CV: %coefficient of variation.**

Sample	N=	Mean	Within-run CV (%)	Between-day CV (%)	Total CV (%)
Control L	80	5.2	1.26	0.51	1.36
Control H	80	9.0	0.85	0.26	1.06
Sample 1	80	5.5	0.73	0.80	1.12
Sample 2	80	11.1	1.11	0.70	1.37
Sample 3	80	12.1	1.14	1.01	1.52

**Table 5 Results from the external site. %CV: %coefficient of variation.**

Sample	N=	Mean	Within-run CV (%)	Between-day CV (%)	Total CV (%)
Control L	80	5.0	1.08	0.75	1.31
Control H	80	8.9	0.65	0.53	0.90
Sample 1	80	5.2	1.18	0.56	1.34
Sample 2	80	8.8	0.82	0.47	1.05

**Table 6 Reproducibility estimated from the results of two sites. %CV: %coefficient of variation.**

Sample	N	Overall mean	Between-site CV (%)	Total CV (%)
Control L	160	5.1	2.55	2.81
Control H	160	9.0	0.56	0.94

### ***12.5 Point of Care (POC) Studies***

External validation of the A1c GEAR System was performed at POC sites to evaluate precision and method comparisons.

In the precision study, three levels of controls were analyzed for 20 days and three levels of patient samples were analyzed for 10 days by POC operators.

**Table 7 Results of precision study at three external sites. %CV: %coefficient of variation.**

Sample	N =	Site	Mean	Within-site CV (%)	Overall mean	Reproducibility Total CV (%)
Control 1	120	1	5.19	2.85%	5.22	2.26%
	120	2	5.21	1.92%		
	120	3	5.27	1.46%		
Control 2	120	1	7.01	2.53%	7.06	2.11%
	120	2	7.07	1.73%		
	120	3	7.10	1.73%		
Control 3	120	1	11.05	3.37%	11.04	2.55%
	120	2	11.09	2.48%		
	120	3	11.00	1.35%		
Sample Low	60	1	5.80	3.14%	5.84	3.12%
	60	2	5.83	3.28%		
	60	3	5.89	2.63%		
Sample Middle	120	1	8.01	3.31%	8.07	4.16%
	128	2	7.87	2.30%		
	120	3	8.34	2.91%		
Sample High	120	1	10.55	3.22%	10.84	5.25%
	128	2	10.59	2.46%		
	120	3	11.38	3.21%		

In the method comparison study, a finger-stick sample and a venous EDTA sample were collected from each donor. The finger-stick blood samples were analyzed with the A1c GEAR System by POC operators and the venous blood samples were analyzed with an ion-exchange HPLC (Tosoh, G8) reference method by qualified laboratory technicians at a reference laboratory.

**Table 8 Linear regression analysis of method comparison study at three external sites. N: number of samples, r: coefficient of correlation.**

Study site	N	Min	Max	Slope (95% confidence interval)	Intercept (95% confidence interval)	r
1	47	4.9	11.9	0.968 (0.941 to 0.994)	0.04 (-0.16 to 0.24)	0.995
2	41	5.4	10.8	0.976 (0.936 to 1.015)	0.12 (-0.15 to 0.40)	0.990
3	46	5.0	9.6	0.989 (0.952 to 1.027)	0.08 (-0.18 to 0.35)	0.990



### **12.6 Interference**

No significant interference was observed up to the following concentrations in both EDTA and NaF whole blood samples, or commercial controls:

- Free - bilirubin 37.0 mg/dl
- Conjugated - bilirubin 40.4 mg/dl
- Rheumatoid factor 550 IU/ml
- Chyle (mixture of lipids) 3120 FTU (Formazine Turbidity Unit)  
includes:
  - Triglycerides 170 mg/dl
  - Phospholipids 182 mg/dl
  - Free fatty acids 124 µEq/dl (approx. 1.24 mmol/l)
- Triglycerides (separate study) 2,000 mg/dl
- Acetaminophen 20.0 mg/dl
- Ibuprofen 50.0 mg/dl
- Glibenclamide 0.2 mg/dl
- Metformin 5.1 mg/dl
- Ascorbic acid 6.0 mg/dl

**NOTE:** It is possible that other substances and/or factors not listed above may interfere with the test and cause false results.

### **12.7 Analytical Specificity**

#### **12.7.1 Hemoglobin (Hb) Variants**

A hemoglobin variant study was performed using commercial samples known to contain Hemoglobin variants C, D, E, S and F. Samples contained both low and high levels of %HbA1c at concentrations from 4.6-11.6%. These variant samples were tested in duplicate using the A1c GEAR System versus a reference method (Primus Ultra Boronate Affinity HPLC). The results indicated samples containing Hemoglobin C were elevated by 24%, samples containing Hemoglobin D were elevated by 16%, samples containing Hemoglobin E were elevated by 13% and samples containing Hemoglobin S were elevated by 14%. Samples containing >10% Hemoglobin F were decreased by 32%. All variants tested were shown to interfere with this device.

#### **12.7.2 Modified Hemoglobin**

The following modified hemoglobin was prepared by incubating with the substance in parentheses and found not to affect the A1c GEAR System for both EDTA and NaF whole blood samples:

- Carbamylated hemoglobin (sodium cyanate, 10 mg/dl)
- Acetylated hemoglobin (acetylsalicylic acid, 200 mg/dl)
- Labile hemoglobin (D-glucose, 2000 mg/dl)



### ***12.8 Limit of Detection***

To estimate the lowest detectable value of %HbA1c for the A1c GEAR System, limit of detection (LOD) studies were performed and LOD was calculated to be 2.6% and LOB was calculated to be 2.3%.

### ***12.9 Stability- Real-Time***

A real time shelf life stability study was performed for MEDIDAS HbA1c using the A1c GEAR analyzer. From these results, it was concluded that the reagent cartridge can be stored for up to one year at 2-8 °C (36-46 °F).

### ***13 Conclusions***

Performance studies were conducted and the data obtained indicate the A1c GEAR System is substantially equivalent to the predicate device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

April 4, 2014

SAKAE CORPORATION  
C/O ERICA AMMIRATI  
575 SHIRLYNN COURT  
LOS ALTOS CA 94022

Re: K130014

Trade/Device Name: A1c Gear System  
Regulation Number: 21 CFR 864.7470  
Regulation Name: Glycosylated hemoglobin assay  
Regulatory Class: II  
Product Code: LCP, JJE  
Dated: March 28, 2014  
Received: March 31, 2014

Dear Ms. Ammirati:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulations (21 CFR Parts 801 and 809), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**Courtney H. Lias -S**

Courtney H. Lias, Ph.D.  
Director  
Division of Chemistry and Toxicology Devices  
Office of In Vitro Diagnostics  
and Radiological Health  
Center for Devices and Radiological Health

Enclosure

**Indications for Use**

510(k) Number (if known)  
K130014

Device Name  
A1c GEAR System

Indications for Use (Describe)

The A1c GEAR System is intended for in vitro diagnostic use only for the quantitative measurement of the percent hemoglobin A1c (%HbA1c) from finger-stick blood or venous whole blood collected in either EDTA or sodium fluoride (NaF) for clinical laboratory and point of care use. The measurement of HbA1c is recommended to monitor long-term glycemic control of persons previously diagnosed with diabetes mellitus. This test is not for screening or diagnosis of diabetes.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.

FOR FDA USE ONLY

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

**Katherine Serrano -S**