



Israel Biomedical Innovations Ltd.  
2 Ha-Eshel St. P.O.BOX 3081  
Caesarea Industrial Park 38900, Israel

JUL 29 2013

**510(K) SUMMARY**

[as required by section 807.92(c)]

**EndoFast Reliant™ SCP and EndoFast Reliant™ LAP**

510(k) Number K\_\_\_\_\_

**Applicant's Name:**

Israel Biomedical Innovations Ltd.  
2 Ha-Eshel St. P.O.BOX 3081  
Caesarea Industrial Park 38900, Israel  
Tel.: +972-4-6277166, Ext. 202  
Fax: +972-4-6277266

**Contact Person:**

Elissa Burg  
VP Quality Assurance and Regulatory Affairs  
Telephone: +972 4 6277166 ext. 107  
Mobile: +972-52 6633572  
Fax: +972 4 6277266  
E-mail: elissa@ibimedical.com

**Prepared On:**

January 7, 2013

**Trade Name:**

EndoFast Reliant™ SCP  
EndoFast Reliant™ LAP

*\*SCP = Sacrocolpopexy*

*\*LAP = Laparoscopic*

**Classification:**

**Name:** Fixation, non-absorbable or absorbable, for pelvic use

**Product Code:** PBQ

**Regulation No:** 884.4530

**Class:** 2

**Classification Panel:** Obstetrics/Gynecology

**Name:** Mesh, Surgical, Polymeric



Israel Biomedical Innovations Ltd.  
2 Ha-Eshel St. P.O.BOX 3081  
Caesarea Industrial Park 38900, Israel

**Name:** Mesh, surgical, synthetic, urogynecologic, for apical vaginal and uterine prolapse, transabdominally placed.

**Product Code:** OTO

**Regulation No:** 878.3300

**Class:** 2

**Classification Panel:** Obstetrics/Gynecology

**Predicate Devices:**

EndoFast Reliant™ System–K080836

Autosuture™ TACKER - K090470

**Device Description:**

The EndoFast Reliant™ SCP (Sacrocolpopexy) is a sterile, single use system which consists of the following parts:

- 3 metal Fixation Devices; each preloaded with a spider fastener composed of metal and includes the MRI Safety Information that is stated in the labeling.
- Plastic Handle; connects to the Fixation Device.
- Surgical mesh composed of monofilament polypropylene material; cut in a rectangular shape of 8x20cm.

The device is introduced through a 5mm trocar in a laparoscopic approach. The mesh is attached to the vaginal vault and the vagina is pulled towards the posterior pelvic bony part. The other end of the mesh is attached to the promontorium, using two fixation points, to restore the normal anatomical position of the uterus and cervix.

The EndoFast Reliant™ LAP (Laparoscopic) is a sterile, single use system which consists of the following parts:

- 3 metal Fixation Devices; each preloaded with a spider fastener composed of metal and include the MRI Safety Information that is stated in the labeling.
- Plastic Handle; connects to the Fixation Device.

The fasteners are introduced through a 5mm trocar in laparoscopic procedures.

**Intended Use:**

The EndoFast Reliant™ SCP is indicated for fixation of surgical mesh to tissue for tissue reinforcement during a laparoscopic Sacrocolpopexy approach.

The EndoFast Reliant™ LAP is indicated for fixation of surgical prosthetic material to tissues for tissue reinforcement during a laparoscopic Sacrocolpopexy approach.

**Technological Characteristics:**

The technological characteristics of the EndoFast Reliant™ SCP and EndoFast Reliant™ LAP are similar to those of the EndoFast Reliant™ System cleared under K080836. The dimensions of the spider fasteners and fixation device in the EndoFast



Israel Biomedical Innovations Ltd.  
2 Ha-Eshel St. P.O.BOX 3081  
Caesarea Industrial Park 38900, Israel

Reliant™ SCP and EndoFast Reliant™ LAP have been decreased to accommodate a laparoscopic approach to be used in Sacrocolpopexy procedures. In addition, the shape of the fixation device in the EndoFast Reliant™ SCP and EndoFast Reliant™ LAP has been changed to accommodate a laparoscopic approach to be used in Sacrocolpopexy procedures.

#### **Performance Testing:**

The appropriate tests to determine substantial equivalence were completed. These include testing in accordance with Guidance for the Preparation of a Premarket Notification Application for a Surgical Mesh (March 22, 1999) and tests designed by the company to address substantial equivalence.

#### **1. The following non-clinical tests have been performed: Performance Tests:**

- a) Functional Characteristics Test to demonstrate appropriate functionality characteristics of EndoFast Reliant SCP after sterilization, Vibration and Drop test and a shelf life study.
  - b) EndoFast Reliant SCP Fastener Support Endurance to Cyclic Loads to validate the endurance of the SCP Spider Fastener in supporting the cyclic load on the surgical mesh.
  - c) SCP Fastener Deployment Reliability and Pull-Out Force Test to verify the deployment reliability and pullout force required for fastener failure, and to verify that the surgical mesh and the artificial tissue have no influence/ effect on the deployment of the Fastener of the SCP Fixation Device
  - d) SCP Fixation Device and Handle Reliability and Force Test to verify the reliability of the SCP Fixation Device and the Handle that is attached to it in deploying the Fasteners without fault and to test the force applied during the deployment.
  - e) Two cadaver studies were performed utilizing the EndoFast Reliant™ SCP.
  - f) Comparative pull-out tests were conducted in order to assess the new products' equivalency with predicate devices.
2. **Gamma Sterilization Validation** for the Surgical Instruments - Using a Radiation dose of 25 kGy in order to assure products' SAL  $10^{-6}$ .
  3. **EtO Sterilization Validation** for the Surgical Mesh - Using Overkill Approach in order to assure products' SAL  $10^{-6}$ .
  4. **Shelf Life Tests** – Using accelerating aging in order to assure the products' adequate performances during their shelf life.
  5. **Biocompatibility** – The products' were evaluated with accordance to ISO 10993-1 and FDA Blue Book Memorandum, G95-1.
  6. **MR compatibility** – In order to assure the safe use of the device at MRI environment.



Israel Biomedical Innovations Ltd.  
2 Ha-Eshel St. P.O.BOX 3081  
Caesarea Industrial Park 38900, Israel

**Conclusion:**

Israel Biomedical Innovations Ltd. believes that, based on the information provided in this submission, the EndoFast Reliant™ SCP and the EndoFast Reliant™ LAP are substantially equivalent to their predicate devices without raising any new safety and/or effectiveness issue.



July 29, 2013

IBI Israel Biomedical Innovations Ltd.  
% Elissa Burg  
VP Quality Assurance and Regulatory Affairs  
2 Ha-Eshel Street, P.O. Box 3081  
CAESAREA INDUSTRIAL PARK 38900  
ISRAEL

Re: K130059  
Trade/Device Name: EndoFast Reliant™ SCP, Endo Fast Reliant™ LAP  
Regulation Number: 21 CFR 884.4530  
Regulation Name: Obstetric-gynecologic specialized manual instrument  
Regulatory Class: Class II  
Product Code: PBQ, OTO  
Dated: June 10, 2013  
Received: June 17, 2013

Dear Elissa Burg,

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must

comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**Benjamin R. Fisher -S**

Benjamin R. Fisher, Ph.D.  
Director  
Division of Reproductive, Gastro-Renal,  
and Urological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure



Israel Biomedical Innovations Ltd.  
2 Ha-Eshel St. P.O.BOX 3081  
Caesarea Industrial Park 38900, Israel

**INDICATIONS FOR USE STATEMENT**

510(k) Number (if known):     K130059    

Device Name:       EndoFast Reliant™ SCP  
\* SCP = Sacrocolpopexy

**Indications for Use:**

The EndoFast Reliant™ SCP is indicated for fixation of surgical mesh to tissue for tissue reinforcement during a laparoscopic Sacrocolpopexy approach.

Prescription Use   X    
(Per 21 CFR 801 Subpart D)

AND/OR

Over-the-Counter Use         
(21 CFR 801 Subpart C)

---

(PLEASE DO NOT WRITE BELOW THIS LINE -CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

**Benjamin R. Fisher -S**

2013.07.29 16:29:29 -04'00'

---

(Division Sign-Off)  
Division of Reproductive, Gastro-Renal, and  
Urological Devices       K130059  
510(k) Number



Israel Biomedical Innovations Ltd.  
2 Ha-Eshel St. P.O. BOX 3081  
Caesarea Industrial Park 38900, Israel

### INDICATIONS FOR USE STATEMENT

510(k) Number (if known): K130059

Device Name: EndoFast Reliant™ LAP

\* LAP = Laparoscopic

#### Indications for Use:

The EndoFast Reliant™ LAP is indicated for fixation of surgical prosthetic material to tissues for tissue reinforcement during a laparoscopic Sacrocolpopexy approach.

Prescription Use X  
(Per 21 CFR 801 Subpart D)

AND/OR

Over-the-Counter Use \_\_\_\_\_  
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE -CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

**Benjamin R. Fisher -S**

2013.07.29 16:30:20 -04'00'

(Division Sign-Off)

Division of Reproductive, Gastro-Renal, and  
Urological Devices

510(k) Number K130059