

K130136

510(k) Summary

Manufacturer: Technomed Europe
Amerikalaan 71
6199 AE Maastricht Airport
The Netherlands

JUL 19 2013

Submitted by: Technomed Europe
Amerikalaan 71
6199 AE Maastricht Airport
The Netherlands
Tel.: (+31) 43-408 6868
Fax: (+31) 43-408 6888

Contact person: Mr. Pierre Vreuls
Manager Regulatory Affairs and Quality Assurance
E-mail: pvreuls@technomed.nl

Date: June 6, 2013

Proprietary Name: Disposable Monopolar and Subdermal Needle Electrodes

Common/usual Name: Needle Electrode

Classification Name: Needle electrodes are classified as class II per 21 CFR section 882.1350 product code GXZ and section 890.1385 product code IKT (diagnostic electromyograph needle electrode).

Legally Marketed Predicate Devices: K990015: Technomed Europe Needle Electrodes
K050325: Technomed Europe Disposable Concentric Probe, Disposable Bipolar Probe, Disposable Monopolar Probe
K062437: Technomed Europe Disposable Hypodermic EMG Needle Electrode
K112034: Carefusion 209, Inc. TECA™ ELITE Disposable Concentric Needles
K072276: Xian Friendship Electronics Co., Ltd., Subdermal Needle Electrodes

Device description: Needle Electrodes are passive devices used for recording, monitoring and stimulation during EEG, EMG, nerve conducting studies and IONM. The Needle Electrodes have

different tip shapes for performing different recording and stimulating procedures.

The Needle Electrodes are disposable and designed for single use and are labeled accordingly. The disposable Needle Electrodes are delivered sterile and can be used after opening of the sterile package.

The needles are comprised of a stainless steel or platinum alloy needle, electrically connected to a lead wire, or with a connector for connection to a separate lead wire. Monopolar needles are coated with a polytetrafluoroethylene coating. In case of a pre-connected lead wire the other end is a regular industry standard DIN 42802 safety connector.

The needles are invasive as they are used subcutaneously. Positioning and use is under supervision of a licensed physician.

Intended Use: Needle Electrodes for Neurological purposes are intended for use with recording, monitoring equipment for the recording of biopotential signals including electroencephalograph (EEG), electromyography (EMG) and nerve potential signals, and are intended for stimulation/recording with stimulation/recording equipment for electromyograph (EMG) and nerve potential signals.

Comparison to predicates:

Characteristic	Predicate device disposable needle electrode EEG / EMG [K990015]	Predicate device Subdermal Needle Electrodes[K072276]	Predicate device disposable concentric needles [K112034]	Disposable monopolar and subdermal needle electrodes[K130136]
Name	Technomed EEG / EMG needle electrodes	Xian Friendship Subdermal Needle Electrodes	Carefusion 209, Inc. TECA™ ELITE Disposable Concentric Needles	Disposable monopolar needle electrode. Disposable subdermal needle electrode
Device class	Class II	Class II	Class II	Same
Product code	IKT and GXZ	GXZ	IKT	IKT and GXZ
Device type	Disposable needle electrode	Disposable subdermal needle electrode	Disposable concentric needle	Disposable monopolar needle electrode. Disposable subdermal needle electrode

Indications for use	Technomed Europe diagnostic needle electrodes are intended to be inserted in the subdermal muscle or nerve tissue only to sense bio-electric. EMG or EEG. signals distally, and will proximal be connected to: Electromyography/electroencephalogram recording equipment	Subdermal Needle Electrodes are intended for use with recording, monitoring and stimulation/recording equipment for the recording of biopotential signals, including electroencephalograph (EEG), electromyograph (EMG) and nerve potential signals and for stimulation during the intraoperative diagnosis of acute dysfunction in corticospinal axonal conduction. The electrodes are sterile and for single patient use only.	Disposable Concentric Needles are intended for use with recording, monitoring and stimulation/recording equipment for the stimulation/recording of biopotential signals Including electromyography (EMG) and nerve potential signals.	Disposable Monopolar and Subdermal Needles are intended for use with recording, monitoring and stimulation/recording equipment for the stimulation/recording of biopotential signals Including electromyography (EMG) and nerve potential signals.
Target population	All patients	All patients	All patients	Same
Anatomical sites	Peripheral nerves and muscles	Peripheral nerves and muscles	Peripheral nerves and muscles	Same
Where used	Hospital	Hospital	Hospital	Same
Device design	Diameters: 0.30mm to 0.60mm Lengths: 25 to 75mm (monopolar), 13 mm (straight subdermal)	Diameter: 0.4mm and 0.6mm Length: 13mm (subdermal) and 23mm (corkscrew)	Diameters: 0.30mm to 0.60mm Lengths: 25mm to 75mm	Diameters: 0.30mm to 0.60mm Lengths: 25 to 75mm (monopolar), 7 to 20 mm (subdermal), 23 mm (corkscrew)
• Tip geometry	Tip geometry: front bevel and pencil tip	Tip geometry: front bevel	Tip geometry: back bevel	Tip geometry: front bevel and pencil tip
Recording / stimulation area	0.42mm ² - 44.5mm ²	16.3mm ² - 44.5 mm ²	0.025mm ² - 0.068mm ²	Same as K990015
Electrode materials	Stainless steel	Stainless steel	Stainless steel, Pt/Ir or W	Stainless steel, Pt/Ir
Coatings	PTFE	n/a	Polyesterimide and low friction lubricant	Same as K990015
Cables	PVC insulated tin plated copper lead wire	PVC insulated tin plated copper lead wire	Detachable reusable cable	Same as K990015
Connectors	DIN 42802 1.5mm touch proof	DIN 42802 Touch proof connector	5 pole DIN	Same as K990015 and K072276
Electrical insulation	Electrical insulation on all surfaces not intended to provide electrical contact with the patient and connection	Electrical insulation on all surfaces not intended to provide electrical contact with the patient and connection	Electrical insulation on all surfaces not intended to provide electrical contact with the patient and connection	Same
Impedance	<200kΩ	Unknown	<300kΩ	Same as K990015

Sterilization method	EO ethylene oxide	EO ethylene oxide	EO ethylene oxide	EO ethylene oxide
Sterility assurance level (SAL)	10 ⁻⁶	Unknown	10 ⁻⁶	Same as K990015
Manufacturers	Technomed Europe	Xian Friendship	CareFusion	Technomed Europe

Biocompatibility: The stainless steel monopolar needle with the polytetrafluoroethylene coating is identical to the Technomed Europe disposable hypodermic needles (K062437) in formulation, processing, sterilization, and no other chemicals have been added (e.g., plasticizers, fillers, color additives, cleaning agents, mold release agents). No further biocompatibility testing with the final finished device was determined to be necessary.

The stainless steel or the platinum/iridium material used in the subdermal needle is identical to that of the Technomed Europe needle electrodes (K990015). No further biocompatibility testing with the final finished device was determined to be necessary.

Performance testing: Electrical safety and compatibility was verified, including impedance testing. Sterilization validation was confirmed to comply equally as for Technomed Europe devices already released to be legally marketed under K050325, and to conform to all requirements.

Conclusion: The comparison to the predicate devices demonstrates that the Needle Electrodes are safe and effective for its intended use and are substantially equivalent to the predicate devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

July 19, 2013

Technomed Europe
c/o Mr. Pierre Vreuls
Amerikalaan 71
6199 AE Maastricht Airport
The Netherlands

Re: K130136

Trade/Device Name: Disposable Monopolar Needle Electrode and Disposable Subdermal
Needle Electrode

Regulation Number: 21 CFR 882.1350

Regulation Name: Needle Electrode

Regulatory Class: Class II

Product Code: GXZ and IKT

Dated: June 6, 2013

Received: June 12, 2013

Dear Mr. Vreuls

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

Page 2 - Mr. Pierre Vreuls

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Victor Krauthamer -S

Victor Krauthamer, Ph.D.
Acting Director
Division of Neurological and Physical
Medicine Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K130136

Device Name: Disposable Monopolar Needle Electrode and Disposable Subdermal Needle Electrode

Indications For Use:

Needle Electrodes for Neurological purposes are intended for use with recording, monitoring equipment for the recording of biopotential signals including electroencephalograph (EEG), electromyograph (EMG) and nerve potential signals, and are intended for stimulation/recording with stimulation/recording equipment for electromyograph (EMG) and nerve potential signals.

Prescription Use x
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use
(21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

<p>Victor Krauthamer -S 2013.07.18 18:26:51 -04'00'</p> <hr/> <p>(Division Sign Off) Division of Neurological and Physical Medicine Devices (DNPMD)</p> <p>510(k) Number <u> K130136 </u></p>
