

**K130536**  
**Traditional 510(k) Summary**

**JUL 16 2013**

**Submitter:** Medtronic Vascular  
37A Cherry Hill Drive  
Danvers, MA 01923-5186

**Contact Person:** Anu Gaur  
Senior Regulatory Affairs Specialist  
Phone: (978) 739-3080  
Fax: (978) 750-8204

**Alternate Contact** Fred Boucher  
Director of Regulatory Affairs  
Phone: (978) 739-3116  
Fax: (978) 750-8204

**Date Prepared:** May 20, 2013

**Trade Name:** Export Advance™ Aspiration Catheter

**Common Name:** Percutaneous Catheter

**Classification Name:** Embolectomy Catheter  
Class II per 21 CFR 870.5150, Product Code DXE

**Predicate Devices:** K120808 - Medtronic Vascular  
Export® AP Catheter.

**Device Description:** The Export Advance™ Aspiration Catheter is a dual lumen catheter used for the aspiration of thrombus and/or debris from a vascular site. The Export Advance™ may also be used for the infusion of diagnostic or therapeutic agents to a desired vascular site.

**Statement of Intended Use:** The Export Advance™ Aspiration Catheter is indicated for:

- Removal/aspiration of embolic material (thrombus/debris) from vessels of the arterial system, and
- To sub selectively infuse/deliver diagnostics or therapeutics agents with or without vessel occlusion.

**Comparison to the predicate devices**

The Export Advance™ Aspiration Catheter represents a series of incremental performance improvements over its predicate device Export® AP Catheter, with primary attributes including, improved deliverability, improved kink resistance, and improved aspiration rate; and including a design feature of a removable preloaded stylet.

**Summary of Technological Characteristics:**

The Export Advance™ includes the following features:

- i. Removable Stylet
- ii. Luer Hub
- iii. Strain Relief
- iv. Inner Liner
- v. Braid Wire
- vi. Proximal Shaft
- vii. Distal Shaft (Dual Lumen/Oversleeve)
- viii. Soft Tip
- ix. Microlumen (wire lumen)
- x. Distal Tip
- xi. Marker band
- xii. Hydrophilic lubricous coating

**Summary of Non-clinical Data:**

Design verification (bench) testing qualification and biocompatibility testing was conducted in accordance with the recommendations presented from the relevant FDA guidance to demonstrate that the subject device Export Advance™ Aspiration Catheter has met the acceptance criteria and performance similar to the predicate device.

**Design Verification Testing:** The design verification (bench) testing was performed based upon the subject device performance specifications. The tests performed for bench testing included:

1. Profile Dimensions (Major & Minor Profile)
2. Guide Wire Lumen ID
3. Working Length
4. Proximal Shaft Tensile
5. Microlumen Tear
6. Tip Tensile
7. Marker Band Tensile
8. Hub Tensile
9. Stylet Hub Tensile
10. Vacuum Integrity

11. Pressure Integrity
12. Air Aspiration
13. Proximal Shaft Crush
14. Proximal Shaft Buckle
15. Evacuation Flow Rate
16. Particle Retrieval
17. 2D Track and Lesion Cross
18. Lubricity & Durability
19. Proximal Shaft Stiffness – Room Temp
20. Proximal Shaft Stiffness – Body Temp
21. Torque Strength
22. Distal Kink
23. Particulate Generation

**Pre-clinical Study (Non-GLP):**

Medtronic Vascular conducted pre-clinical *in vivo* (non-GLP) studies for design evaluation on Export Advance. These *in vivo* studies results provided confirmatory evidence that the Export Advance design is suitable to meet the incremental performance attributes as compared with its predicate device, Export AP Catheter, and related comparative evaluations to support the substantial equivalence.

**Biocompatibility Testing (GLP):** Pursuant to the ISO 10993-1:2009/AC: 2010- *Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process*; and 21 CFR 58 -*Good Laboratory Practice for Nonclinical Laboratory Studies*; Medtronic Vascular has concluded full biocompatibility testing on the subject device.

1. Cytotoxicity
2. *In vitro* Hemolysis
3. Systemic Toxicity
4. Sensitization
5. C3a Complement Activation
6. Sc5b9 Complement Activation
7. ISO Intracutaneous Reactivity
8. USP Material Mediated Pyrogen Study
9. *In vivo* Thromboresistance

No new safety or effectiveness issues were raised during the testing. The bench testing qualification and biocompatibility testing demonstrated that the subject device Export Advance™ Aspiration Catheter is safe, effective, and substantially equivalent to the predicate device.

**Summary of Clinical  
Data:**

No clinical investigation has been performed on the subject device Export Advance™ catheter.

**Conclusion from Data:**

Medtronic Vascular has demonstrated that the subject device Export Advance™ Aspiration Catheter is substantially equivalent to the predicate device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

July 16, 2013

Medtronic Inc.  
C/O Colleen Mullins  
37a Cherry Hill Drive  
Danvers, MA 01923 US

Re: K130536  
Trade/Device Name: Export Advance™ Aspiration Catheter  
Regulation Number: 21 CFR 870.5150  
Regulation Name: Embolectomy Catheter  
Regulatory Class: Class II  
Product Code: DXE  
Dated: June 5, 2013  
Received: June 6, 2013

Dear Ms. Mullins:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act

or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



for  
Bram D. Zuckerman, M.D.  
Director  
Division of Cardiovascular Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number:     K130536    

Device Name: Medtronic Export Advance™ Aspiration Catheter

### Indications for Use:

The Export Advance™ Aspiration Catheter is indicated for:

- Removal/aspiration of embolic material (thrombus/debris) from vessels of the arterial system, and
- To sub selectively infuse/deliver diagnostics or therapeutics agents with or without vessel occlusion.

Prescription Use   X    
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE  
IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)



Page 1 of \_\_