

R131169

DEC 24 2013



5. 510(k) SUMMARY

510(k) Summary

This 510(k) Summary is submitted in accordance with the requirements of 21 CFR Part 807, Section 92.

5.1 Applicant:

Meril Life Sciences Private Limited
Bilakhia house, Survey No. 135/139
Muktanand Marg,
Chala, Vapi
Gujarat
396 191
INDIA

5.2 Contact Person:

Utpal Thakor
Director
Meril Life Sciences Private Limited
Bilakhia house, Survey No. 135/139
Muktanand Marg,
Chala, Vapi
Gujarat
396 191
INDIA
Phone: +91 - 2603052100
Mobile: +91 - 9601260306
Fax: +91 260-3052125
E mail: utpal.thakor@merillife.com

5.3 Date prepared: April 11, 2013



5.4 Device information:

Proprietary Name:	Mozec™ – Rx PTCA Balloon Dilatation Catheter
Common / Usual Name:	Rapid Exchange PTCA Balloon Dilatation Catheter
Regulation name:	Percutaneous Transluminal Coronary Angioplasty (PTCA) Catheter – 21CFR 870.5100 (a)
Product Code:	LOX
Device Class:	Class II

5.5 Predicate device:

1. **Sprinter Legend – Rx PTCA Balloon Dilatation Catheter (P790017)**
2. **Apex Monorail Balloon Catheter (P860019)**
3. **Tamarin blue® PTCA Rx Dilatation Catheter (K112735)**

5.6 Device description:

Mozec™ - Rx PTCA Balloon Dilatation Catheter is a sterile, single use rapid exchange catheter consisting of a semi-compliant balloon, a soft tip, a dual lumen distal shaft and a single lumen proximal shaft. Radiopaque platinum/iridium marker(s) contained in the balloon segment facilitate fluoroscopic visualization of the proximal and distal ends of the Mozec™ balloon's working length (2.00 mm to 4.50 mm sizes) or the center of the balloon's working length (1.50 mm size). The Mozec™ balloons are offered in diameters ranging from 1.50mm to 4.50mm and lengths varying from 9mm to 41mm.

5.7 Indication for use:

The Mozec™ - Rx PTCA Balloon Dilatation Catheter is indicated for balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis for the purpose of improving myocardial perfusion. The Mozec™ - Rx PTCA Balloon Dilatation Catheter (balloon models 2.25mm to 4.50mm) is also indicated for post-delivery expansion of balloon expandable stents.



5.8 Comparison of Technological characteristics:

The Mozec™ – Rx PTCA Balloon Dilatation Catheter is similar to the predicate devices with respect to intended use, device design, materials, rated burst pressure, and method of sterilization.

5.9 Non clinical Performance data:

To ensure that the device design and construction are suitable for the intended use, the Mozec™ – Rx PTCA Balloon Dilatation Catheter was subjected to the performance testing recommended in the *Guidance for Industry and FDA Staff- Class II Special Controls Guidance Document for Certain Percutaneous Transluminal Coronary Angioplasty (PTCA) Catheters* (September 8, 2010).

The safety and effectiveness of the Mozec™ – Rx PTCA Balloon Dilatation Catheter has been evaluated in the following non Clinical tests;

- Dimensional verification
- Balloon preparation, deployment and retraction
- Balloon rated burst pressure
- Balloon compliance
- Balloon fatigue
- Balloon inflation and deflation time
- Catheter bond strength(s)
- Tip pull test
- Flexibility and kink test
- Radiopacity
- Coating integrity
- Particulate evaluation
- Balloon rated burst pressure (In stent)
- Balloon fatigue (In Stent)
- Biocompatibility testing in compliance with the ISO 10993-1
 - Cytotoxicity
 - Sensitization
 - Irritation / Intracutaneous reactivity
 - Systemic Toxicity



- Hemolysis
- C3a complement activation
- SC5b-9 complement activation
- In vivo thromboresistance
- Material mediated pyrogenicity

5.10 Conclusion

The MozecTM – Rx PTCA Balloon Dilatation Catheter met all predetermined acceptance criteria as specified by applicable standards, FDA guidance documents and test protocols. No safety and efficacy issues were raised during the testing program.

The MozecTM – Rx PTCA Balloon Dilatation Catheter is similar to the predicate devices with respect to intended use, device design, materials, rated burst pressure, and method of sterilization. Therefore, Meril Life Sciences Pvt. Ltd. believes the MozecTM – Rx PTCA Balloon Dilatation Catheter is considered substantially equivalent to the predicate devices.



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

December 24, 2013

Meril Life Sciences Private Limited
C/O Judith Danielson
Senior Regulatory Consultant
CardioMed Device Consultants, LLC
5523 Research Park Drive, Suite 360
Baltimore, MD 21228

Re: K131169
Trade/Device Name: Mozec™ – Rx PTCA Balloon Dilatation Catheter
Regulation Number: 21 CFR 870.5100
Regulation Name: Percutaneous Transluminal Coronary Angioplasty (PTCA) Catheter
Regulatory Class: Class II
Product Code: LOX
Dated: December 6, 2013
Received: December 9, 2013

Dear Ms. Danielson,

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



for Bram D. Zuckerman, M.D.

Director
Division of Cardiovascular Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure



4. INDICATIONS FOR USE STATEMENT

Indications for Use

510(k) Number (if known): K131169

Device Name: Mozec™ – Rx PTCA Balloon Dilatation Catheter

Indications for Use: The Mozec™ - Rx PTCA Balloon Dilatation Catheter is indicated for balloon dilatation of the stenotic portion of a coronary artery or bypass graft stenosis for the purpose of improving myocardial perfusion. The Mozec™ - Rx PTCA Balloon Dilatation Catheter (balloon models 2.25mm to 4.50mm) is also indicated for post-delivery expansion of balloon expandable stents.

Prescription Use X

AND/OR

Over-The-Counter Use

(Part 21 CFR 801 Subpart D)

(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NECESSARY)

Concurrence of CDRH, Office of Device Evaluation (ODI) *[Signature]* Date: 2013.12.24
09:02:28
-05'00"