

K131351

510(k) SUMMARY

510(k) Owner:	<p>Alfa Wassermann Diagnostic Technologies, LLC 4 Henderson Drive West Caldwell, NJ 07006</p> <p>Contact: Hkatz@AlfaWassermannUS.com Hyman Katz, Ph.D. Phone: 973-852-0158 Fax: 973-852-0237</p> <p style="text-align: right;">AUG 15 2013</p>
Date Summary Prepared:	July 3, 2013
Device:	<p>Trade Name: ACE Alkaline Phosphatase Reagent Classification: Class 2 Common/Classification Name: Nitrophenylphosphate, Alkaline Phosphatase Or Isoenzymes (21 C. F.R. § 862.1050) Product Code CJE</p> <p>Trade Name: ACE Amylase Reagent Classification: Class 2 Common/Classification Name: Saccharogenic, Amylase (21 C. F.R. § 862.1070) Product Code CIJ</p> <p>Trade Name: ACE ALT Reagent Classification: Class 1 Common/Classification Name: NADH Oxidation/NAD Reduction, ALT/SGPT (21 C.F.R. § 862.1030) Product Code CKA</p> <p>Trade Name: ACE AST Reagent Classification: Class 2 Common/Classification Name: NADH Oxidation/NAD Reduction, AST/SGOT (21 C.F.R. § 862.1100) Product Code CIT</p>
Predicate Devices:	<p>Manufacturer for reagent system predicates: Alfa Wassermann ACE and ACE Axcel Clinical Chemistry Systems and ACE Reagents (K113253, K931786, K930104, K113436, K113382)</p>

<p>Device Descriptions:</p>	<p>In the ACE Alkaline Phosphatase Reagent assay, alkaline phosphatase catalyzes the hydrolysis of colorless p-nitrophenyl phosphate to p-nitrophenol and inorganic phosphate. In an alkaline solution (pH 10.5), p-nitrophenol is in the phenoxide form and has a strong absorbance at 408 nm. The rate of increase in absorbance, monitored bichromatically at 408 nm/486 nm, is directly proportional to the alkaline phosphatase activity in the sample.</p> <p>In the ACE Amylase Reagent assay, α-amylase hydrolyzes the 2-chloro-p-nitrophenyl-α-D-maltotrioxide substrate to release 2-chloro-p-nitrophenol and form 2-chloro-p-nitrophenyl-α-D-maltoside, maltotriose and glucose. The rate of increase in absorbance, monitored bichromatically at 408 nm/ 647 nm, is directly proportional to the α-amylase activity in the sample.</p> <p>In the ACE ALT Reagent assay, alanine aminotransferase converts the L-alanine and α-ketoglutarate substrates in the reagent to L-glutamate and pyruvate, respectively. Lactate dehydrogenase (LDH) catalyzes the oxidation of the reduced cofactor to the cofactor. The rate of conversion of the reduced cofactor to the cofactor can be determined by monitoring the decrease in absorbance bichromatically at 340 nm/647 nm. This rate of conversion from the reduced cofactor to the cofactor is a function of the activity of ALT in the sample.</p> <p>In the ACE AST Reagent assay, aspartate aminotransferase converts the L-aspartate and α-ketoglutarate in the reagent to oxaloacetate and L-glutamate, respectively. The oxaloacetate undergoes reduction, with concurrent oxidation of NADH to NAD⁺ in the malate dehydrogenase-catalyzed indicator reaction. NADH absorbs strongly at 340 nm, whereas NAD⁺ does not. Therefore, the rate of conversion of NADH to NAD⁺ can be determined by monitoring the decrease in absorbance bichromatically at 340 nm/647 nm. This rate of conversion from NADH to NAD⁺ is a function of the activity of AST in the sample. Lactate dehydrogenase is added to prevent interference from endogenous pyruvate, which is normally present in blood.</p>
<p>Intended Use:</p>	<p>Indications for Use:</p> <p>The ACE Alkaline Phosphatase Reagent is intended for the quantitative determination of alkaline phosphatase activity in serum and lithium heparin plasma using the ACE, ACE <i>Alera</i>, and ACE Axcel Clinical Chemistry Systems. Measurements of alkaline phosphatase are used in the diagnosis and treatment of liver, bone, parathyroid, and intestinal diseases. This test is intended for use in clinical laboratories and physician office laboratories. For <i>in vitro</i> diagnostic use only.</p> <p>The ACE Amylase Reagent is intended for the quantitative determination of α-amylase activity in serum and lithium heparin plasma using the ACE, ACE <i>Alera</i>, and ACE Axcel Clinical Chemistry Systems. Amylase measurements are used primarily for the diagnosis and treatment of pancreatitis (inflammation of the pancreas). This test is intended for use in clinical laboratories and physician office laboratories. For <i>in vitro</i> diagnostic use only.</p>

Intended Use:	<p>The ACE ALT Reagent is intended for the quantitative determination of alanine aminotransferase activity in serum and lithium heparin plasma using the ACE, ACE <i>Alera</i>, and ACE Axcel Clinical Chemistry Systems. Alanine aminotransferase measurements are used in the diagnosis and treatment of certain liver diseases (e.g., viral hepatitis and cirrhosis) and heart diseases. This test is intended for use in clinical laboratories and physician office laboratories. For <i>in vitro</i> diagnostic use only.</p> <p>The ACE AST Reagent is intended for the quantitative determination of aspartate aminotransferase activity in serum and lithium heparin plasma using the ACE, ACE <i>Alera</i>, and ACE Axcel Clinical Chemistry Systems. Measurements of aspartate aminotransferase are used in the diagnosis and treatment of certain types of liver and heart disease. This test is intended for use in clinical laboratories and physician office laboratories. For <i>in vitro</i> diagnostic use only.</p>
Technological Characteristics:	<p>The ACE Alkaline Phosphatase Reagent is composed of two reagent bottles (Buffer and Substrate Reagent). The reagents contain AMP Buffer (pH 10.45), magnesium acetate, and p-nitrophenyl phosphate.</p> <p>The ACE Amylase Reagent is composed of a single reagent bottle. The reagents contain 2-chloro-p-nitrophenyl-α-D-maltotrioxide, sodium chloride, calcium acetate, potassium thiocyanate, and MES buffer (pH 6.0).</p> <p>The ACE ALT Reagent consists of two reagent bottles (Substrate and Coenzyme). The reagents contain L-alanine, α-ketoglutarate, nicotinamide adenine dinucleotide-reduced (NADH), lactate dehydrogenase, and Tris buffer.</p> <p>The ACE AST Reagent consists of two reagent bottles (Substrate and Coenzyme). The reagents contain L-aspartate, α-ketoglutarate, nicotinamide adenine dinucleotide-reduced (NADH), malate dehydrogenase, lactate dehydrogenase, and Tris buffer.</p>

Device Comparison with Predicate	<u>Comparison of similarities and differences with predicate device</u>		
	ACE Alkaline Phosphatase Reagent		
	ALP	Candidate Device	Predicate Device K931786 (ACE ALP)
	Intended Use/ Indications for Use	The ACE Alkaline Phosphatase Reagent is intended for the quantitative determination of alkaline phosphatase activity.	Same
	Platforms	ACE, ACE <i>Alera</i> [®] , and ACE Axcel Clinical Chemistry Systems	ACE Clinical Chemistry System
	Method	Photometric	Same
	Calibration Stability	Not a calibrated test	Same
	On-Board Stability	SA2002: 20 Days RX2002: 7 Days	Same
	Sample Type	Serum and lithium heparin plasma	Serum
	Sample Volume	4 µL	Same
	Reaction Volume	169 µL	Same
	Expected Values	44 - 147 U/L	Same
	Measuring Range	9 - 1400 U/L	Same
	Sample Stability	Serum ALP is stable for 7 days at 4-8°C and for 2 months at -20°C	Same
	ACE Amylase Reagent		
	AMYLASE	Candidate Device	Predicate Device K931786 (ACE Amylase)
	Intended Use/ Indications for Use	The ACE Amylase Reagent is intended for the quantitative determination of α -amylase activity.	Same
	Platforms	ACE, ACE <i>Alera</i> [®] , and ACE Axcel Clinical Chemistry Systems	ACE Clinical Chemistry System
	Method	Photometric	Same
	Calibration Stability	Not a calibrated test	Same
	On-Board Stability	30 Days	Same
	Sample Type	Serum and lithium heparin plasma	Serum
	Sample Volume	3 µL	Same
Reaction Volume	168 µL	Same	
Expected Values	20 - 104 U/L	Same	
Measuring Range	9 - 1900 U/L	Same	
Sample Stability	Serum amylase is stable for 7 days at room temperature (18-26°C) and at 4°C for one month. Recommended storage is at 4°C.	Same	

Device Comparison with Predicate	ACE ALT Reagent		
	ALT	Candidate Device	Predicate Device K930104 (ACE ALT)
	Intended Use/ Indications for Use	The ACE ALT Reagent is intended for the quantitative determination of alanine aminotransferase activity concentration.	Same
	Platforms	ACE, ACE <i>Alera</i> [®] , and ACE Axcel Clinical Chemistry Systems	ACE Clinical Chemistry System
	Method	Photometric	Same
	Calibration Stability	Not a calibrated test	Same
	On-Board Stability	30 Days	Same
	Sample Type	Serum and lithium heparin plasma	Serum
	Sample Volume	13 µL	Same
	Reaction Volume	185 µL	Same
	Expected Values	5 - 30 U/L	Same
	Measuring Range	4 - 480 U/L	Same
	Sample Stability	Specimen is stable for 7 days at 4-8°C and -20°C.	Same
	ACE AST Reagent		
	AST	Candidate Device	Predicate Device K930104 (ACE AST)
	Intended Use/ Indications for Use	The ACE AST Reagent is intended for the quantitative determination of aspartate aminotransferase activity.	Same
	Platforms	ACE, ACE <i>Alera</i> [®] , and ACE Axcel Clinical Chemistry Systems	ACE Clinical Chemistry System
	Method	Photometric	Same
	Calibration Stability	Not a calibrated test	Same
	On-Board Stability	30 Days	Same
Sample Type	Serum and lithium heparin plasma	Serum	
Sample Volume	13 µL	Same	
Reaction Volume	185 µL	Same	
Expected Values	7 - 31 U/L	Same	
Measuring Range	4 - 450 U/L	Same	
Sample Stability	Specimen activity is stable for 28 days at 4°C and at least one year at -20°C	Same	

Performance
Data:

In-House
Precision –
Serum vs.
Plasma

Performance data for the Alfa Wassermann ACE Reagents run on the Alfa Wassermann ACE, ACE *Alera* and ACE Axcel Clinical Chemistry Systems

In-House Precision: Serum vs. Plasma – ACE Alkaline Phosphatase Reagent

Precision (SD, %CV)									
ALP (n=20)	ACE			<i>Alera</i>			Axcel		
	Mean (U/L)	Within- Run	Total	Mean (U/L)	Within- Run	Total	Mean (U/L)	Within- Run	Total
Serum Low	92	1.2 , 1.3%	1.9 , 2.0%	91	0.8 , 0.9%	2.0 , 2.2%	90	1.2 , 1.3%	2.1 , 2.3%
Serum Mid	649	6.3 , 1.0%	8.6 , 1.3%	642	6.4 , 1.0%	6.4 , 1.0%	645	6.8 , 1.0%	6.9 , 1.1%
Serum High	1198	20.2 , 1.7%	20.8 , 1.7%	1190	5.6 , 0.5%	8.7 , 0.7%	1194	6.1 , 0.5%	7.9 , 0.7%
Plasma Low	76	1.9 , 2.5%	2.8 , 3.7%	75	0.8 , 1.1%	3.4 , 4.6%	74	1.1 , 1.5%	3.7 , 5.1%
Plasma Mid	614	5.8 , 0.9%	24.4 , 4.0%	609	5.1 , 0.8%	20.2 , 3.3%	613	3.4 , 0.6%	20.5 , 3.3%
Plasma High	1163	6.8 , 0.6%	33.5 , 2.9%	1149	5.9 , 0.5%	32.9 , 2.9%	1155	7.6 , 0.7%	35.7 , 3.1%

In-House Precision: Serum vs. Plasma – ACE Amylase Reagent

Precision (SD, %CV)									
AMY (n=20)	ACE			<i>Alera</i>			Axcel		
	Mean (U/L)	Within- Run	Total	Mean (U/L)	Within- Run	Total	Mean (U/L)	Within- Run	Total
Serum Low	46	0.7 , 1.5%	0.9 , 2.1%	46	1.0 , 2.1%	1.5 , 3.3%	46	1.2 , 2.6%	1.6 , 3.4%
Serum Mid	830	9.2 , 1.1%	11.5 , 1.4%	825	5.1 , 0.6%	11.9 , 1.4%	826	9.2 , 1.1%	11.9 , 1.4%
Serum High	1597	15.0 , 0.9%	19.8 , 1.2%	1577	9.9 , 0.6%	29.4 , 1.9%	1586	8.5 , 0.5%	18.6 , 1.2%
Plasma Low	41	0.7 , 1.8%	1.3 , 3.2%	42	0.5 , 1.3%	1.5 , 3.5%	41	1.0 , 2.5%	2.5 , 6.0%
Plasma Mid	806	8.1 , 1.0%	12.5 , 1.5%	801	4.4 , 0.5%	13.1 , 1.6%	805	5.3 , 0.7%	12.2 , 1.5%
Plasma High	1604	15.4 , 1.0%	29.3 , 1.8%	1596	18.8 , 1.2%	33.8 , 2.1%	1604	21.6 , 1.3%	25.0 , 1.6%

Performance
Data:

In-House
Precision –
Serum vs.
Plasma

In-House Precision: Serum vs. Plasma – ACE ALT Reagent

Precision (SD, %CV)									
ALT (n=24)	ACE			Alera			Axcel		
	Mean (U/L)	Within- Run	Total	Mean (U/L)	Within- Run	Total	Mean (U/L)	Within- Run	Total
Serum Low	36	0.9 , 2.6%	1.4 , 4.0%	36	0.6 , 1.7%	1.2 , 3.3%	37	1.1 , 3.1%	1.4 , 3.7%
Serum Mid	114	1.3 , 1.2%	1.8 , 1.6%	114	1.3 , 1.1%	2.2 , 2.0%	115	1.4 , 1.2%	2.3 , 2.0%
Serum High	216	2.9 , 1.4%	4.0 , 1.9%	216	1.8 , 0.8%	3.5 , 1.6%	218	1.4 , 0.6%	3.2 , 1.5%
Plasma Low	32	0.9 , 2.7%	1.3 , 4.0%	32	0.8 , 2.6%	1.5 , 4.7%	34	1.0 , 2.9%	1.4 , 4.1%
Plasma Mid	112	1.1 , 1.0%	1.5 , 1.4%	112	0.9 , 0.8%	1.1 , 1.0%	113	1.0 , 0.9%	1.5 , 1.4%
Plasma High	219	1.1 , 0.5%	1.6 , 0.7%	219	2.1 , 0.9%	3.0 , 1.4%	222	1.7 , 0.7%	2.7 , 1.2%

In-House Precision: Serum vs. Plasma – ACE AST Reagent

Precision (SD, %CV)									
AST (n=24)	ACE			Alera			Axcel		
	Mean (U/L)	Within- Run	Total	Mean (U/L)	Within- Run	Total	Mean (U/L)	Within- Run	Total
Serum Low	26	0.7 , 2.6%	1.2 , 4.5%	25	0.7 , 2.7%	0.8 , 3.3%	25	1.1 , 4.3%	1.4 , 5.7%
Serum Mid	155	1.3 , 0.8%	2.7 , 1.8%	155	1.0 , 0.7%	2.8 , 1.8%	155	0.7 , 0.5%	2.3 , 1.5%
Serum High	301	3.5 , 1.2%	4.8 , 1.6%	302	2.2 , 0.7%	4.0 , 1.3%	304	2.7 , 0.9%	4.1 , 1.4%
Plasma Low	26	0.8 , 3.2%	1.4 , 5.4%	26	0.6 , 2.5%	1.0 , 3.7%	26	0.9 , 3.5%	1.0 , 4.0%
Plasma Mid	157	1.6 , 1.0%	2.3 , 1.4%	157	1.5 , 0.9%	1.7 , 1.1%	158	1.4 , 0.9%	2.0 , 1.3%
Plasma High	304	3.5 , 1.2%	4.5 , 1.5%	303	2.6 , 0.9%	3.9 , 1.3%	305	3.2 , 1.1%	4.3 , 1.4%

Performance Data: In-House Matrix Comparison – Serum vs. Plasma	Performance data for the Alfa Wassermann ACE Reagents run on the Alfa Wassermann ACE, ACE <i>Alera</i> and ACE Axcel Clinical Chemistry Systems	
	<u>In-House Matrix Comparison: Serum vs. Plasma – ACE ALP Reagent</u>	
	System	Range
	Results - Serum vs. Plasma	
ACE	9 - 1274 U/L	Slope: 0.998
108 pairs		Intercept: -8.3
		Correlation: 0.9980
		Std. Error Est: 13.5
		Confidence Interval Slope: 0.986 to 1.010
	Confidence Interval Intercept: -11.5 to -5.1	
ACE <i>Alera</i>	9 - 1202 U/L	Slope: 0.983
108 pairs		Intercept: -6.4
		Correlation: 0.9952
		Std. Error Est: 20.2
		Confidence Interval Slope: 0.965 to 1.002
	Confidence Interval Intercept: -11.2 to -1.6	
ACE Axcel	11 - 1222 U/L	Slope: 1.017
62 pairs		Intercept: -6.5
		Correlation: 0.9982
		Std. Error Est: 14.5
		Confidence Interval Slope: 1.001 to 1.033
	Confidence Interval Intercept: -11.1 to -1.8	
<u>In-House Matrix Comparison: Serum vs. Plasma – ACE Amylase Reagent</u>		
System	Range	Results - Serum vs. Plasma
ACE	11 - 1766 U/L	Slope: 0.977
104 pairs		Intercept: 1.7
		Correlation: 0.9995
		Std. Error Est: 8.8
		Confidence Interval Slope: 0.970 to 0.983
	Confidence Interval Intercept: -0.2 to 3.6	
ACE <i>Alera</i>	11 - 1703 U/L	Slope: 0.979
101 pairs		Intercept: 0.9
		Correlation: 0.9994
		Std. Error Est: 9.0
		Confidence Interval Slope: 0.972 to 0.986
	Confidence Interval Intercept: -1.0 to 2.9	
ACE Axcel	10 - 1890 U/L	Slope: 0.994
52 pairs		Intercept: -1.76
		Correlation: 0.9996
		Std. Error Est: 11.54
		Confidence Interval Slope: 0.986 to 1.002
	Confidence Interval Intercept: -5.33 to 1.80	

Performance Data: In-House Matrix Comparison – Serum vs. Plasma	<u>In-House Matrix Comparison: Serum vs. Plasma – ACE ALT Reagent</u>		
	System	Range	Results - Serum vs. Plasma
	ACE	4 - 460 U/L	Slope: 1.003
	54 pairs		Intercept: -3.6
			Correlation: 0.9994
			Std. Error Est: 2.8
			Confidence Interval Slope: 0.994 to 1.013
		Confidence Interval Intercept: -4.5 to -2.8	
	ACE Alera	5 - 463 U/L	Slope: 1.000
	52 pairs		Intercept: -3.6
Correlation: 0.9986			
Std. Error Est: 4.3			
Confidence Interval Slope: 0.985 to 1.015			
	Confidence Interval Intercept: -4.9 to -2.2		
ACE Axcel	6 - 469 U/L	Slope: 0.985	
56 pairs		Intercept: -3.35	
		Correlation: 0.9993	
		Std. Error Est: 2.94	
		Confidence Interval Slope: 0.976 to 0.995	
	Confidence Interval Intercept: -4.24 to -2.47		
<u>In-House Matrix Comparison: Serum vs. Plasma – ACE AST Reagent</u>			
System	Range	Results - Serum vs. Plasma	
ACE	4 - 404 U/L	Slope: 0.981	
53 pairs		Intercept: 0.7	
		Correlation: 0.9992	
		Std. Error Est: 3.2	
		Confidence Interval Slope: 0.971 to 0.992	
	Confidence Interval Intercept: -0.3 to 1.7		
ACE Alera	4 - 396 U/L	Slope: 0.999	
52 pairs		Intercept: -0.6	
		Correlation: 0.9989	
		Std. Error Est: 3.9	
		Confidence Interval Slope: 0.985 to 1.012	
	Confidence Interval Intercept: -1.9 to 0.6		
ACE Axcel	4 - 408 U/L	Slope: 1.001	
52 pairs		Intercept: 0.22	
		Correlation: 0.9989	
		Std. Error Est: 3.97	
		Confidence Interval Slope: 0.987 to 1.014	
	Confidence Interval Intercept: -1.04 to 1.49		

Performance
Data:

Precision -
POL

POL – Precision for ACE and ACE Alera Clinical Chemistry Systems

(Note: Please refer to previously cleared submissions k113436 (ALP and Amylase) and k113382 (ALT and AST) for ACE Axcel POL data)

ALP			ACE Results		Mean	ACE Alera Results	
Lab	Sample	Mean	U/L SD, %CV			Within-Run	Total
			Within-Run	Total	Within-Run		
In-House	1	59	1.2	1.4	60	1.1	1.3
			2.0%	2.5%		1.8%	2.1%
POL 1	1	55	0.7	2.3	56	0.8	1.7
			1.4%	4.1%		1.4%	3.0%
POL 2	1	54	0.9	1.6	59	1.2	2.1
			1.6%	2.9%		2.0%	3.5%
POL 3	1	57	1.0	1.5	56	1.4	3.0
			1.8%	2.6%		2.5%	5.4%
In-House	2	648	5.5	6.7	653	4.5	7.0
			0.8%	1.0%		0.7%	1.1%
POL 1	2	632	6.9	13.4	626	7.4	16.7
			1.1%	2.1%		1.2%	2.7%
POL 2	2	631	6.2	11.9	659	4.2	18.0
			1.0%	1.9%		0.6%	2.7%
POL 3	2	642	3.7	8.6	640	5.9	21.9
			0.6%	1.3%		0.9%	3.4%
In-House	3	1191	7.0	10.7	1192	9.4	13.4
			0.6%	0.9%		0.8%	1.1%
POL 1	3	1145	10.3	23.0	1135	19.0	25.0
			0.9%	2.0%		1.7%	2.2%
POL 2	3	1159	12.8	17.2	1209	9.6	29.3
			1.1%	1.5%		0.8%	2.4%
POL 3	3	1185	6.2	7.3	1165	6.6	37.2
			0.5%	0.6%		0.6%	3.2%

Performance
Data at POL:Precision -
POL**POL – Precision for ACE and ACE Alera Clinical Chemistry Systems**

AMY			ACE Results		Mean	ACE Alera Results	
Lab	Sample	Mean	U/L SD, %CV			Within-Run	Total
			Within-Run	Total	Within-Run		
In-House	1	39	1.0	1.3	39	0.8	1.4
			2.7%	3.2%		2.1%	3.5%
POL 1	1	37	0.5	1.2	38	0.9	1.7
			1.5%	3.2%		2.4%	4.4%
POL 2	1	39	0.7	0.9	40	1.3	1.3
			1.9%	2.4%		3.2%	3.2%
POL 3	1	39	0.6	1.1	39	1.0	1.1
			1.6%	2.7%		2.5%	2.8%
In-House	2	741	7.1	8.6	747	4.4	7.3
			1.0%	1.2%		0.6%	1.0%
POL 1	2	725	4.6	11.5	723	4.7	7.0
			0.6%	1.6%		0.6%	1.0%
POL 2	2	727	7.4	8.2	770	4.3	6.1
			1.0%	1.1%		0.6%	0.8%
POL 3	2	737	9.6	11.3	747	5.8	7.4
			1.3%	1.5%		0.8%	1.0%
In-House	3	1429	7.5	13.1	1437	11.6	12.8
			0.5%	0.9%		0.8%	0.9%
POL 1	3	1389	18.1	27.8	1388	19.5	21.6
			1.3%	2.0%		1.4%	1.6%
POL 2	3	1401	22.4	23.8	1500	10.3	11.7
			1.6%	1.7%		0.7%	0.8%
POL 3	3	1410	12.9	15.7	1435	8.4	14.4
			0.9%	1.1%		0.6%	1.0%

Performance
Data at POL:

Precision -
POL

POL – Precision for ACE and ACE Alera Clinical Chemistry Systems

ALT			ACE Results		Mean	ACE Alera Results	
Lab	Sample	Mean	U/L SD, %CV			Within-Run	Total
			Within-Run	Total	Within-Run		
In-House	1	33	0.9	1.3	32	1.1	1.2
			2.7%	3.9%		3.5%	3.9%
POL 1	1	26	0.5	1.7	28	1.4	2.4
			1.9%	6.4%		5.1%	8.4%
POL 2	1	27	0.9	1.8	26	0.9	2.1
			3.2%	6.7%		3.6%	8.2%
POL 3	1	30	1.0	1.7	29	2.1	2.4
			3.4%	5.7%		7.5%	8.4%
In-House	2	191	4.8	4.8	190	4.0	4.1
			2.5%	2.5%		2.1%	2.1%
POL 1	2	189	1.2	1.6	193	1.7	2.0
			0.6%	0.8%		0.9%	1.0%
POL 2	2	192	2.3	5.9	194	2.2	2.5
			1.2%	3.1%		1.1%	1.3%
POL 3	2	188	2.9	3.9	195	3.5	5.0
			1.5%	2.1%		1.8%	2.6%
In-House	3	309	2.2	3.7	307	4.0	4.1
			0.7%	1.2%		1.3%	1.3%
POL 1	3	304	4.0	4.4	309	3.1	3.6
			1.3%	1.4%		1.0%	1.2%
POL 2	3	311	2.9	10.2	314	2.3	2.6
			0.9%	3.3%		0.7%	0.8%
POL 3	3	303	3.4	5.0	310	8.5	9.1
			1.1%	1.6%		2.8%	3.0%

Performance
Data at POL:
Precision -
POL

POL – Precision for ACE and ACE Alera Clinical Chemistry Systems

AST			ACE Results		Mean	ACE Alera Results	
Lab	Sample	Mean	U/L SD, %CV			Within-Run	Total
			Within-Run	Total	Within-Run		
In-House	1	28	1.4	2.1	28	1.2	1.7
			5.1%	7.4%		4.5%	6.0%
POL 1	1	23	0.9	1.4	26	1.0	1.3
			3.7%	6.0%		3.8%	5.1%
POL 2	1	28	1.5	1.7	27	2.0	2.4
			5.6%	6.0%		7.5%	8.9%
POL 3	1	26	0.7	2.0	29	2.9	2.9
			2.7%	7.7%		9.9%	9.9%
In-House	2	223	6.4	6.6	222	4.8	7.1
			2.9%	3.0%		2.2%	3.2%
POL 1	2	220	2.2	4.0	220	2.0	3.2
			1.0%	1.8%		0.9%	1.5%
POL 2	2	227	2.0	4.8	233	3.8	5.0
			0.9%	2.1%		1.6%	2.1%
POL 3	2	223	2.6	4.7	229	5.3	7.2
			1.2%	2.1%		2.3%	3.2%
In-House	3	410	6.3	6.8	406	2.9	6.8
			1.5%	1.7%		0.7%	1.7%
POL 1	3	416	5.0	5.5	416	7.8	9.2
			1.2%	1.3%		1.9%	2.2%
POL 2	3	420	9.0	10.2	428	5.1	5.6
			2.2%	2.4%		1.2%	1.3%
POL 3	3	407	3.9	7.3	417	8.2	12.1
			0.9%	1.8%		2.0%	2.9%

Performance
Data:
Method
Comparison -
POL on ACE

POL – Method Comparison for ACE Clinical Chemistry System

Reagent	Statistic	In-House ACE (x)	In-House ACE (x)	In-House ACE(x)
		vs. POL 1 ACE (y)	vs. POL 2 ACE (y)	vs. POL 3 ACE (y)
ALP	n	49	50	50
	Range	58 to 1199	58 to 1199	58 to 1199
	Regression	$y = 0.989x - 9.5$	$y = 0.977x - 7.9$	$y = 0.982x - 2.8$
	Correlation	0.9987	0.9997	0.9995
	Std. Error Est.	12.4	5.9	7.4
	CI Slope	0.975 to 1.004	0.970 to 0.984	0.973 to 0.990
	CI Intercept	-13.8 to -5.1	-9.9 to -5.8	-5.4 to -0.2
	AMY	n	51	51
Range	28 to 1732	28 to 1732	28 to 1732	
Regression	$y = 0.970x + 1.5$	$y = 0.973x + 3.6$	$y = 0.974x + 3.9$	
Correlation	0.9995	0.9998	0.9998	
Std. Error Est.	11.1	7.5	7.4	
CI Slope	0.962 to 0.979	0.967 to 0.978	0.968 to 0.980	
CI Intercept	-2.1 to 5.2	1.1 to 6.0	1.5 to 6.3	
ALT	n	44	47	49
	Range	6 to 442	6 to 442	6 to 442
	Regression	$y = 1.006x - 4.7$	$y = 1.021x - 4.0$	$y = 0.982x - 2.3$
	Correlation	0.9978	0.9993	0.9981
	Std. Error Est.	6.5	3.7	5.7
	CI Slope	0.985 to 1.026	1.010 to 1.033	0.964 to 0.999
	CI Intercept	-6.9 to -2.4	-5.2 to -2.8	-4.1 to -0.4
	AST	n	50	50
Range		6 to 413	6 to 413	6 to 413
Regression		$y = 0.999x - 0.6$	$y = 1.019x + 2.4$	$y = 0.992x + 0.6$
Correlation		0.9993	0.9989	0.9994
Std. Error Est.		3.9	4.9	3.6
CI Slope		0.988 to 1.010	1.005 to 1.033	0.982 to 1.003
CI Intercept		-1.9 to 0.6	0.9 to 4.0	-0.6 to 1.7

Performance
Data at POL:

Method
Comparison -
POL on ACE
Alera

POL – Method Comparison for ACE *Alera* Clinical Chemistry System

Reagent	Statistic	In-House ACE (x)	In-House ACE (x)	In-House ACE (x)
		vs. POL 1 <i>Alera</i> (y)	vs. POL 2 <i>Alera</i> (y)	vs. POL 3 <i>Alera</i> (y)
ALP	n	50	50	50
	Range	58 to 1199	58 to 1199	58 to 1199
	Regression	$y = 0.997x - 4.6$	$y = 1.029x - 4.1$	$y = 1.010x - 6.6$
	Correlation	0.9992	0.9991	0.9986
	Std. Error Est.	10.1	10.8	13.0
	CI Slope	0.985 to 1.008	1.016 to 1.041	0.995 to 1.025
	CI Intercept	-8.1 to -1.1	-7.9 to -0.4	-11.2 to -2.1
AMY	n	51	51	51
	Range	28 to 1732	28 to 1732	28 to 1732
	Regression	$y = 0.960x + 3.0$	$y = 1.010x + 5.8$	$y = 0.990x + 3.7$
	Correlation	0.9991	0.9995	0.9995
	Std. Error Est.	15.1	11.7	11.3
	CI Slope	0.948 to 0.971	1.001 to 1.018	0.981 to 0.999
	CI Intercept	-1.9 to 7.9	2.0 to 9.6	0.0 to 7.4
ALT	n	50	48	50
	Range	6 to 442	6 to 442	6 to 442
	Regression	$y = 1.019x - 0.5$	$y = 1.012x - 3.5$	$y = 0.970x + 2.4$
	Correlation	0.9986	0.9985	0.9977
	Std. Error Est.	5.1	5.3	6.1
	CI Slope	1.003 to 1.035	0.995 to 1.028	0.951 to 0.990
	CI Intercept	-2.1 to 1.1	-5.3 to -1.8	0.5 to 4.4
AST	n	50	50	50
	Range	6 to 413	6 to 413	6 to 413
	Regression	$y = 1.028x + 1.4$	$y = 1.040x + 0.5$	$y = 1.004x + 1.8$
	Correlation	0.9995	0.9992	0.9995
	Std. Error Est.	3.5	4.3	3.3
	CI Slope	1.018 to 1.037	1.027 to 1.052	0.994 to 1.013
	CI Intercept	0.3 to 2.5	-0.8 to 1.9	0.8 to 2.9

Performance
Data:
ACE Alera

**Performance data for the Alfa Wassermann ACE Reagents run on the
Alfa Wassermann ACE *Alera* Clinical Chemistry Systems**

Detection Limits - ACE *Alera* Clinical Chemistry System

<i>ACE Alera</i>	ALP	Amylase	ALT	AST
LoB (U/L)	2.8	0.2	1.6	2.2
LoD (U/L)	4.8	0.9	3.1	3.3
LoQ (U/L)	4.8	5.6	4.1	3.3

Linearity - ACE *Alera* Clinical Chemistry System

ACE Reagents	Low level tested	Upper level tested	Linear to:	Linear Regression Equation	R ²
ALP	4.0 U/L	1401 U/L	1400 U/L	$y = 0.998x - 0.5$	0.9993
Amylase	4.0 U/L	2012 U/L	1900 U/L	$y = 1.013x + 0.2$	0.9974
ALT	3.1 U/L	504 U/L	480 U/L	$y = 1.007x - 0.17$	0.9992
AST	3.0 U/L	491 U/L	450 U/L	$y = 1.013x + 0.24$	0.9992

Performance
Data:

ACE Alera

Interferences - ACE Alera Clinical Chemistry System

Interferents on ACE Alera	No Significant Interference at or below:			
	ALP	Amylase	ALT	AST
Icterus	70.6 mg/dL	30.0 mg/dL	50 mg/dL	50 mg/dL
Hemolysis	62.5 mg/dL	62.5 mg/dL	500 mg/dL	62.5 mg/dL
Lipemia	1000 mg/dL	1000 mg/dL	419 mg/dL	439 mg/dL
Ascorbic Acid	6 mg/dL	6 mg/dL	6 mg/dL	6 mg/dL

Precision - ACE Alera Clinical Chemistry System

on ACE Alera		Precision (SD, %CV)		
		Mean	Within-Run	Total
ALP U/L	Low	43	0.9 , 2.2%	1.9 , 4.5%
	Mid	164	2.8 , 1.7%	5.5 , 3.4%
	High	339	5.3 , 1.6%	8.9 , 2.6%
AMYLASE U/L	Low	45	0.9 , 2.1%	1.5 , 3.3%
	Mid	139	1.7 , 1.2%	3.9 , 2.8%
	High	311	5.0 , 1.6%	7.7 , 2.5%
ALT U/L	Low	37	1.2 , 3.3%	1.8 , 4.9%
	Mid	181	2.2 , 1.2%	3.2 , 1.8%
	High	320	5.2 , 1.6%	5.2 , 1.6%
AST U/L	Low	34	1.1 , 3.3%	1.5 , 4.4%
	Mid	175	3.5 , 2.0%	3.6 , 2.1%
	High	312	4.8 , 1.5%	4.9 , 1.6%

Performance Data: <i>ACE Alera</i>	<u>Method Comparison - ACE Alera Clinical Chemistry System</u>				
	In-House ACE (x) vs. In-House ACE Alera (y)				
		ALP	Amylase	ALT	AST
	n	60	57	50	50
	Range (U/L)	9 to 1199	10 to 1732	6 to 442	6 to 413
	Slope	0.991	0.995	0.988	1.006
	Intercept	0.5	2.9	1.3	1.5
	Correlation Coefficient	0.9999	0.9998	0.9999	0.9998
	Std. Error	3.5	6.6	1.6	1.9
	CI Slope	0.987 to 0.995	0.990 to 1.000	0.983 to 0.992	1.001 to 1.012
CI Intercept	-0.5 to 1.6	0.9 to 4.9	0.8 to 1.8	0.9 to 2.2	
Conclusions:	Based on the foregoing data, the device is safe and effective for use in clinical laboratories and physician office laboratories. These data indicate substantial equivalence for lithium heparin plasma sample collection tubes to the predicate device's use of serum sample collection tubes.				



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

August 15, 2013

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

Alfa Wassermann Diagnostic Technologies, LLC
C/O Hyman Katz, Ph.D.
4 Henderson Drive
WEST CALDWELL NJ 07006

Re: K131351

Trade/Device Name: ACE Alkaline Phosphatase Reagent
ACE Amylase Reagent
ACE ALT Reagent
ACE AST Reagent

Regulation Number: 21 CFR 862.1050

Regulation Name: Alkaline phosphatase or isoenzymes test system

Regulatory Class: II

Product Code: CJE, CIJ, CKA, CIT

Dated: July 8, 2013

Received: July 9, 2013

Dear Dr. Katz:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Carol C. Benson -S for

Courtney H. Lias, Ph.D.
Director
Division of Chemistry and Toxicology Devices
Office of In Vitro Diagnostics
and Radiological Health
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): k131351

Device Name: ACE Alkaline Phosphatase Reagent

Indications for Use: The ACE Alkaline Phosphatase Reagent is intended for the quantitative determination of alkaline phosphatase activity in serum and lithium heparin plasma using the ACE, ACE *Alera*, and ACE Axcel Clinical Chemistry Systems. Measurements of alkaline phosphatase are used in the diagnosis and treatment of liver, bone, parathyroid, and intestinal diseases. This test is intended for use in clinical laboratories and physician office laboratories. For *in vitro* diagnostic use only.

Device Name: ACE Amylase Reagent

Indications for Use: The ACE Amylase Reagent is intended for the quantitative determination of α -amylase activity in serum and lithium heparin plasma using the ACE, ACE *Alera*, and ACE Axcel Clinical Chemistry Systems. Amylase measurements are used primarily for the diagnosis and treatment of pancreatitis (inflammation of the pancreas). This test is intended for use in clinical laboratories and physician office laboratories. For *in vitro* diagnostic use only.

Prescription Use X
(21 CFR Part 801 Subpart D)

AND/OR

Over-The-Counter Use.
(21 CFR Part 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE; CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Devices or Radiological Health (OIR)

Ruth A. Chesler -S

Division Sign-Off
Office of In Vitro Devices or Radiological Health
510(k) k131351

Indications for Use

510(k) Number (if known): k131351

Device Name: ACE ALT Reagent

Indications for Use: The ACE ALT Reagent is intended for the quantitative determination of alanine aminotransferase activity in serum and lithium heparin plasma using the ACE, ACE *Alera*, and ACE Axcel Clinical Chemistry Systems. Alanine aminotransferase measurements are used in the diagnosis and treatment of certain liver diseases (e.g., viral hepatitis and cirrhosis) and heart diseases. This test is intended for use in clinical laboratories and physician office laboratories. For *in vitro* diagnostic use only.

Device Name: ACE AST Reagent

Indications for Use: The ACE AST Reagent is intended for the quantitative determination of aspartate aminotransferase activity in serum and lithium heparin plasma using the ACE, ACE *Alera*, and ACE Axcel Clinical Chemistry Systems. Measurements of aspartate aminotransferase are used in the diagnosis and treatment of certain types of liver and heart disease. This test is intended for use in clinical laboratories and physician office laboratories. For *in vitro* diagnostic use only.

Prescription Use X
(21 CFR Part 801 Subpart D)

AND/OR

Over-The-Counter Use.
(21 CFR Part 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE; CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of In Vitro Devices or Radiological Health (OIR)

Ruth A. Chesler -S

Division Sign-Off
Office of In Vitro Devices or Radiological Health
510(k) k131351