

## 510(k) Summary of Safety and Effectiveness

This summary of 510(k) safety and effectiveness information is being submitted in accordance with the requirements of 21 CFR 807.92.

The assigned 510(k) number is:     K132059    

### Submitter

Advanced Instrumentations, Inc.  
6800 N.W. 77<sup>th</sup> Court  
Miami, FL 33166  
Telephone: 305-477-6331  
Fax: 305-477-5351

Registration # 1066270

### Official correspondent:

Jorge Millan, PhD  
Email: [jmillan@hiatec.org](mailto:jmillan@hiatec.org)  
601 West 20 St  
Hialeah, FL 33010  
Phone : (305) 925-1260

### Date Prepared:

October 24, 2013

OCT 25 2013

### Device name and classification:

- **Device Name:** DUS-5000 Diagnostic Ultrasound System
- **Classification Name:** 892.1550 System, Imaging, Pulsed Doppler Ultrasonic  
**Product code:** IYN  
892.1560 Ultrasonic, Pulsed echo, Imaging  
**Product code:** IYO  
892.1570 Transducer, Ultrasonic, Diagnostic  
**Product code:** ITX

- **Regulatory Class:** Class II

**Predicate Device:**

Model U50 Diagnostic Ultrasound System, K123249 Manufacturer: EDAN Instruments

**Device Description:**

The DUS-5000 Diagnostic Ultrasound System, which applies advanced technologies such as Phased Inversion Harmonic Compound Imaging (eHCI), Multi-Beam-Forming (mBeam), Speckle Resistance Imaging (eSRI), and Spatial Compounding Imaging, etc. Various image parameter adjustments, 12.1 inch LCD and diverse probes are configured to provide clear and stable images.

Its function is to acquire and display Ultrasound images in B-mode, M-mode, PW-mode, Color-mode, PDI/DPDI mode. This system provides a series of probes that include linear array, convex array, micro-convex array with a frequency range of approximately 2.5 MHz to 11 MHz.

**Intended Use:**

The diagnostic ultrasound system (DUS-5000) is applicable for adults, pregnant women, pediatric patients' ultrasound evaluation in hospitals and clinics. It is intended for use in abdominal, obstetrics, gynecology, pediatric, small parts, urology, peripheral vascular, musculoskeletal (conventional and superficial), transvaginal and cardiac clinical applications, by or on the order of a physician \*or similarly qualified health care professional.

**Effectiveness and Safety Contraindications:**

**Clinical Test**

Clinical testing is not required

**Non-clinical test:**

DUS-5000 Diagnostic Ultrasound System

---

The following safety standards are conducted on the subject device:

1. IEC 60601 -1 Electrical Safety
2. IEC 60601-1-2 Electromagnetic Compatibility
3. Acoustic output testing as per the guideline "Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Transducers" dated September 9, 2008.
4. ISO 10993-1, ISO 10993-5 and ISO 10993-10

**Comparison to the predicate device:**

The subject device has identical technology characteristics, same intended use, same product design, materials and manufacturing process, as well as same performance effectiveness, performance safety and the same needle-guide bracket material, property and sterilization methods as the predicate device.

**Substantially Equivalent Determination:**

This premarket notification submission demonstrates that DUS-5000 Diagnostic Ultrasound System is substantially equivalent to the predicate device.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

October 25, 2013

Advanced Instrumentations, Inc.  
% Jorge Millan, Ph.D.  
Official Correspondent  
Hialeah Technology Center, Inc.  
601 West 20 Street  
HIALEAH FL 33010

Re: K132059

Trade/Device Name: DUS-5000 Diagnostic Ultrasound System  
Regulation Number: 21 CFR 892.1550  
Regulation Name: Ultrasonic pulsed doppler imaging system  
Regulatory Class: II  
Product Code: IYN, IYO, ITX  
Dated: October 8, 2013  
Received: October 10, 2013

Dear Dr. Millan:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

This determination of substantial equivalence applies to the following transducers intended for use with the DUS-5000 Diagnostic Ultrasound System, as described in your premarket notification:

Transducer Model Number

C352UB  
L1042UB  
L742UB

E612UB  
C612UB  
C6152UB

C422UB  
L552UB

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



for

Janine M. Morris  
Director, Division of Radiological Health  
Office of In Vitro Diagnostics  
and Radiological Health  
Center for Devices and Radiological Health

Enclosure

Advanced Instrumentations  
DUS-5000 Diagnostic Ultrasound System  
510K Submission

## Indications for Use

510(k) Number (if known): K132059

Device Name: DUS-5000 Diagnostic Ultrasound System

### Intended Use:

The diagnostic ultrasound system (DUS-5000) is applicable for adults, pregnant women, pediatric patients' ultrasound evaluation in hospitals and clinics. It is intended for use in abdominal, obstetrics, gynecology, pediatric, small parts, urology, peripheral vascular, musculoskeletal (conventional and superficial), transvaginal and cardiac clinical applications, by or on the order of a physician or similarly qualified health care professional.

Prescription Use  Or Over the Counter Use   
(21 CFR Part 801 Subpart D) (21 CFR Part 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH; Office of In Vitro Diagnostics and Radiological Health (OIR)



---

(Division Sign-Off)  
Division of Radiological Health  
Office of In Vitro Diagnostics and Radiological Health  
510(k)     K132059

## Diagnostic Ultrasound Indications for Use Form

### DUS-5000 Diagnostic Ultrasound System

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

General	Clinical Application	Mode of Operation						
		B	M	PW	CW	Color	Combined (Specify) <sup>[1]</sup>	Other (Specify) <sup>[2][3]</sup>
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal / Obstetrics	P	P	P		P	P	P
	Abdominal	P	P	P		P	P	P
	Intra-operative (Specify)							
	Intra-operative (Neuro logical)							
	Laparoscopic							
	Pediatric	P	P	P		P	P	P
	Small Organ (Specify) *	P	P	P		P	P	P
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal	P	P	P		P	P	P
	Trans-urethral							
	Musculo-skeletal (Conventional)	P	P	P		P	P	P
	Musculo-skeletal (Superficial)	P	P	P		P	P	P
	Intravascular							
Other (Specify) **	P	P	P		P	P	P	
Cardiac	Adult Cardiac	P	P	P		P	P	P
	Pediatric Cardiac	P	P	P		P	P	P
	Intravascular (Cardiac)							
	Trans-esoph. (Cardiac)							
	Intra- cardiac							
Peripheral vascular	Peripheral vascular	P	P	P		P	P	P
	Other (Specify)							

N = new indication; P = previously cleared by FDA; E = added under this appendix PDI= Power Doppler Imaging

Additional comments: Combined mode: B+M, B+PW, B+Color, B+PDI/DPDI, B+Color+PW, B+PDI/DPDI +PW

Note \* Small Organ includes Thyroid, Testes, Breast | Cleared applications under K123249

\*\* Other use includes Urology, Kidney, Gynecology

[1]: PDI: Power Doppler Imaging, DPDI: Directional Power Doppler Imaging

[2]: Biopsy Guidance

[3]: Harmonic Imaging. This feature does not use contrast agent.

(PLEASE DO NOT WRITE BELOW THIS LINE: - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH: Office of In Vitro Diagnostics and Radiological Health (OIR)

Prescription Use (Per 21 CFR 801.109)

## Diagnostic Ultrasound Indications for Use Form

### DUS-5000 with C352UB

#### Transducer

Intended Use: Diagnostic ultrasound imaging or fluid flow analysis of the human body as follows:

General	Clinical Application	Mode of Operation						
		B	M	PW	CW	Color	Combined (Specify) <sup>[1]</sup>	Other (Specify) <sup>[2][3]</sup>
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal / Obstetrics	P	P	P		P	P	P
	Abdominal	P	P	P		P	P	P
	Intra-operative (Specify)							
	Intra-operative (Neuro logical)							
	Laparoscopic							
	Pediatric							
	Small Organ (Specify) *							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Musculo-skeletal (Conventional)							
	Musculo-skeletal (Superficial)							
	Intravascular							
Other (Specify) **	P	P	P		P	P	P	
Cardiac	Adult Cardiac							
	Pediatric Cardiac							
	Intravascular (Cardiac)							
	Trans-esoph. (Cardiac)							
	Intra- cardiac							
Peripheral vascular	Peripheral vascular							
	Other (Specify)							

N = new indication; P = previously cleared by FDA; E = added under this appendix PDI= Power Doppler Imaging

Additional comments: Combined mode: B+M, B+PW, B+Color, B+PDI/DPDI, B+Color+PW, B+PDI/DPDI +PW

Note \* Small Organ includes Thyroid, Testes, Breast |Cleared applications under K123249

\*\* Other use includes Urology, Kidney, Gynecology

[1]: PDI: Power Doppler Imaging ,DPDI: Directional Power Doppler Imaging

[2]: Biopsy Guidance

[3]: Harmonic Imaging. This feature does not use contrast agent.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH; Office of In Vitro Diagnostics and Radiological Health (OIR)

Prescription Use (Per 21 CFR 801.109)

## Diagnostic Ultrasound Indications for Use Form

### DUS-5000 with L1042UB

#### Transducer

General	Clinical Application	Mode of Operation						
		B	M	PW	CW	Color	Combined (Specify) <sup>[1]</sup>	Other (Specify) <sup>[2][3]</sup>
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal / Obstetrics							
	Abdominal							
	Intra-operative (Specify)							
	Intra-operative (Neuro logical)							
	Laparoscopic							
	Pediatric							
	Small Organ (Specify) *	P	P	P		P	P	P
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Musculo-skeletal(Conventional)	P	P	P		P	P	P
	Musculo-skeletal (Superficial)	P	P	P		P	P	P
	Intravascular							
Other (Specify) **								
Cardiac	Adult Cardiac							
	Pediatric Cardiac							
	Intravascular(Cardiac)							
	Trans-esoph.(Cardiac)							
	Intra- cardiac							
Peripheral vascular	Peripheral vascular	P	P	P		P	P	P
	Other (Specify)							

N = new indication; P = previously cleared by FDA; E = added under this appendix PDI= Power Doppler Imaging

Additional comments: Combined mode: B+M, B+PW, B+Color, B+PDI/DPDI, B+Color+PW, B+PDI/DPDI +PW

Note \* Small Organ includes Thyroid, Testes, Breast (Cleared applications under K123249)

\*\* Other use includes Urology, Kidney, Gynecology

[1]: PDI; Power Doppler Imaging, DPDI: Directional Power Doppler Imaging

[2]: Biopsy Guidance

[3]: Harmonic Imaging, This feature does not use contrast agent.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH; Office of In Vitro Diagnostics and Radiological Health (OIR)

Prescription Use (Per 21 CFR 801.109)

**Diagnostic Ultrasound Indications for Use Form**

**DUS-5000 with L742UB  
Transducer**

General	Clinical Application	Mode of Operation						
		B	M	PW	CW	Color	Combined (Specify) <sup>[1]</sup>	Other (Specify) <sup>[2][3]</sup>
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal / Obstetrics							
	Abdominal							
	Intra-operative (Specify)							
	Intra-operative (Neuro logical)							
	Laparoscopic							
	Pediatric							
	Small Organ (Specify) *	P	P	P		P	P	P
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Musculo-skeletal (Conventional)	P	P	P		P	P	P
	Musculo-skeletal (Superficial)	P	P	P		P	P	P
	Intravascular							
Other (Specify) **								
Cardiac	Adult Cardiac							
	Pediatric Cardiac							
	Intravascular(Cardiac)							
	Trans-esoph.(Cardiac)							
	Intra- cardiac							
Peripheral vascular	Peripheral vascular	P	P	P		P	P	P
	Other (Specify)							

N = new indication; P = previously cleared by FDA; E = added under this appendix PDI= Power Doppler Imaging

Additional comments: Combined mode: B+M, B+PW, B+Color, B+PDI/DPDI, B+Color+PW, B+PDI/DPDI +PW

Note \* Small Organ includes Thyroid, Testes, Breast (Cleared applications under K123249)

\*\* Other use includes Urology, Kidney, Gynecology

[1]: PDI: Power Doppler Imaging, DPDI: Directional Power Doppler Imaging

[2]: Biopsy Guidance

[3]: Harmonic Imaging, This feature does not use contrast agent.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH: Office of In Vitro Diagnostics and Radiological Health (OIR)

Prescription Use (Per 21 CFR 801.109)

**Diagnostic Ultrasound Indications for Use Form**

**DUS-5000 with E612UB**  
**Transducer**

General	Clinical Application	Mode of Operation						
		B	M	PW	CW	Color	Combined (Specify) <sup>[1]</sup>	Other (Specify) <sup>[2][3]</sup>
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal / Obstetrics							
	Abdominal							
	Intra-operative (Specify)							
	Intra-operative (Neuro logical)							
	Laparoscopic							
	Pediatric							
	Small Organ (Specify) *							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal	P	P	P		P	P	P
	Trans-urethral							
	Musculo-skeletal(Conventional)							
	Musculo-skeletal (Superficial)							
	Intravascular							
Other (Specify) **								
Cardiac	Adult Cardiac							
	Pediatric Cardiac							
	Intravascular(Cardiac)							
	Trans-esoph.(Cardiac)							
	Intra- cardiac							
Peripheral vascular	Peripheral vascular							
	Other (Specify)							

N = new indication; P = previously cleared by FDA; E = added under this appendix PDI= Power Doppler Imaging

Additional comments: Combined mode: B+M, B+PW, B+Color, B+PDI/DPDI, B+Color+PW, B+PDI/DPDI +PW

Note \* Small Organ includes Thyroid, Testes, Breast |Cleared applications under K123249

\*\* Other use includes Urology, Kidney, Gynecology

[1]: PDI: Power Doppler Imaging, DPDI: Directional Power Doppler Imaging

[2]: Biopsy Guidance

[3]: Harmonic Imaging, This feature does not use contrast agent.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH; Office of In Vitro Diagnostics and Radiological Health (OIR)

Prescription Use (Per 21 CFR 801.109)

**Diagnostic Ultrasound Indications for Use Form**

**DUS-5000 with C612UB**  
**Transducer**

General	Clinical Application	Mode of Operation						
		B	M	PW	CW	Color	Combined (Specify) <sup>[1]</sup>	Other (Specify) <sup>[2][3]</sup>
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal / Obstetrics							
	Abdominal							
	Intra-operative (Specify)							
	Intra-operative (Neuro logical)							
	Laparoscopic							
	Pediatric	P	P	P		P	P	P
	Small Organ (Specify) *							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Musculo-skeletal(Conventional)							
	Musculo-skeletal (Superficial)							
	Intravascular							
Other (Specify) **								
Cardiac	Adult Cardiac							
	Pediatric Cardiac	P	P	P		P	P	P
	Intravascular(Cardiac)							
	Trans-esoph.(Cardiac)							
	Intra- cardiac							
Peripheral vascular	Peripheral vascular							
	Other (Specify)							

N = new indication; P = previously cleared by FDA; E = added under this appendix PDI= Power Doppler Imaging

Additional comments: Combined mode: B+M, B+PW, B+Color, B+PDI/DPDI, B+Color+PW, B+PDI/DPDI +PW

Note \* Small Organ includes Thyroid, Testes, Breast. (Cleared applications under K123249)

\*\* Other use includes Urology, Kidney, Gynecology

[1]: PDI: Power Doppler Imaging ,DPDI: Directional Power Doppler Imaging

[2]: Biopsy Guidance

[3]: Harmonic Imaging, This feature does not use contrast agent.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH: Office of In Vitro Diagnostics and Radiological Health (OIR)

Prescription Use (Per 21 CFR 801.109)

## Diagnostic Ultrasound Indications for Use Form

### DUS-5000 with C6152UB Transducer

General	Clinical Application	Mode of Operation						
		B	M	PW	CW	Color	Combined (Specify) <sup>[1]</sup>	Other (Specify) <sup>[2][3]</sup>
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal / Obstetrics							
	Abdominal							
	Intra-operative (Specify)							
	Intra-operative (Neuro logical)							
	Laparoscopic							
	Pediatric	P	P	P		P	P	P
	Small Organ (Specify) *							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Musculo-skeletal(Conventional)							
	Musculo-skeletal (Superficial)							
	Intravascular							
	Other (Specify) **							
Cardiac	Adult Cardiac							
	Pediatric Cardiac	P	P	P		P	P	P
	Intravascular(Cardiac)							
	Trans-esoph.(Cardiac)							
	Intra- cardiac							
Peripheral vascular	Peripheral vascular							
	Other (Specify)							

N = new indication; P = previously cleared by FDA; E = added under this appendix PDI= Power Doppler Imaging

Additional comments: Combined mode: B+M, B+PW, B+Color, B+PDI/DPDI, B+Color+PW, B+PDI/DPDI +PW

Note \* Small Organ includes Thyroid, Testes, Breast |Cleared applications under K123249

\*\* Other use includes Urology, Kidney, Gynecology

[1]: PDI: Power Doppler Imaging ,DPDI: Directional Power Doppler Imaging

[2]: Biopsy Guidance

[3]: Harmonic Imaging. This feature does not use contrast agent.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH; Office of In Vitro Diagnostics and Radiological Health (OIR)

Prescription Use (Per 21 CFR 801.109)

## Diagnostic Ultrasound Indications for Use Form

DUS-5000 with C422UB

Transducer

General	Clinical Application	Mode of Operation						
		B	M	PW	CW	Color	Combined (Specify) <sup>[1]</sup>	Other (Specify) <sup>[2][3]</sup>
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal / Obstetrics							
	Abdominal	P	P	P		P	P	P
	Intra-operative (Specify)							
	Intra-operative (Neuro logical)							
	Laparoscopic							
	Pediatric							
	Small Organ (Specify) *							
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Musculo-skeletal(Conventional)							
	Musculo-skeletal (Superficial)							
	Intravascular							
Other (Specify) **								
Cardiac	Adult Cardiac	P	P	P		P	P	P
	Pediatric Cardiac							
	Intravascular(Cardiac)							
	Trans-esoph.(Cardiac)							
	Intra- cardiac							
Peripheral vascular	Peripheral vascular							
	Other (Specify)							

N = new indication; P = previously cleared by FDA; E = added under this appendix PDI= Power Doppler Imaging

Additional comments: Combined mode: B+M, B+PW, B+Color, B+PDI/DPDI, B+Color+PW, B+PDI/DPDI +PW

Note \* Small Organ includes Thyroid, Testes, Breast (Cleared applications under K123249)

\*\* Other use includes Urology, Kidney, Gynecology

[1]: PDI: Power Doppler Imaging, DPDI: Directional Power Doppler Imaging

[2]: Biopsy Guidance

[3]: Harmonic Imaging, This feature does not use contrast agent.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH; Office of In Vitro Diagnostics and Radiological Health (OIR)

Prescription Use (Per 21 CFR 801.109)

## Diagnostic Ultrasound Indications for Use Form

### DUS-5000 with L552UB Transducer

General	Clinical Application	Mode of Operation						
		B	M	PW	CW	Color	Combined (Specify) <sup>[1]</sup>	Other (Specify) <sup>[2][3]</sup>
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal / Obstetrics							
	Abdominal							
	Intra-operative (Specify)							
	Intra-operative (Neuro logical)							
	Laparoscopic							
	Pediatric	P	P	P		P	P	P
	Small Organ (Specify) *	P	P	P		P	P	P
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Musculo-skeletal (Conventional)	P	P	P		P	P	P
	Musculo-skeletal (Superficial)	P	P	P		P	P	P
	Intravascular							
Other (Specify) **								
Cardiac	Adult Cardiac							
	Pediatric Cardiac							
	Intravascular (Cardiac)							
	Trans-esoph. (Cardiac)							
	Intra- cardiac							
Peripheral vascular	Peripheral vascular	P	P	P		P	P	P
	Other (Specify)							

N = new indication; P = previously cleared by FDA; E = added under this appendix PDI= Power Doppler Imaging

Additional comments: Combined mode: B+M, B+PW, B+Color, B+PDI/DPDI, B+Color+PW, B+PDI/DPDI +PW

Note \* Small Organ includes Thyroid, Testes, Breast |Cleared applications under K123249

\*\* Other use includes Urology, Kidney, Gynecology

[1]: PDI: Power Doppler Imaging, DPDI: Directional Power Doppler Imaging

[2]: Biopsy Guidance

[3]: Harmonic Imaging, This feature does not use contrast agent.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH; Office of In Vitro Diagnostics and Radiological Health (OIR)

Prescription Use (Per 21 CFR 801.109)

### Diagnostic Ultrasound Indications for Use Form

#### DUS-5000 with L552UB Transducer

General	Clinical Application	Mode of Operation						
		B	M	PW	CW	Color	Combined (Specify) <sup>[1]</sup>	Other (Specify) <sup>[2][3]</sup>
Ophthalmic	Ophthalmic							
Fetal Imaging & Other	Fetal / Obstetrics							
	Abdominal							
	Intra-operative (Specify)							
	Intra-operative (Neuro logical)							
	Laparoscopic							
	Pediatric	P	P	P		P	P	P
	Small Organ (Specify) *	P	P	P		P	P	P
	Neonatal Cephalic							
	Adult Cephalic							
	Trans-rectal							
	Trans-vaginal							
	Trans-urethral							
	Musculo-skeletal(Conventional)	P	P	P		P	P	P
	Musculo-skeletal (Superficial)	P	P	P		P	P	P
	Intravascular							
Other (Specify) **								
Cardiac	Adult Cardiac							
	Pediatric Cardiac							
	Intravascular(Cardiac)							
	Trans-esoph.(Cardiac)							
	Intra- cardiac							
Peripheral vascular	Peripheral vascular	P	P	P		P	P	P
	Other (Specify)							

N = new indication; P = previously cleared by FDA; E = added under this appendix PDI= Power Doppler Imaging

Additional comments: Combined mode: B+M, B+PW, B+Color, B+PDI/DPDI, B+Color+PW, B+PDI/DPDI +PW

Note \* Small Organ includes Thyroid, Testes, Breast |Cleared applications under K123249

\*\* Other use includes Urology, Kidney, Gynecology

[1]: PDI: Power Doppler Imaging ,DPDI: Directional Power Doppler Imaging

[2]: Biopsy Guidance

[3]: Harmonic Imaging, This feature does not use contrast agent.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRI: Office of In Vitro Diagnostics and Radiological Health (OIR)

Prescription Use (Per 21 CFR 801.109)