

PREMARKET NOTIFICATION 510(k) SUMMARY
As required by 21 CFR §807.92(c)

Submitter

DEC 12 2013

510(k) Owner: MedApps, Inc., DBA Alere Connect
Owner / Operator: 10027842
Registration: 3005916763
Address: 8767 E. Via De Ventura, Suite 300, Scottsdale, AZ 85258
Telephone: 480-305-6323
Fax Number: 480-305-6328
Contact Person: Kent Dicks
Contact Person Title: CEO
Date Prepared: September 3, 2013

Device Information

Trade Name: MedApps 2.0 - Remote Patient Monitoring System
Common Name: Remote Patient Monitoring System
Classification Status: Class II per regulations 870.2910
Classification Name: Transmitters and Receivers, Physiological Signal,
Radiofrequency (21 CFR 870.2910, Product Code DRG)

A. LEGALLY MARKETED PREDICATE DEVICE

Legally marketed predicate devices are:

K040966 Carematix Modified System
K083862 MedApps 2.0 - Remote Patient Monitoring System
K124000 MedApps 2.0 - Remote Patient Monitoring System (with MobileLink added)

B. INDICATIONS FOR USE

The MedApps (Alere Connect) 2.0 - Remote Patient Monitoring System consists of 1) a cellular communication hub (MedApps' HealthPAL, Alere™ MobileLink or Alere™ HomeLink) an over-the-counter device that resides with the end-user (patient), which connects to commercially available FDA cleared accessory devices, specifically glucose meters, scales, blood pressure monitors, pulse oximeters, and PT/INR monitors (excluding HomeLink) and 2) web-based health data management application (MedApps' HealthCOM), that provides access to collected data stored on a secure host server system.

MedApps / Alere Connect Remote Patient Monitoring devices receive and store measurements collected from the described accessory devices, either wirelessly using short-range radio protocols (e.g. Bluetooth, Zigbee, WiFi, Bluetooth Low Energy (BLE), Fitlinxx Radios) or tethered via cable (e.g. USB, serial, etc.). Regardless of connectivity mode, the MedApps / Alere Connect monitoring devices do not alter the indications for use of the described peripheral accessory health devices.

MedApps / Alere Connect devices indicate successful or failed data reception and transmission with visual and audio feedback using a combination of any of the

following: Display, LED Lights, verbal messages, and / or audio tones / chimes. MedApps / Alere Connect devices store collected data and forward / transmit to server for access in HealthCOM via commercially available, FCC compliant, wireless telecommunication protocols (including but not limited to cellular GSM, CDMA and WiMax).

Healthcare professionals, clinicians and other authorized personnel can review the transmitted information within the MedApps HealthCOM system, where they can review collected readings, establish parameters to indicate readings exceptions to set thresholds, or trigger Interactive Voice Response (IVR) messages to the patient remotely to issue information such as reminders (e.g. "We haven't received readings from you today, please take and send your readings") or possibly educational information for conditions such as diabetes, hypertension, CHF, etc. Additionally, HealthCOM can port collected data to the healthcare providers' clinical back-end system(s) of choice.

The MedApps /Alere Connect 2.0 Remote Patient Monitoring System is not intended for diagnosis or as a substitute for medical care, nor is it intended to provide real-time / time-critical data. The device is contraindicated for patients requiring direct medical supervision or emergency intervention.

C. MedApps 2.0 SYSTEM DESCRITPION

The MedApps 2.0 - Remote Patient Monitoring System consists of:

(1) HealthPAL hardware:

The physical component of the HealthPAL is an electronic device contained in a plastic enclosure with an OLED screen, built-in M2M cellular chip, speaker, smart cable connection, smart cables, wireless module, LED lights to indicate activity, timer button to assist patients with their reading schedule (i.e. remind them to take their reading in X minutes), last reading button, volume up and down buttons. The HealthPAL Model 105 contains a GSM cellular module while the HealthPAL Model 106 contains a CDMA cellular module.

(2) HealthHUB hardware / software:

The HealthHUB hardware is an extension of the HealthPAL functionality. HealthHUB acts as a "docking" station for the HealthPAL in order to act as a conduit for the AC power adaptor connecting the electrical wall outlet to the HealthPAL providing power and battery charging capability. The Hubs also provide additional connections to off the shelf Glucose Meters, Scales, Blood Pressure Monitors and Pulse Oximeters, via smart cables (per validated in HealthPAL software). The HealthHUB model MA200 allows for multiple wired connections for accessory devices. HealthHUB Model 205 is specific for the HealthPAL MA105, and the HealthHUB 206 is specific for the HealthPAL MA 106 with both Hubs having one wired connector.

(3) HealthPAL firmware / software:

The firmware captures data from commercially available health monitors, and stores and transmits the information to the HealthCOM server, via the embedded communication chip / platform.

The firmware allows HealthPAL to receive information via wire or via embedded wireless module from accessory medical devices that are compatibly wireless enabled, which have been paired to the MedApps HealthPAL.

The firmware has many additional functions including:

- Download of user profiles from the server to configure HealthPAL remotely.
- HealthPAL has audio capability to deliver verbal announcement of readings and acknowledgment of data transmission from all connected accessory medical devices, time settings, volume control, educational content and reminders, in any language that is loaded to the device.
- Timer capability, activated by the user to provide assistance with adhering to a reading schedule (reminders to take readings within a set timeframe).
- OLED screen displays information regarding the HealthPAL's status including battery level, volume level, data transmission status, transmission pending indicator, activity icons / messages and other information to provide ease of use and promote patient adherence; as well as information received from accessory medical devices, such as the type of device, measurement, date and time of the last reading collected.
- Battery charging, isolation circuits, and interfaces to individual accessory medical devices / protocols via the smart cable.

(4) MobileLink hardware / software:

AC020 MobileLink is a modified MA105 HealthPAL device. The physical component of the MobileLink is an electronic device contained in a plastic enclosure with built-in M2M cellular chip, speaker, standard USB cable and USB Smart Cable connection, OLED screen to review the reading, and LED lights to indicate activity regarding the receiving and transmitting of collected data.

Like the HealthPAL, MobileLink's firmware / software captures, data from commercially available retail health monitors, and stores and transmits information to the HealthCOM server, via the embedded communication chip / platform.

The firmware allows MobileLink to receive information via wire, either standard USB or with a customized USB Smart Cable, from accessory medical devices.

The firmware has many additional functions including:

- Download of user profiles from the server to configure MobileLink remotely.
- MobileLink's Audio feature uses audio tones to indicate acknowledgment of collected readings from all connected accessory medical devices as well as reading transmission via the cellular network.
- MobileLink's visual user interface utilizes an OLED display to display collected readings from attached accessory medical devices as well as reading transmission acknowledgements. MobileLink's visual interface also contains a LED light to show power and provide reading request indication capability.

(5) HomeLink hardware / software:

AC200 HomeLink is a modified MA105 HealthPAL device. The physical component of the HomeLink is an electronic device contained in a plastic enclosure with built-in M2M cellular chip, speaker, standard USB cable and USB Smart Cable connection, touch screen to review the reading and allow patients with ability to have "IVR"

(questions and answers) functionality that was formally part of the HealthCOM IVR section, and wireless module.

Like the HealthPAL or MobileLink, HomeLink's firmware / software captures biometric data from FDA cleared commercially available accessory devices, and stores and transmits information to the HealthCOM server, via the embedded communication chip / platform.

The firmware allows HomeLink to receive information via wire or via embedded wireless module from accessory medical devices that are compatibly wireless enabled, which have been paired to the HomeLink device.

The firmware has many additional functions including:

- Download of user profiles from the server to configure HomeLink remotely.
- HomeLink's has audio capability to deliver verbal announcement of readings and acknowledgment of data transmission from all connected accessory medical devices, volume control, user feedback, in any language that is loaded to the device.
- HomeLink's visual user interface utilizes a touch screen display to display collected readings from attached accessory medical devices as well as other user interactions. One such patient HomeLink user interface includes the MedApps 2.0 System IVR functionality, which originally was executed from the HealthCOM webportal by professionals to patient's phones. The updated version allows end users (patients) to provide answers on the HomeLink device based on professional predefined questions that are reviewed in HealthCOM. This allows an improved user interaction of the previous IVR procedure. HomeLink screen also displays information regarding the device's status including volume level, data transmission status, preferred language, transmission pending indicator, activity icons / messages and other information to provide ease of use and promote patient adherence; as well as information received from accessory medical devices, such as the type of device, measurement, date and time of the last reading collected.

(6) MedApps HealthCOM software application:

The HealthCOM software application allows caregivers access to review patient data collected from accessory medical devices using MedApps hardware on the secure HealthCOM website. HealthCOM software allows professional caregivers to set patient readings.

HealthCOM software also allows the patient to establish an account and to direct / authorize their data to be directed to an outside, validated Personal Health Record (PHR), Electronic Health Record or Medical Record (EHR or EMR).

(7) MedApps IVR software application from HealthCOM:

The IVR (Interactive Voice Response) software application provides the ability to contact the patient remotely, by phone (designated in the user profile), and executes an pre-approved ("canned") scripts to deliver pre-approved ("canned") reminder messages ("Your nurse would like to talk to you, can I connect you now", "We haven't received a reading from you today, please send us your readings"), educational content and gather survey information.

In addition, the MedApps IVR application will send out Email, SMS / Text Messages, Paging, IM and other forms of communications in order to contact patients or caregivers. This will include reminders and alerts, based on clinically defined parameters / thresholds established in HealthCOM by the professional care provider.

The original IVR system has moved to the HomeLink device for easier interaction with the patient compared to the need for the patient to receive an IVR phone call on their personal phone, which includes professionally defined questions that the patient can answer and both the questions and answers will be recorded in the HealthCOM system for professional review.

(8) Accessory Device Descriptions:

OEM Accessory Device	Description	Integrated OEM FDA cleared accessory devices			FDA Clearance number
		Manufacturer	Model	Alere Connect Hub Device	
Blood Pressure Monitor	Off the shelf device that measures blood pressure, pulse of an individual	A&D Medical	UA-767PBT	HealthPAL & HomeLink	K012472
		A&D Medical	UA-767PC	MobileLink	
Scale	Off the shelf device that measures weight of an individual	A&D Medical	UC-321PBT	HealthPAL & HomeLink	Class 1 no 510k number
			UC-321PL	MobileLink	
Glucose Meters	Off the shelf devices that measure glucose level of blood for an individual.	LifeScan (J&J) OneTouch	Ultra and Ultra2	HealthPAL, MobileLink & HomeLink	K053529
		Abbott	Freestyle Lite		K070850
		Abbott	Freestyle Freedom Lite		K051839
		Bayer	Contour		K091820
		Nipro	TrueBalance		K090495
		Nipro	TrueResult		K080641
Pulse Oximeter	Prescription Use Device that measures oxygen concentration of an individual	Nonin	9560	HealthPAL & HomeLink	K081285
		Oxypulse	OxiPro-6	MobileLink	K100815
PT/INR	Prescription Use Device that measures prothrombin time of whole blood for monitoring of oral anticoagulation therapy.	Alere (formally Hemosense)	InRatio2	HealthPAL & MobileLink	K021923

D. TECHNOLOGICAL CHARACTERISTICS SUMMARY – as required by 807.92(a)(6)

Feature	Carematix Modified System K040966	MedApps Submission K083862	MedApps 2.0 Submission K124000	MedApps 2.0 Submission (This submission)
Indications of Use	Physiological monitoring system that collects, accumulates and transmits patient vital signs and other data from a patient who may be remote from the healthcare practitioner to the practitioner.	The MedApps 2.0 - Remote Patient Monitoring System consists of a patient device, MedApps HealthPAL and HealthCOM software monitoring application.	Same as MedApps 2.0 with the exception of updates to include MobileLink device and PT/INR monitor as an accessory device	Same as MedApps 2.0 System with the addition of the HomeLink device (similar to HealthPAL and MobileLink Hub devices)
Intended Use	Telemedicine System	Telemedicine System	Same as Carematix and MedApps 2.0	Same as Carematix and MedApps 2.0
Intended Users	Home users and Healthcare providers	Same	Same	Same
Site of Use	Home & Clinic	Same Home (HealthPAL) Clinic (HealthCOM)	Same Home (HealthPAL/ MobileLink); Clinic (HealthCOM)	Same Home (HealthPAL/ MobileLink/ HomeLink); Clinic (HealthCOM)
Data Collection Software	Proprietary Software	MedApps Proprietary Software	MedApps Proprietary Software	MedApps Proprietary Software
Data Collection Software Functionality	Transmit data from Sensor devices to Central Database	Same	Same	Same
Communication method of hub with Central Server	Via modem over telephone line	Via Embedded Cellular Technology	Via Embedded Cellular Technology	Via Embedded Cellular Technology
Types of sensors which can be interfaced (wired or wirelessly) to receiver hub	Medical Devices designed for Home use: Glucose, Scale Blood Pressure Pulse Ox, FEX/PEF,PT/INR Temperature	Medical Devices designed for Home use: Glucose, Scale Blood Pressure Pulse Ox	Same as MedApps (Glucose, Scale, BP, Pulse Ox with PT/INR) Same as CareMatix (PT/INR)	Same as MedApps (Glucose, Scale, BP, Pulse Ox with PT/INR) HomeLink does not interface with PT/INR
Maximum number and type of measurement devices that can be connected to the devices	Determined by vital sign devices that are designed for Home use, and have a data port. (Wireless or Wired)	Same	Same	Same
Maximum data throughput under worst case conditions	Multiple readings are stored on the medical devices and act as a backup if data needs to be re-sent to the server	Same	Same	Same

Feature	Carematix Modified System K040966	MedApps 2.0 Submission K083862	MedApps 2.0 Submission K124000	MedApps 2.0 Submission (This Submission)
Time Delay in the processing of data collected and transmitted	Readings stored in the medical devices can be sent up to the server when the connection is restored.	Same	Same	Same
Implementation method of collecting data from sensors	Modify OTS sensors with previous 510k approval by adding communications interface without altering sensors	Short range radio system using Wireless (Bluetooth) and Wired (tethered) cables.	Short range radio system using Wired (tethered)	Short range radio system using Wireless (Bluetooth) and Wired (tethered) cables.
Sensor Software	Sensor Software unchanged	Same	Same	Same
Connectivity	Wired or wireless to hub	Short range radio system using Bluetooth and Wired (tethered) cables.	Short range radio system using Wired (tethered)	Short range radio system using Bluetooth and Wired (tethered) cables.
Communication method of hub with devices	Wireless RF protocol	Short range radio system using Wireless (Bluetooth) and Wired (tethered) cables.	Short range radio system using Wired (tethered)	Short range radio system using Wireless (Bluetooth) and Wired (tethered) cables.
Communication s Protocol	Proprietary	Wireless (Bluetooth) V2.0 and Wired (Tethered)	Wired (Tethered)	Wireless (Bluetooth) V2.0 and Wired (Tethered)
Communication Frequency	915 MHz FCC assigned channel	(HealthPAL) Bluetooth : 2.402 to 2.480 GHz GSM: 850 / 900 / 1800 / 1950 Mhz	(HealthPAL or MobileLink) GSM: 850 / 900 / 1800 / 1950 Mhz	(HomeLink) Bluetooth : 2.402 to 2.480 GHz GSM: 850 / 900 / 1800 / 1950 Mhz
Power Source	Wall power plug (120 VAC/50-60) and Batteries in Device	Wall power plug (120 VAC/50-60) or Rechargeable Battery (HealthPAL)	Wall power plug (120 VAC/50-60) or Rechargeable Battery (HealthPAL or MobileLink)	Wall power plug (120 VAC/50-60) (HomeLink)
Display	On devices and hub, and monitors connected to central server	Same	Same	Same
Communication with Patients	On screen display	On screen display of Readings, Voice Output and Interactive Voice Response (IVR)	On screen display with audio tones instead of voice.	On screen display of Readings, Voice Output and Interactive Voice Response (IVR) through HomeLink
Use of Thresholds / Algorithms for determining how Thresholds are set and changed	Thresholds are set by Healthcare professionals in Server Software	Same	Same	Same

Information presented to the user, if it is different from that presented by the measurement devices	On screen display	Audio/visual reading feedback on screen and by speaker; and Interactive Voice Response (IVR) System for patient contact	Visual reading feedback on screen and audio tone by speaker; and Interactive Voice Response (IVR) System for patient contact	Audio/visual reading feedback on screen and by speaker; and Interactive Voice Response (IVR) System for patient contact is through HomeLink device
Messages and Instructions that can be sent to the User.	On screen display	On screen display of Readings, Voice Output and Interactive Voice Response (IVR)	On screen display of Readings, Voice Output and Interactive Voice Response (IVR)	On screen display of Readings, Voice Output and Interactive Voice Response (IVR) functionality uses HomeLink device

Below is a Technological Characteristics Summary Comparison between MA105 HealthPAL, AC020 MobileLink, and AC200 HomeLink MedApps 2.0 System devices:

Feature	MA105 HealthPAL	AC020 MobileLink	AC200 HomeLink
Indications of Use	Enables healthcare providers to monitor and manage biometric patient data collected remotely	Same	Same
Intended Use	Telemedicine System	Same	Same
Intended Users	Home users and patients outside of the clinical setting, as well as Healthcare providers for HealthCOM	Same	Same
Site of Use	Remote setting (e.g. Home / Work), Clinic	Same	Same
Data Collection Software (firmware)	MedApps Proprietary Software	Same	Same
Data Collection Software Functionality	Transmit data from Sensor devices to Central Database	Same	Same
Communication method of device hub with Central Server	Via Embedded Cellular Technology (GSM or CDMA)	Same	Same
Types of sensors which can be interfaced (wired or wirelessly) to receiver hub	Medical Devices designed for Home use: Glucose, Scale, Blood Pressure, Pulse Ox, PT/INR	Medical Devices designed for Home use: Glucose, Scale, Blood Pressure, Pulse Ox, PT/INR	Medical Devices designed for Home use: Glucose, Scale, Blood Pressure, Pulse Ox
Transmission	Transmits information to the MedApps secure host server called "HealthCOM"	Same	Same
Implementation method of collecting data from sensors and general Connectivity	Short range radio system using Wireless (Bluetooth) and Wired / tethered (Smart Cables).	MobileLink uses wired / tethered connection (USB, Smart Cables)	Short range radio system using Wireless (Bluetooth) and Wired / tethered (Smart Cables).
Communication method of hub with devices	Short range radio system using Wireless (Bluetooth) and Wired / tethered (Smart Cables).	MobileLink uses wired / tethered connection (USB, Smart Cables)	Short range radio system using Wireless (Bluetooth) and Wired / tethered (Smart Cables).

Communication Frequency	Bluetooth : 2.402 to 2.480 GHz GSM: 850 / 900 / 1800 / 1900 Mhz	No Bluetooth capability GSM: 850 / 900 / 1800 / 1900 Mhz	Bluetooth : 2.402 to 2.480 GHz GSM: 850 / 900 / 1800 / 1900 Mhz
Power Source	AC adaptor Wall power plug (120 VAC/50-60) and Rechargeable Batteries in Device	AC Adaptor that is 60601-1 3 rd Edition compatible with Lithium battery only used for soft shut down functionality and not powering device.	AC Adaptor that is 60601-1 3 rd Edition compatible for powering device.
Device Communication with Patients	On screen display and audio voice feedback	On screen display with audio tones instead of voice.	Touch screen display with audio voice feedback.
Applicable Certification Testing	Safety 60601-1, EMC/EMI/FCC (60601-1-2), ESD & Radiated Immunity, FCC Bluetooth, (PTCRB), CTIA (battery), ETSI	Safety 60601-1 3 rd Edition, EMC/EMI/FCC (60601-1-2), ESD & Radiated Immunity, PTCRB, ETSI (Declaration of Conformity)	Safety 60601-1 3 rd Edition, EMC/EMI/FCC (60601-1-2), ESD & Radiated Immunity, PTCRB ETSI (See Declaration of Conformity)

Technological Characteristics Summary Comparison between MA105 HealthPAL, AC020 MobileLink, and AC200 HomeLink devices as part of the MedApps 2.0 System focusing on the HomeLink modification to other two devices:

Feature	MA105 HealthPAL	AC020 MobileLink	AC200 HomeLink	Modification Description
Device Dimensions	4.2"x 2.1" x 0.9", weighs 3 oz.	4.5" x 1.9" x 0.75", weighs 2.8 oz.	5.1" x 7.9" x 0.7", weighs 17 oz.	Modified to increase the size of the device
User Interface	1.8" Full Color Passive Matrix OLED and buttons	1.54" single color OLED and buttons	7" Full Color Touch Screen	Modified to a touch screen instead of buttons
Voice Interface	Voice and tones	Tones only	Voice and tones	Same as HealthPAL
Microprocessors	C8051F344 (Silicon Labs); 8051 architecture	PIC32MX975F512L (Microchip); MIPS architecture	PIC32MX975F512L (Microchip); MIPS architecture	Same as MobileLink
Memory Devices	4GB MicroSD Card; SDHC	4GB or 8GB MicroSD Card; SDHC	4GB or 8GB MicroSD; SDHC	Same as MobileLink
Energy Sources	Main Power Supply: SHWA12A-S022(MEPOS)6061-1 3 rd Ed.[external] Switching Power Supplies: TPS61080,BQ24014, TPS63001(TI); LTC3525(Linear)[internal]	Main Power Supply: 18ur-09-2000(Group West) Power Supply components: BQ24014 (TI), TPS61080(TI); LTC3122(Liner); A4490(Allegro) [Internal]	Main Power Supply: 18ur-09-2000 (Group West) Power Supply components: BQ24014(TI); LTC2427,LT1617, LT1930(Linear); XC9119(Torex); A4490(Allegro)[Internal]	Power Supply same as MobileLink; More internal power supply components compared to other two devices
Safety Features	ESD Protection: TPD6E001 (TI)	ESD Protection: TPD6E001 (TI)	ESD Protection: TPD6E001 (TI)	Same as both devices
Communications to Server	Cell Module: GSM0308-10(Enfora);	Cell Module: SL8081 (Sierra Wireless);	Cell Module: SL8081 (Sierra Wireless);	Same as MobileLink

Communications to Accessory Devices	Bluetooth: WT11-A-Ai4(BlueGiga); Serial Communication through Smart Cable (30 pin)	Serial Communication through Smart Cable (USB)	Bluetooth: WT11-A-Ai4(BlueGiga); BR-LE4.0-S2A(Blue Radios); FitLinxx; Serial Communication through Smart Cable (USB)	Bluetooth same as HealthPAL; Added Bluetooth low energy and Fitlinxx
Programming Language	C Programming Language	C Programming Language	C Programming Language	Same as HealthPAL and MobileLink
Firmware Architecture	Single threaded	Multi-threaded (RTOS)	Multi-threaded (RTOS)	Same as MobileLink
IVR functionality communication	Professional use HealthCOM; End Users use their phone	Professional use HealthCOM; End Users use their phone	Professional use HealthCOM; End Users use HomeLink	HomeLink is used for IVR functionality instead of the End User's phone

Data Collection:

The one predicate device and the MedApps 2.0 System solution connect to medical devices (designed for home use) by either wired (cable) connection or wireless (HealthPAL/HomeLink- Bluetooth). The data is collected from the devices and sent to a secure central server using various communication methods.

Telecommunication Platform to Central Server:

Carematrix uses modem off telephone line; MedApps 2.0 System devices uses embedded Machine to Machine (M2M) module to transmit data via cellular connectivity.

Patient Feedback Technology:

The predicate device and the MedApps solution allow data and messages to be displayed on a screen (for the HealthPAL/HomeLink) for the patient to read and acknowledge. MobileLink also uses audio and visual acknowledgement / feedback. The MedApps 2.0 System solution also uses an Interactive Voice Response (IVR) system in order to communicate with the patient and ask questions, gather survey information, or issue reminders. The HomeLink includes the user interface (IVR process) as part of the user interaction instead of using a patient's phone.

Backend Data Storage:

All systems (predicate device and the MedApps 2.0 system solution), provide a backend system that allows data to be stored, and healthcare professionals to access and monitor collected patient data.

Types of sensors which can be interfaced (wired or wirelessly) to receiver hub:

CareMatix Modified System and the proposed MedApps 2.0 System both include PT/INR monitor as an FDA cleared accessory device to the receiver hub. The CareMatix uses wireless connectivity to the hub whereas the MedApps 2.0 System uses both wired and wireless connectivity to the hub. CareMatix System 510k (K040966) clearance summary letter includes the predicate Avid Care (K011779 and K010029) telemedicine system that connects to PT/INR monitor as an accessory device. The HomeLink includes all of the original MedApps 2.0 System (K083862) accessory devices without interfacing with the PT/INR monitor at this time.

E. NON-CLINICAL PERFORMANCE DATA TESTING AND REVIEW – as required by 807.92(b)(1)

Non-Clinical Testing

The submitted 2.0 System has undergone MedApps' design control verification and validation testing. MedApps 2.0 validation testing include testing of all executable code and functionality and confirmation that all identified risks have been adequately addressed by software functionality, the user interface, documentation or user SOP.

MedApps 2.0 System verification and validation activities as part of the design control process include testing of all Design Specifications (Design Control Inputs) based on risk analysis, certification standards, and Verification plans. MedApps Product Verification and Release Plan execution on HealthPAL, MobileLink, and HomeLink devices ensures each medical device works with each applicable type of user accessory medical device (glucose, blood pressure monitor, scale, pulse oximeter and PT/INR (excluding HomeLink) as part of the MedApps 2.0 System including integration to HealthCOM backend software application. The output of these design control verification analysis documents **MedApps 2.0 - Remote Patient Monitoring System** shall meet its requirements and design specifications as intended. Additionally, the HomeLink question and answer functionality, which moved to the device instead of the patient's personal phone as part of the MedApps 2.0 System was included in the verification testing to ensure this user interface is verified.

Lastly, all relevant certification testing such as EMC (60601-1-2) and Safety (60601-1) are described in MedApps' Declaration of Conformity.

F. SUBSTANTIAL EQUIVALENT

The MedApps 2.0 Remote Patient Monitoring System is substantially equivalent to the predicate devices in terms of data collection software functionality, operating system for the patient device, communication method of patient device with central server, types of sensors which can be integrated to the patient medical device, implementation methods of collecting data from sensors, sensor software, connectivity, communication protocol, power source and general display method.

The HomeLink (this submission) is substantially equivalent to the HealthPAL (described in 510(k) K083862) and MobileLink (K124000) as all three devices, as part of the MedApps 2.0 System, connect to commercially available Glucose Meters, Scales, Blood Pressure Monitors, Pulse Oximeters, and PT/INR; data is collected, stored and transmitted using off-the-shelf FCC approved wireless / cellular connectivity. All devices provide audio and visual feedback / acknowledgement that readings have been collected and transmitted to MedApps' secure host server called "HealthCOM". The HomeLink (this submission) is substantially equivalent to the HealthPAL (described in 510(k) K083862) as both devices, as part of the MedApps 2.0 System, connect to commercially available Glucose Meters, Scales, Blood Pressure Monitors, Pulse Oximeters; data is collected, stored and transmitted using off-the-shelf FCC approved wireless / cellular connectivity. Both provide audio (voice) and visual feedback / acknowledgement that readings have been collected and transmitted to MedApps' secure host server called "HealthCOM". Both devices also use wireless connectivity (Bluetooth) for accessory devices that have this capability.

G. SAFETY AND EFFICACY

The MedApps 2.0 Remote Patient Monitoring System does not rely on an assessment of clinical performance data. The device conforms to FDA's recognized consensus standards and relies on its conformity to demonstrate its safety and efficacy. The device does not introduce any new questions concerning the safety or efficacy and is, therefore, substantially equivalent to the predicate devices.



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

December 12, 2013

Medapps Inc., Dba Alere Connect
Kent Dicks
8767 E. Via De Ventura, Ste 300
Scottsdale, AZ 85258 US

Re: K132803
Trade/Device Name: Medapps 2.0 Remote Patient Monitoring System, Healthpal,
Healthcom, Mobile Link
Regulation Number: 21 CFR 870.2910
Regulation Name: Medapps Remote Patient Monitoring System
Regulatory Class: Class II
Product Code: DRG, NBW
Dated: 11/12/2013
Received: 11/11/2013

Dear Kent Dicks:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Owen P. Faris -S

for

Bram D. Zuckerman, M.D.
Director
Division of Cardiovascular Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number: K132803

Preparation Date: September 3, 2013

Device Name: **MedApps (Alere Connect) 2.0 - Remote Patient Monitoring System**

Indications For Use:

The MedApps (Alere Connect) 2.0 - Remote Patient Monitoring System consists of 1) a **cellular communication hub** (MedApps' HealthPAL, Alere™ MobileLink or Alere™ HomeLink) an over-the-counter device that resides with the end-user (patient), which connects to commercially available FDA cleared accessory devices, specifically glucose meters, scales, blood pressure monitors, pulse oximeters, and PT/INR monitors (excluding HomeLink) and 2) web-based health data management application (MedApps' HealthCOM), that provides access to collected data stored on a secure host server system.

MedApps / Alere Connect Remote Patient Monitoring devices receive and store measurements collected from the described accessory devices, either wirelessly using short-range radio protocols (e.g. Bluetooth, Zigbee, WiFi, Bluetooth Low Energy (BLE), Fitlinxx Radios) or tethered via cable (e.g. USB, serial, etc). Regardless of connectivity mode, the MedApps / Alere Connect monitoring devices do not alter the indications for use of the described peripheral accessory health devices.

MedApps / Alere Connect devices indicate successful or failed data reception and transmission with visual and audio feedback using a combination of any of the following: Display, LED Lights, verbal messages, and / or audio tones / chimes. MedApps / Alere Connect devices store collected data and forward / transmit to server for access in HealthCOM via commercially available, FCC compliant, wireless telecommunication protocols (including but not limited to cellular GSM, CDMA and WiMax).

Healthcare professionals, clinicians and other authorized personnel can review the transmitted information within the MedApps HealthCOM system, where they can review collected readings, establish parameters to indicate readings exceptions to set thresholds, or trigger Interactive Voice Response (IVR) messages to the patient remotely to issue information such as reminders (e.g. "We haven't received readings from you today, please take and send your readings") or possibly educational information for conditions such as diabetes, hypertension, CHF, etc. Additionally, HealthCOM can port collected data to the healthcare providers' clinical back-end system(s) of choice.

The MedApps /Alere Connect 2.0 Remote Patient Monitoring System is not intended for diagnosis or as a substitute for medical care, nor is it intended to provide real-time / time-critical data. The device is contraindicated for patients requiring direct medical supervision or emergency intervention.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use _____ OR Over-The-Counter Use X
(Per 21 CFR 801.109)

Page Exhibit01-1  Digitally signed by Owen P. Faris -S Date: 2013.12.12 14:51:39 -05'00'