

**NxStage Medical, Inc.
NxStage System One with NxView
510(k) Premarket Notification**

MAR - 3 2014

This 510(k) Summary is being submitted in accordance with the requirements of the Safe Medical Device Act (SMDA) of 1990. The content contained in this 510(k) summary has been provided in conformance with 21 CFR §807.92

A. Date Prepared: November 18, 2013

B. Submitter's Information:

Name: NxStage Medical, Inc.

Address: 350 Merrimack Street
Lawrence, MA 01843

**FDA Establishment
Owner/Operator
Number:** 9045797

Contact Person: Mary Lou Stroumbos
Director, Regulatory Affairs

Phone: (978) 687-4872
Fax: (978) 687-4750

Manufacturer: NxStage Medical, Inc.
350 Merrimack Street
Lawrence, MA 01843

**FDA Establishment
Registration Number:** 3003464075

Sterilization Site: Steris Isomedix (NxStage Cartridge
Express)
1000 S. Sarah Place
Ontario, CA 91761

**NxStage Medical, Inc.
NxStage System One with NxView
510(k) Premarket Notification**

C. Device Name:

Trade/Proprietary Name:	NxStage System One with NxView
Common/Usual Name:	Hemodialysis System
Classification Name:	High Permeability Hemodialysis System
Regulation Number:	876.5860
Product Code:	78 KDI
Device Classification:	Class II
Device Panel:	Gastroenterology/Urology

D. Substantial Equivalence:

The NxStage System One with NxView has the same intended use and utilizes the same fundamental technology as the predicate NxStage System One. The NxStage System One with NxView has been compared to the legally marketed predicate devices as cleared through K122051 (April 23, 2013) and K040074 (April 8, 2004) and was found to be substantially equivalent.

E. Device Description/Indications for Use:

The NxStage System One is comprised of the NxStage Cyclor, an electromechanical control unit; the NxStage Cartridge, a sterile, single-use extracorporeal blood and fluid management circuit (with or without a pre-attached high permeability filter) that mounts integrally within the NxStage Cyclor; and NxView, a flat panel touch screen interface which is mounted on top of the cyclor and provides online instructions for use, summarized system information, and remote viewing of treatment information. The combined system is designed to deliver hemofiltration, hemodialysis and/or ultrafiltration in an acute or chronic care facility. The NxStage System One with NxView is also indicated for Therapeutic Plasma Exchange in a clinical environment.

**NxStage Medical, Inc.
NxStage System One with NxView
510(k) Premarket Notification**

Indications for use:

NxStage System One:

The NxStage System One is indicated for the treatment of acute and chronic renal failure or fluid overload using hemofiltration, hemodialysis, and/or ultrafiltration, in an acute or chronic care facility.

The NxStage System One is also indicated for Therapeutic Plasma Exchange in a clinical environment.

All treatments must be administered under physician's prescription, and must be observed by a trained and qualified person, considered to be competent in the use of this device by the prescribing physician.

NxView:

NxView is a computer-based touch screen user interface that provides on-line instructions for use, summarized system information and remote access.

NxView is contraindicated as the sole method of monitoring a patient during treatment.

F. Technological Characteristics:

The proposed device has the same technological characteristics and is similar in design and configuration as compared to the predicate device. The proposed device is designed with similar software, components and features also used in the predicate device.

Table 1 Device Technological Characteristics Comparison Table		
Parameter	Proposed Device NxStage System One with NxView	Predicate Device NxStage System One (K122051)
<i>Intended Use</i>		
<i>Hemodialysis</i>	Yes	Yes
<i>Hemofiltration</i>	Yes	Yes
<i>Ultrafiltration</i>	Yes	Yes
<i>Technology / Components:</i>		
<i>Pumps</i>	Same	4 peristaltic pumps
<i>Valves (clamps)</i>	Same Same	2 solenoid actuated pinch clamps 8 cam driven pinch clamps

NxStage Medical, Inc.
NxStage System One with NxView
510(k) Premarket Notification

Table 1 Device Technological Characteristics Comparison Table										
Parameter	Proposed Device NxStage System One with NxView	Predicate Device NxStage System One (K122051)								
<i>Air / fluid detectors</i>	Same	3 ultrasonic air/ fluid detectors								
<i>Blood leak detector</i>	Same	1 optical blood leak detector								
<i>Pressure transducers</i>	Same	5 electronic pressure transducers								
<i>Temperature sensors</i>	Same	1 electronic temperature sensor								
<i>Flow Rates:</i>										
<i>Blood</i>	Same	10-600 ml/min								
<i>Prescription Fluid /Dialysate Flow</i>	0-12000 ml/hr (NX1000-5)	0 -12000 ml/hr (NX1000-1 & NX1000-2) 0-18000 ml/hr (NX1000-3)								
<i>Ultrafiltration</i>	Same	0-2400 ml/hr								
<i>Transmembrane Pressure Monitoring Specification</i>	Same	Yes								
<i>Venous Pressure Monitor</i>	Same	0 to 400 mmHg								
<i>Effluent fluid Pressure Monitor</i>	Same	0 to 500 mmHg								
<i>Air Detector</i>	Same	Reduction of detector signal lasting 6 ms minimum (Approximates a 60 micro liter bubble at 400 mmHg venous pressure and 600 ml/min blood flow)								
<i>Blood Leak Detector</i>	15% reduction in detector signal lasting 20 seconds minimum. Signal reduction % based on a 0.35 ml/min leak rate of 32 Hct blood.	15% reduction in detector signal lasting 20 seconds minimum. Signal reduction % based on a 0.45 ml/min leak rate of 32 Hct blood.								
<i>Effluent Volume Accuracy</i>	Same	Greater of 300 ml/ 12 hr or 3% of exchange volume (For software versions 4.7 and below) For software versions 4.8 and higher: <table border="1" style="margin-left: auto; margin-right: auto;"> <thead> <tr> <th style="text-align: center;">Therapy Fluid Flow Rate L/hr</th> <th colspan="2" style="text-align: center;">Specification greater of</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">> 3</td> <td rowspan="2" style="text-align: center;">+ 5% UF* or</td> <td style="text-align: center;">±100 ml/hr*</td> </tr> <tr> <td style="text-align: center;">≤ 3</td> <td style="text-align: center;">± 25 ml/hr*</td> </tr> </tbody> </table>	Therapy Fluid Flow Rate L/hr	Specification greater of		> 3	+ 5% UF* or	±100 ml/hr*	≤ 3	± 25 ml/hr*
Therapy Fluid Flow Rate L/hr	Specification greater of									
> 3	+ 5% UF* or	±100 ml/hr*								
≤ 3		± 25 ml/hr*								
<i>IV Prescription Fluid</i>	Same	Off-line, sterile- physician prescribed, indicated for infusion								

**NxStage Medical, Inc.
NxStage System One with NxView
510(k) Premarket Notification**

Table 1 Device Technological Characteristics Comparison Table		
Parameter	Proposed Device NxStage System One with NxView	Predicate Device NxStage System One (K122051)
<i>Dialysate</i>	Same	Dialysate available as pre-packaged pre-mixed sterile fluids or via the PureFlow SL (K043436 K060296, K090919 & K111174)
<i>Compatible Blood Tubing Set</i>	Same	NxStage Standard Cartridge
<i>Software</i>	Software version 4.9	Software version 4.8
<i>NxView</i>	NxView touch screen interface included as standard with NxStage System One	OneView (K040074) touch screen interface is an optional accessory

G. Summary of Non-Clinical Test/Performance Testing - Bench

NxStage believes that the information and data provided in this submission clearly describes the proposed device and demonstrates that the device is adequately designed for the labeled indications for use. Performance, verification and validation testing was conducted to characterize performance of the proposed device and the predetermined acceptance criteria was met. Results of this testing have documented that the proposed NxStage System One with NxView is substantially equivalent to the predicate devices and is suitable for the labeled indications for use.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

March 3, 2014

NxStage Medical, Inc.
Mary Lou Stroumbos
Director, Regulatory Affairs
350 Merrimack Street
Lawrence, MA 01843

Re: K133547
Trade/Device Name: NxStage[®] System One™ with NxView™
Regulation Number: 21 CFR§ 876.5860
Regulation Name: High permeability hemodialysis system
Regulatory Class: II
Product Code: KDI
Dated: January 24, 2014
Received: January 27, 2014

Dear Mary Lou Stroumbos,

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies.

You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOffices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Benjamin R. Fisher -S

Benjamin R. Fisher, Ph.D.
Director
Division of Reproductive, Gastro-Renal,
and Urological Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

INDICATIONS FOR USE

510(k) Number (if known): K133547

Device Name: NxStage[®] System One[™] with NxView[™]

Indications for Use: **NxStage System One:**
The NxStage System One is indicated for the treatment of acute and chronic renal failure, or fluid overload using hemofiltration, hemodialysis, and/or ultrafiltration, in an acute or chronic care facility.

The NxStage System One is also indicated for Therapeutic Plasma Exchange in a clinical environment.

All treatments must be administered under physician's prescription, and must be observed by a trained and qualified person, considered to be competent in the use of this device by the prescribing physician.

NxView:
NxView is a computer-based touch screen user interface that provides on-line instructions for use, summarized system information and remote access.

NxView is contraindicated as the sole method of monitoring a patient during treatment.

Prescription Use X
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use _____
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page 1 of 1

Benjamin R-Fisher -S
2014.03.02 18:23:28 -05'00'