

510(k) Summary

1. Submitter Identification

510(k) Submitter	GRANDWAY TECHNOLOGY (SHENZHEN)
Address	Block 6 and 7, Zhu Keng Industrial Zone District, Shenzhen, Guang Dong, People's Republic of China
Phone Number	(00852)-2851-6789
Fax Number	(00852)-2851-6278
Contact Person	Mr. Patrick Chow
Date of Submission	19 th November, 2013

2. Device Identification

Trade Name	Digital Automatic Wrist Blood Pressure Series [Model No.: MD22xy]
	x --- The first character (2 & 3) refers to the device type, i.e. wrist or cabinet (housing).
	y --- The second character (0, 1, 2, 3) refers to the device change revision of device. The number refers to those device changes that do not affect the safety or efficacy of the device. The number refers to conformity test results of EMC and performance, i.e. IEC60601-1 and IEC60601-2-3.
Common Name	Non-invasive Blood Pressure Measurement Device

4. Device Description

Digital Automatic Wrist Blood Pressure Monitor WBPM22 Series is a non-invasive blood pressure measurement system for use by medical professional or patient to measure the systolic and diastolic blood pressure, and pulse rate measurement and then display the readings on a digital panel.

The device utilizes the oscillometric methodology, in which an inflatable cuff is placed around the Wrist of an individual, for blood pressure measurement. The device detects your blood's movement through your brachial artery and converts it into a digital reading.

The table below illustrate the feature presence in Digital Automatic Wrist Blood Pressure Monitor WBPM22 Series.

Model	Blood Pressure Measurement	Pulse Rate Measurement	WHO Classification	Irregular Heartbeat
MD2220	✓	✓	✓	✓
MD2222	✓	✓	✓	✓
MD2223	✓	✓	✓	✓
MD2230	✓	✓	✓	✓
MD2231	✓	✓	✓	✓
MD2232	✓	✓	✓	✓

5. Indication for Use

6. Comparison of Technological Characteristics between Predicate Devices

Digital Automatic Wrist Blood Pressure Monitor WBPM22 Series is a predicate device, WBPM15 Series (K120673) in the device comparison.

Comparison between WBPM22 Series and predicate device		
Item	Predicate Device	WBPM22 Series
Indication for Use	Digital Automatic Wrist Blood Pressure Monitor WBPM15 Series is for use by medical professional or at home. The WBPM15 Series is intended to measure the systolic and diastolic blood pressure, and pulse rate (heartbeat rate) of an individual by using a non-invasive technique, in which an inflatable cuff is wrapped around the wrist of an individual. The inflatable cuff circumference is limited to 13.5cm - 19.5cm.	Digital Automatic Wrist Blood Pressure Monitor WBPM22 Series is for use by medical professional or at home. The WBPM22 Series is intended to measure the systolic and diastolic blood pressure, and pulse rate of an individual by using a non-invasive technique, in which an inflatable cuff is wrapped around the wrist.
Measurement Method	Non-invasive, Oscillometric	Non-invasive, Oscillometric
Irregular Heartbeat Detection (IHB)	Yes	Yes
Patient Population	Age 16 or above	Age 16 or above
Blood Pressure Measurement Range	Cuff Pressure: 0 - 300 mmHg Systolic Pressure: 50 - 250 mmHg Diastolic Pressure: 30 - 200 mmHg	Cuff Pressure: 0 - 300 mmHg Systolic Pressure: 50 - 250 mmHg Diastolic Pressure: 30 - 200 mmHg
Number of User	2 independent users	2 independent users
Memory Space	2 users × 120 memory space	2 users × 60 memory space 2 users × 120 memory space
Blood Pressure Measurement Accuracy	± 3 mmHg or 2% of reading	± 3 mmHg or 2% of reading
Pulse Rate	30 - 180 beats/min	30 - 180 beats/min

Comparison between WBPM22 Series and predicate device			
Item	Predicate Device	WBPM22 Series	Comment
Storage and Transportation Condition	Temperature: -20 - 60 °C Humidity: 10 - 95 % R.H. max Atmospheric Pressure: 700 - 1060 kPa	Temperature: -20 - 60 °C Humidity: 10 - 95 % R.H. max Atmospheric Pressure: 700 - 1060 kPa	Identical
Material	Resistances, capacitance, transistors, amplifiers, pressure sensor, CPU, PCB, cuff ABS button, ABS cabinet, batteries and packaging	Resistances, capacitance, transistors, amplifiers, pressure sensor, CPU, PCB, cuff ABS button, ABS cabinet, batteries and packaging	Identical
Compatibility with Environment and Other Devices	No influence with environment and other device	No influence with environment and other device	Identical
Applicable Standard	<ul style="list-style-type: none"> ◇ EN 1060-1:1995/A2:2009 ◇ EN 1060-3:1997/A2:2009 ◇ IEC 60601-1:2005+ CORR.1(2006)+CORR. 2 (2007) ◇ EN 60601-1-2:2007 ◇ FCC Part 15 ◇ ISO 10993-5:2009 ◇ ISO 10993-10:2002 + A1:2006 ◇ EN 60601-1-4:2007 ◇ ANSI/AAMI SP-10:2002 	<ul style="list-style-type: none"> ◇ EN 1060-1:1995+A2:2009 ◇ EN 1060-3:1997+A2:2009 ◇ IEC 60601-1:2012 ◇ EN 60601-1-2:2007 ◇ FCC Part 15 Subpart B ◇ ISO 10993-5:2009 ◇ ISO 10993-10:2010 ◇ IEC 62304:2006 ◇ IEC 81060-2:2009 	Equivalent, change upon FDA's recognized standard

Digital Automatic Wrist Blood Pressure Monitor WBPM22 Series is a non-invasive measuring device and utilizes the oscillometric methodology to measure the blood pressure reading. The key components of device are: a pressure sensor, a electric valve and an electronic control module together with an electric pump. The electric pump inflate (and deflate) the inflatable cuff automatically according to our designed architecture. The predicate device adopts exactly same methodology and key components for measuring blood pressure.

7. Clinical and Non-clinical Tests

Clinical Test Summary

Testing to insure clinical accuracy of the device in accordance with ISO 81060-2:2009 as documented in Clinical Test report.

One hundred patients (49 males and 51 females) were invited for the study. Standard auscultation method was used as the reference blood pressure monitor measuring in the left wrist. Blood pressure measurements were repeated alternatively with the device and auscultation in the same arm according to the sequence in ISO 81060-2:2009.

Non-Clinical Test Summary

Digital Automatic Wrist Blood Pressure Monitor WBPM22 Series underwent non-clinical tests to show that all requirement specifications and standards are met. The tests includes the follows:

- ❖ EN 1060-1:1995+A2:2009
- ❖ EN 1060-3:1997+A2:2009
- ❖ IEC 60601-1:2012
- ❖ EN 60601-1-2:2007
- ❖ FCC Part 15 Subpart B
- ❖ ISO 10993-5:2009
- ❖ ISO 10993-10:2010
- ❖ IEC 62304:2006

As all of the clinical and non-clinical testing performed on Digital Automatic Wrist Blood Pressure Monitor WBPM22 Series are same as the predicate device, a comparison testing is conducted to show the performance of Digital Automatic Wrist Blood Pressure Monitor WBPM22 Series is equivalent to the predicate device.

8. Conclusion

Digital Automatic Wrist Blood Pressure Monitor WBPM22 Series has the same design and same technological characteristics as the predicate device, Digital Automatic Wrist Blood Pressure Monitor WBPM15 Series (K120673). Moreover both devices have undergone testing has demonstrated that no differences in the technology or design are questioning on safety or effectiveness to be raised. Thus, Digital Automatic Wrist Blood Pressure Monitor WBPM22 Series is substantially equivalent to the predicate device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

March 5, 2014

Grandway Technology (Shenzhen) Limited
Patrick Chow
Zhu Keng Industrial Zone
Ping Shan, Long Gang District
Shenshen, Guang Dong, 518118 CH

Re: K133618

Trade/Device Name: Digital Automatic Wrist Blood Pressure M

Regulation Number: 21 CFR 870.1130

Regulation Name: Non-Invasive Blood Pressure Measurement S

Regulatory Class: Class II

Product Code: DXN

Dated: February 12, 2014

Received: November 25, 2013

Dear Patrick Chow,

We have reviewed your Section 510(k) premarket notification of intent to
referenced above and have determined the device is substantially equivalent
for use stated in the enclosure) to legally marketed predicate devices mar
commerce prior to May 28, 1976, the enactment date of the Medical De

Page 2 - Patrick Chow

found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In a future letter, we will publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination letter indicates that FDA has made a determination that your device complies with other applicable Federal or any Federal statutes and regulations administered by other Federal agencies, in addition to those with all the Act's requirements, including, but not limited to: registration (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting and record-keeping for device-related adverse events) (21 CFR 803); good manufacturing practice requirements for medical devices (the quality systems (QS) regulation (21 CFR Part 820); and if applicable, radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000.10-1000.15.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 807.97), please contact the Division of Small Manufacturers, International and Consumer Assistance, toll free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Also, please note the regulation entitled, "Misbranding by reference to product labeling" (21 CFR Part 807.97). For questions regarding the reporting of adverse events (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance, toll free number (800) 638 2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,



Section 5 Indication for Use Statement

DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration Indications for Use	Form Approved: OMB No. 0910-0120 Expiration Date: December 31, 2013 See PRA Statement on last page.
---	---

510(k) Number (if known)

Device Name
Digital Automatic Wrist Blood Pressure Monitor WBPM22 Series

Indications for Use (Describe)
Digital Automatic Wrist Blood Pressure Monitor WBPM22 Series is for use by medical professional or home user. The WBPM22 Series is intended to measure the systolic and diastolic blood pressure, and pulse rate of an adult individual by using a non-invasive technique, in which an inflatable cuff is wrapped around the wrist.


Type of Use (Select one or both, as applicable)

- Prescription Use (Part 21 CFR 801 Subpart D)
 Over-The-Counter Use (21 CFR 801 Subpart C)

PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON A SEPARATE PAGE IF NEEDED.

FOR FDA USE ONLY

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)



Date: 2014.03.05

16:08:17 -05'00'

for Bram Zuckerman

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 70 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."