

JAN 24 2014

Siemens Healthcare Diagnostics, Inc.
Dimension Vista Creatinine (CRE2) Method Premarket Notification

1. 510(k) Summary

1.1 Description

Dimension Vista® System Creatinine (CRE2) Flex® reagent cartridge

This summary of 510(k) safety and effectiveness information is submitted in accordance with the requirements of SMDA 1990 and 21 CFR §807.92.

1.2. Assigned 510(k) number

The assigned 510(k) number is:

K133728

1.3 Applicant and Date

Applicant: Laura J. Duggan, Ph. D., RAC
Siemens Healthcare Diagnostics Inc.
P.O. Box 6101
Newark, DE 19714-6101
Email: laura.j.duggan@siemens.com
Phone: 302-631-7654
Fax: 302-631-6299

Date: December 5, 2013

1.4 Proprietary and Established Names

Dimension Vista® Creatinine (CRE2) Flex® reagent cartridge

1.5 Common Name

Creatinine

1.6 Regulatory Information

Dimension Vista® Creatinine (CRE2) Flex® reagent cartridge

The CRE2 method is an *in vitro* diagnostic test for the quantitative measurement of creatinine in human serum, plasma, and urine on the Dimension Vista® System. Creatinine measurements are used in the diagnosis and treatment of certain renal disease, in monitoring renal dialysis, and as a calculation basis for measuring other urine analytes.

Siemens Healthcare Diagnostics, Inc.
Dimension Vista Creatinine (CRE2) Method Premarket Notification

Classification Name:	Creatinine test system
Regulation Section:	21CFR862.1225 – Creatinine test system
Classification:	Class II
Product Code:	CGX
Panel:	Clinical Chemistry

1.7 Predicate Device

The predicate device used to demonstrate substantial equivalence to the Dimension Vista® Creatinine (CRE2) Flex® reagent cartridge is the Dimension Vista® Creatinine (CREA) Flex® reagent cartridge previously cleared under K061238.

1.8 Device Description / Test Principle

The Dimension Vista® Creatinine (CRE2) Flex® reagent cartridge uses a modified kinetic Jaffe technique. In the presence of a strong base such as sodium hydroxide, picrate reacts with creatinine to form a red chromophore. The rate of increasing absorbance at 510 nm due to the formation of this chromophore is directly proportional to the creatinine concentration in the sample and is measured using a bichromatic (510, 577 nm) rate technique.

1.9 Intended Use

The CRE2 method is an *in vitro* diagnostic test for the quantitative measurement of creatinine in human serum, plasma, and urine on the Dimension Vista® System. Creatinine measurements are used in the diagnosis and treatment of certain renal disease, in monitoring renal dialysis, and as a calculation basis for measuring other urine analytes.

1.10 Indication(s) for Use

The CRE2 method is an *in vitro* diagnostic test for the quantitative measurement of creatinine in human serum, plasma, and urine on the Dimension Vista® System. Creatinine measurements are used in the diagnosis and treatment of certain renal disease, in monitoring renal dialysis, and as a calculation basis for measuring other urine analytes.

1.11 Substantial Equivalence Information

Both the Dimension Vista® Creatinine (CRE2) Flex® reagent cartridge and the predicate Dimension Vista® Creatinine (CREA) Flex® reagent cartridge employ prepackaged reagents for use on automated clinical chemistry test systems. A comparison of the similarities and differences between the devices is provided in the following tables:

Similarities for the Dimension Vista Creatinine (CRE2) Flex® reagent cartridge:

	New Device	Predicate Device
Feature	Dimension Vista® System Creatinine (CRE2) Flex® reagent cartridge	Dimension Vista® System Creatinine (CREA) Flex® reagent cartridge (K061238)
Intended Use	The CRE2 method is an <i>in vitro</i> diagnostic test for the quantitative measurement of creatinine in human serum, plasma, and urine on the Dimension Vista® System. Creatinine measurements are used in the diagnosis and treatment of certain renal disease, in monitoring renal dialysis, and as a calculation basis for measuring other urine analytes.	The CREA method is an <i>in vitro</i> diagnostic test for the quantitative measurement of creatinine in human serum, plasma and urine on the Dimension Vista® System.
Device Technology (detection)	Modified Jaffe Methodology (creatinine alkaline picrate) with photometric detection	Same
Calibration Interval	90 days - same reagent lot	Same
Limit of Blank/ Analytical Sensitivity	Limit of Blank: 0.05 mg/dL	Analytical Sensitivity : 0.1 mg/dL
Calibration	Chem 1 Calibrator (K061838) 2 levels (n=5)	Same
Traceability & Standardization	NIST SRM 914a (IDMS assigned crystalline creatinine standard prepared by the National Institute of Standards and Technology)	Same

Differences for the Dimension Vista Creatinine (CRE2) Flex® reagent cartridge:

	New Device	Predicate Device
Feature	Dimension Vista® System Creatinine (CRE2) Flex® reagent cartridge	Dimension Vista® System Creatinine (CREA) Flex® reagent cartridge (K061238)
Measuring Range (serum)	0.15 - 20.0 mg/dL	0.10 – 20.00 mg/dL
Measuring Range (urine)	5.00 – 300 mg/dL	0.10 – 200.00 mg/dL
Reportable Range (serum)	0.15 – 40.0 mg/dL	0.10 – 40.00 mg/dL
Reportable Range (urine)	5.00 – 900 mg/dL	0.10 – 200.00 mg/dL (no specific dilution factor recommended)
Sample Volume	13.2 µL	8 µL
Reagents	Reagent 2 = Lithium Picrate (125 mM) Reagent 1 = Sodium Hydroxide (1000 mM)	Reagent 1 = Lithium Picrate (125 mM) Reagent 2 = Sodium Hydroxide (2000 mM) and potassium ferricyanide
Reagent Volumes	Volume of Reagent 1 used = 30.6 µL Volume of Reagent 2 used = 16.5 µL	Volume of Reagent 1 used = 29.6 µL Volume of Reagent 2 used = 7.2 µL
Detection Conditions	Wavelength = 510 and 577 nm Type of Measurement = Bichromatic rate	Wavelength = 510 and 600 nm Type of Measurement = Bichromatic rate
Expected Values	Serum and Plasma Males: 0.700 – 1.30 mg/dL Females: 0.550 – 1.02 mg/dL Urine Males: 0.95 – 2.49 g/24 hr	Serum and Plasma 0.6–1.3 mg/dL Urine 0.6–2.5 g/24hr

	Females: 0.60 – 1.80 g/24 hr	
Interferences	<p>No significant interference at a Creatinine concentration of 1.5 mg/dL in serum from:</p> <p>Hemoglobin at 1000 mg/dL Bilirubin (conjugated) at 20 mg/dL, Lipemia (Intralipid) at 1000 mg/dL</p> <p>Bilirubin (unconjugated) at 20 mg/dL</p>	<p>No significant interference at a Creatinine concentration of 1.3 mg/dL in serum from:</p> <p>Hemoglobin at 500 mg/dL, Bilirubin (conjugated) at 20 mg/dL, Lipemia (Intralipid) at 3000 mg/dL</p> <p>No significant interference at a Creatinine concentration of 1.1 mg/dL from:</p> <p>Bilirubin (unconjugated) at 10 mg/dL</p>

1.12 Standard/Guidance Document Reference

- Stability Testing of In Vitro Diagnostic Reagents (CEN 13640)
- Interference Testing in Clinical Chemistry; Approved Guideline (EP07-A2)
- Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline (EP09-A2-IR)
- Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline (EP05-A2)
- Evaluation of the Linearity of Quantitative Measurement (EP-06-A)
- Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline (EP17-A2)
- Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline- Third Edition (C28-A3c)
- In Vitro Diagnostic Creatinine Test System - Guidance for Industry - July 2, 1998
- In Vitro Diagnostic Devices: Guidance for the Preparation of 510(k) Submissions – Jan. 1997
- Format for Traditional and Abbreviated 510(k)'s – Guidance for Industry and Staff – Nov. 17, 2005
- Guidance for Industry and FDA Staff: Administrative Procedures for CLIA Categorization – May 7, 2008
- eCopy Program for Medical Device Submissions - Guidance for Industry and Food and Drug Administration Staff - October 10, 2013
- Refuse to Accept Policy for 510(k)s - Guidance for Industry and Food and Drug Administration Staff – December 31, 2012

1.13 Performance Characteristics

The following data represent typical method performance. These data were collected on the Dimension Vista[®] 500 System.

1.13.1 Method Comparison

The predicate device selected for the method comparison was the Dimension Vista[®] Creatinine (CREA) Flex[®] reagent cartridge cleared under K061238. Remnant de-identified serum samples were tested. No patient history information was obtained on these samples. Inclusion/exclusion data criteria are not applicable. The study included both native and spiked samples.

These studies were conducted internally by Siemens Healthcare Diagnostic Inc. R&D organization personnel. The personnel conducting the study were laboratory technicians with training similar to personnel who would conduct the tests in a hospital laboratory setting. They were trained on the operation of both the device and the predicate device. A split sample method comparison, following EP09-A2, demonstrated good agreement between the Dimension Vista[®] Creatinine (CRE2) Flex[®] reagent cartridge and the predicate Dimension Vista[®] Creatinine (CREA) Flex[®] reagent cartridge with serum patient samples.

One hundred forty serum patient samples across the assay range were tested on the Dimension Vista[®] Creatinine (CRE2) Flex[®] reagent cartridge and the Dimension Vista[®] Creatinine (CREA) Flex[®] reagent cartridge. In a second study, 112 urine samples were tested on the Dimension Vista[®] Creatinine (CRE2) Flex[®] reagent cartridge and the Dimension Vista[®] Creatinine (CREA) Flex[®] reagent cartridge. The results across the full assay range were analyzed by linear regression. Although the samples were tested in duplicate, only the first result was used for the analysis.

Comparative Method	Range (mg/dL)	Slope	Intercept (mg/dL)	Correlation Coefficient	n	Sample type
Vista CREA Assay	0.38 – 18.93	1.02	-0.11	1.000	140	serum
Vista CREA Assay	8.39 – 299.60	1.05	-4.29	0.999	112	urine

The model equation for the regression statistics is: [results for Dimension Vista[®] Creatinine (CRE2)] = slope x [comparative method results] + intercept.

An additional study was completed comparing the Dimension Vista[®] Creatinine (CRE2) Flex[®] reagent cartridge with the IDMS reference method. Remnant de-identified serum samples were tested. No patient history information was obtained on these samples. Inclusion/exclusion data criteria are not applicable. All of the samples were native.

Forty eight patient samples were tested on the Dimension Vista[®] Creatinine (CRE2) Flex[®] reagent cartridge and the IDMS reference method. The results were analyzed by linear

regression. Although the samples were tested in duplicate, only the first result was used for the analysis.

Comparative Method	Range (mg/dL)	Slope	Intercept (mg/dL)	Correlation Coefficient	n	Sample type
IDMS Reference Method	0.18– 6.32	1.06	-0.03	0.997	48	serum

The model equation for the regression statistics is: [results for Dimension Vista[®] Creatinine (CRE2)] = slope x [comparative method results] + intercept.

1.13.2 Serum Plasma Equivalency

Serum and lithium heparin plasma equivalency was demonstrated for the Dimension Vista[®] Creatinine (CRE2) Flex[®] reagent cartridge. Fifty six matched serum and lithium heparin plasma samples were tested using the Dimension Vista[®] Creatinine (CRE2) Flex[®] reagent cartridge. The table below summarizes the linear regression statistics.

Serum vs.	Slope	Intercept	Correlation Coefficient (r)	Range	n
Lithium Heparin Plasma	1.03	-0.002	0.998	0.418 – 17.9	56

One replicate of each sample was processed. All samples in the study were fresh and never frozen. The eight spiked sample sets were prepared by spiking equal amounts of purified creatinine into the matched serum and lithium heparin plasma samples.

1.13.3 Precision

Precision testing was performed in accordance with CLSI EP05-A2 Evaluation of Precision Performance of Quantitative Measurement Methods: Approved Guideline – Second Edition. Samples consisted of two (2) serum pools, three (3) levels of BioRad Multiquel material, two (levels) of BioRad Liquicheck material and two (2) urine pools. Testing was performed over twenty (20) days, two (2) separate runs with two test samples for each test material. Analysis of variance (ANOVA) was used to evaluate the data consistent with the recommendations of EP05-A2. The data are summarized in the following table:

Siemens Healthcare Diagnostics, Inc.
 Dimension Vista Creatinine (CRE2) Method Premarket Notification

		Repeatability			Within-Lab	
Sample		Mean (mg/dL)	SD	%CV	SD	%CV
Serum	Serum Pool 1	1.30	0.03	2.2	0.04	2.7
	Serum Pool 2	16.8	0.09	0.6	0.11	0.7
	BioRad Multiquel Level 1	0.530	0.02	2.9	0.02	4.0
	BioRad Multiquel Level 2	1.87	0.02	1.1	0.05	2.5
	BioRad Multiquel Level 3	7.23	0.05	0.7	0.08	1.1
Urine	Urine Pool 1	37.2	0.91	2.4	1.33	3.6
	Urine Pool 2	254	1.88	0.7	2.58	1.0
	BioRad Liquicheck Level 1	62.0	0.75	1.2	1.38	2.2
	BioRad Liquicheck Level 2	145	1.53	1.1	2.93	2.0

BioRad® is a registered trademark of Bio-Rad Laboratories, Irvine, CA 92618, USA.
 Multiquel® is a registered trademark of Bio-Rad Laboratories, Irvine, CA 92618, USA.
 Liquicheck™ is a trademark of Bio-Rad Laboratories, Irvine, CA 92618, USA.

1.13.4 Limit of Blank and Limit of Detection

The Limit of Blank (LoB) and Limit of Detection (LoD) were evaluated in accordance with CLSI EP17-A2 Protocols for Determination of Limits of Detection and Limits of Quantitation: Approved Guideline.

Dimension Vista[®] CRE2 Limit of Detection Results with Serum		
Limit	Protocol	Value
LoB	4 samples with no analyte were tested (N=5) for 3 days, one run per day, 2 reagent lots,	0.030 mg/dL
LoD	4 low patient serum samples were tested (N=5) for 3 days, one run per day, 2 reagent lots	0.056 mg/dL

Dimension Vista[®] CRE2 Limit of Detection Results with Urine		
Limit	Protocol	Value
LoB	4 samples with no analyte were tested (N=5) for 3 days, one run per day, 2 reagent lots,	0.438 mg/dL
LoD	4 low patient urine samples were tested (N=5) for 3 days, one run per day, 2 reagent lots	0.828 mg/dL

The nonparametric approach described in EP-17A2 was followed to determine the Limit of Detection.

LoB = Mean of Blank Measurement + 1.645 x Standard Deviation of Blank Measurements

LoD = Limit of Blank + C_pSD_s

- C_p is a correction factor for the 95% CI normal variate to account for bias in the SD_s estimate.
- SD_s is an estimate of method imprecision pooled from replicates of the low analyte samples

Siemens Healthcare Diagnostics, Inc.
Dimension Vista Creatinine (CRE2) Method Premarket Notification

The LoB was determined to be 0.030 mg/dL with serum samples and is consistent with the claim of 0.05 mg/dL. With urine samples, the LoB was determined to be 0.438 mg/dL and is consistent with the claim of 1.0 mg/dL.

The LoD was determined to be 0.056 mg/dL with serum samples and is consistent with the claim of 0.1 mg/dL. With urine samples, the LoD was determined to be 0.828 mg/dL and is consistent with the claim of 2.0 mg/dL.

1.13.5 Limit of Quantitation

The Limit of Quantitation (LoQ) for CRE2 for serum and plasma is 0.15 mg/dL and based on allowable total error of 0.15 mg/dL, determined consistent with CLSI Guideline EP17-A2. Total Error is calculated using: $TE = bias + 2 * SD$. The reference values are traceable to Isotope Dilution Mass Spectrometry (IDMS). The established allowable total error for creatinine published by CLIA¹, CAP², AAB³, NYS⁴ and WLSH⁵ is +/- 0.3 mg/dL or +/- 15% (greater). The LoQ of 0.15 mg/dL based on allowable total error of 0.15 mg/dL meets these recommendations.

CLIA¹ - CLIA '88 Proficiency Testing Limits, U.S. Federal Register.

CAP² - College of American Pathologists Participant Summary, April 2004.

*AAB³ - American Association of Bioanalysts Table of Grading Limits (undated, approx Oct 2005).
For details, visit www.aab.org/pts/qrdlim.htm*

NYS⁴ - Wadsworth Center Clinical Laboratory Evaluation Program, Guide to Program Requirements and Services, New York State Department of Health (undated, approx June, 2004). For details, visit www.wadsworth.org/labcert.

WLSH⁵ - WSLH Proficiency Testing is a national, full-service PT program located at the Wisconsin State Laboratory of Hygiene on the campus of the University of Wisconsin - Madison. (undated, approx June, 2004). For more information, visit www.slh.wisc.edu/pt.

The Limit of Quantitation (LoQ) for CRE2 for urine is 5.00 mg/dL and based on allowable total error of 3.00 mg/dL, determined consistent with CLSI Guideline EP17-A2. Total Error is calculated using: $TE = bias + 2 * SD$. The reference values are traceable to Isotope Dilution Mass Spectrometry (IDMS). There is no recognized allowable total error limit for urine creatinine measurements.

1.13.6 Linearity (Measurement Range)

Linearity was evaluated for serum and urine samples by using 12 equally spaced samples which spanned the assay range. Each was prepared by mixing high and low creatinine concentration samples across the measurement range as described in CLSI Evaluation of the Linearity of Quantitative Measurement Procedure (EP06-A). The low sample represents an artificial sample matrix prepared with bovine serum albumin.

Regression Statistics

Range of samples	Slope	Intercept	Correlation Coefficient	N
serum				
0.00 – 21.9 mg/dL	0.996	-0.05	1.0	12
urine				
1.71 – 338 mg/dL	0.995	0.47	1.0	12

1.14 Analytical Specificity

12.5.1 Non-interfering Substances

CLSI EP7A2 was followed for the interference testing. The interference study was conducted using a "paired difference worst case scenario" approach where these compounds were spiked into fresh sample pools containing either low or high levels of creatinine analyte in both serum pools and urine pools.

Siemens Healthcare Diagnostics, Inc.
Dimension Vista Creatinine (CRE2) Method Premarket Notification

DV CRE2 Serum/Plasma Substance		DV CRE2					
Substance	Concentration of Substance (mg/dL)	LOW pool (~ 1.5mg/dL)			HIGH pool (~ 5.0 mg/dL)		
		Mean	Control	% Diff	mean	control	% Diff
Acetaminophen	20	1.41	1.40	1%	5.03	5.00	1%
Acetoacetate	20	1.35	1.38	-3%	4.87	4.87	0%
Amikacin	8	1.40	1.40	0%	5.05	5.00	1%
Ampicillin	5.3	1.43	1.40	2%	5.06	5.00	1%
Ascorbic Acid	6	1.41	1.40	0%	4.99	5.00	0%
Caffeine	6	1.40	1.42	-2%	4.97	4.96	0%
Carbamazepine	3	1.43	1.42	1%	5.01	5.06	-1%
Cephalexin	25	1.35	1.38	-3%	4.86	4.87	0%
Cephapirin	25	1.43	1.39	3%	4.99	4.97	0%
Cephadrine	25	1.44	1.39	3%	5.00	4.97	1%
Chloramphenicol	5	1.43	1.41	1%	5.01	5.03	0%
Chlordiazepoxide	1	1.40	1.43	-2%	5.01	5.02	0%
Chlorpromazine	0.2	1.40	1.40	0%	5.00	4.99	0%
Cholesterol Supertrate	503	1.34	1.30	3%	4.58	4.57	0%
Cimetidine	2	1.40	1.42	-1%	4.97	4.96	0%
Dextran 40	6000	1.56	1.51	3%	5.13	5.23	-2%
Diazepam	0.51	1.44	1.43	1%	5.04	5.02	0%
Digoxin	6.1 ng/mL	1.41	1.42	0%	4.95	4.93	0%
EDTA	200	1.37	1.39	-2%	4.93	4.91	0%
Erythromycin	6	1.44	1.41	2%	5.00	5.03	-1%
Ethanol	400	1.41	1.40	0%	4.98	5.00	0%
Ethosuximide	25	1.40	1.42	-1%	4.97	4.96	0%
Furosemide	6	1.43	1.41	1%	5.03	5.03	0%
Gentamicin	1	1.42	1.40	1%	5.00	5.00	0%
Heparin (196 Units/mg)	3 U/mL	1.40	1.42	-1%	4.94	4.99	-1%
Ibuprofen	50	1.43	1.42	1%	5.01	5.06	-1%
Immunoglobulin G (IgG)	5000	1.56	1.51	3%	5.14	5.23	-2%
Isopropanol	1.0 g/dL	1.40	1.39	0%	4.91	4.91	0%
Lidocaine	1.2	1.43	1.41	1%	5.08	5.03	1%
Lithium	2.2	1.40	1.40	0%	5.00	4.99	0%
Nicotine	0.1	1.44	1.42	1%	5.05	5.06	0%
Nortriptyline	1000 ng/mL	1.39	1.39	-1%	4.96	4.91	1%
Penicillin G (1654 Units/mg)	25 U/mL	1.41	1.40	1%	4.98	4.99	0%
Pentobarbital	8	1.43	1.42	0%	4.99	4.98	0%
Phenobarbital	10	1.44	1.42	1%	5.00	4.98	0%
Phenytoin	5	1.42	1.43	-1%	5.02	5.02	0%
Potassium oxalate	500 mg/dL	1.40	1.38	1%	4.86	4.87	0%
Primidone	4	1.42	1.42	0%	5.04	5.06	0%
Propoxyphene	0.16	1.42	1.43	0%	5.02	5.02	0%
Protein, Albumin	6000	1.56	1.51	3%	4.88	5.345	-9%
Protein, Total	12g/dL	1.57	1.48	6%	5.13	5.22	-2%
Salicylic acid	60	1.42	1.41	1%	5.03	5.03	0%
Sodium fluoride	400	1.30	1.30	0%	4.62	4.67	-1%
Theophylline	4	1.39	1.42	-2%	5.00	4.96	1%
Urea	500	1.43	1.40	2%	5.04	4.99	1%
Uric acid	20	1.43	1.41	1%	4.97	4.99	0%
Valproic acid	50	1.41	1.42	-1%	4.97	4.96	0%
Vancomycin	10	1.41	1.40	1%	4.99	5.00	0%

Siemens Healthcare Diagnostics, Inc.
 Dimension Vista Creatinine (CRE2) Method Premarket Notification

DV CRE2 Urine Method Specific		DV CRE2					
Substance	Concentration of Substance	LOW pool (~40 mg/dL)			HIGH pool (~175 mg/dL)		
		Mean	Control	%Diff	mean	control	%Diff
50% Acetic Acid	25ml/24 hr collection	40.75	40.25	1%	180.0	178.8	1%
6N Hydrochloric Acid	0.60%	41.30	41.35	0%	178.8	182.5	-2%
6N Nitric Acid	0.60%	41.18	41.35	0%	181.3	182.5	-1%
Acetone	100 mg/dL	42.93	41.98	2%	184.5	186.8	-1%
Albumin	0.5 g/dL	42.08	42.98	-2%	183	185.8	-1%
Ascorbic Acid	1.5 g/dL	38.65	40.48	-5%	173.5	176.0	-1%
Bilirubin (conjugated)	2 mg/dL	42.08	41.75	1%	182.0	179.8	1%
Boric Acid	1% w/v	41.83	41.75	0%	185.8	184.0	1%
Ethanol	1 g/dL	41.98	41.65	1%	183.8	179.8	2%
Gamma Globulin	0.5 g/dL	42.75	41.88	2%	184.3	185.5	-1%
Glucose	2 g/dL	42.65	41.88	2%	184.8	185.5	0%
Hemoglobin	115 mg/dL	41.78	41.58	0%	183.8	183.8	0%
Oxalic Acid	0.1 g/dL	40.28	40.48	0%	178.3	176.0	1%
Sodium Carbonate	5g/24 hr collection	41.08	40.65	1%	182.0	183.0	-1%
Sodium Fluoride	1% w/v	41.05	41.75	-2%	187.3	184.0	2%

12.5.2 Interfering Substances

The Dimension Vista® Creatinine (CRE2) Flex® reagent cartridge was evaluated for interference according to CLSI EP7-A2. Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference. Dilution studies were conducted to determine the level at which the spiked substance no longer displayed significant interference. These studies were conducted at two creatinine analyte concentrations, if both sample pools show significant interference. This study was conducted as needed for both serum pools and urine pools.

Siemens Healthcare Diagnostics, Inc.
Dimension Vista Creatinine (CRE2) Method Premarket Notification

Substance	Concentration of Substance (mg/dL)	Mean Test CRE2	Mean Control	% Diff
Acetone	150	5.77	4.87	18%
Acetone	75	5.32	4.87	9%
Acetone	37.5	5.15	4.87	6%
Acetone	37.5	1.59	1.39	15%
Acetone	18.75	1.50	1.39	8%
Cefoxitin	5	5.29	4.97	6%
Cefoxitin	2.5	1.58	1.40	13%
Cefoxitin	1.25	1.49	1.40	6%
Cephalothin	15	1.53	1.39	10.1%
Cephalothin	15	5.09	4.97	2%
Cephalothin	13.75	1.49	1.40	6%
Glucose	400mg/dL	1.54	1.37	12%
Glucose	300mg/dL	1.51	1.37	10%
Glucose	200mg/dL	1.45	1.37	6%
Glucose	600mg/dL	5.15	4.91	5%
Pyruvate	10.5 mg/dL	5.67	4.87	16%
Pyruvate	5.26 mg/dL	5.31	4.92	8%
Pyruvate	1.32 mg/dL	1.49	1.38	8%
Triglycerides	3000	4.93	4.39	12%
Triglycerides	2500	4.97	4.57	9%
Triglycerides	2000	5.01	4.65	8%
Triglycerides	2000	1.49	1.32	13%
Triglycerides	1500	1.48	1.35	10%
Triglycerides	1000	1.50	1.39	8%

12.5.3 Hemolysis, Icterus, Lipemia (HIL) Interference

The Dimension Vista[®] Creatinine (CRE2) Flex reagent cartridge was evaluated for interference according to CLSI EP7-A2. Bias is the difference in the results between the control sample (without the interferent) and the test sample (contains the interferent) expressed in percent. Bias exceeding 10% is considered interference.

Siemens Healthcare Diagnostics, Inc.
Dimension Vista Creatinine (CRE2) Method Premarket Notification

DV CRE2		HIL Interference		
Substance	Concentration of Substance (mg/dL)	Mean Test CRE2 Result (mg/dL)	Mean Control CRE2 Result (mg/dL)	% Diff
Intralipid 20%	2000mg/dL	5.27	4.56	15.5%
Intralipid 20%	1500mg/dL	5.15	4.71	9.5%
Intralipid 20%	1500mg/dL	1.51	1.31	15.1%
Intralipid 20%	1000mg/dL	1.45	1.36	7.0%
Hemoglobin	1000mg/dL	1.45	1.35	7.0%
Hemoglobin	1000mg/dL	4.87	4.72	3.2%
Bilirubin (unconj)	40mg/dL	0.98	1.40	-29.6%
Bilirubin (unconj)	20mg/dL	1.44	1.45	-0.2%
Bilirubin (unconj)	40mg/dL	4.61	4.93	-6.6%
Bilirubin (conj)	40mg/dL	1.16	1.38	-16.1%
Bilirubin (conj)	20mg/dL	1.38	1.44	-4.5%
Bilirubin (conj)	40mg/dL	4.72	4.94	-4.5%
Bilirubin (conj)	20mg/dL	4.96	5.04	-1.5%

1.15 Expected Values

Serum and Plasma

Transference was used to validate established reference ranges for serum and plasma creatinine following CLSI C28-A3c.

Serum samples from forty normal healthy adult female donors and forty-three normal healthy adult male donors were tested N=1 using the Dimension Vista[®] Creatinine (CRE2) Flex[®] reagent cartridge. The range was considered validated if ≤10% of the samples fell outside of established range.

Gender	Males	Females
Range	0.700 -1.30 mg/dL	0.550 -1.02 mg/dL
# Samples	53	40
# Sample outside Range	5	2
% Samples outside range	9.4%	5.0%

These data support the published ranges for Males (0.7 to 1.3 mg/dL (Tietz 1999) and Females (0.55 to 1.02 mg/dL (Clin Chem 54:3 (2008))).

Urine

These data support the published ranges for Males (0.7 to 1.3 mg/dL (Tietz 1999) and Females (0.55 to 1.02 mg/dL (Clin Chem 54:3 (2008))).

Similarly published reference ranges were validated for urine specimens using samples from twenty normal healthy adult female donors and twenty-two normal healthy adult male donors following CLSI C28-A3c. Samples were processed N=1 using the Dimension Vista® Creatinine (CRE2) Flex® reagent cartridge. The range was considered validated if ≤10% of the samples fell outside of established range.

Gender	Males	Females
Reference Range	0.95-2.49 g/day	0.60-1.80 g/day
# Samples	22	20
# Sample outside Range	2	1
% Samples outside range	9.1%	5.0%

These data support the published ranges for Males (0.95 to 2.49 g/day (Clin Chem Acta 344 (2004) 137-148) and Females (0.6 to 1.8 g/day (Tietz 1999))).

1.16 Results Outside of the Measuring (Reportable) Range

The instrument generates a flag which states "Above Assay Range" or "Below Assay Range" when it obtains a result outside of the measuring range. Serum and plasma samples with results in excess of 20 mg/dL and urine samples above 300 mg/dL are reported as "above Assay Range" and should be repeated on dilution.

Autodilution for Serum and Plasma Samples: The autodilute sample volume is 6.6 µL and recommended dilution factor 1:2 for serum and plasma. This extends the serum and plasma reportable range to 40 mg/dL.

Autodilution (AD) for Urine Samples: Autodilutions for urines are a two-step process. Urine samples automatically get a 1:15 AUD in the aliquot tray followed by a second 1:3 dilution in the aliquot tray making a final 1:45 dilution. This extends the urine reportable range to 900 mg/dL [79560 µmol/L]. Refer to your Dimension Vista® Operator's Guide.

These protocols were tested with samples at the upper end of the undiluted range. The mean percent recovery was calculated and met the acceptance criteria of 90 - 110% recovery.

DV CRE2 Reportable Range Dilution of above assay range results

Serum

Serum results > 20 mg/dL can be automatically re-run using 1/2 sample (6.6 μ L vs 13.2 μ L).
This extends the serum reportable range to 40.0 mg/dL.

Sample Volume (μ L)	Mean of N=5 reps
13.2 μ L (normal volume)	16.5
6.6 μ L (auto-dilute volume)	16.0
% Bias	-3.0%

Urine

Urine results > 300 mg/dL can be automatically diluted 3x with water and re-run.
This extends the urine reportable range to 900 mg/dL.

	Mean of N=5 reps
undiluted (normal)	268
auto-dilute 3x with water	245
% Bias	-8.6%

Conclusion

Results support serum extended range up to 40.0 mg/dL, and urine extended range up to 900 mg/dL.

1.17 Conclusion

The Dimension Vista[®] Creatinine (CRE2) Flex[®] reagent cartridge (K1033A) is substantially equivalent in principle and performance to the Dimension Vista[®] Creatinine (CREA) Flex[®] reagent cartridge cleared under K061238. Comparative testing described in the submission demonstrates substantially equivalent performance.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

January 24, 2014

SIEMENS HEALTHCARE DIAGNOSTICS INC.
LAURA DUGGAN, PH.D.
REGULATORY TECHNICAL SPECIALIST
P.O. BOX 6101
NEWARK DE 19714-6101

Re: K133728

Trade/Device Name: Dimension Vista® Creatinine (CRE2) Flex® Reagent Cartridge
Regulation Number: 21 CFR 862.1225
Regulation Name: Creatinine test system
Regulatory Class: II
Product Code: CGX
Dated: December 05, 2013
Received: December 6, 2013

Dear Dr. Duggan:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Parts 801 and 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulations (21 CFR Parts 801 and 809), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

 Courtney H. Lias -S

Courtney H. Lias, Ph.D.
Director
Division of Chemistry and Toxicology Devices
Office of In Vitro Diagnostics
and Radiological Health
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

k133728

Device Name

Dimension Vista Creatinine (CRE2) Flex Reagent Cartridge

Indications for Use (Describe)

The CRE2 method is an in vitro diagnostic test for the quantitative measurement of creatinine in human serum, plasma, and urine on the Dimension Vista® System. Creatinine measurements are used in the diagnosis and treatment of certain renal disease, in monitoring renal dialysis, and as a calculation basis for measuring other urine analytes.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.

FOR FDA USE ONLY

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

Ruth A. Chesler -S