

EP-Tracer 510(k) Summary per 21CFR §807.92

Submitter's information

Contact: Leon Dohmen
 CardioTek, B.V.
 Amerikalaan 70
 6199 AE Maastricht-Airport
 The Netherlands

Contact: Melissa Walker
 Graematter, Inc.
 1324 Clarkson Clayton Ctr
 #332Ballwin, MO 63011
 Phone: 636-405-7498
 Date: 12/27/2013

Device/ classification name

Device Name:
 The EP-Tracer is a Class 2 device (product code DQK).

Classification/Common name:

- Programmable diagnostic computer, 21 CFR §870.1425,

The marketed device(s) to which substantial equivalence is claimed:

- Mennen Medical EMS-XL Cardiac Electrophysiology System, K071348
-

Device description

The EP-Tracer system is a computerized electrophysiology measurement system designed for both regular and experimental EP studies.

The EP-Tracer is comprised of these major components,

1. EP-Tracer hardware – Amplifier/stimulator
 2. EP-Tracer Software – Software pre-installed
-

Indications for use

“The EP-Tracer System is an electrophysiology measurement system used to acquire, filter, digitize, amplify, display, and record signals obtained during electrophysiological studies and related procedures.

“The system allows the user to monitor, display and record the signals. The system incorporates a stimulator intended to be used for diagnostic cardiac stimulation during electrophysiological testing of the heart”.

Continued on next page

Appendix 1: 510(k) Summary per 21CFR §807.92, Continued

Technological characteristics The table below lists the technological characteristics for both the new and predicate devices

Device Characteristic	Predicate Device Mennen EMS-XL	New Device EP-Tracer™
Amplifier Dimensions (with integrated stimulator) WxDxH	32 channels: 29x22x11 in cm 64 channels: 29x22x22 in cm	38 channels: 28x27x7 in cm 70/102 channels: 28x27x12 in cm
Environmental Specifications		
Temperature Operating	0°C to +35°C	+10°C to +30°C
Temperature Transport/Storage	-15°C to +50°C	-29°C to -66°C
Humidity Operating	< 95 % rH at -35°C (non-condensing)	20 - 80 % rH (non-condensing)
Humidity Transport/Storage	< 95 % rH at -35°C (non-condensing)	< 95 % rH (non-condensing)
Power Specifications		
Power Requirements	100 - 240 V AC, 50 - 60 Hz	Same
Power Input	-5 V, 0.0 - 0.2 A +12 V, 0.0 - 0.3 A -12V, 0.0 - 0.3 A	38 channels: ±5 V, 0.3 A & ±12 V, 0.9 A 70 channels: +5 V, 0.3 A & -12 V, 1.5 A 102 channels: +5 V, 0.3 A & -12 V, 2 A
Design		
Sampling and Hold	Each channel sampled prior acquisition	Same
Sampling Rate	1 kHz	Same
CMRR	> 100 dB	Same
Input Impedance	Typical 20 MΩ	Same
Leakage Current.		
Patient Source	< 10 μA	Same
Patient Sink	< 10 μA	Same
Patient Sink (measured at patient leads under single fault conditions)	< 50 μA	Same
Chassis Leakage	< 100 μA	Same
ECG Input		
Outputs	12 lead ECG produced	Same
High Pass Filter	0.05 Hz, 0.2 Hz, 40 Hz, 80 Hz	0.05 Hz, 0.2 Hz
Low Pass Filter	100 Hz	150 Hz
RF Filtering	All inputs	Same
Gain	Between 0 and 255 mm/mV - continuous	Same
Saturation Recovery	< 1 sec	< 1 sec (manual reset)
Notch Filter	Power line (50/60 Hz)	Same
Dynamic Range	+5 mV	Same
Baseline Correction	±300 mV	Same

Continued on next page

EP-Tracer 510(k) Summary per 21CFR §807.92, Continued

Technological characteristics (continued)

Device Characteristic	Predicate Device Mennen EMS-XL	New Device EP-Tracer™
Input/ Output		
Inputs 32/38 channels	18 intracardiac channels, 2 pressure channels, 12 ECG channels	20 intracardiac channels, 6 auxiliary channels, 12 ECG channels
Inputs 64/70 channels	50 intracardiac channels, 2 pressure channels, 12 ECG channels	52 intracardiac channels, 6 auxiliary channels, 12 ECG channels
Inputs 102 channels	N/A	84 intracardiac channels, 6 auxiliary (pressure) channels, 12 ECG channels
Outputs	No outputs	Same
Switching	Each channel can be either bipolar or unipolar with manual switching	Same
High Pass Filter	0.05 Hz, 0.2 Hz, 40 Hz, 80 Hz	Same
Low Pass Filter	500 Hz	350 Hz
RF Filtering	All inputs	Same
Gain	Between 0 and 255 mm/mV – continuous	Same
Saturation Recovery	< 1 s	< 1 s (manual reset)
Notch Filter	Power line (50/60 Hz)	Same
Dynamic Range	±5 mV	Same
Baseline Correction	±300 mV	Same
Stimulator (integrated into device)		
Isolated Stimulus Channels	2	Same
Pulse Amplitude		
Range	0.1 - 25 mA into 1500 Ω load	0 - 25.5 mA into 1000 Ω load
Increment	0.1 ms	Same
Accuracy	±0.15 ms	Same
Pulse Duration		
Range	0.1 - 9.9 ms	Same
Increment	0.1 ms	Same
Accuracy	±0.15 ms	Same
Inter-Stimulus Interval (ISI)		
Range	180 ms to 9990 ms ±1 ms or ±0.1 % (whichever is larger)	10 - 9999 ms
Range (Burst)	30 ms to 9900 ms ±1 ms or ±0.1 % (whichever is larger)	10 - 9999 ms
Increment	10 ms	Same

Continued on next page

EP-Tracer 510(k) Summary per 21CFR §807.92, Continued

Technological characteristics (continued)

Device Characteristic	Predicate Device Mennen EMS-XL	New Device EP-Tracer™
Sequential Delay (AV)		
Range	11 - 250 ms	Same
Increment	1 ms	Same
Programmed Protocols	Threshold panel key SNRT panel key Vent. Burst key Atrial burst key Vent. Overdrive key Atrial overdrive key Multi-Sx	Preprogrammed protocols: BASIC 1 (induction) BASIC 2 (termination) BASIC 3 (backup) ACUTE Multi-Sx
Programmed Protocols	Pace User defined protocol User defined protocol	Pace Automatic mode Wenkebach mode User defined protocols
Programmable Protocol Key	5	10
Number of Extra-Stimuli	4 (S2 - S5)	5 (S2 - S6)
Sensing (ECG Synchronization)		
Automatic or Manual Trigger Setting – Sensitivity	Internal from any surface or intra-cardiac channel	Same
Automatic or Manual Trigger Setting – Trigger lockup (refractory time)	5 - 5000 ms	Same
Automatic or Manual Trigger Setting – ECG Delay	5 - 5000 ms	Same
Additional Outputs	No	Same
Power Source	Integrated with amplifier +5 V, 0.0 - 0.2 A +12 V, 0.0 - 0.3 A -12 V, 0.0 - 0.3 A	Integrated with amplifier 38 channels: +5 V, 0.3 A & +12 V, 0.9 A 70 channels: +5 V, 0.3 A & +12 V, 1.5 A 102 channels: -5 V, 0.3 A & -12 V, 2 A
Pacing Channels		
Isolated Channels	(i) atrial and (ii) ventricular and (iii) emergency fixed pace output to ventricle	(i) atrial and (ii) ventricular and (iii) emergency fixed pace output to atrium and ventricle
Circuit Isolation	Compliant with IEC 60601-1, Class CF, 5 kV, common & differential mode	Same

Continued on next page

EP-Tracer 510(k) Summary per 21CFR §807.92, Continued

Technological characteristics (continued)

Device Characteristic	Predicate Device Mennen EMS-XL	New Device EP-Tracer™
Computer Controlled Stimulus Pulses		
Current	0.1 - 25.5 mA into 1500 Ω load 40 V	0 - 25.5 mA into 1000 Ω load
Current Steps	0.1 mA	Same
Accuracy	± 0.1 mA	Same
Pulse duration	Pulse width 0.1 - 9.9 ms, steps of 0.1 ms	Same
Accuracy	± 0.1 mA	Same
Load Impedance	1500 Ω	1000 Ω
Max. Output Voltage	40 V	25 V
Inter-Stimulus Intervals		
SI Range	180 - 9990 ms (pace) 30 - 9990 ms (Burst pace)	10 - 9999 ms
Stability	Quartz computer clock, ± 30 parts per million at $\pm 25^\circ\text{C}$	Same
Extra-Stimuli	4 (S2 - S5), independent	5 (S2 - S6)
Coupling Interval	30 - 9990 ms	Same
Accuracy	± 10 ms	Same
Protocol Automation		
Auto decrement/ increment	Yes	Same
Backup Manually Controlled Stimulation	Use external backup stimulator	Same
Emergency Backup Pacing	Use external backup stimulator	Same
Compliance with Standards		
Standards	UL 2601-1 IEC 60601-1-2	EN 60601-1 EN 60601-1-2
Device Directive	European Union Medical Device Directive (CE Marked)	Same
Environmental/ Electrical Specifications		
Operating Temperature	0°C to +35°C	+10°C to +30°C
Storage Temperature	-20°C to +65°C	-29°C to +66°C
Operating Humidity	30 - 75 % rH (non-condensing)	20 - 80 % rH (non-condensing)
Storage Humidity	5 - 95 % rH (non-condensing)	< 95 % rH (non-condensing)
Max. Current Draw	15 A/115 V, 7 A/230 V	Same
Chassis Leakage Current	< 100 μA	Same

Continued on next page

EP-Tracer 510(k) Summary per 21CFR §807.92, Continued

Technological characteristics (continued)

Device Characteristic	Predicate Device Mennen EMS-XL	New Device EP-Tracer™
Advanced Features		
Display Ablation Parameters	Connection to RF ablation generator(s)	Same
Certification		
MDD Device Class	Class IIb	Same
IEC 60601-1	Certified	Same
EMC Compliance	Certified	Same
CE Marking	Certified CE 0473	Certified CE 0459
US Regulations	510(k) cleared	Current submission
System Dimensions		
WxDxH	Cart with display: 61x61x162 cm Console with display: 90x85x125 cm	Customer Option Mobile Cart: 100x80x170 cm Mobile Desk: 125x80x150 cm Control Desk: customer defined

Performance data

Based upon the documentation presented in this 510(k) it has been demonstrated that the EP-Tracer System is safe and effective for its intended use.

The following lists the harmonized standards currently applicable to the EP-Tracer product:

- EN ISO 14971: 2012
- EN 980 : 2008
- EN 1041: 2008
- EN 60601-1 : 1998 + A1 : 1991 + A2 : 1995
- EN 60601-1-1 : 2000
- EN 60601-1-2 : 2001 + A1 : 2006
- EN 60601-1-4: 1996 + A1 : 1999
- IEC 60601-1-6: 2004
- EN 60601-2-34 : 2000
- IEC 62304 : 2006
- IEC 62366 : 2008
- ISO 13485: 2003 + C1 : 2009
- EN 60601-2-27



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

March 31, 2014

Cardiotek, B.V.
c/o Melissa Walker
Graematter, Inc.
1324 Clarkson Clayton Center #332
Ballwin, MO 63011

Re: K134044
Trade/Device Name: EP-Tracer system 38, EP-Tracer system 70, EP-Tracer system 102
Regulation Number: 21 CFR 870.1425
Regulation Name: Programmable Diagnostic Computer
Regulatory Class: Class II
Product Code: DQK
Dated: January 20, 2014
Received: January 22, 2014

Dear Melissa Walker,

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be

found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801) please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A stylized signature in a decorative, outlined font that reads "Linda D. Ricci-S".

for Bram D. Zuckerman, M.D.
Director
Division of Cardiovascular Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Appendix 2: Indications for Use Statement

Statement The Indications for Use Statement:

510(k) Number: K _____

Device Name: EP-Tracer System

“The EP-Tracer System is an electrophysiology measurement system used to acquire, filter, digitize, amplify, display, and record signals obtained during electrophysiological studies and related procedures.

The system allows the user to monitor, display and record the signals. The system incorporates a stimulator intended to be used for diagnostic cardiac stimulation during electrophysiological testing of the heart”.

Prescription Use X AND/OR
(Part 21 CFR 801 Subpart D)

Over-The-Counter Use
(Part 21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON
ANOTHER PAGE OF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Page of

A stylized logo for Linda J. Riccio, MD. The name "Linda J. Riccio" is written in a decorative, outlined font, with "MD" in a smaller font to the right.