



MAR 27 2014

**Special 510(k) Summary  
as required by 21 CFR 807.92(a)  
K140106**

A ) Submitted by: Renovis Surgical Technologies, Inc.  
1901 W. Lugonia Ave, Ste 340  
Redlands, CA 92374  
Phone: 909-557-2360  
Fax: 909-307-8571

Official Contact: Anthony DeBenedictis  
Vice President of Quality Assurance

Consultant: Sharyn Orton, Ph.D.  
MEDIcept, Inc.  
200 Homer Ave  
Ashland, MA 01721

Date: February 28, 2014

B) Device Name: Intervertebral Fusion Device With Bone Graft, Lumbar  
Common Name: Intervertebral body fusion device  
Proprietary Name: S128 Anterior Lumbar Interbody Fusion (ALIF) System  
Device Class: Class II – 888.3080  
Regulation and Product code: 888.3080, OVD - Intervertebral body fusion device  
Classification panel: Orthopedic

C) Predicates:

- K131122 Renovis S128 Anterior Lumbar Interbody Fusion (ALIF) System

D) Device Description:

The Renovis S128 Anterior Lumbar Interbody Fusion (ALIF) System is cleared under K131122. The S128 ALIF implants (cages) are to be used with the bone screws and anterior cover plate assembly and requires no additional supplementary fixation systems. The Renovis S128 ALIF System implants are available in a variety of sizes (widths, height, depths, and bone screw sizes; see below) to suit the individual pathology and anatomical conditions of the patient. The implants are manufactured from PEEK or additively manufactured and machined Titanium. The bone screws and cover plate

assembly are both manufactured from Titanium alloy. The PEEK markers are manufactured from Tantalum.

This Special 510(k) Premarket Notification is submitted for the additional offering of gamma sterilized S128 implants (PEEK cages; titanium cages, screws and cover plates).

#### E) Indications For Use:

The Renovis S128 Anterior Lumbar Interbody Fusion (ALIF) System is indicated for intervertebral body fusion procedures in skeletally mature patients with degenerative disc disease (DDD) of the lumbar spine at one or two contiguous levels from L2-S1. Degenerative disc disease is defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies. These DDD patients may have up to Grade 1 spondylolisthesis or retrolisthesis at the involved level(s). Renovis S128 ALIF System implants are to be used with autogenous bone graft.

Patients should be skeletally mature and have at least six months of non-operative treatment.

The Renovis S128 ALIF System is a stand-alone device and is intended to be used with the cover plate and screws provided and requires no additional supplementary fixation. The anterior cover plate must be utilized whenever the device is implanted using the bone screws provided. Should the physician choose to use fewer than the four screws provided, additional supplemental fixation cleared by the FDA for use in the lumbar spine must be used.

#### F) Substantial Equivalence Comparison and Discussion

	Renovis S128 ALIF System	Renovis S128 ALIF System K131122
Product code	OVD	
Implant Material (cages, cover plate, and screws)	PEEK (ZA-500) per ASTM F2026 Titanium alloy Ti-6Al-4V per ASTM F136	
PEEK marker material	Tantalum per ASTM F 560-08	
Dimensions (mm)		
A/P	26, 28, 30	
M/L	30, 34, 38	
H	11 - 21	
Lordosis	7°, 12°	
Number of screws	4	
Screw Diameter (mm)	4.5, 5.	
Screw Length (mm)	20, 25, 30, 35	
Cover plate (mm)	8.3 x 22	
Provided sterile?	Yes - gamma	No

### *Conclusion*

Based upon the same intended use, design, function, technology, and materials, the Renovis S128 ALIF System is substantially equivalent to the predicate devices and does not raise new issues of safety or effectiveness.

#### G) Performance Data

Implants are sterilized by  $^{60}\text{Co}$  Gamma irradiation validated to a sterility assurance level (SAL) of  $10^{-6}$  by selecting and substantiating a 25 kGy dose by the VDmax<sup>25</sup> method, according to ISO 11137-1. Titanium alloy components are not affected by gamma sterilization and/or aging. PEEK components were tested to an average dose of 200kGy and underwent accelerated aging to simulated 10 or more years. After aging the samples were tested, and the results did not show any significant difference between untreated PEEK and gamma treated and aged PEEK.

Packaging has been validated to maintain sterility for 3 years in compliance with ISO 11607-2 demonstrates compliance with accelerated aging simulation per ASTM F1980-7 and real time aging; and performance following distribution per ISTA 2A.

### *Conclusion*

Gamma sterilization does not have a negative effect on Renovis S128 ALIF System implants.

#### H) Compliance with Design Controls

The results of design validation support that the Renovis S128 ALIF System is substantially equivalent to the predicate device and the offering of gamma irradiated implant components does not raise new issues of safety or effectiveness.

#### I) Compliance with Consensus Standards and FDA Guidance

##### **Standards - Renovis complies with:**

- ASTM F2026 Standard Specification for PEEK Polymers for Surgical Implant Applications
- ASTM F136 Standard Specification for Wrought Titanium-6Aluminum-4Vanadium ELI (Extra Low Interstitial) Alloy for Surgical Implant Applications, and
- ASTM F 560-08, Standard Specification for Unalloyed Tantalum for Surgical Implant Applications
- ISO 11137-1: 2006: Sterilization of health care products -- Radiation -- Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices
- ANSI/AAMI/ISO 11137-2:2006 Sterilization of health care products -- Radiation -- Part 2: Establishing the sterilization dose
- ISO11607-2: 2006 Packaging for terminally sterilized medical devices -- Part 2: Validation requirements for forming, sealing and assembly processes
- ASTM F1980-07: 2007 Standard Guide for Accelerated Aging of Sterile Barrier Systems for Medical Devices

- ASTM D4169:2009 Standard Practice for Performance Testing of Shipping Containers and Systems
- ISTA 2A, 2011 Partial-Simulation Performance Test Procedure: Packaged Products 150lb (68 kg) or Less

**FDA guidance:**

- Updated 510(k) Sterility Review Guidance K90-1; Guidance for Industry and FDA, August 2002
- Guidance for Industry and FDA Staff - Class II Special Controls Guidance Document: Intervertebral Body Fusion Device, October 2007
- Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: FDA Review Guidance, April 1996
- Guidance for Industry and FDA Staff Spinal System 510(k)'s, May 2004 (for labeling language)



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

March 27, 2014

Renovis Surgical Technologies, Incorporated  
% Sharyn Orton, Ph.D.  
MEDIcept, Incorporated  
200 Homer Avenue  
Ashland, Massachusetts 01721

Re: K140106

Trade/Device Name: Renovis S128 Anterior Lumbar Interbody Fusion (ALIF) System  
Regulation Number: 21 CFR 888.3080  
Regulation Name: Intervertebral body fusion device  
Regulatory Class: Class II  
Product Code: OVD  
Dated: March 5, 2014  
Received: March 7, 2014

Dear Dr. Orton:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

Page 2 – Sharyn Orton, Ph.D.

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Ronald P. Jean -S for

Mark N. Melkerson  
Director  
Division of Orthopedic Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**Indications for Use Statement**

510(k) Number (if known): K140106

Device Name: Renovis S128 Anterior Lumbar Interbody Fusion (ALIF) System

**Indications for Use:**

The Renovis S128 Anterior Lumbar Interbody Fusion (ALIF) System is indicated for intervertebral body fusion procedures in skeletally mature patients with degenerative disc disease (DDD) of the lumbar spine at one or two contiguous levels from L2-S1. Degenerative disc disease is defined as discogenic back pain with degeneration of the disc confirmed by history and radiographic studies. These DDD patients may have up to Grade I spondylolisthesis or retrolisthesis at the involved level(s). Renovis S128 ALIF System implants are to be used with autogenous bone graft.

Patients should be skeletally mature and have at least six months of non-operative treatment.

The Renovis S128 ALIF System is a stand-alone device and is intended to be used with the cover plate and screws provided and requires no additional supplementary fixation. The anterior cover plate must be utilized whenever the device is implanted using the bone screws provided. Should the physician choose to use fewer than the four screws provided, additional supplemental fixation cleared by the FDA for use in the lumbar spine must be used.

Prescription Use   X   AND/OR Over-the-Counter Use \_\_\_\_\_  
(21 CFR 801, Subpart D) (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NECESSARY)

---

*Concurrence of CDRH, Office of Device Evaluation (ODE)*

**Anton E. Dmitriev, PhD**

**Division of Orthopedic Devices**