

3M Health Care  
3M Infection Prevention

3M Center  
2510 Conway, Bldg 275-05-W-06  
St. Paul, MN 55144  
651 733 1110

k140250

**JUL 25 2014**



**510(k) SUMMARY**  
**Sponsor Information:**

Applicant Name: 3M Company  
3M Health Care  
Infection Prevention Division  
3M Center  
2510 Conway Ave., Bldg 275-5W-06  
St. Paul, MN 55144

Contact Person: Kristin Totushek  
Regulatory Affairs

Phone Number: 651-736-6117  
Fax Number: 651-737-5320  
e-mail: ktotushek@mmm.com

Date Prepared: July 23, 2014

**Device Name and Classification:**

Trade Name: 3M™ Steri-Drape™ CHG  
Antimicrobial Incise Drape

Common or Usual Name: Antimicrobial Incise Drape

Classification Name: Surgical Drape and Drape Accessories

Product Code: KXX per 21CFR 878.4370

Predicate Device: 3M™ Ioban™ 2 Antimicrobial Incise Drape  
(K801550)

**Description of Device:**

The 3M™ Steri-Drape™ CHG Antimicrobial Incise Drape is a sterile adhesive film that is applied pre-operatively to skin. The film adheres and conforms to body contours, providing a sterile surface on top of the skin prior to surgery. The surgeon incises through the adhesive film and skin, ensuring a sterile film barrier all the way to the edge of the incision. The CHG contained in the adhesive provides antimicrobial activity to reduce the risk of microbial contamination in the surgical wound.

3M Health Care

**3M Infection Prevention**

3M Center  
2510 Conway, Bldg 275-05-W-06  
St. Paul, MN 55144  
651 733 1110

The 3M™ Steri-Drape™ CHG Antimicrobial Incise Drape consists of a transparent, conformable, breathable fluid barrier film laminated to a lightly-tinted, acrylic pressure-sensitive adhesive that contains 2% w/w Chlorhexidine Gluconate (CHG). The adhesive drape is delivered on a silicone-coated release liner that is discarded after application.

*In vitro* testing demonstrates that 3M™ Steri-Drape™ CHG Antimicrobial Incise Drape has broad spectrum antimicrobial activity.

**Indications for Use:**

3M™ Steri-Drape™ CHG Antimicrobial Incise Drape is indicated for use as an incise drape with continuous antimicrobial activity based on *in vitro* time kill data out to 90 minutes. It is intended for external use only.

Catalog #	Color	Size	
8840EZ	Light pink	13 in x 13 in	34 cm x 35 cm
8848EZ	Light pink	22 in x 23 in	56 cm x 60 cm
8850EZ	Light pink	22 in x 17 in	56 cm x 45 cm
8851EZ	Light pink	22 in x 33 in	56 cm x 85 cm

**Substantial Equivalence:**

3M™ Steri-Drape™ CHG Antimicrobial Incise Drape is substantially equivalent to the predicate device, 3M™ Ioban™ 2 Antimicrobial Incise Drape cleared under K801550.

3M™ Steri-Drape™ CHG Antimicrobial Incise Drape is composed of the same or similar components and has the same or similar performance, intended use and indications for use as the predicate device.

Based on the similarities to the predicate device, the minor differences do not present any new safety or effectiveness issues and the device is substantially equivalent to the predicate.

**Comparative Data for Determining Substantial Equivalence of New Device to Predicate Device:**

Characteristic	3M™ Steri-Drape™ CHG Antimicrobial Incise Drape	3M™ Ioban™ 2 Antimicrobial Incise Drape	Comparison
Intended Use	3M™ Steri-Drape™ CHG Antimicrobial	3M™ Ioban™ 2 Antimicrobial Incise	Same

Characteristic	3M™ Steri-Drape™ CHG Antimicrobial Incise Drape	3M™ Ioban™ 2 Antimicrobial Incise Drape	Comparison
	Incise Drape is indicated for use as an incise drape with continuous antimicrobial activity. It is intended for external use only.	Drape is indicated for use as an incise drape with continuous antimicrobial activity. It is intended for external use only.	
Design/Materials	Antimicrobial impregnated adhesive coated on breathable film with silicone-coated release liner.	Antimicrobial impregnated adhesive coated on breathable film with silicone-coated release liner.	Same
Antimicrobial	Chlorhexidine Gluconate (CHG)	Iodine	Both are well characterized, broad spectrum, antimicrobials. Safety and efficacy testing demonstrates substantial equivalence.
Sold Sterile	Yes	Yes	Same
Sterilization Method	Ethylene Oxide sterilized, 10 <sup>-6</sup> SAL per ISO 11135	Gamma Sterilized, 10 <sup>-6</sup> SAL per ISO 11137	Different sterilization method but same sterility assurance level (SAL). Iodine is not compatible with Ethylene Oxide. CHG is not compatible with gamma irradiation. Both devices are sterilized to 10 <sup>-6</sup> SAL per respective ISO standards.
Packaging	Film/Tyvek Pouch	Foil Pouch	Foil laminate is used to protect Ioban from EO exposure if sterilized in a kit. Film/Tyvek is used to enable sterilization of CHG drape with EO.
Sterility Assurance Level	SAL 10 <sup>-6</sup>	SAL 10 <sup>-6</sup>	Same

Characteristic	3M™ Steri-Drape™ CHG Antimicrobial Incise Drape	3M™ Ioban™ 2 Antimicrobial Incise Drape	Comparison
Principals of Operation	The 3M™ Steri-Drape™ CHG Antimicrobial Incise Drape is a sterile adhesive film that is applied pre-operatively to skin. The film adheres and conforms to body contours, providing a sterile surface on top of the skin prior to surgery. The surgeon incises through the adhesive film and skin, ensuring a sterile film barrier all the way to the edge of the incision. The CHG contained in the adhesive provides additional antimicrobial activity to reduce the risk of microbial contamination of the surgical wound.	The 3M™ Ioban™ 2 Antimicrobial Incise Drape is a sterile adhesive film that is applied pre-operatively to skin. The film adheres and conforms to body contours, providing a sterile surface on top of the skin prior to surgery. The surgeon incises through the adhesive film and skin, ensuring a sterile film barrier all the way to the edge of the incision. The iodine contained in the adhesive provides additional antimicrobial activity to reduce the risk of microbial contamination of the surgical wound.	Same
Biocompatibility	Not cytotoxic, slight irritant, not a potential skin sensitizer	Not cytotoxic, slight irritant, not a potential skin sensitizer	Same

### Bench Testing Performance Data:

Bench testing was performed to demonstrate substantial equivalence for this submission, as recommended by the FDA Guidance on Premarket Notification [510(k)] Submissions for Surgical Gowns and Surgical Drapes and the additional recommendations from FDA Draft Guidance on Premarket Notification [510(k)]

3M Health Care

**3M Infection Prevention**

3M Center

2510 Conway. Bldg 275-05-W-06

St. Paul, MN 55144

651 733 1110

Submissions for Medical Devices that Include Antimicrobial Agents. Data submitted includes:

- Barrier performance per AAMI PB70:2012 using ASTM F1670 "Resistance to Synthetic Blood Penetration"
- Tensile and Elongation
- Flammability testing was performed per 16 CFR 1610
- Moisture Vapor Transmission Rate (MVTR)
- *in vitro* Antimicrobial Activity Spectrum and Minimum Effective Concentration (MEC) per ASTM E2315
- CHG Release Kinetics

Characteristic	Test Method Description	Device Component Tested	Data Summary		Comparison
			Test Article: 3M™ Steri-Drape™ CHG Anti-microbial Incise Drape	Predicate: 3M™ Ioban™ 2 Anti-microbial Incise Drape	
Tensile and Elongation	Based on PSTC-131, ASTM D882, ASTM D3759	Finished Drape Incise Area, sample cut to size	Tensile > 1365 g/25mm and Elongation >500%	Tensile > 1365 g/25mm and Elongation >500%	Same
MVTR	Based on ASTM E96 and Payne Cup Method	Finished Drape sample, cut to size	> 400 g/m <sup>2</sup> /24 hrs	> 400 g/m <sup>2</sup> /24 hrs	Same
Flammability	16 CFR 1610	Finished Drape Sample	Class I (Normal)	Class I (Normal)	Same
Barrier Performance	AAMI PB70:2012	Finished Drape Sample, cut to size	Level 4	Level 4	Same
Antimicrobial performance	ASTM E2315	Finished Drape Sample, cut to size	Broad Spectrum, > 4 log reduction at 90 min	Broad Spectrum, > 4 log reduction at 90 min	Same

**In vivo Animal Efficacy Data**

*In vivo* animal testing was performed to demonstrate antimicrobial efficacy compared with a placebo incise drape.

Characteristic	Description of Test Method	Device Component Tested		Data Summary	Conclusion
		Test Article	Placebo		
Antimicrobial efficacy study- in vivo porcine model	Wound contamination rates of test versus placebo drape in an in vivo porcine surgical model	Finished Drape sample, cut to size.	Placebo drape sample, cut to size.	T <sub>0</sub> Log <sub>10</sub> CFU/sample CHG drape lower than placebo (observed difference 0.35 log <sub>10</sub> , P=0.026)	Test article reduces microbial contamination of surgical wound compared to the placebo
				T <sub>4</sub> Log <sub>10</sub> CFU/sample CHG drape lower than placebo (observed difference 0.68 log <sub>10</sub> , P=0.001)	
				Difference in change in log <sub>10</sub> counts/sample T <sub>0</sub> to T <sub>4</sub> , observed difference 0.34 log <sub>10</sub> , P=0.14	
				Log <sub>10</sub> CFU/cm <sup>2</sup> under drape lower than placebo (observed difference 0.70 log <sub>10</sub> , P=0.0002)	

### Clinical Performance Data:

In house clinical panels were performed to demonstrate good adhesion to skin. Outside clinical adhesion studies were conducted to demonstrate safety and adhesion to various skin types.

Safety data were collected in these studies to demonstrate substantial equivalence for this submission to the predicate, 3M™ Ioban™ 2 Antimicrobial Incise Drape.

Test	Device Component Tested	Testing Summary	Summary	Comparison
Adhesion to Skin – in house clinical panels	Finished Drape Samples, cut to size	Adhesion comparison of 3M™ Steri-Drape™ CHG Antimicrobial Incise Drape and 3M™ Ioban™ 2 Antimicrobial Incise Drape on skin prepped with DuraPrep Solution or a leading PVP-I prep solution during a wet challenge wear study	No statistical difference in skin adhesion over either DuraPrep Solution or PVP-I prep solutions	Same performance as predicate
Adhesion to skin – in house clinical panels	Finished Drape Samples, cut to size	Adhesion comparison of 3M™ Steri-Drape™ CHG Antimicrobial Incise Drape and 3M Ioban 2 Antimicrobial Incise Drape on skin during a dry challenge wear study	Removal force of test article is higher than predicate, no difference in skin reactions	Similar performance to predicate
Adhesion to skin – external clinical study	Finished Drape Samples, cut to size	Evaluation of skin condition and removal force post- wear on subjects with very dry and potentially fragile skin types	No adverse events associated with test article post-wear; removal force slightly higher for test article	Similar performance to predicate

3M Health Care

**3M Infection Prevention**

3M Center

2510 Conway, Bldg 275-05-W-06

St. Paul, MN 55144

651 733 1110

**Performance, Safety and Efficacy Conclusion:**

Based on bench, animal and clinical studies, 3M™ Steri-Drape™ CHG Antimicrobial Incise Drape is as safe, as effective and is demonstrated to be substantially equivalent to the predicate device, 3M™ Ioban™ 2 Antimicrobial Incise Drape.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

July 25, 2014

3M Health Care  
Infection Prevention Division  
Ms. Kristin Totushek  
Regulatory Affairs Specialist  
3M Center, 2510 Conway Ave. Bldg 275-5W-06  
St. Paul, MN 55144

Re: K140250  
Trade/Device Name: 3M™ Steri-Drape™ CHG Antimicrobial Incise Drape  
Regulation Number: 21 CFR 878.4370  
Regulation Name: Surgical Drape Accessories  
Regulatory Class: Class II  
Product Code: KXX  
Dated: June 23, 2014  
Received: June 24, 2014

Dear Ms. Totushek,

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801 and Part 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and Part 809), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638 2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

*Tejashri Purohit-Sheth, M.D.*      Tejashri Purohit-Sheth, M.D.  
Clinical Deputy Director  
DAGRID/ODE/CDRH    FOR

Erin I. Keith, M.S.  
Director  
Division of Anesthesiology, General  
Hospital,  
Respiratory, Infection Control and  
Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**Indications for Use**

510(k) Number (if known)  
K140250

Device Name  
3M™ Steri-Drape™ CHG Antimicrobial Incise Drape

Indications for Use (Describe)

3M™ Steri-Drape™ CHG Antimicrobial Incise Drape is indicated for use as an incise drape with continuous antimicrobial activity based on in vitro time kill data out to 90 minutes. It is intended for external use only.

Catalog #	Color	Size
8840EZ	Light pink	13 in x 13 in 34 cm x 35 cm
8848EZ	Light pink	22 in x 23 in 56 cm x 60 cm
8850EZ	Light pink	22 in x 17 in 56 cm x 45 cm
8851EZ	Light pink	22 in x 33 in 56 cm x 85 cm

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.**

**FOR FDA USE ONLY**

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

**Sreekanth Gutala -S**

Digitally signed by Sreekanth Gutala -S  
DN: c=US, o=U.S. Government, ou=HHS, ou=FDA,  
ou=People, 0.9.2342.19200300.100.1.1=2000540490,  
cn=Sreekanth Gutala -S  
Date: 2014.07.24 11:04:33 -04'00'

---

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*\*An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number.\**