



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WO66-G609
Silver Spring, MD 20993-0002

May 11, 2015

Well Brain International, Ltd.
Victor K. Wai
Managing Director
Room 03, 14/F, Fook Yip Bldg.
53-57 Kwai Fung Crescent
Kwai Chung, N. T. Hong Kong, China

Re: K142055
Trade/Device Name: GYMFORM® ABS&CORE, Model: VDPGYCSET0042
Regulation Number: 21 CFR 890.5850
Regulation Name: Powered Muscle Stimulator
Regulatory Class: Class II
Product Code: NGX
Dated: April 3, 2015
Received: April 10, 2015

Dear Mr. Wai,

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the [Federal Register](#).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-

related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Felipe Aguel -S

for Carlos L. Peña, PhD, MS
Director
Division of Neurological
and Physical Medicine Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K142055

Device Name

GYMFORM® ABS&CORE, Model: VDPGYCSET0042

Indications for Use (Describe)

GYMFORM®ABS& CORE is intended to stimulate healthy muscles in order to improve or facilitate muscle performance. The ABS& CORE may be considered a technique or method for muscle training. 2-area belt is intended for use on the muscles in abdomen or lower back separately. Mini belt is intended for use on the muscles in arms, legs, thighs or buttocks areas separately.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

Chapter 5. 510(k) Summary

510(k) Summary

This summary of 510(K) safety and effectiveness information is being submitted in accordance with the requirement of 21 CFR 807.92.

1. Date of the summary prepared: May 8, 2015

2. Submitter's Information

- ◆ 510(k) Owner's Name: Well Brain International Ltd.
- ◆ Establishment Registration Number: 3004950644
- ◆ Address: Room 03, 14/F, Fook Yip Bldg., 53-57 Kwai Fung Crescent, Kwai Chung, N. T. Hong Kong, China
- ◆ Phone: (852) 2619-0833
- ◆ Fax: (852) 2429-0960
- ◆ Contact Person: Victor K Wai
- ◆ Email: victor@wellbrain-intl.com

3. Subject Device Information

- ◆ Trade Name: GYMFORM[®] ABS & CORE, Model: VDPGYCSET0042
- ◆ Common Name: Powered muscle stimulator
- ◆ Classification name: Stimulator, Muscle, Powered, For muscle conditioning
- ◆ Review Panel: Physical Medicine
- ◆ Product Code: NGX
- ◆ Regulation Class: II
- ◆ Regulation Number: 890.5850

4. Predicate Device Information

Predicate Device	Predicate Device I	Predicate Device II
Sponsor	Well Brain International Ltd.	Contour Technology

Sponsor: Well Brain International Ltd.

Subject Device: GYMFORM® ABS& CORE, Model: VDPGYCSET0042

File No.: 510(k) submission report (V1.0), Chapter 5

Predicate Device	Predicate Device I	Predicate Device II
Device Name	Gymform® ABS-A-ROUND, model: VDPGYCIND0016	Contour Technology Muscle Stimulator
510(k) Number	K130074	K111476
Product Code	NGX	NGX
Regulation Number	21 CFR 890.5850	21 CFR 890.5850
Regulation Class	II	II

5. Device Description

GYMFORM® ABS&CORE is a two channels battery operated muscle stimulation system specifically designed to stimulation the muscles.

The 2-area belt is intended for use on the muscles in abdomen or lower back separately.

Mini belt is intended for use on the muscles in arms, legs, thighs and buttocks areas separately.

It is comprised of a console for signal generation, two belts (2-area belt and mini belt) for fixation, and electrode pads for signal connection to skin. The electrode pads are replaceable.

Power is derived from 2 “AAA” batteries located in a compartment protected by a removable battery cover.

There is no current passed from side to side. The user cannot access the wiring or connectors within the belt.

The stimulator sends gentle electrical current to targeted muscle group through the electrodes placed on the skin. The parameters of the device are controlled by the buttons. Its intensity level can be adjustable by user.

6. Intended Use / Indications for Use

GYMFORM® ABS& CORE is intended to stimulate healthy muscles in order to improve or facilitate muscle performance. The ABS& CORE may be considered a technique or method for muscle training.

2-area belt is intended for use on the muscles in abdomen or lower back separately.

Mini belt is intended for use on the muscles in arms, legs, thighs or buttocks areas separately.

7. Test Summary

ABS & CORE has been evaluated the safety and performance by lab bench testing as following:

- ◆ Electrical safety test according to IEC 60601-1 and IEC 60601-2-10 standards

Sponsor: Well Brain International Ltd.

Subject Device: GYMFORM® ABS & CORE, Model: VDPGYCSET0042

File No.: 510(k) submission report (V1.0), Chapter 5

- ◆ Electromagnetic compatibility test according to IEC 60601-1-2 standard
- ◆ Software verification and validation test according to the requirements of the FDA “Guidance for Pre Market Submissions and for Software Contained in Medical Devices”
- ◆ Waveform test report to verify the output specifications of the device according to IEC 60601-2-10 and Guidance for Powered Muscle Stimulator.
- ◆ Dispersion test and Shelf test according to ASTM F 1980-07, Standard Guide for Accelerated Aging of Sterile Barrier Systems for Medical Devices and Guidance: Shelf Life of Medical Device.
- ◆ Adhesion test according to Section 5.4 of AAMI EC 12_2000_(R) 2010 Disposable ECG electrodes.

8. Comparison to predicate device and conclusion

The technological characteristics, features, specifications, materials, mode of operation, and intended use of GYMFORM® ABS & CORE is substantially equivalent to the predicate devices quoted above.

The differences between the subject device and predicate devices do not raise new issues of safety or effectiveness.

Elements of Comparison	Subject Device	Predicate Device I	Predicate Device II	Remark
Basic Unit Characteristics				
Device Name and Model	GYMFORM® ABS & CORE, Model: VDPGYCSET0042	Gymform® ABS-A-ROUND, model: VDPGYCIND0016	Contour Technology Muscle Stimulator	--
510 (K) Number	Applying	K130074	K111476	--
Intended Use	GYMFORM® ABS & CORE is intended to stimulate healthy muscles in order to improve or facilitate muscle performance. The ABS & CORE may be considered a technique or method for muscle training. 2-area belt is	Intended Use / Indications for Use: The GYMFORM® ABS-A-ROUND is intended to stimulate healthy muscles in order to improve or facilitate muscle performance. The ABS-A-ROUND	The Contour Technology Muscle Stimulator is intended to stimulate healthy muscles in order to improve or facilitate muscle performance. The Contour Technology Muscle Stimulator may therefore be considered a	SE

Sponsor: Well Brain International Ltd.

Subject Device: GYMFORM® ABS& CORE, Model: VDPGYCSET0042

File No.: 510(k) submission report (V1.0), Chapter 5

Elements of Comparison		Subject Device	Predicate Device I	Predicate Device II	Remark
		intended for use on the muscles in abdomen or lower back separately. Mini belt is intended for use on the muscles in arms, legs, thighs or buttocks areas separately.	may be considered a technique or method for muscle training. The 3-area belt is intended for use on the muscles in abdomen, left waist and right waist alternately. The Mini belt is intended for use on the muscles in arms, legs (lower extremities), thighs and buttocks areas separately.	technique or method for muscle training. The Contour Technology Muscle Stimulator Ab Belt accessory is intended for use on abdominal muscles only for strengthening and toning of abdominal muscles. The Contour Technology Muscle Stimulator BackPad accessory is intended for use on the lower back muscles only.	
Stimulated muscles		Abdomen, lower back, arms, legs, thighs and buttocks	Abdomen, lower back, arms, legs, thighs and buttocks	Abdomen, lower back	SE
Power Sources		2 x 1.5V AAA batteries	3 x 1.5V AAA batteries	4 x 1.5V AAA batteries	SE Note 1
Method of Line Current Isolation		Battery Supply N/A	Battery Supply N/A	Battery Supply N/A	SE
Patient Leakage Current	Normal Condition	10µA	5.8µA	--	SE Note 1
	Single Fault Condition	50µA	8.5µA	--	SE Note 1
Number of Modes		6	6	--	SE

Sponsor: Well Brain International Ltd.

Subject Device: GYMFORM[®] ABS & CORE, Model: VDPGYCSET0042

File No.: 510(k) submission report (V1.0), Chapter 5

Elements of Comparison		Subject Device	Predicate Device I	Predicate Device II	Remark
Number of Channels		2	3	--	SE Note 1
-Synchronous or Alternating		Alternating	Alternating	--	SE
-Method of Channel Isolation		Press MODE button for 3 seconds	Press FLR button	--	SE Note 1
Number Intensity Level		31 steps	99 steps	--	SE Note 1
Regulated Current or Regulated Voltage		Regulated Voltage	Regulated Voltage	--	SE
Software/Firmware/Microprocessor control		Yes	Yes	--	SE
Automatic Overload Trip		No	No	--	SE
Automatic No-load Trip		Yes.	Yes	--	SE
Automatic Shut Off		Yes.	Yes	--	SE
Patient Override Control		Yes	Yes	--	SE
Indicator Display	On/Off Status	Yes	Yes	--	SE
	Low Battery	Yes	Yes	--	SE
	Voltage/Current Level	Yes	No	--	SE
Timer Range		Default time is 10 minutes	Default time is 19 minutes	--	SE Note 1
Console weight		50g(Without batteries)	70g (Without batteries)	--	SE Note 1
Accessories weight		Big belt: 150g Small belt: 65g Electrode pad: 15g Carry bag: 60g	3-area belt: 310g Mini belt: 25g Electrode pad (big): 40g Electrode pad	--	SE Note 1

Sponsor: Well Brain International Ltd.

Subject Device: GYMFORM® ABS & CORE, Model: VDPGYCSET0042

File No.: 510(k) submission report (V1.0), Chapter 5

Elements of Comparison	Subject Device	Predicate Device I	Predicate Device II	Remark
		(small): 26g		
Console dimensions	100 mm(L) × 68 mm (W) × 24.5 mm (H)	91.8 mm(L) × 25.5 mm (W) × 82 mm (H)	--	SE Note 1
Electrode pad dimension	40 mm (L) × 70 mm (W) × 3mm (H)	Small: 33.0 cm ² Big: 34.5 cm ²	40.5 cm ²	SE Note 1
Housing Materials and Construction	Console: ABS plastic Belt: Polyester Electrode pads: Glycerine, Polyacrylic acid, Water and Salt	Console: ABS plastic Belt: Polyester Electrode pads: Glycerine, Polyacrylic acid, Water and Salt	--	SE
Output Specification				
Waveform	Symmetrical	Symmetrical	--	SE
Shape	Rectangular	Rectangular	--	SE
Maximum Output Voltage ($V_{\text{peak-to-peak}}$) (+/- 10%)	132V @ 500Ω 138V @ 2kΩ 140V @ 10kΩ	108V @ 500Ω 124V @ 2kΩ 126V @ 10kΩ	--	SE Note 2
Maximum Current Density ($I_{\text{peak-to-peak}}$)	264mA @ 500Ω 69mA @ 2kΩ 14mA @ 10kΩ	216mA @ 500Ω 62mA @ 2kΩ 12.6mA @ 10kΩ	--	SE Note 2
Frequency range	2 Hz, 10 Hz, 50 Hz, 90 Hz, 120 Hz	2 Hz, 10 Hz, 50 Hz, 90 Hz, 120 Hz	1 to 120 Hz	SE Note 2
Pulse width range	108μs / 124μs	100 μs / 120 μs	--	SE Note 2

Sponsor: Well Brain International Ltd.

Subject Device: GYMFORM[®] ABS& CORE, Model: VDPGYCSET0042

File No.: 510(k) submission report (V1.0), Chapter 5

Elements of Comparison	Subject Device	Predicate Device I	Predicate Device II	Remark
Pulse duration(Only changes with the mode)	Mode 1: 500ms; Mode 2: 11.1ms; Mode 3: 8.3ms; Mode 4: Front 90s: 11.1ms; Back 90s: 8.3ms; Mode 5: 100ms; Mode 6: 20ms	Mode 1: 500ms; Mode 2: 11.1ms; Mode 3: 8.33ms; Mode 4: Front 90s: 11.1ms; Back 90s: 8.33ms; Mode 5: 100ms; Mode 6: 20ms	340 μ s	SE Note 2
Phase duration(Only changes with the mode)	Mode 1: Not applicable(Continuous pulse) Mode 2: 2s Mode 3: 2s Mode 4: 2s Mode 5: 10s Mode 6: 16s	Mode 1: Not applicable(Continuous pulse) Mode 2: 2s Mode 3: 2s Mode 4: 2s Mode 5: 10s Mode 6: 16s	--	SE
Net Charge	19.2 μ C@ 500 Ω	15.7 μ C @ 500 Ω	--	SE Note 2
Maximum Phase Charge	16.4 μ C@ 500 Ω	13.0 μ C@ 500 Ω	--	SE Note 2
Maximum Current Density	0.082 mA/cm ² @ 500 Ω	0.057 mA/cm ² @ 500 Ω	0.55 mA/cm ² @ 500 Ω	SE Note 2
Maximum Power Density	94.8 μ W/cm ² @ 500 Ω	53.8 μ W/cm ² @ 500 Ω	--	SE Note 2
ON time	0.5s	1s	--	SE Note 2
OFF time	0.5s	1s	--	SE Note 2
Contraction and Relaxation time	Adjustable, due to different modes.	Adjustable, due to different modes.	--	SE
Burst Mode				
Pulse per burst	1~397	1~397	--	SE

Sponsor: Well Brain International Ltd.

Subject Device: GYMFORM® ABS& CORE, Model: VDPGYCSET0042

File No.: 510(k) submission report (V1.0), Chapter 5

Elements of Comparison	Subject Device	Predicate Device I	Predicate Device II	Remark
Bursts per second	0.125~1	0.125~1	--	SE
Burst Duration(s)	1~8	1~8	--	SE
Duty Cycle	0.02%~1.28%	0.02%~1.19%	--	SE Note 2
Additional Features				
Environment for operating	Temperature: 5 ~ 40° C Humidity: 20 ~65% RH	Temperature: 5 ~ 45° C Humidity: 20 ~65% RH	--	SE
Environment for storage	Temperature: 0 ~40° C Humidity: 10 ~90% RH	Temperature: 5 ~ 45° C Humidity: 20 ~65% RH	--	SE
Biocompatibility	All user directly contacting materials are compliance with ISO10993-5 and ISO10993-10 requirements.	All user directly contacting materials are compliance with ISO10993-5 and ISO10993-10 requirements.	All user directly contacting materials are compliance with ISO10993-5 and ISO10993-10 requirements.	SE
Electrical Safety	Comply with IEC 60601-1 and IEC 60601-2-10	Comply with IEC 60601-1 and IEC 60601-2-10	Comply with IEC 60601-1 and IEC 60601-2-10	SE
EMC	Comply with IEC 60601-1-2	Comply with IEC 60601-1-2	Comply with IEC 60601-1-2	SE

Comparison in Detail(s):

Note 1:

Although the power sources, patient leakage current in normal condition and single fault condition, number of channels, synchronous or alternating, method of channel isolation, number intensity level, timer range console weight, accessories weight, console dimensions, electrode pad dimension of subject device are a

Sponsor: Well Brain International Ltd.

Subject Device: GYMFORM[®] ABS& CORE, Model: VDPGYCSET0042

File No.: 510(k) submission report (V1.0), Chapter 5

little different from the predicate devices, they are all compliant with requirements of IEC 60601-1, IEC 60601-1-2 and Guidance for Powered Muscle Stimulator. So the differences of the function specifications will not raise any safety or effectiveness issue.

Note 2:

Although the maximum output voltage, maximum current density, frequency range, pulse width range, pulse duration, contraction time, net charge, maximum phase charge, maximum current density, maximum power density, on time, off time, contraction time, relaxation time and duty cycle of subject device are a little different from the predicate devices, they are all compliant with the requirements of IEC 60601-1, IEC 60601-2-10, and Guidance for Powered Muscle Stimulator. So the differences of function specification will not raise any safety or effectiveness issue.

Final Conclusion:

The subject device “GYMFORM[®] ABS& CORE” is Substantial Equivalence to the predicate device.