



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

November 25, 2014

Prismatik Dentalcraft, Incorporated  
Mr. Brandon Shepard  
Regulatory Affairs & Quality Assurance Specialist  
2212 Dupont Drive, Suite P  
Irvine, CA 92612

Re: K142118  
Trade/Device Name: Inclusive® Titanium Abutments compatible with: Straumann Bone Level, Nobel Biocare Branemark, and Nobel Biocare NobelActive Implants.  
Regulation Number: 21 CFR 872.3630  
Regulation Name: Endosseous Dental Implant Abutment  
Regulatory Class: II  
Product Code: NHA  
Dated: October 24, 2014  
Received: October 28, 2014

Dear Mr. Shepard,

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

A handwritten signature in black ink that reads "Susan Runno DDS, MA". The signature is written in a cursive style. A faint, semi-transparent "FDA" watermark is visible behind the signature.

Erin Keith  
Acting Director  
Division of Anesthesiology,  
General Hospital, Respiratory, Infection  
Control, and Dental Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

004

**Indications for Use Statement**DEPARTMENT OF HEALTH AND HUMAN SERVICES  
Food and Drug Administration

Form Approved: OMB No. 0910-0120

Expiration Date: January 31, 2017

See PRA Statement below.

**Indications for Use**510(k) Number (if known)  
K142118

## Device Name

Inclusive® Titanium Abutments compatible with: Straumann Bone Level, Nobel Biocare Branemark, and Nobel Biocare NobelActive Implants.

## Indications for Use (Describe)

Inclusive® Titanium Abutments are premanufactured prosthetic components directly connected to endosseous dental implants and are intended for use as an aid in prosthetic rehabilitation.

Inclusive® Titanium Abutments are compatible with:

- Straumann: Bone Level NC and RC implant sizes
- Nobel Biocare: Branemark RP size implant
- Nobel Biocare NobelActive NP and RP internal connection implants

## Type of Use (Select one or both, as applicable)

 Prescription Use (Part 21 CFR 801 Subpart D) Over-The-Counter Use (21 CFR 801 Subpart C)**PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.****FOR FDA USE ONLY**

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

005

## 510(k) Summary

[As Required by 21 CFR 807.92]

This 510(k) summary of safety and effectiveness information is being submitted in accordance with the requirements 21 CFR 807.92.

### A. SUBMITTER INFORMATION

**Company Name:** PRISMATIK DENTALCRAFT, INC.

**Company Address:** 2212 Dupont Dr., Suite P  
Irvine, CA 92612

**Company Phone / Fax:** (949) 225-1269 / (978) 313-0850

**Contact Person:** Primary Contact:  
Brandon Shepard, (949) 225-1243

Secondary Contact:  
Marilyn Pourazar, (949) 225-1269

**Date Summary Prepared:** November 24, 2014

### B. DEVICE IDENTIFICATION

**Trade/Proprietary Name:** Inclusive<sup>®</sup> Titanium Abutments  
compatible with: Straumann Bone Level, Nobel  
Biocare Branemark, and Nobel Biocare  
NobelActive Implants.

**Common Name:** Endosseous Dental Implant Abutment

**Regulation Number:** 872.3630

**Product Code:** NHA

**Device Class:** 2

**Review Panel:** Dental

---

**C. IDENTIFICATION OF PREDICATE DEVICE**

Trade/Proprietary Name:

- Inclusive Titanium Abutment Blanks (K083192)

**D. DEVICE DESCRIPTION**

Inclusive Titanium Abutments are endosseous implant abutments which are placed into the dental implant to provide support for a prosthetic restoration. The abutment is placed over the implant shoulder and is fastened into the implant with a screw. Abutments are made of titanium grade 23 (Ti-6AL-4V ELI) and meet ASTM F-136 Standard.

**E. INDICATIONS FOR USE**

Inclusive<sup>®</sup> Titanium Abutments are premanufactured prosthetic components directly connected to endosseous dental implants and are intended for use as an aid in prosthetic rehabilitation.

Inclusive<sup>®</sup> Titanium Abutments are compatible with:

- Straumann: Bone Level NC and RC implant sizes
- Nobel Biocare: Branemark RP size implant
- Nobel Biocare NobelActive NP and RP internal connection implants

**F. NON-CLINICAL TESTING**

Non-clinical test data was used to evaluate the device's safety and effectiveness, and determine substantial equivalence with predicate devices.

Non-clinical testing was performed in accordance with FDA Guidance "*Class II Special Controls Guidance Document: Root-form Endosseous Dental Implants and Endosseous Dental Implant Abutments*" and it consisted of Reliability Calculation, Fatigue Testing and Static Load Failure Testing of finished assembled implant/abutment systems.

In addition, sterilization validation information and a recommended sterilization method based on ANSI/AAMI ST79 and ISO 17665-1 is provided in the Information for Use.

The testing performed demonstrated implant to abutment compatibility and established that the predicate device is as safe, as effective, and performs as well or better than the predicate device(s).

**G. SUBSTANTIAL EQUIVALENCE**

Inclusive<sup>®</sup> Titanium Abutments compatible with Straumann Bone Level, Nobel Biocare Branemark, and Nobel Biocare NobelActive Implants are substantially equivalent to the Inclusive Titanium Abutment Blanks (K083192) identified in Section C above. They are substantially equivalent in intended use, materials, design and performance.

*(See Comparison Tables below).*

### Comparison of Predicate and Proposed Devices: Straumann Bone Level Platform

	Predicate Device	Proposed Device	
	Inclusive Titanium Abutment Blanks (K083192)	Inclusive Titanium Abutments compatible with: Straumann Bone Level Implants	Similarities and Differences
Manufacturer	Prismatik Dentalcraft, Inc.	Prismatik Dentalcraft, Inc.	N/A
Dimensions of Abutment	Cylindrical Diameter: 9.4mm Connection Length: 2.735mm	Cylindrical Diameter: 9.4mm Connection Length: 2.735mm	Same
Abutment Screw Manufacturer	OEM ( <i>Straumann</i> )	Inclusive ( <i>Prismatik</i> )	Equivalent screw from different manufacturer
Dimensions of Abutment Screw	Length 7.85mm; M1.6 Thread	Length 7.85mm; M1.6 Thread	Same
Indications for Use	The Inclusive Abutment is intended to be used in conjunction with endosseous implants in the maxillary and/or mandibular arch to provide support for crowns, bridges or overdenture prostheses. The prosthesis can be cement retained to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.	Inclusive Titanium Abutments are premanufactured prosthetic components directly connected to endosseous dental implants and are intended for use as an aid in prosthetic rehabilitation.	Same Intended Use. Different wording in the Indications for Use
Platform Compatibility	Straumann Bone Level NC and RC	Straumann Bone Level NC and RC	Same
Connection	Slot	Slot	Same
Design/Construction	Machined	Machined	Same
Anatomical Site	Oral Cavity	Oral Cavity	Same
Abutment Angle	0°-20°	0°-30°	Similar; increased angulation range
Implant Seat	Taper	Taper	Same
Screw Seat	Taper	Taper	Same
Material	Titanium Alloy	Titanium Alloy	Same

### Comparison of Predicate and Proposed Devices: Nobel Biocare Branemark Platform

	Predicate Device	Proposed Device	
	Inclusive Titanium Abutment Blanks (K083192)	Inclusive Titanium Abutments compatible with: Nobel Biocare Branemark Implants	Similarities and Differences
Manufacturer	Prismatik Dentalcraft, Inc.	Prismatik Dentalcraft, Inc.	N/A
Dimensions of Abutment	Cylindrical Diameter: 9.4mm Across Flats Length: 2.71mm	Cylindrical Diameter: 9.4mm Across Flats Length: 2.71mm	Same
Abutment Screw Manufacturer	OEM ( <i>Nobel Biocare</i> )	Inclusive ( <i>Prismatik</i> )	Equivalent screw from different manufacturer
Dimensions of Abutment Screw	Length 7.3mm; M2.0 Thread	Length 7.3mm; M2.0 Thread	Same
Indications for Use	The Inclusive Abutment is intended to be used in conjunction with endosseous implants in the maxillary and/or mandibular arch to provide support for crowns, bridges or overdenture prostheses. The prosthesis can be cement retained to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.	Inclusive Titanium Abutments are premanufactured prosthetic components directly connected to endosseous dental implants and are intended for use as an aid in prosthetic rehabilitation.	Same Intended Use. Different Indications for Use statement.
Platform Compatibility	Nobel Biocare Branemark RP	Nobel Biocare Branemark RP	Same
Connection	External Hex	External Hex	Same
Design/Construction	Machined	Machined	Same
Anatomical Site	Oral Cavity	Oral Cavity	Same
Abutment Angle	0°-20°	0°-30°	Similar; increased angulation range
Implant Seat	Flat	Flat	Same
Screw Seat	Flat	Flat	Same
Material	Titanium Alloy	Titanium Alloy	Same Titanium material

### Comparison of Predicate and Proposed Devices: Nobel Biocare NobelActive Platform

	<b>Predicate Device</b>	<b>Proposed Device</b>	
	<b>Inclusive Titanium Abutment Blanks (K083192)</b>	<b>Inclusive Titanium Abutments compatible with: Nobel Biocare NobelActive Implants</b>	<b>Similarities and Differences</b>
Manufacturer	Prismatik Dentalcraft, Inc.	Prismatik Dentalcraft, Inc.	N/A
Dimensions of Abutment	Cylindrical Diameter: 9.4mm Across Flats Length: 2.235mm	Cylindrical Diameter: 9.4mm Across Flats Length: 2.235mm	Same
Abutment Screw Manufacturer	OEM ( <i>Nobel Biocare</i> )	Inclusive ( <i>Prismatik</i> )	Equivalent screw from different manufacturer
Dimensions of Abutment Screw	Length 7.3mm; M1.6 Thread	Length 7.3mm; M1.6 Thread	Same
Indications for Use	The Inclusive Abutment is intended to be used in conjunction with endosseous implants in the maxillary and/or mandibular arch to provide support for crowns, bridges or overdenture prostheses. The prosthesis can be cement retained to the abutment. The abutment screw is intended to secure the abutment to the endosseous implant.	Inclusive Titanium Abutments are premanufactured prosthetic components directly connected to endosseous dental implants and are intended for use as an aid in prosthetic rehabilitation.	Same Intended Use. Different Indications for Use statement.
Platform Compatibility	Nobel Biocare NobelActive NP, RP	Nobel Biocare NobelActive NP, RP	Same
Connection	Internal Hex	Internal Hex	Same
Design/Construction	Machined	Machined	Same
Anatomical Site	Oral Cavity	Oral Cavity	Same
Abutment Angle	0°-20°	0°-30°	Similar; increased angulation range
Implant Seat	Taper	Taper	Same
Screw Seat	Taper	Taper	Same
Material	Titanium Alloy	Titanium Alloy	Same