



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

May 26, 2015

J.C.E.C. Company, Inc.  
c/o Ms. Elizabeth N. Dupras, RAC  
B&H Consulting Services, Inc.  
50 Division Street, Suite 206  
Somerville, NJ 08876

Re: K142549

Trade/Device Name: Oral7™ Moisturizing Gel  
Oral7™ Moisturizing Mouthwash  
Oral7™ Moisturizing Mouth Spray

Regulation Number: Unclassified  
Regulation Name: Saliva, Artificial  
Product Code: LFD  
Dated: April 24, 2015  
Received: April 27, 2015

Dear Ms. Dupras:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

 Tina  
Kiang -S

for Erin I. Keith, M.S.  
Director  
Division of Anesthesiology, General Hospital,  
Respiratory, Infection Control and  
Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**Indications for Use**

510(k) Number (if known)  
K142549

Device Name  
Oral7™ Moisturizing Gel, Oral7™ Moisturizing Mouthwash and Oral7™ Moisturizing Mouth Spray

Indications for Use (Describe)  
Relieve symptoms of dry mouth, refresh, moisturize, clean and soothe oral irritation, and lubricate oral dryness.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)       Over-The-Counter Use (21 CFR 801 Subpart C)

**PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.**

**FOR FDA USE ONLY**

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 5. 510(k) SUMMARY (K142549)

### 5.1. Submitter Information

- **Sponsor**

J.C.E.C. Company, Inc. (doing business as Fountain Health Care, LLC)  
12 Frieda Lane  
Kendall Park, NJ 08824

John Canvin  
President

Phone: (908) 420-3759  
email: canvinj@comcast.net

- **Regulatory Agent**

B&H Consulting Services, Inc.  
50 Division Street  
Suite 206  
Somerville, NJ 08876

Primary Contact:

Elizabeth N. Dupras, RAC  
Phone: 908-704-1691, ext. 288  
Fax: 908-704-1693  
email: edupras@bhconsultingservices.com

- **Date of Summary**

21 May 2015

### 5.2. Device Name

The proposed Oral7™ devices are outlined in Table 5-1.

**Table 5-1: Proposed Oral7™ Device Names**

Proprietary Name	Common/Usual Name	Classification Name (Product Code)
Oral7™ Moisturizing Gel	Saliva, Artificial	Saliva, Artificial (LFD)
Oral7™ Moisturizing Mouthwash	Saliva, Artificial	Saliva, Artificial (LFD)
Oral7™ Moisturizing Mouth Spray	Saliva, Artificial	Saliva, Artificial (LFD)

Regulatory Class: Unclassified  
Product Code: LFD

---

### 5.3. Equivalent Device Identification

The primary predicate device for Oral7™ Moisturizing Gel is Oral Balance Gel (K061331); the primary predicate device for Oral7™ Moisturizing Mouthwash and Oral7™ Moisturizing Mouth Spray is Biotène® Dry Mouth Oral Rinse (K101477).

The formulation of Oral7™ Moisturizing Mouthwash is identical to the formulation of Oral7™ Moisturizing Mouth Spray. The products have the same intended use, technical characteristics and physiological purpose; the only difference is in the mode of delivery. Therefore, Biotène® Moisturizing Mouth Spray (K103745) is a reference predicate for Oral7™ Moisturizing Mouth Spray.

The predicate devices are outlined in Table 5-2.

**Table 5-2: Predicate Devices**

Product	Notification Holder	510(k) Number
<i>Primary Predicates</i>		
Oral Balance Gel	Laclede, Inc.	K061331
Biotène® Dry Mouth Oral Rinse	GlaxoSmithKline Consumer Healthcare	K101477
<i>Reference Predicate</i>		
Biotène® Moisturizing Mouth Spray	GlaxoSmithKline Consumer Healthcare	K103745

### 5.4. Device Descriptions

Oral7™ products are specially formulated saliva substitutes which contain moisturizers, humectants, proteins and salivary enzymes that collectively have lubricating, moisturizing, soothing and refreshing properties to relieve and treat the symptoms of dry mouth. The Oral7™ products are listed in Table 5-3.

**Table 5-3: Oral7™ Product Presentations**

Product	Presentation
Oral7™ Moisturizing Gel	1.6 fluid ounces (48 mL) in plastic aluminum barrier laminate tube with membrane
Oral7™ Moisturizing Mouthwash	8.5 and 16.9 fluid ounces (250 and 500 mL, respectively) in white, polyethylene terephthalate bottles
Oral7™ Moisturizing Mouth Spray	1.5 fluid ounces (44 mL) in white, polyethylene terephthalate bottles

## 5.5. Intended Use

The intended use of the devices is symptomatic treatment of xerostomia.

*Indications for Use:* Relieve symptoms of dry mouth, refresh, moisturize, clean and soothe oral irritation, and lubricate oral dryness.

## 5.6. Comparison Table

A comparison of the Oral7™ products and the primary and reference predicate devices is provided in Table 5-4.

**Table 5-4: Comparison of Oral7™ Products to the Primary and Reference Predicate Devices**

Parameter	Primary Predicate Devices				Reference Predicate Device	
	Gel		Mouthwash and Mouth Spray/ Oral Rinse		Mouth Spray	
	Oral7™	Oral Balance (K061331)	Oral7™	Biotène® (K101477)	Oral7™	Biotène® (K103745)
Intended Use	Symptomatic treatment of xerostomia	Symptomatic treatment of xerostomia	Symptomatic treatment of xerostomia	Symptomatic treatment of xerostomia	Symptomatic treatment of xerostomia	Symptomatic treatment of xerostomia
Method of Use	Ready to use gel	Ready to use gel	Ready to use liquid	Ready to use liquid	Ready to use spray	Ready to use spray
Applications/Day	As needed	As needed	As needed	As needed	As needed	As needed
Disease State	Xerostomia	Xerostomia	Xerostomia	Xerostomia	Xerostomia	Xerostomia
Area of Use	Oral cavity	Oral cavity	Oral cavity	Oral cavity	Oral cavity	Oral cavity
Type of Product	Gel	Gel	Liquid	Liquid	Liquid	Liquid
Presentation	Non-sterile	Non-sterile	Non-sterile	Non-sterile	Non-sterile	Non-sterile

## 5.7. Summary of Testing

### 5.7.1. Nonclinical Testing

Biocompatibility testing for the final finished devices was conducted according to ISO 10993-1. Testing for cytotoxicity (ISO 10993-5), sensitization and irritation (ISO 10993—10) were completed. The cytotoxicity testing showed mild reactivity, which was not unexpected due to the known sensitivity/nature of the test (i.e., use of bare cells) and the presence of lysozyme, lactoperoxidase, mentha piperita and potassium thiocyanate in the formulations. The mild reactivity of the formulations supports the use in the proposed indication and does not impact safe use of the formulations. The sensitization and oral mucosal irritation testing demonstrated that the formulations are non-sensitizing and non-irritating.

### 5.7.2. Clinical Testing

The submission does not contain clinical testing.

### **5.7.3. Bench Testing**

Stability testing has been conducted on three (3) batches each of Oral7™ Moisturizing Gel and Oral7™ Moisturizing Mouthwash. Oral7™ Moisturizing Mouth Spray has the identical formulation as Oral7™ Moisturizing Mouthwash. Both formulations are filled into polyethylene terephthalate bottles; therefore, stability data for the Oral7™ Moisturizing Mouthwash support the shelf-life of Oral7™ Moisturizing Mouth Spray.

The data demonstrate that there are no significant changes in any of the product characteristics over the three year period of study for Oral7™ products. Available data support a three (3) year shelf-life for the Oral7™ products.

### **5.8. Conclusion**

The subject Oral7™ products have the same materials (i.e., moisturizers, humectants, proteins and salivary enzymes) as the primary predicate devices. All ingredients in the Oral7™ products are commonly used for their intended functions at equivalent or higher levels than in the proposed formulation.

In addition, the Oral7™ products have the same intended use and the same technological properties as the primary and reference predicate devices. The Oral7™ products were shown to be safe for the intended use in tests for cytotoxicity, sensitization and irritation.

The materials and manufacturing/processing of the Oral7™ products are the same as the primary predicate devices; therefore, there is no impact on safety and effectiveness of the subject devices.

Given the similarities in chemical composition, available delivery forms and physical properties between the Oral7™ products and the Biotène primary and reference predicate devices, Oral7™ Moisturizing Gel, Oral7™ Moisturizing Mouthwash and Oral7™ Moisturizing Mouth Spray are as safe and as effective as the primary and reference predicate devices.