



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center – WO66-G609
Silver Spring, MD 20993-0002

January 16, 2015

Ellipse Technologies, Incorporated
Ms. Rebecca Shelburne
Regulatory Affairs Project Manager
13900 Alton Parkway, Suite 123
Irvine, California 92618

Re: K142599
Trade/Device Name: PRECICE Trauma Nail System
Regulation Number: 21 CFR 888.3020
Regulation Name: Intramedullary fixation rod
Regulatory Class: Class II
Product Code: HSB
Dated: December 22, 2014
Received: December 23, 2014

Dear Ms. Shelburne:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to

<http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Industry and Consumer Education at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address

<http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Lori A. Wiggins -S

for

Mark N. Melkerson

Director

Division of Orthopedic Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES
Food and Drug Administration

Form Approved: OMB No. 0910-0120
Expiration Date: January 31, 2017
See PRA Statement on last page.

Indications for Use

510(k) Number (if known)
K142599

Device Name
PRECICE Trauma Nail System

Indications for Use (Describe)
The Ellipse PRECICE Trauma Nail System is indicated for open and closed fracture fixation, pseudoarthrosis, or mal-unions and non-unions of long bones.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D) Over-The-Counter Use (21 CFR 801 Subpart C)

PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.

FOR FDA USE ONLY

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

“An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number.”



13900 Alton Parkway, Suite 123
 Irvine, California 92618 USA
 Office: 949-837-3600
 Fax: 949-837-3664
 www.ellipse-tech.com

**Ellipse PRECICE Trauma Nail System
 510(k) Summary – K TBD
 September 2014**

1. Company: Ellipse Technologies, Incorporated
 13900 Alton Parkway, Suite 123
 Irvine, CA 92618

Contact: Rebecca Shelburne
 Regulatory Affairs Specialist
 Phone: (949) 837-3600 x227
 Fax: (949) 837-3664

Date of Submission: September 12, 2014

2. Proprietary Trade Name: Ellipse PRECICE Trauma Nail System

3. Classification Name: Intramedullary Fixation Rod (21 CFR 888.3020)

4. Product Code: HSB (Rod, Fixation, Intramedullary and Accessories)

5. Product Description: The Ellipse PRECICE Trauma Nail System is composed of the PRECICE Trauma Nail (supplied sterile), locking screws, surgical instruments and an external remote controller (ERC). The Nail is available in various diameters, lengths and screw hole configurations to accommodate a variety of patient anatomies. The locking screws are also available in a variety of diameters and lengths. The PRECICE Trauma Nail is supplied sterile by gamma radiation while the locking screws and accessories are supplied non-sterile and must be sterilized prior to use. The Nail contains an enclosed rare earth magnet, telescoping lead screw/nut assembly, and planetary gearing. The Nail is supplied pre-distracted by 10 mm to allow for compression fracture reduction techniques.

6. Indications: The Ellipse PRECICE Trauma Nail System is indicated for open and closed fracture fixation, pseudoarthrosis, or mal-unions and non-unions of long bones.

7. Substantial equivalence: A detailed comparison to the predicate device demonstrates that the Ellipse PRECICE Trauma Nail System is substantially equivalent to the following 510(k) cleared device:

Trade Name:	PRECICE® Trauma Nail System
Common Name:	Intramedullary Fixation Rod
510(k) Clearance Number:	K141447



13900 Alton Parkway, Suite 123
Irvine, California 92618 USA
Office: 949-837-3600
Fax: 949-837-3664
www.ellipse-tech.com

Substantial equivalence is based on same intended use, identical technological characteristics and principles of operation.

The Ellipse PRECICE Trauma Nail System and the predicate device have the same intended use. Specifically, the PRECICE Trauma Nail System and the predicate are both designed to treat long bones using compression fracture reduction techniques. The PRECICE Trauma Nail System subject of this premarket notification is indicated for open and closed fracture fixation, pseudoarthrosis, or mal-unions and non-unions of long bones, specifically.

The Ellipse PRECICE Trauma Nail has the same technological characteristics and principles of operation as that of the predicate. The design of both devices are identical. Both devices are inserted into the intramedullary canal of the long bone and secured with locking screws. Both devices are adjusted non-invasively by the Ellipse External Remote Controller (ERC).

There are no changes to the design of the PRECICE Trauma Nail being made as a result of this submission, therefore all testing that was performed on the predicate PRECICE Trauma Nail System is applicable. These tests include mechanical testing according to the methods outlined in the standard ASTM F1264-03, validation of the gamma radiation sterilization cycle in accordance with the VD_{max}^{25} methodology as given in ANSI/AAMI/ISO 11137-2 to verify that the gamma radiation sterilization process provides a sterility assurance level of 10^{-6} , and design verification and validation, shelf life testing for the packaging after accelerated aging, O-ring seal performance testing and biocompatibility in accordance with ISO 10993-1 for the intended use of the device.

There are no changes to the design of the ERC being made as a result of this submission, therefore all testing that was performed on the predicate PRECICE Trauma Nail System for the ERC are applicable.

Conclusions can be drawn from these tests that the PRECICE Trauma Nail System is substantially equivalent to the predicate device.